




*** * Φ3.5 空心钉安装器械包 91449901-91826501**

Φ3.5 Cannulated Screw Instrument Set



91440401-91042001 空心丝锥
Hollow Tap,Φ3.5mm



91121003-91600001 导针
Guide Pin,Φ1.0×150mm



91150402-91050001 限位器扳手
Limiter Wrench,SW2.5



91150810-91010001 钻头限位器
Limiter For Drill Bit, Φ2.5mm



91440301-91061001 钻套
Drill Sleeve,Φ2.5mm



91142501-91042001 空心六角起子
Hollow Hexagonal Screwdriver,SW2.5



91440201-91001001 Φ1.0导针钻套
Drill Sleeve For Guide Pin,Φ1.0mm



91132501-91081001 SW2.5六角起子
Hexagonal Screwdriver,SW2.5



91440501-91042001 空心埋头钻
Hollow Countersink Drill



91440601-91021001 打滑取出器
Slippage Extractor,SW2.5



91150502-91041001 导针测深器
Depth Gauge For Guide Pin



91150703-91610001 清洁导针
Cleaning Guide Pin,Φ1.0mm



91150906-91001001 持钉套
Screw Holding Sleeve



91150302-91090001 螺钉夹持钳
Screw Holding Forceps

91440001-91066001 空箱
Box