




*** * Φ3.5 无头加压空心钉安装器械包 91629902-91863401**

Φ3.5 Headless Compression Cannulated Screw Instrument Set




91621101-91006001 托盘
Tray



91620101-91042001 双头导向器
Double-End Guider,Φ1.0/Φ2.1mm



91121001-91600001 导针
Guide Pin,Φ1.0×120mm



91121002-91010001 螺纹导针
Threaded Guide Pin,Φ1.0×120mm



91620201-91081001 空心钻(D型快装式)
Hollow Drill Bit (Quick Coupling),Φ2.1mm



91620401-91081001 沉孔空心钻(D型快装式)
Hollow Countersink Drill (Quick Coupling)



91150401-91050001 限位器扳手
Limiter Wrench,SW2.0



91142301-91022001 空心梅花起子(D型快装式)
Hollow Torx Screwdriver(Quick Coupling),T2.3



91132303-91061001 梅花起子(D型快装式)
Torx Screwdriver(Quick Coupling),T2.3



91150110-91892006 快换手柄 (D型快装)
Quick Coupling Handle



91150805-91010001 限位器
Limiter,Φ2.1mm



91150501-91001001 导针测深器
Depth Gauge For Guide Pin



91150301-91090001 螺钉夹持镊子
Screw Holding Forceps



91150204-91012001 空心快换接头 (D型快装)
Hollow Quick Coupling Joint



91150701-91610001 清洁导针
Cleaning Guide Pin,Φ1.0mm

91620002-91003102 空箱
Box