




* * $\Phi 4.0$ 空心钉安装器械包 91439901-91026401

$\Phi 4.0$ Cannulated Screw Instrument Set



91121201-91600001 导针
Guide Pin, $\Phi 1.2 \times 150\text{mm}$

91103001-91041001 空心钻
Hollow Drill Bit, $\Phi 3.0 \times 150\text{mm}$

91430301-91061001 钻套
Drill Sleeve, $\Phi 3.0\text{mm}$

91121202-91610001 螺纹导针
Threaded Guide Pin, $\Phi 1.2 \times 150\text{mm}$

91430401-91042001 空心丝锥
Hollow Tap, $\Phi 4.0\text{mm}$

91430201-91001001 导针钻套
Drill Sleeve For Guide Pin, $\Phi 1.2\text{mm}$

91150402-91050001 限位器扳手
Limiter Wrench, SW2.5

91150906-91001001 持钉套
Screw Holding Sleeve

91142501-91042001 空心六角起子
Hollow Hexagonal Screwdriver, SW2.5

91150811-91010001 钻头限位器
Limiter For Drill Bit, $\Phi 3.0\text{mm}$

91150302-91090001 螺钉夹持镊
Screw Holding Forceps

91132501-91081001 SW2.5六角起子
Hexagonal Screwdriver, SW2.5

91440501-91042001 空心埋头钻
Hollow Countersink Drill

91440601-91021001 打滑取出器
Slippage Extractor, SW2.5

91150502-91041001 导针测深器
Depth Gauge For Guide Pin

91150704-91610001 清洁导针
Cleaning Guide Pin, $\Phi 1.2\text{mm}$

91430001-91066001 空箱
Box