




\* \*  $\phi$  6.5无头加压空心钉安装器械包 91649901-91004801

## $\phi$ 6.5 Headless Compression Cannulated Screw Instrument Set



91640501-91006001 托盘  
Instrument box PPSU




91122507-91600001 导针  
Guide Pin,  $\phi$  2.5 x 180



91150102-91042001 空心T型快换手柄(三角快装)  
Trial Template



91150202-91012001 空心快换接头(三角快装)  
Hollow Quick Coupling Joint



91122508-91010001 螺纹导针  
Threaded Guide Pin,  $\phi$  2.5 x 180



91150125-91042001 空心快装手柄(三角快装)  
Trial Template



91150301-91090001 螺钉夹持钳  
Screw Holding Forceps



91640101-91042001 双头导向器  
Double-End Guider,  $\phi$  4.5



91150403-91050001 限位器扳手  
Limiter Wrench, SW3.0



91150819-91010001 限位器  
Limiter, SW3.0



91150709-91610001 清洁导针  
Clean the needle,  $\phi$  2.5 x 200



91134003-91061001 六角起子(三角快装)  
Trial Template, SW4.0



91640301-91042001 空心埋头钻(三角快装)  
Hollow Countersink Drill



91150510-91001001 导针测深器  
Depth Gauge for Guide Pin,  $\phi$  2.5 x 180



91144002-91022001 空心六角起子(三角快装)  
Trial Template,  $\phi$  2.5/SW4.0



91650401-91081001 空心钻头  
Drill Bit,  $\phi$  4.5/  $\phi$  2.5

91640001-91003101 空箱  
Box

G