






* * $\Phi 7.3$ 空心钉安装器械包 91419901

$\Phi 7.3$ Cannulated Screw Instrument Set



91122504-91010001 导针
Guide Pin, $\Phi 2.5 \times 250\text{mm}$



91122505-91610001 螺纹导针
Threaded Guide Pin, $\Phi 2.5 \times 250\text{mm}$



91105002-91031001 空心钻
Hollow Drill Bit, $\Phi 5.0 \times 250\text{mm}$



91150812-91010001 钻头限位器
Limiter For Drill Bit, $\Phi 5.0\text{mm}$



91150506-91021001 导针测深器
Depth Gauge For Guide Pin, $\Phi 7.3\text{mm}$



91410201-91001001 导针钻套
Drill Sleeve For Guide Pin, $\Phi 2.5\text{mm}$



91410301-91001001 扩张器
Expander



91410401-91021001 钻套
Drill Sleeve, $\Phi 5.0\text{mm}$



91410501-91021001 扶持器
Support

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object. 20220629



* * $\Phi 7.3$ 空心钉安装器械包 91419901

$\Phi 7.3$ Cannulated Screw Instrument Set



91410601-91063001 多孔导向器
Multi Holes Guider



91410701-91042001 两孔导向器
Two Holes Guider



91410801-91060001 开口扳手
Opening Wrench,SW10



91150403-91050001 限位器扳手
Limiter Wrench,SW3.0



91410901-91042001 空心丝锥
Hollow Tap, $\Phi 7.3\text{mm}$



91144001-91042001 空心六角起子
Hollow Hexagonal Screwdriver,SW4.0



91154001-91001001 持钉套
Screw Holding Sleeve



91134001-91081001 六角起子
Hexagonal Screwdriver,SW4.0




91411001-91003001 空心套筒
Hollow Sleeve,SW6.0



91411101-91042001 空心埋头钻
Hollow Countersink Drill



91411201-91021001 打滑取出器
Slippage Extractor,SW4.0



91150706-91610001 清洁导针
Cleaning Guide Pin, $\Phi 2.5 \times 250\text{mm}$

91410001-91066001 空箱
Box