

# Femoral Intramedullary Nail System

Femoral Intramedullary Nail T8194/T8110/T8111



| Product Name                | Material       | Diameter | Length              | Product No.       |
|-----------------------------|----------------|----------|---------------------|-------------------|
| Femoral Intramedullary Nail | Titanium Alloy | Φ9.4     | 320/340/360/380     | T8194320-T8194380 |
|                             |                | Φ10      | 320/340/360/380/400 | T8110320-T8110400 |
|                             |                | Φ11      | 360/380/400         | T8111360-T8111400 |

Femoral Retrograde Intramedullary Nail T8290/T8210/T8211



| Product Name                           | Material       | Diameter | Length                      | Product No.       |
|----------------------------------------|----------------|----------|-----------------------------|-------------------|
| Femoral Retrograde Intramedullary Nail | Titanium Alloy | Φ9       | 260/280/300                 | T8290260-T8290300 |
|                                        |                | Φ10      | 240/260/280/300/320/340/360 | T8210240-T8210360 |
|                                        |                | Φ11      | 300/320/340/360             | T8211300-T8211360 |

Φ4.5 Full Thread Locking Screw T852



| Product Name              | Material       | Diameter | Length                     | Product No.       |
|---------------------------|----------------|----------|----------------------------|-------------------|
| Full Thread Locking Screw | Titanium Alloy | Φ4.5     | 30                         | T8524030          |
|                           |                |          | 32                         | T8524032          |
|                           |                |          | 35                         | T8524035          |
|                           |                |          | 38                         | T8524038          |
|                           |                |          | 40/45/50/55/60/65/70/75/80 | T8524040-T8524080 |

Φ6.5 Full Thread Locking Screw T853




| Product Name              | Material       | Diameter | Length                  | Product No.       |
|---------------------------|----------------|----------|-------------------------|-------------------|
| Full Thread Locking Screw | Titanium Alloy | Φ6.5     | 60/65/70/75/80/85/90/95 | T8536560-T8536595 |

Femoral Nail Instrument Set 92269901



Femoral Nail Aiming Tools

Femoral Retrograde Nail Instrument Set 92249901



Femoral Retrograde Nail Aiming Tools