




Craniomaxillofacial Fixation Series

B

Screw Features Self
Retaining, Easy to Operate.


Φ1.5mm

Self Tapping Blade


Anti-slip Cross Head



MCF Cross Screwdriver Blade
40mm 91130206



MCF Straight Plate T705A T707A




Product Name	Material	Thickness	Hole x Length	Product No.
MCF Straight Plate	Pure Titanium	0.6mm	2×15	T705A0200
			4×20	T705A0400
			6×29	T705A0600
			8×38	T705A0800
		0.8mm	2×16	T707A0200
			4×21	T707A0400
			6×30	T707A0600
			8×39	T707A0800

MCF Straight Plate T705A T707A




Product Name	Material	Thickness	Hole x Length	Product No.
MCF Straight Plate	Pure Titanium	0.6mm	12×53	T705A1200
		0.8mm	12×54	T707A1200

MCF Oblique L Plate (Left/Right) T705B T707B




Product Name	Material	Thickness	Hole x Length	Product No.
MCF Oblique L Plate (Left/Right)	Pure Titanium	0.6mm	4×17L	T705B0401
			4×17R	T705B0402
		0.8mm	4×17L	T707B0401
			4×17R	T707B0402

MCF Oblique L Plate (Left/Right) T705B T707B




Product Name	Material	Thickness	Hole x Length	Product No.
MCF Oblique L Plate (Left/Right)	Pure Titanium	0.6mm	5×21L	T705B0501
			5×21R	T705B0502
		0.8mm	5×22L	T707B0501
			5×22R	T707B0502

MCF Oblique L Plate (Left/Right) T705B T707B



Product Name	Material	Thickness	Hole x Length	Product No.
MCF Oblique L Plate (Left/Right)	Pure Titanium	0.6mm	7×26L	T705B0701
			7×26R	T705B0702
		0.8mm	7×27L	T707B0701
			7×27R	T707B0702

MCF T Plate T705C T707C




Product Name	Material	Thickness	Hole x Length	Product No.
MCF T Plate	Pure Titanium	0.6mm	7×25	T705C0700
		0.8mm	7×26	T707C0700

MCF Y Plate T705D T707D




Product Name	Material	Thickness	Hole x Length	Product No.
MCF Y Plate	Pure Titanium	0.6mm	5×20	T705D0500
		0.8mm	5×21	T707D0500

MCF Arc-Shaped Plate T705E T707E




Product Name	Material	Thickness	Hole x Length	Product No.
MCF Arc-Shaped Plate	Pure Titanium	0.6mm	6×30	T705E0600
			8×39	T705E0800
		0.8mm	6×30	T707E0600
			8×39	T707E0800

MCF Double Y Plate T705F T707F




Product Name	Material	Thickness	Hole x Length	Product No.
MCF Double Y Plate	Pure Titanium	0.6mm	6×20	T705F0600
		0.8mm	6×21	T707F0600

MCF X Plate T705G T707G




Product Name	Material	Thickness	Hole x Length	Product No.
MCF X Plate	Pure Titanium	0.6mm	4×10	T705G0400
		0.8mm	4×11	T707G0400

MCF Anatomic L Plate Orthognathic T707H




MCF Oblique L Plate Orthognathic T707J



MCF Double Bend Jaw Plate Orthognathic T710C




MCF Straight Sagittal Split Plate Orthognathic T710B



MCF Curved Sagittal Split Plate Orthognathic T710A



MCF Orbital Floor Anatomical Reconstruction Plate T705H



Product Name	Material	Thickness	Hole	Bridge Length	Product No.
MCF Anatomic L Plate Orthognathic	Pure Titanium	0.8mm	6L	4	T707H6041
			6R	4	T707H6042
			6L	8	T707H6081
			6R	8	T707H6082
			6L	12	T707H6121
MCF Oblique L Plate Orthognathic	Pure Titanium	0.8mm	6R	12	T707H6122
			4L	2	T707J4021
			4R	2	T707J4022
			4L	6	T707J4061
			4R	6	T707J4062
			4L	11	T707J4111
			4R	11	T707J4112
			6L	4	T707J6041
			6R	4	T707J6042
			6L	8	T707J6081
			6R	8	T707J6082
			6L	13	T707J6131
MCF Double Bend Jaw Plate Orthognathic	Pure Titanium	0.8mm	5	4	T710C0504
				6	T710C0506
				8	T710C0508
				10	T710C0510

Product Name	Material	Thickness	Hole	Bridge Length	Product No.
MCF Straight Sagittal Split Plate Orthognathic	Pure Titanium	1.0mm	4	6	T710B0406
				8	T710B0408
				10	T710B0410
				12	T710B0412
MCF Curved Sagittal Split Plate Orthognathic	Pure Titanium	1.0mm	6	6	T710A0660
				8	T710A0860
				10	T710A1060
				12	T710A1260

Product Name	Material	Thickness	Size	Bridge Length	Product No.
MCF Orbital Floor Anatomical Reconstruction Plate	Pure Titanium	0.6mm	Large L	42x37	T705H001L
			Large R	42x37	T705H001R
			Small L	42x36	T705H002L
			Small R	42x36	T705H002R

MCF 0.6mm Maxillofacial Plate Series



Compatible Screws

MCF $\Phi 1.5$ SDST Screw



MCF $\Phi 1.5$ ST Screw



MCF $\Phi 1.7$ Correction Screw



MCF 0.8mm Maxillofacial Plate Series



Compatible Screws

MCF $\Phi 2.0$ SDST Screw



MCF $\Phi 2.0$ ST Screw



MCF 0.8 Standard Orthognathic Plate Series



MCF 1.0 Standard Orthognathic Plate Series



MCF Orbital Floor Anatomical Reconstruction Plate



Compatible Screws

MCF $\Phi 2.0$ SDST Screw



MCF $\Phi 2.0$ ST Screw



MCF $\Phi 1.5$ Self Drilling Self Tapping Screw 32325



Product Name	Material	Diameter	Length	Product No.
MCF $\Phi 1.5$ Self Drilling Self Tapping Screw	Titanium Alloy	$\Phi 1.5$	4	32325004
			5	32325005
			6	32325006
			7	32325007

MCF $\Phi 1.5$ Self Tapping Screw 32321



Product Name	Material	Diameter	Length	Product No.
MCF $\Phi 1.5$ Self Tapping Screw	Titanium Alloy	$\Phi 1.5$	4	32321004
			5	32321005
			6	32321006
			7	32321007

MCF $\Phi 1.7$ Self Tapping Correction Screw 32322



Product Name	Material	Diameter	Length	Product No.
MCF $\Phi 1.7$ Self Tapping Correction Screw	Titanium Alloy	$\Phi 1.7$	5	32322005
			6	32322006
			7	32322007
			8	32322008

MCF $\Phi 2.0$ Self Drilling Self Tapping Screw 32326




Product Name	Material	Diameter	Length	Product No.
MCF $\Phi 2.0$ Self Drilling Self Tapping Screw	Titanium Alloy	$\Phi 2.0$	5	32326005
			6	32326006
			7	32326007
			8	32326008

MCF $\Phi 2.0$ Self Tapping Screw 32323




Product Name	Material	Diameter	Length	Product No.
MCF $\Phi 2.0$ Self Tapping Screw	Titanium Alloy	$\Phi 2.0$	5	32323005
			6	32323006
			7	32323007
			8	32323008

MCF Plate Instrument Set 91279901




MCF Straight Plate (Sterile Kit)



Compatible Tool

MCF All-In-One Screwdriver Cross 138mm



MCF Straight Plate T705A T707A




Product Name	Material	Thickness	Hole x Length	Product No.
MCF Straight Plate	Pure Titanium	0.6	2×15	T705A0200

MCF Φ1.5 Self Drilling Self Tapping Screw 32325




Product Name	Material	Diameter	Length	Product No.
MCF Φ1.5 Self Drilling Self Tapping Screw	Titanium Alloy	Φ1.5	4	32325004
			5	32325005
			6	32325006
			7	32325007

Digital 3D Titanium Mesh Customization




CT.MRI.DSA Data of the Patient



Reconstruction Image




Titanium Mesh Customized Design



Digital 3D Titanium Mesh Shaping Completed

B

Titanium Mesh T708




Product Name	Material	Thickness	Specification	Product No.
Titanium Mesh	Pure Titanium	0.8	90×90	T708W0090
			80×100	T708W0100
			100×119	T708W0119
			112×132	T708W0132
			145×150	T708W0150
			142×158	T708W0158

MCF Φ2.0 Self Drilling Self Tapping Screw 32326




Product Name	Material	Diameter	Length	Product No.
MCF Φ2.0 Self Drilling Self Tapping Screw	Titanium Alloy	Φ2.0	5	32326005
			6	32326006
			7	32326007
			8	32326008

B



Product Code Directive



T710N0411



T710N0822



T710N1235



T710N1646

B

MCF Φ2.0 Self Drilling Self Tapping Screw 32326



Product Name	Material	Diameter	Length	Product No.
MCF Φ2.0 Self Drilling Self Tapping Screw	Titanium Alloy	Φ2.0	5	32326005
			6	32326006
			7	32326007
			8	32326008

Compatible Tools

MCF All-In-One Screwdriver Cross 138mm 91130221



Product Features: Self Tapping Self Drilling, Self Retaining

Metal Bone Fracture Screw Cross T5022



Product Name	Material	Diameter	Length	Product No.
Metal Bone Fracture Screw Cross	Titanium Alloy	Φ2.0	5	T5022105
			6	T5022106
			7	T5022107
			8	T5022108

Compatible Tools

Limited Drill Bit Φ1.5×8mm 91111501



Crosshead Screwdriver 91130201



Product Features: Non Self Tapping Self Drilling, Not Self Retaining