




Common Plate Screw Series

Phalangeal/Metacarpal Plate Series

Straight Micro Plate (Mid) T7510



T Micro Plate (Mid) T7520



L Micro Plate (Mid)(Left/Right) T7530



Straight Micro Plate (Large) T7610



T Micro Plate (Large) T7620



L Micro Plate (Large)(Left/Right) T7630



Small T Reconstruction Plate T2010



Product Name	Material	Holes	Length	Product No.
Straight Micro Plate (Mid)	Pure Titanium	4	23	T7510004
		5	29	T7510005
		6	35	T7510006
		7	41	T7510007
T Micro Plate (Mid)		4	18	T7520004
		5	24	T7520005
		6	30	T7520006
		7	36	T7520007
L Micro Plate (Mid)(Left/Right)		4	18	T7530004L/R
		5	24	T7530005L/R
		6	30	T7530006L/R
Note: Use HA2.7 Screw				

Product Name	Material	Holes	Length	Product No.
Straight Micro Plate (Large)	Pure Titanium	4	31	T7610004
		5	39	T7610005
		6	47	T7610006
		8	63	T7610008
T Micro Plate (Large)		5	32	T7620005
		6	40	T7620006
		7	48	T7620007
		8	56	T7620008
L Micro Plate (Large)(Left/Right)		5	31	T7630005L/R
		6	39	T7630006L/R
		7	47	T7630007L/R
Note: Use HA2.5 Screw				

Product Name	Material	Holes	Length	Product No.
Small T Reconstruction Plate	Pure Titanium	6	43	T2010006
		7	54	T2010007
		8	65	T2010008
Note: Use HB4.0 Screw on Head, HA3.5 Screw on the rest (Plate shown in figure: 7 Holes)				

Curved Pelvic Plate T3820



Product Name	Material	Holes	Length	Product No.
Curved Pelvic Plate	Pure Titanium	5	60	T3820005
		6	72	T3820006
		8	96	T3820008
		10	120	T3820010
		12	144	T3820012
Note: Use HB4.0 Screw (Plate shown in figure: 10 Holes)				

Posterior Wall Acetabular Plate (Left/Right) T3821



Product Name	Material	Holes	Product No.
Posterior Wall Acetabular Plate (Left/Right)	Pure Titanium	9	T3821009L/R
Note: Use HB4.0 Screw			
Product with Patent			

Pubis Lineae Actuate Plate (Left/Right) T3840




Product Name	Material	Holes	Product No.
Pubis Lineae Actuate Plate (Left/Right)	Pure Titanium	12	T3840012L/R
Note: Use HB4.0 Screw			
Product with Patent			

Symphysis Iliopectineal Plate (Left/Right) T3850



Product Name	Material	Holes	Product No.
Symphysis Iliopectineal Plate (Left/Right)	Pure Titanium	13	T3850013L/R
Note: Use HB4.0 Screw			
Product with Patent			

Anterior Articulation Sacroiliac Plate T3880



Product Name	Material	Holes	Product No.
Anterior Articulation Sacroiliac Plate	Pure Titanium	4	T3880004
Note: Use HB4.0 Screw			

Symphysis Pubis Plate T3851



Product Name	Material	Holes	Product No.
Symphysis Pubis Plate	Pure Titanium	8	T3851008
Note: Use HB4.0 Screw			

Posterior Wall Hook Acetabular Plate T3821



Product Name	Material	Holes	Product No.
Posterior Wall Hook Acetabular Plate	Pure Titanium	3	T3821003
Note: Use HB4.0 Screw			

Product with Patent

PM Pelvic Reconstruction Plate T1330



Product Name	Material	Holes	Length	Product No.
PM Pelvic Reconstruction Plate	Pure Titanium	8	115	T1330008
		12	175	T1330012
Note: Use HB4.0 Screw (Plate shown in figure: 6 Holes)				

Restricted Compression Plate (Upper Limb)
T1110



Product Name	Material	Holes	Length	Product No.
Restricted Compression Plate (Upper Limb)	Pure Titanium	6	77	T1110006
		7	89	T1110007
		8	101	T1110008
		10	125	T1110010
		12	149	T1110012
Note: Use HA3.5 Screw (Plate shown in figure: 6 Holes)				

Restricted Compression Plate (Tibia)
T1120



Product Name	Material	Holes	Length	Product No.
Restricted Compression Plate (Tibia)	Pure Titanium	7	130	T1120007
		8	146	T1120008
		10	178	T1120010
		12	210	T1120012
Note: Use HA4.5 Screw (Plate shown in figure: 8 Holes)				

LC-DCP Ulna Radius Plate T1220



Product Name	Material	Holes	Length	Product No.
LC-DCP Ulna Radius Plate	Pure Titanium	5	62	T1220005
		6	74	T1220006
		7	86	T1220007
		8	98	T1220008
		10	122	T1220010
Note: Use HA3.5 Screw (Plate shown in figure: 6 Holes)				

One-third Tubular Plate T1410



Product Name	Material	Holes	Length	Product No.
One-third Tubular Plate	Pure Titanium	5	61	T1410005
		6	73	T1410006
		7	85	T1410007
		8	97	T1410008
		10	121	T1410010
Note: Use HA3.5 Screw (Plate shown in figure: 6 Holes)				

Reconstruction Plate T1310



Product Name	Material	Holes	Length	Product No.
Reconstruction Plate	Pure Titanium	5	65	T1310005
		6	78	T1310006
		8	104	T1310008
		10	130	T1310010
		12	156	T1310012
Note: Use HA3.5, HB4.0 Screw (Plate shown in figure: 6 Holes)				

Metallic Bone Fracture Screw Series

Ø1.5 Metal Bone Fracture Screw (Triangle)
T54115



Head Type	Material	Diameter	Length	Product No.
Triangle	Titanium Alloy	1.5	6/8/10/12/14	T5411506- T5411514

Ø2.0 Metal Bone Fracture Screw (Triangle)
T54120



Head Type	Material	Diameter	Length	Product No.
Triangle	Titanium Alloy	2.0	6/8/10/12/14/16/ 18/20/22	T5412006- T5412022

Ø2.5 Metal Bone Fracture Screw (Torx Type)
T54125



Head Type	Material	Diameter	Length	Product No.
Torx	Titanium Alloy	2.5	6/8/10/12/14/16/ 18/20/22/24/26	T5412506- T5412526

Ø2.7 Metal Bone Fracture Screw (Torx Type)
T54127



Head Type	Material	Diameter	Length	Product No.
Torx Type	Titanium Alloy	2.7	6/8/10/12/14/16/ 18/20/22/24/26/ 28/30	T5412706- T5412730

Metal Bone Fracture Screw (Torx Type)
(HA Full Thread) T54135



Head Type	Material	Diameter	Length	Product No.
Torx Type	Titanium Alloy	3.5	12/14/16/18/20/ 22/24/26/28/30/ 32/34/36/38/40/ 42/44	T5413512- T5413544

HA4.5 Metal Bone Fracture Screw (Torx Type)
(HA Full Thread) T54345



Head Type	Material	Diameter	Length	Product No.
Torx Type	Titanium Alloy	4.5	20/22/24/26/28/ 30/32/34/36/38/ 40/42/44/46/48/ 50	T5434520- T5434550
			55/60	T5434555- T5434560

HB4.0 Metal Bone Fracture Screw (Torx Type)
(HB Full Thread) T545



Head Type	Thread Type	Material	Diameter	Length	Product No.
Torx Type	HB Full Thread	Titanium Alloy	4.0	12/14/16/18/20/ 22/24/26/28/30/ 32/34/36/38/40/ 42/44/46/48/50/ 52/54/56/58/60	T5454012- T5454060

HB6.5 Metal Bone Fracture Screw (Torx Type)
(HB Full Thread) T549



Head Type	Thread Type	Material	Diameter	Length	Product No.
Torx Type	HB Full Thread	Titanium Alloy	6.5	30/35/40/45/ 50/55/60/65/ 70/75/80/85	T5496530- T5496585

HB6.5 Metal Bone Fracture Screw (Torx Type)
(HB Half Thread) T548



Head Type	Thread Type	Material	Diameter	Length	Product No.
Torx Type	HB Half Thread	Titanium Alloy	6.5	50/55/60/65/70/ 75/80/85/90	T5486550- T5486590