



C



QML Femoral Reverse Instrument Set 92279901



QML Multi-Dimensional Femoral Reverse Intramedullary Nail 63301/63302/63303



| Product Name                                              | Material       | Diameter | Length | Product No. |
|-----------------------------------------------------------|----------------|----------|--------|-------------|
| QML Multi-Dimensional Femoral Reverse Intramedullary Nail | Titanium Alloy | 9.4      | 180    | 63301180    |
|                                                           |                |          | 220    | 63301220    |
|                                                           |                |          | 260    | 63301260    |
|                                                           |                |          | 300    | 63301300    |
|                                                           |                |          | 320    | 63301320    |
|                                                           |                |          | 340    | 63301340    |
|                                                           |                |          | 360    | 63301360    |
|                                                           |                |          | 380    | 63301380    |
|                                                           |                | 10       | 400    | 63301400    |
|                                                           |                |          | 180    | 63302180    |
|                                                           |                |          | 220    | 63302220    |
|                                                           |                |          | 260    | 63302260    |
|                                                           |                |          | 300    | 63302300    |
|                                                           |                |          | 320    | 63302320    |
|                                                           |                |          | 340    | 63302340    |
|                                                           |                |          | 360    | 63302360    |
|                                                           |                |          | 380    | 63302380    |
|                                                           |                |          | 400    | 63302400    |
|                                                           |                | 11       | 180    | 63303180    |
|                                                           |                |          | 220    | 63303220    |
|                                                           |                |          | 260    | 63303260    |
|                                                           |                |          | 300    | 63303300    |
|                                                           |                |          | 320    | 63303320    |
|                                                           |                |          | 340    | 63303340    |
|                                                           |                |          | 360    | 63303360    |
|                                                           |                |          | 380    | 63303380    |
|                                                           |                |          | 400    | 63303400    |

QML  $\Phi$ 5.0 Full Thread Locking Screw 63207



| Product Name                             | Material       | Diameter | Length                       | Product No.           |
|------------------------------------------|----------------|----------|------------------------------|-----------------------|
| QML $\Phi$ 5.0 Full Thread Locking Screw | Titanium Alloy | 5.0      | 30-85<br>(in 5mm increments) | 63207030<br>-63207085 |

$\Phi$ 4.5 Retaining Double Thread Locking Screw 63009




| Product Name                                     | Material       | Diameter | Length                       | Product No.           |
|--------------------------------------------------|----------------|----------|------------------------------|-----------------------|
| $\Phi$ 4.5 Retaining Double Thread Locking Screw | Titanium Alloy | 4.5      | 30-85<br>(in 5mm increments) | 63009030<br>-63009085 |

QML Locking End Cap 63304



| Product Name        | Material       | Length | Product No. |
|---------------------|----------------|--------|-------------|
| QML Locking End Cap | Titanium Alloy | 0      | 63304100    |

QML End Cap 63305



| Product Name | Material       | Length | Product No. |
|--------------|----------------|--------|-------------|
| QML End Cap  | Titanium Alloy | 0      | 63305100    |
|              |                | 5      | 63305105    |
|              |                | 10     | 63305110    |
|              |                | 15     | 63305115    |