

Tibial Intramedullary Nail Aiming Tools



Tibial Intramedullary Nail T8384/T8390/T8310



Product Name	Material	Diameter	Length	Product No.
Tibial Intramedullary Nail	Titanium Alloy	Φ8.4	260/280/300/320	T8384260-T8384320
		Φ9.0	280/300/320/340	T8390280-T8390340
		Φ10	300/320/340	T8310300-T8310340

Φ4.5 Full Thread Locking Screw T852



Product Name	Material	Diameter	Length	Product No.
Full Thread Locking Screw	Titanium Alloy	Φ4.5	30	T8524030
			32	T8524032
			35	T8524035
			38	T8524038
			40/45/50/55/60/65/70/75/80	T8524040-T8524080

Tibial Nail Instrument Set 92329901

