



Digital 3D Titanium Mesh Customization




CT.MRI.DSA Data of the Patient



Reconstruction Image




Titanium Mesh Customized Design



Digital 3D Titanium Mesh Shaping Completed

B

Titanium Mesh T708



Product Name	Material	Thickness	Specification	Product No.
Titanium Mesh	Pure Titanium	0.8	90×90	T708W0090
			80×100	T708W0100
			100×119	T708W0119
			112×132	T708W0132
			145×150	T708W0150
			142×158	T708W0158

MCF Φ2.0 Self Drilling Self Tapping Screw 32326



Product Name	Material	Diameter	Length	Product No.
MCF Φ2.0 Self Drilling Self Tapping Screw	Titanium Alloy	Φ2.0	5	32326005
			6	32326006
			7	32326007
			8	32326008