



## WX Locking Plate Screw Series

G

## WX (Variable Angle) Locking Plate Series

VA Distal Medial Radius Locking Plate  
(Left/Right) 23467



| Product Name                                                                                         | Material      | Holes | Length | Product No. |
|------------------------------------------------------------------------------------------------------|---------------|-------|--------|-------------|
| VA Distal Medial Radius Locking Plate (Left/Right)                                                   | Pure Titanium | 6L    | 50     | 23467061    |
|                                                                                                      |               | 6R    | 50     | 23467062    |
| Note: Use Φ2.7 Variable Angle Locking Screw (Anodized), HA2.7 Screw (Plate shown in figure: 6 Holes) |               |       |        |             |

VA Distal Lateral Radius Locking Plate 23468



| Product Name                                                                                         | Material      | Holes | Length | Product No. |
|------------------------------------------------------------------------------------------------------|---------------|-------|--------|-------------|
| VA Distal Lateral Radius Locking Plate                                                               | Pure Titanium | 5     | 50     | 23468050    |
|                                                                                                      |               | 6     | 55     | 23468060    |
| Note: Use Φ2.7 Variable Angle Locking Screw (Anodized), HA2.7 Screw (Plate shown in figure: 6 Holes) |               |       |        |             |

VA 2.7 Radius Head Locking Plate (Left/Right)  
23475



| Product Name                                                                                                              | Material      | Holes | Length | Product No. |
|---------------------------------------------------------------------------------------------------------------------------|---------------|-------|--------|-------------|
| VA 2.7 Radius Head Locking Plate (Left/Right)                                                                             | Pure Titanium | 3L    | 40     | 23475031    |
|                                                                                                                           |               | 3R    | 40     | 23475032    |
|                                                                                                                           |               | 4L    | 48     | 23475041    |
|                                                                                                                           |               | 4R    | 48     | 23475042    |
|                                                                                                                           |               | 5L    | 56     | 23475051    |
|                                                                                                                           |               | 5R    | 56     | 23475052    |
| Note: Use Φ2.7 Variable Angle Locking Screw (Anodized), Φ2.7 Locking Screw (Anodized)<br>(Plate shown in figure: 5 Holes) |               |       |        |             |

VA S Clavicle Locking Plate (Low Stress)  
(Anodized) 23477



| Product Name                                                                                                                        | Material      | Holes    | Length | Product No. |
|-------------------------------------------------------------------------------------------------------------------------------------|---------------|----------|--------|-------------|
| VA S Clavicle Locking Plate<br>(Low Stress)(Anodized)                                                                               | Pure Titanium | 7L       | 90     | 23477071    |
|                                                                                                                                     |               | 7R       | 90     | 23477072    |
|                                                                                                                                     |               | 8L       | 102    | 23477081    |
|                                                                                                                                     |               | 8R       | 102    | 23477082    |
|                                                                                                                                     |               | 9L       | 114    | 23477091    |
|                                                                                                                                     |               | 9R       | 114    | 23477092    |
|                                                                                                                                     |               | 10L      | 126    | 23477101    |
|                                                                                                                                     |               | 10R      | 126    | 23477102    |
|                                                                                                                                     |               | 12L      | 150    | 23477121    |
| 12R                                                                                                                                 | 150           | 23477122 |        |             |
| Note: Use Φ3.5 Variable Angle Locking Screw (Anodized), Φ3.5 Locking Screw (Anodized), HA3.5 Screw (Plate shown in figure: 9 Holes) |               |          |        |             |

VA Clavicle Stem Locking Plate (Left/Right)  
23472/23464



| Product Name                                                                                                                               | Material      | Holes     | Length | Product No. |
|--------------------------------------------------------------------------------------------------------------------------------------------|---------------|-----------|--------|-------------|
| VA Clavicle Stem Locking Plate (Left/Right)                                                                                                | Pure Titanium | 7L(Short) | 90     | 23472071    |
|                                                                                                                                            |               | 7R(Short) | 90     | 23472072    |
|                                                                                                                                            |               | 7L        | 109    | 23464071    |
|                                                                                                                                            |               | 7R        | 109    | 23464072    |
|                                                                                                                                            |               | 8L        | 122    | 23464081    |
|                                                                                                                                            |               | 8R        | 122    | 23464082    |
|                                                                                                                                            |               | 9L        | 136    | 23464091    |
|                                                                                                                                            |               | 9R        | 136    | 23464092    |
|                                                                                                                                            |               | 11L       | 162    | 23464111    |
|                                                                                                                                            |               | 11R       | 162    | 23464112    |
| Note: Use Φ3.5 Variable Angle Locking Screw(Anodized), Φ3.5 Locking Screw(Anodized), HA3.5 Screw(Anodized) (Plate shown in figure: 9holes) |               |           |        |             |

VA Clavicle Hook Locking Plate (Left/Right)  
23470



| Product Name                                                                                                                                        | Material      | Height | Holes | Length | Product No |
|-----------------------------------------------------------------------------------------------------------------------------------------------------|---------------|--------|-------|--------|------------|
| VA Clavicle Hook Locking Plate (Left/Right)                                                                                                         | Pure Titanium | 14     | 8L    | 54     | 23470081   |
|                                                                                                                                                     |               |        | 8R    | 54     | 23470082   |
|                                                                                                                                                     |               |        | 9L    | 66     | 23470091   |
|                                                                                                                                                     |               |        | 9R    | 66     | 23470092   |
|                                                                                                                                                     |               |        | 11L   | 92     | 23470111   |
|                                                                                                                                                     |               |        | 11R   | 92     | 23470112   |
| Note: Use Φ2.7 Variable Angle Locking Screw (Anodized), Φ2.7/3.5 Locking Screw (Anodized), HA3.5 Screw (Anodized) (Plate shown in figure: 11 Holes) |               |        |       |        |            |

VA Clavicle Hook Locking Plate (Left/Right)  
23476



| Product Name                                                                                                                                                         | Material      | Height | Holes | Length | Product No |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|--------|-------|--------|------------|
| VA Clavicle Hook<br>Locking Plate<br>(Left/Right)                                                                                                                    | Pure Titanium | 17     | 8L    | 54     | 23476081   |
|                                                                                                                                                                      |               |        | 8R    | 54     | 23476082   |
|                                                                                                                                                                      |               |        | 9L    | 66     | 23476091   |
|                                                                                                                                                                      |               |        | 9R    | 66     | 23476092   |
|                                                                                                                                                                      |               |        | 11L   | 92     | 23476111   |
|                                                                                                                                                                      |               |        | 11R   | 92     | 23476112   |
| Note: Used For $\Phi 2.7$ Variable Angle Locking Screw (Anodized), $\Phi 2.7/3.5$ Locking Screw (Anodized), HA3.5 Screw (Anodized) (Plate shown in figure: 11 Holes) |               |        |       |        |            |

VA Distal Clavicle Locking Plate (Low Stress)  
(Left/Right) 23478



| Product Name                                                                                                                                       | Material      | Holes | Length | Product No. |
|----------------------------------------------------------------------------------------------------------------------------------------------------|---------------|-------|--------|-------------|
| VA Distal Clavicle Locking Plate<br>(Low Stress)(Left/Right)                                                                                       | Pure Titanium | 4L    | 82     | 23478041    |
|                                                                                                                                                    |               | 4R    | 82     | 23478042    |
|                                                                                                                                                    |               | 5L    | 96     | 23478051    |
|                                                                                                                                                    |               | 5R    | 96     | 23478052    |
|                                                                                                                                                    |               | 6L    | 110    | 23478061    |
|                                                                                                                                                    |               | 6R    | 110    | 23478062    |
|                                                                                                                                                    |               | 7L    | 123    | 23478071    |
|                                                                                                                                                    |               | 7R    | 123    | 23478072    |
|                                                                                                                                                    |               | 8L    | 136    | 23478081    |
|                                                                                                                                                    |               | 8R    | 136    | 23478082    |
| Note: Use Φ2.7 Variable Angle Locking Screw (Anodized), Φ2.7/3.5 Locking Screw (Anodized), HA3.5 Screw (Anodized) (Plate shown in figure: 6 Holes) |               |       |        |             |

VA Proximal Humerus Locking Plate 23465



| Product Name                                                                                                                                  | Material      | Holes | Length | Product No |
|-----------------------------------------------------------------------------------------------------------------------------------------------|---------------|-------|--------|------------|
| VA Proximal Humerus Locking Plate                                                                                                             | Pure Titanium | 3     | 91     | 23465030   |
|                                                                                                                                               |               | 5     | 115    | 23465050   |
| Note: Use Φ3.5 Variable Angle Locking Screw (Anodized), Φ3.5 Locking Screw (Anodized), HA3.5 Screw(Anodized) (Plate shown in figure: 3 Holes) |               |       |        |            |

VA Proximal Humerus Locking Plate (Long)  
23466



| Product Name                                                                                                                                   | Material      | Holes | Length | Product No. |
|------------------------------------------------------------------------------------------------------------------------------------------------|---------------|-------|--------|-------------|
| VA Proximal Humerus Locking Plate(Long)                                                                                                        | Pure Titanium | 5     | 144    | 23466050    |
|                                                                                                                                                |               | 7     | 180    | 23466070    |
|                                                                                                                                                |               | 9     | 216    | 23466090    |
| Note: Use Φ3.5 Variable Angle Locking Screw (Anodized), Φ3.5 Locking Screw (Anodized), HA3.5 Screw (Anodized) (Plate shown in figure: 7 Holes) |               |       |        |             |

VA Distal Humerus Subcondyle Locking Plate (Left/Right) 23471



| Product Name                                                                                                                                       | Material      | Holes | Length | Product No |
|----------------------------------------------------------------------------------------------------------------------------------------------------|---------------|-------|--------|------------|
| VA Distal Humerus Subcondyle Locking Plate (Left/Right)                                                                                            | Pure Titanium | 12L   | 142    | 23471121   |
|                                                                                                                                                    |               | 12R   | 142    | 23471122   |
|                                                                                                                                                    |               | 15L   | 188    | 23471151   |
|                                                                                                                                                    |               | 15R   | 188    | 23471152   |
| Note: Use Φ2.7 Variable Angle Locking Screw(Anodized), Φ2.7/3.5 Locking Screw (Anodized), HA3.5 Screw (Anodized) (Plate shown in figure: 12 Holes) |               |       |        |            |

VA Proximal Ulna Locking Plate (Left/Right) 23473



| Product Name                                                                                                                                 | Material      | Holes | Length | Product No. |
|----------------------------------------------------------------------------------------------------------------------------------------------|---------------|-------|--------|-------------|
| VA Proximal Ulna Locking Plate (Left/Right)                                                                                                  | Pure Titanium | 4L    | 112    | 23473041    |
|                                                                                                                                              |               | 4R    | 112    | 23473042    |
|                                                                                                                                              |               | 6L    | 140    | 23473061    |
|                                                                                                                                              |               | 6R    | 140    | 23473062    |
|                                                                                                                                              |               | 8L    | 168    | 23473081    |
|                                                                                                                                              |               | 8R    | 168    | 23473082    |
| Note: Use Φ3.5 Variable Angle Locking Screw(Anodized), Φ3.5 Locking Screw (Anodized), HA3.5 Screw(Anodized) (Plate shown in figure: 6 Holes) |               |       |        |             |

VA Proximal Ulna Locking Compression Hook Plate 23474



| Product Name                                                                                                                                               | Material      | Holes | Length | Product No. |
|------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|-------|--------|-------------|
| VA Proximal Ulna Locking Compression Hook Plate                                                                                                            | Pure Titanium | 5     | 64     | 23474050    |
| Note: Use $\Phi 3.5$ Variable Angle Locking Screw (Anodized), $\Phi 3.5$ Locking Screw (Anodized), HA3.5 Screw (Anodized) (Plate shown in figure: 5 Holes) |               |       |        |             |

VA Distal Fibula Reconstruction Locking Plate (Left/Right) T4415



| Product Name                                                                                                                       | Material      | Holes | Length | Product No. |
|------------------------------------------------------------------------------------------------------------------------------------|---------------|-------|--------|-------------|
| VA Distal Fibula Reconstruction Locking Plate (Left/Right)                                                                         | Pure Titanium | 6L    | 88     | T4415006L   |
|                                                                                                                                    |               | 6R    | 88     | T4415006R   |
|                                                                                                                                    |               | 8L    | 108    | T4415008L   |
|                                                                                                                                    |               | 8R    | 108    | T4415008R   |
| Note: Use Φ2.7 Variable Angle Locking Screw(Anodized), Φ2.7 Locking Screw (Anodized), HA2.7 Screw (Plate shown in figure: 6 Holes) |               |       |        |             |

VA Distal Lateral Fibula Locking Plate-II (Left/Right) 23563



| Product Name                                                                                                                                             | Material      | Holes    | Length | Product No |
|----------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|----------|--------|------------|
| VA Distal Lateral Fibula Locking Plate-II(Left,Right)                                                                                                    | Pure Titanium | 3L       | 86     | 23563031   |
|                                                                                                                                                          |               | 3R       | 86     | 23563032   |
|                                                                                                                                                          |               | 4L       | 99     | 23563041   |
|                                                                                                                                                          |               | 4R       | 99     | 23563042   |
|                                                                                                                                                          |               | 5L       | 112    | 23563051   |
|                                                                                                                                                          |               | 5R       | 112    | 23563052   |
|                                                                                                                                                          |               | 6L       | 125    | 23563061   |
|                                                                                                                                                          |               | 6R       | 125    | 23563062   |
|                                                                                                                                                          |               | 7L       | 138    | 23563071   |
| 7R                                                                                                                                                       | 138           | 23563072 |        |            |
| Note: Use Φ2.7 Variable Angle Locking Screw (Anodized), Φ2.7/3.5 Locking Screw (Anodized), HA2.7/HA3.5 Screw (Anodized) (Plate shown in figure: 5 Holes) |               |          |        |            |

VA Combo Minimally Invasive Calcaneus Locking Plate (Max)(Left/Right) 23559



| Product Name                                                                         | Material      | Holes | Product No. |
|--------------------------------------------------------------------------------------|---------------|-------|-------------|
| VA Combo Minimally Invasive Calcaneus Locking Plate (Max)(Left/Right)                | Pure Titanium | 8L    | 23559081    |
|                                                                                      |               | 8R    | 23559082    |
| Note: Use Φ3.5 Variable Angle Locking Screw(Anodized), Φ3.5 Locking Screw (Anodized) |               |       |             |

VA Combo Minimally Invasive Calcaneus Locking Plate (Mini)(Left/Right) 23560



| Product Name                                                                         | Material      | Holes | Product No. |
|--------------------------------------------------------------------------------------|---------------|-------|-------------|
| VA Combo Minimally Invasive Calcaneus Locking Plate(Mini)(Left/Right)                | Pure Titanium | 8L    | 23560081    |
|                                                                                      |               | 8R    | 23560082    |
| Note: Use Φ3.5 Variable Angle Locking Screw(Anodized), Φ3.5 Locking Screw (Anodized) |               |       |             |

VA Calcaneus Anatomical Locking Plate (Mini)(Left/Right) 23564



| Product Name                                                                         | Material      | Holes | Length | Product No |
|--------------------------------------------------------------------------------------|---------------|-------|--------|------------|
| VA Calcaneus Anatomical Locking Plate (Mini)(Left/Right)                             | Pure Titanium | 11L   | 58     | 23564111   |
|                                                                                      |               | 11R   | 58     | 23564112   |
| Note: Use Φ3.5 Variable Angle Locking Screw (Anodized), Φ3.5 Locking Screw(Anodized) |               |       |        |            |

VA Calcaneus Anatomical Locking Plate (Mid)(Left/Right) 23565



| Product Name                                                                        | Material      | Holes | Length | Product No |
|-------------------------------------------------------------------------------------|---------------|-------|--------|------------|
| VA Calcaneus Anatomical Locking Plate (Mid)(Left/Right)                             | Pure Titanium | 11L   | 64     | 23565111   |
|                                                                                     |               | 11R   | 64     | 23565112   |
| Note: Use Φ3.5 Variable Angle Locking Screw(Anodized), Φ3.5 Locking Screw(Anodized) |               |       |        |            |

VA Calcaneus Anatomical Locking Plate (Max)(Left/Right) 23566



| Product Name                                                                                    | Material      | Holes | Length | Product No |
|-------------------------------------------------------------------------------------------------|---------------|-------|--------|------------|
| VA Calcaneus Anatomical Locking Plate (Max)(Left/Right)                                         | Pure Titanium | 11L   | 71     | 23566111   |
|                                                                                                 |               | 11R   | 71     | 23566112   |
| Note: Use $\Phi$ 3.5 Variable Angle Locking Screw(Anodized), $\Phi$ 3.5 Locking Screw(Anodized) |               |       |        |            |

VA Pelvis Reconstruction Locking Compression Plate (Anodized) 23569



| Product Name                                                                                                                          | Material      | Holes | Length | Product No |
|---------------------------------------------------------------------------------------------------------------------------------------|---------------|-------|--------|------------|
| VA Pelvis Reconstruction Locking<br>Compression Plate<br>(Anodized)                                                                   | Pure Titanium | 5     | 62     | 23569050   |
|                                                                                                                                       |               | 6     | 74     | 23569060   |
|                                                                                                                                       |               | 8     | 98     | 23569080   |
|                                                                                                                                       |               | 10    | 122    | 23569100   |
|                                                                                                                                       |               | 12    | 146    | 23569120   |
|                                                                                                                                       |               | 14    | 170    | 23569140   |
| Note: Use $\Phi 3.5$ Variable Angle Locking Screw(Anodized), $\Phi 3.5$ Locking Screw<br>(Anodized) (Plate shown in figure: 10 Holes) |               |       |        |            |

VA 3.5 Proximal Medial Tibia Locking Plate  
(Left/Right) 23562



| Product Name                                                                                                                                   | Material      | Holes | Length | Product No |
|------------------------------------------------------------------------------------------------------------------------------------------------|---------------|-------|--------|------------|
| VA 3.5 Proximal Medial Tibial Locking Plate (Left/Right)                                                                                       | Pure Titanium | 5L    | 106    | 23562051   |
|                                                                                                                                                |               | 5R    | 106    | 23562052   |
|                                                                                                                                                |               | 7L    | 134    | 23562071   |
|                                                                                                                                                |               | 7R    | 134    | 23562072   |
|                                                                                                                                                |               | 9L    | 162    | 23562091   |
|                                                                                                                                                |               | 9R    | 162    | 23562092   |
|                                                                                                                                                |               | 11L   | 190    | 23562111   |
|                                                                                                                                                |               | 11R   | 190    | 23562112   |
| Note: Use Φ3.5 Variable Angle Locking Screw (Anodized), Φ3.5 Locking Screw (Anodized), HA3.5 Screw (Anodized) (Plate shown in figure: 9 Holes) |               |       |        |            |

VA KSL Proximal Lateral Tibia Locking Plate  
(Left/Right) 23558



| Product Name                                                                                                                                  | Material      | Holes | Length | Product No |
|-----------------------------------------------------------------------------------------------------------------------------------------------|---------------|-------|--------|------------|
| VA KSL Proximal Lateral Tibial Locking Plate (Left/Right)                                                                                     | Pure Titanium | 5L    | 98     | 23558051   |
|                                                                                                                                               |               | 5R    | 98     | 23558052   |
|                                                                                                                                               |               | 7L    | 130    | 23558071   |
|                                                                                                                                               |               | 7R    | 130    | 23558072   |
|                                                                                                                                               |               | 9L    | 162    | 23558091   |
|                                                                                                                                               |               | 9R    | 162    | 23558092   |
|                                                                                                                                               |               | 11L   | 194    | 23558111   |
|                                                                                                                                               |               | 11R   | 194    | 23558112   |
| Note: Use Φ3.5 Variable Angle Locking Screw (Anodized), Φ3.5 Locking Screw (Anodized), HA3.5 Screw(Anodized) (Plate shown in figure: 7 Holes) |               |       |        |            |

VA Distal Anterolateral Tibia Locking Plate  
(Left/Right) 23561



| Product Name                                                                                                                                                   | Material      | Holes | Length | Product No. |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|-------|--------|-------------|
| VA Distal Anterolateral Tibial Locking Plate (Left/Right)                                                                                                      | Pure Titanium | 5L    | 123    | 23561051    |
|                                                                                                                                                                |               | 5R    | 123    | 23561052    |
|                                                                                                                                                                |               | 7L    | 155    | 23561071    |
|                                                                                                                                                                |               | 7R    | 155    | 23561072    |
|                                                                                                                                                                |               | 9L    | 187    | 23561091    |
|                                                                                                                                                                |               | 9R    | 187    | 23561092    |
|                                                                                                                                                                |               | 11L   | 219    | 23561111    |
|                                                                                                                                                                |               | 11R   | 219    | 23561112    |
| Note: Use $\Phi$ 2.7 Variable Angle Locking Screw (Anodized), $\Phi$ 2.7/3.5 Locking Screw (Anodized), HA3.5 Screw (Anodized) (Plate shown in figure: 7 Holes) |               |       |        |             |

VA Distal Anteromedial Tibia Locking Plate  
(Left/Right) 23567



| Product Name                                                                                                                                  | Material      | Holes | Length | Product No |
|-----------------------------------------------------------------------------------------------------------------------------------------------|---------------|-------|--------|------------|
| VA Distal Anteromedial Tibial Locking Plate (Left/Right)                                                                                      | Pure Titanium | 4L    | 115    | 23567041   |
|                                                                                                                                               |               | 4R    | 115    | 23567042   |
|                                                                                                                                               |               | 6L    | 143    | 23567061   |
|                                                                                                                                               |               | 6R    | 143    | 23567062   |
|                                                                                                                                               |               | 8L    | 171    | 23567081   |
|                                                                                                                                               |               | 8R    | 171    | 23567082   |
|                                                                                                                                               |               | 10L   | 199    | 23567101   |
|                                                                                                                                               |               | 10R   | 199    | 23567102   |
|                                                                                                                                               |               | 12L   | 227    | 23567121   |
|                                                                                                                                               |               | 12R   | 227    | 23567122   |
| Note: Use Φ3.5 Variable Angle Locking Screw (Anodized), Φ3.5 Locking Screw (Anodized), HA3.5 Screw (Anodized)(Plate shown in figure: 8 Holes) |               |       |        |            |

VA Distal Medial Tibia Locking Plate  
(Left/Right) 23568



| Product Name                                                                                                                                 | Material      | Holes | Length | Product No. |
|----------------------------------------------------------------------------------------------------------------------------------------------|---------------|-------|--------|-------------|
| VA Distal Medial Tibial Locking Plate (Left/Right)                                                                                           | Pure Titanium | 4L    | 115    | 23568041    |
|                                                                                                                                              |               | 4R    | 115    | 23568042    |
|                                                                                                                                              |               | 6L    | 143    | 23568061    |
|                                                                                                                                              |               | 6R    | 143    | 23568062    |
|                                                                                                                                              |               | 8L    | 171    | 23568081    |
|                                                                                                                                              |               | 8R    | 171    | 23568082    |
|                                                                                                                                              |               | 10L   | 199    | 23568101    |
|                                                                                                                                              |               | 10R   | 199    | 23568102    |
|                                                                                                                                              |               | 12L   | 227    | 23568121    |
|                                                                                                                                              |               | 12R   | 227    | 23568122    |
| Note: Use Φ3.5 Variable Angle Locking Screw(Anodized), Φ3.5 Locking Screw (Anodized), HA3.5 Screw(Anodized) (Plate shown in figure: 8 Holes) |               |       |        |             |

## WX (Variable Angle) Locking Plate Series

VA Ulna Radius Locking Plate 23405



| Product Name                                                                                                                                  | Material      | Holes | Length | Product No. |
|-----------------------------------------------------------------------------------------------------------------------------------------------|---------------|-------|--------|-------------|
| VA Ulna Radius Locking Plate                                                                                                                  | Pure Titanium | 5     | 90     | 23405050    |
|                                                                                                                                               |               | 6     | 106    | 23405060    |
|                                                                                                                                               |               | 7     | 122    | 23405070    |
|                                                                                                                                               |               | 8     | 138    | 23405080    |
|                                                                                                                                               |               | 10    | 154    | 23405100    |
|                                                                                                                                               |               | 12    | 170    | 23405120    |
| Note: Use Φ3.5 Variable Angle Locking Screw(Anodized), Φ3.5 Locking Screw (Anodized), HA3.5 Screw (Anodized) (Plate shown in figure: 7 Holes) |               |       |        |             |

VA One-third Tubular Locking Plate 23406



| Product Name                                                                                                                                  | Material      | Holes | Length | Product No |
|-----------------------------------------------------------------------------------------------------------------------------------------------|---------------|-------|--------|------------|
| VA One-third Tubular Locking Plate                                                                                                            | Pure Titanium | 6     | 85     | 23406060   |
|                                                                                                                                               |               | 8     | 113    | 23406080   |
|                                                                                                                                               |               | 10    | 141    | 23406100   |
| Note: Use Φ3.5 Variable Angle Locking Screw(Anodized), Φ3.5 Locking Screw (Anodized), HA3.5 Screw (Anodized) (Plate shown in figure: 6 Holes) |               |       |        |            |

VA Reconstruction Locking Plate 23401



| Product Name                                                                                                                                 | Material      | Holes | Length | Product No. |
|----------------------------------------------------------------------------------------------------------------------------------------------|---------------|-------|--------|-------------|
| VA Reconstruction Locking Plate                                                                                                              | Pure Titanium | 5     | 87     | 23401050    |
|                                                                                                                                              |               | 6     | 103    | 23401060    |
|                                                                                                                                              |               | 7     | 119    | 23401070    |
|                                                                                                                                              |               | 8     | 135    | 23401080    |
|                                                                                                                                              |               | 9     | 151    | 23401090    |
|                                                                                                                                              |               | 10    | 167    | 23401100    |
|                                                                                                                                              |               | 12    | 199    | 23401120    |
|                                                                                                                                              |               | 14    | 231    | 23401140    |
| Note: Use Φ3.5 Variable Angle Locking Screw(Anodized), Φ3.5 Locking Screw (Anodized), HA3.5 Screw(Anodized) (Plate shown in figure: 7 Holes) |               |       |        |             |

VA Pelvis Reconstruction Locking Compression Plate (Straight) (Anodized) 23047



| Product Name                                                                                                                                 | Material      | Holes | Length | Product No. |
|----------------------------------------------------------------------------------------------------------------------------------------------|---------------|-------|--------|-------------|
| VA Pelvic Reconstruction Locking Compression Plate (Streight) (Anodized)                                                                     | Pure Titanium | 6     | 76     | 23047060    |
|                                                                                                                                              |               | 8     | 102    | 23047080    |
|                                                                                                                                              |               | 10    | 128    | 23047100    |
|                                                                                                                                              |               | 12    | 154    | 23047120    |
|                                                                                                                                              |               | 14    | 180    | 23047140    |
|                                                                                                                                              |               | 16    | 206    | 23047160    |
|                                                                                                                                              |               | 18    | 232    | 23047180    |
|                                                                                                                                              |               | 20    | 258    | 23047200    |
| Note: Use Φ3.5 Variable Angle Locking Screw(Anodized), Φ3.5 Locking Screw (Anodized), HA3.5 Screw(Anodized) (Plate shown in figure: 8 Holes) |               |       |        |             |

VA Humerus Locking Plate 23402



| Product Name                                                                                                                                 | Material      | Holes | Length | Product No. |
|----------------------------------------------------------------------------------------------------------------------------------------------|---------------|-------|--------|-------------|
| VA Humerus Locking Plate                                                                                                                     | Pure Titanium | 6     | 115    | 23402060    |
|                                                                                                                                              |               | 7     | 131    | 23402070    |
|                                                                                                                                              |               | 8     | 147    | 23402080    |
|                                                                                                                                              |               | 9     | 163    | 23402090    |
|                                                                                                                                              |               | 10    | 179    | 23402100    |
|                                                                                                                                              |               | 12    | 211    | 23402120    |
|                                                                                                                                              |               | 14    | 243    | 23402140    |
| Note: Use Φ3.5 Variable Angle Locking Screw(Anodized), Φ3.5 Locking Screw (Anodized), HA3.5 Screw(Anodized) (Plate shown in figure: 8 Holes) |               |       |        |             |