





* * 颈椎前路钢板器械包 93239901

Anterior Cervical Plate Instrument Set

A



93230101-93030001 定位针
Located Pin



93230201-93004001 定位钉植入器
Insert Device For Located Screw



93230301-93061001 快装开孔器 (D型)
Awl, $\Phi 2.0$



93230401-93021001 快装钻头 (D型)
Drill Bit, $\Phi 2.5 \times 13\text{mm}$



93230402-93021002 快装钻头 (D型)
Drill Bit, $\Phi 2.5 \times 14\text{mm}$



93230403-93021003 快装钻头 (D型)
Drill Bit, $\Phi 2.5 \times 15\text{mm}$



93230404-93021004 快装钻头 (D型)
Drill Bit, $\Phi 2.5 \times 16\text{mm}$



93230501-93041001 快装丝锥 (D型)
Tap, $\Phi 4.0$



93230601-93062001 单孔导向器
Single Hole Guider

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object. 20220311



* * 颈椎前路钢板器械包 93239901

Anterior Cervical Plate Instrument Set



93230701-93063001 两孔导向器
Two Holes Guider



93650201-93060001 撑开钉
Distraction Screw

**



93230901-93003001 撑开钉起子
Distraction Screw Screwdriver, SW4.0



93231001-93061101 撑开器
Distraction Device



93231101-93063001 持钉套
Screw Holding Sleeve

**



93231201-93061001 快装内六角起子 (D型)
Inner Hexagonal Screwdriver, SW2.5

**



93231301-93061001 快装一字起子 (D型)
Straight Screwdriver



93231401-93002001 持板器
Plate Holder



93231501-93048001 弯板钳
Plate Bending Clamp



93231601-93042001 快装转换接头 (D型)
Exchanging Joint



93231701-93042003 快装手柄 (D型)
Quick Coupling Handle



93231801-93027101 钉板盒
Screw&Plate Box

93210001-93007001 空箱
Box

A