




*** * 股骨带锁髓内钉安装器械包 92269901**

Femoral Intramedullary Nail Instrument Set



92260101-92002001 髓腔铰刀
Medullary Canal Reamer,Φ9mm



92260201-92002002 髓腔铰刀
Medullary Canal Reamer,Φ10.4mm



92320401-92002003 髓腔铰刀
Medullary Canal Reamer,Φ11mm



92260401-92002004 髓腔铰刀
Medullary Canal Reamer,Φ12mm



92260501-92002005 髓腔铰刀
Medullary Canal Reamer,Φ12.5mm



92320501-92002001 钻套
Drill Sleeve,Φ4.0mm



92320601-92002002 钻套
Drill Sleeve,Φ5.2mm



92320701-92002003 外套管
Outer Sleeve



92320801-92080001 软组织分离器
Soft Tissue Separator

Our company have final right to modify all the specification in this catalogue, If there is any inconformitybetween contant, model, color, please refer to the object.20220629



*** * 股骨带锁髓内钉安装器械包 92269901**

Femoral Intramedullary Nail Instrument Set



91114001-91040001 钻头
Drill Bit, $\Phi 4.0 \times 300\text{mm}$



91150818-91030001 限位器
Limiter, $\Phi 4.0/\text{SW}3.0$



91115201-91040001 钻头
Drill Bit, $\Phi 5.2 \times 300\text{mm}$



92321001-92001001 T型平头钻
T Flat Drill Bit, $\Phi 5.2 \times 220\text{mm}$



92322401-92001001 T型六角起子
T Hexagonal Wrench, SW3.5



91133504-91081001 封帽起子
Sealing Cap Screwdriver, SW3.5



91150908-91001001 持钉套
Screw Holding Sleeve



92261501-92002001 测深器
Depth Gauge



92261601-92080001 定位架锁轮
Locking Nut For Located Bracket



92261701-92080001 导杆锁轮
Locking Nut For Guide Rod



92261801-92001001 U型卡块
U Plug



91150404-91060001 锁轮扳手
Locking Nut Wrench, SW5.0



91150403-91060001 限位器扳手
Limiter Wrench, SW3.0



92261901-92004001 瞄准器手把
Diopter Handle



92262001-92002001 连接螺杆
Connecting Rod



* * 股骨带锁髓内钉安装器械包 92269901

Femoral Intramedullary Nail Instrument Set



92262101-92003001 远端定位架
Distal Located Bracket



92262201-92003001 导杆
Guide Rod



92262301-92021001 定位杆
Located Rod



92262401-9042001 打入器手柄
Impactor Handle



92262501-92002001 打拔器
Inserter- Extractor



92262601-92002001 万向接头
Universal Joint, M8




92101402-92040001 开口扳手
Opening Wrench, SW12



92100706 -92002001 开孔器
Awl



92322201-92001001 直角定位杆
Right-Angled Located Rod



92100202-92010001 导针
Guide Pin, Φ2.5×660mm



92100402-92003001 快装手柄
Quick Coupling Handle



* *
91124001-91010001 导针
Guide Pin, Φ4.0×300mm

92260001-91000101 空箱
Box