


\* \* JZWA脊柱微创固定系统安装器械包 93519901-93054311

## JZWA Minimally Invasive Spinal Fixation System Instrument Set

A



\*\*\*\*\*  
93510101-93061001 螺纹导针  
Threaded Guide Pin,  $\Phi 1.5 \times 450\text{mm}$



93510201-93067001 预紧起子  
Preloaded Screwdriver



93510301-93025001 六角起子  
Hexagonal Screwdriver, SW2.5



\*\*\*  
93510401-93093001 顶丝扳手\*2  
Nut Wrench, T5.3



93510501-93054001 开口器  
Awl



93510601-93093001 空心丝锥  
Hollow Tap,  $\Phi 5.0\text{mm}$



93510701-93093001 空心丝锥  
Hollow Tap,  $\Phi 6.0\text{mm}$




93510801-93093001 空心丝锥  
Hollow Tap,  $\Phi 7.0\text{mm}$



93510901-93065001 钉尾折断器  
Screw Cutter

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object.20220507



\* \* JZWA 脊柱微创内固定系统安装器械包 93519901-93054311

JZWA Minimally Invasive Spinal Fixation System Instrument Set

A



\*\*  
93511001-93825001 扩张套筒(大)  
Expanding Sleeve(Large)



\*\*  
93511101-93294001 扩张套筒(中)  
Expanding Sleeve(Medium)



\*\*  
93511201-93474001 扩张套筒(小)  
Expanding Sleeve(Small)



93511301-93009001 对抗扳手\*2  
Counter Wrench



93511401-93009001 折断对抗器  
Break-off Counter Device



\*\*  
93511501-93067001 置棒器  
Rod Placing Device



93511601-93060101 测量器  
Measurement Device



93511701-93009101 撑开钳  
Distraction Clamp



93511801-93009101 加压钳  
Compression Clamp



93511901-93003101 弯棒钳  
Rod Bending Clamp



\*\*  
93512001-93002101 棘轮手柄  
Ratchet Handle



\*\*  
93512101-93002101 扭力手柄  
Torque Handle, 12.5Nm




\*\*\*\*  
93512201-93081001 稳定套  
Stable Sleeve



\*\*  
93512302-93026001 U型钉上钉器  
Screwdriver For U Pedicle Screw, sw4.0

93510001-93001101 空箱  
Box

**JZWA Minimally Invasive Spinal Fixation System Instrument Set****A****选配件**93512501-93003001 空心丝锥定制款 (选配)  
Hollow Tap,  $\Phi 4.5\text{mm}$ 93512501-93003002 空心丝锥定制款 (选配)  
Hollow Tap,  $\Phi 5.0\text{mm}$ 93512501-93003003 空心丝锥定制款 (选配)  
Hollow Tap,  $\Phi 6.0\text{mm}$ 93512501-93003004 空心丝锥定制款 (选配)  
Hollow Tap,  $\Phi 7.0\text{mm}$ 93512901-93005001 U型螺纹上钉器 (选配)  
SW4.0\*85mm91352902 U型螺纹上钉器 (选配)  
SW4.0\*110mm93311205-93042001 顶丝预紧起子 T型 (选配)  
T5.393327105-93023001 万向钉调钉器 (选配)  
Universal Screw Adjuster93512401-93082001 穿刺针 (选配)  
Screwdriver For U Pedicle Screw93513201-93053001 JZWA 持钉外套管  
Nail-holding outer cannula, 30mm93513301-93053001 JZWA 持杆内套管  
Nail-holding inner cannula, 195mm93513401-93061001 顶丝起子 T5.3 加强型  
Reinforced Set Screw Driver, T5.3\*240mm93513101-93042001 加压撑开辅助器  
Compression Spreader Auxiliary



## D10 脊柱脊柱系统工具包

### JZWA Minimally Invasive Spinal Fixation System Instrument Set

A



D1050102- D10 短尾万向钉上钉器



D1050103-91005001 D10 短尾固定钉上钉器  
Screw Screwdriver



D1050104-91005001 D10 短尾万向钉上钉器 (梅花)  
Short tail multiaxial pedicle screw screwdriver



D1050101-91005001 D10 短尾固定钉上钉器 (梅花)  $\phi 4.5$  D1050201 D10 脊柱钉丝锥  
Screw Screwdriver



$\phi 5.0$  D1050202 D10 脊柱钉丝锥



$\phi 5.5$  D1050203 D10 脊柱钉丝锥



$\phi 6.0$  D1050204 D10 脊柱钉丝锥



$\phi 6.5$  D1050205 D10 脊柱钉丝锥