






*** * KTO 股骨带锁髓内钉安装器械包 92259901**

KTO Femoral Intramedullary Nail Instrument Set



92250101-92002001 锁钉钻套
Drill Sleeve For Locking Screw, Φ4.0mm



92250201-92002001 导针钻套
Drill Sleeve For Guide Pin, Φ2.5mm



92250301-92002001 锁钉套管
Locking Screw Sleeve



92250401-92080001 软组织分离器(长)
Soft Tissue Separator (Long)



92250501-92002001 定位杆钻套
Drill Sleeve For Located Rod, Φ5.2mm




92250601-92002001 定位杆套管
Located Rod Sleeve



92250701-92080001 软组织分离器(短)
Soft Tissue Separator, (Short)



91122502-91020001 螺纹导针
Threaded Guide Pin, Φ2.5×320mm



91122503-91600001 导针
Guide Pin, Φ2.5×320mm



* * KTO 股骨带锁髓内钉安装器械包 92259901

KTO Femoral Intramedullary Nail Instrument Set



91123201-91010001 螺纹导针
Threaded Guide Pin, Φ3.2×400mm



92250801-92061001 导针测深器
Depth Gauge Of Guide Pin, Φ2.5×320mm



92250901-92002001 测深器
Depth Gauge



91114001-91040001 远端骨钻
Distal Bone Drill, Φ4.0×300mm



91150801-91030001 限位器
Limiter, Φ4.0/SW3.5



**

91124001-91010001 导针
Guide Pin, Φ4.0×300mm



92251101-92044001 空心台阶钻
Hollow Step Drill Bit, Φ6.4mm



92251201-92061001 限位器
Limiter



91115201-91040001 定位杆骨钻
Bone Drill For Located Rod, Φ5.2×300mm



92251301-92001001 定位杆平头钻
Flat Drill Bit For Located Rod, Φ5.2mm



92223801-92001001 T型六角起子
T Hexagonal Screwdriver, SW3.5



92251501-92005001 T型股骨颈钉/尾帽扳手
T Femoral Neck Screw / End Cap Wrench, SW4.0/M2.5



92251601-92001001 连接螺杆扳手
Connecting Rod Wrench, SW6.5



92251701-92001001 六角起子(快装式)
Hexagonal Screwdriver (Quick Coupling), SW3.5



92251801-92001001 股骨颈钉/尾帽扳手(快装式)
Femoris Neck Screw / End Cap Wrench, SW4.0



* * KTO 股骨带锁髓内钉安装器械包 92259901

KTO Femoral Intramedullary Nail Instrument Set



92251901-92021001 T型连接螺杆万向扳手
T Universal Wrench For Connecting Rod,SW6.5



92252001-92027001 软杆尾帽万向扳手
Universal Wrench For Soft Rod End Cap,SW4.0



92252101-92030001 扳手
Wrench,Φ4.0mm




92252201-92048001 近端导杆
Proximal Guide Rod



92252301-92048001 远端导杆
Distal Guide Rod



92252401-92027001 股骨近端钉远端导杆
Distal Guide Rod For Proximal Femoral Screw



92252501-92046001 远端定位架
Distal Located Bracket



**

92252601-92080001 导杆锁轮
Locking Nut For Guide Rod



92252602-92080001 导杆连接锁轮
Locking Nut For Connecting Guide Rod



92252701-92001001 定位卡块
Located Plug



92252801-92069001 手把
Handle



92252901-92061001 主钉拔出连接杆
Connecting Rod For Extracting Nail,M8




**

92253001-92080001 连接螺杆
Connecting Rod



92253101-92021001 T型定位杆
T Located Rod



92253201-92001001 显影量尺
Developing Measure Gauge



*** * KTO 股骨带锁髓内钉安装器械包 92259901**

KTO Femoral Intramedullary Nail Instrument Set




92253301-92021001 直角定位杆
Right-Angled Located Rod



92253401-92065001 近端空心钻
Proximal Hollow Drill, $\Phi 13.8\text{mm}$



92100801-92025101 软扩刀杆
Soft Expansion Tool Bar



92100901-92005001 软扩刀盒
Soft Reamer Head Box



91109801-91021001 软扩刀头
Soft Reamer Head, $\Phi 8.5\text{mm}$



91109802-91021002 软扩刀头
Soft Reamer Head, $\Phi 9.0\text{mm}$



91109803-91021003 软扩刀头
Soft Reamer Head, $\Phi 9.5\text{mm}$



91109804-91021004 软扩刀头
Soft Reamer Head, $\Phi 10\text{mm}$



91109805-91021005 软扩刀头
Soft Reamer Head, $\Phi 10.5\text{mm}$



91109806-91021006 软扩刀头
Soft Reamer Head, $\Phi 11\text{mm}$



91109807-91021007 软扩刀头
Soft Reamer Head, $\Phi 11.5\text{mm}$



91109808-91021008 软扩刀头
Soft Reamer Head, $\Phi 12\text{mm}$



91109809-91021009 软扩刀头
Soft Reamer Head, $\Phi 12.5\text{mm}$



91109810-91021010 软扩刀头
Soft Reamer Head, $\Phi 13\text{mm}$



92100702-92023001 开孔器
Awl




* * KTO 股骨带锁髓内钉安装器械包 92259901

KTO Femoral Intramedullary Nail Instrument Set



92100301-92004001 复位棒(导针引入器)
Reduction Rod(Guide Pin Guider)



92100201-92061001 球头导针
Ball Head Guide Pin,Φ2.5×800



92100101-92066001 导针把持器
Guide Pin Holder



92253501-92001001 打拔器螺帽
Inserter-Extractor Nut



92253601-92062001 打拔锤
Inserter-Extractor Hammer



92253701-92063001 打拔器杆
Insertor-extractor Rod



92101001-92042001 保护套筒
Protective Sleeve



92101101-92064001 蜂窝式导针套筒
Honeycomb Guide Pin Sleeve,Φ3.2mm



92101201-92064001 套筒手柄
Sleeve Handle



92100401-92003001 快装手柄
Quick Coupling Handle



92100501-92062001 转换接头
Exchanging Joint

92250001-92063101 空箱
Box