





* * KTO 滑动加压伽玛钉安装器械包 92229901

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set



91123201-91020001 螺纹导针
Threaded Guide Pin, $\Phi 3.2 \times 400\text{mm}$




**

91124001-91010001 导针
Guide Pin, $\Phi 4.0 \times 300\text{mm}$



92230101-92061001 导针测深器
Depth Gauge For Guide Pin



92220201-92001001 显影量尺
Radiographic Ruler



92100701-92023001 开孔器
Awl



92220401-92042001 保护套筒
Protection Sleeve, $\Phi 16.5\text{mm}$



92220501-92081001 导针套管
Guide Pin Sleeve, $\Phi 3.2\text{mm}$



92220601-92065001 近端空心钻
Proximal Cannulated Drill Bit, $\Phi 16.5\text{mm}$




92100101-92066001 导针把持器
Guide Pin Gripper

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object.20220620



* * KTO 滑动加压伽玛钉安装器械包 92229901

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set



92100201-92061001 球头导针
Ball-tipped Guide Pin, $\Phi 2.5 \times 800\text{mm}$



92100301-92004001 复位棒(导针引入器)
Reduction Tool



92221001-92069001 手把
Handle



92221101-92048001 近端导杆
Proximal Aiming Arm



92221201 连接螺杆
Connecting Rod, SW7.0



92231001-92001001 连接螺杆扳手
Connecting Rod Wrench, SW7.0
曾用编码 92210901



92221401-92006001 远端导杆
Distal Aiming Arm



92231201-92080001 导杆锁轮
Locking Bolt For Aiming Arm



92221601-92065001 近端套管
Proximal Sleeve



92221701-92042001 近端皮质套管
Proximal Cortical Sleeve



92231501-92080001 近端套管针
Proximal Trocar



92221901-92061001 支持加压螺母
Buttress Compression Nut



92222001-92061001 子钉开口钻
Entwy Drill For Sub Screw, $\Phi 7.3\text{mm}$



92222101-92002001 子钉钻头
Drill Bit For Sub Screw, $\Phi 6.5\text{mm}$



92222201-92003001 子钉孔定位器
Locator For Sub Screw Hole, 65mm



* * KTO 滑动加压伽玛钉安装器械包 92229901

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set



92222301-92065001 母钉皮质骨空心钻
Cortical Bone Cannulated Bit For Parent Screw, Φ10.5mm



92222401-92046001 母钉限位空心钻
Cannulated Drill Bit For Parent Screw, Φ10.3mm



92222501-92001001 限位器
Limiter



92222601-92046001 母钉空心丝锥
Cannulated Tap For Parent Screw



92232401-92021001 母钉持杆
Holder Parent Screw



92232501-92042001 母钉中空持杆
Cannulated Holder Parent Screw



92222901-92046001 母钉上钉器
Parent Screw Screwdriver



92232701-92081001 母钉加压螺母
Compression Nut For Parent Screw



92223101-92002001 母钉导针套管
Guiding Pin Sleeve For Parent Screw, Φ3.2mm



92232901-92084001 子钉上钉器
Sub Screw Screwdriver, SW4.0/M2.5



91114001-91040001 骨钻
Bone Drill Bit, Φ4.0×300mm



91150801-91030001 限位器
Limiter, Φ4.0/SW3.5



**
92233001-92002001 远端锁钉套管
Distal Locking Screw Sleeve



**
92233101-92002001 钻套
Drill Bit Sleeve, Φ4.0mm



92233201-92080001 套管针
Trocar, Φ4.0mm



* * KTO 滑动加压伽玛钉安装器械包 92229901

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set



92233301-92002001 测深器
Depth Gauge



92233401-92001001 锁钉自持杆
Self Holding Rod For Locking Screw,SW3.5



92233501-92005001 T型六角扳手
Hex T Wrench , SW3.5



92233601-92001001 六角起子(快装式)
Hex Screwdriver(Quick Coupling),SW4.0



92233701-92041001 锁钉扳手
Locking Screw Wrench,SW4.0/M2.5



92233801-92001001 连接螺栓/尾帽扳手
Wrench For Connecting Screw/End Cap,SW7.0



92233901-92027001 软杆尾帽万向扳手
Hexible End Cap Wrench,SW7.0



92100401-92003001 快装手柄
Quick Coupling Handle



92100501-92062001 转换接头
Adapter



92224301-92062001 打入取出连接杆
Inserter For Main Nail, M12



92224401-92062001 主钉取出杆
Extractor Rod For Main Nail,M12



92100601-92027001 滑锤
Mallet



92101401-92090001 组合扳手
Combination Wrench



92235401-92020001 清洁导针
Cleaning Guide Pin,Φ3.0mm
曾用编码 92213401




92100801-92025101 软扩刀杆
Soft Reamer Shaft



* * KTO 滑动加压伽玛钉安装器械包 92229901

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set



92100901-92005001 软扩刀盒
Soft Reamer Head Box



91109801-91021001 软扩刀头
Soft Reamer Head,Φ8.5mm



91109802-91021002 软扩刀头
Soft Reamer Head,Φ9.0mm



91109803-91021003 软扩刀头
Soft Reamer Head,Φ9.5mm



91109804-91021004 软扩刀头
Soft Reamer Head,Φ10mm



91109805-91021005 软扩刀头
Soft Reamer Head,Φ10.5mm



91109806-91021006 软扩刀头
Soft Reamer Head,Φ11mm



91109807-91021007 软扩刀头
Soft Reamer Head,Φ11.5mm



91109808-91021008 软扩刀头
Soft Reamer Head,Φ12mm



91109809-91021009 软扩刀头
Soft Reamer Head,Φ12.5mm



91109810-91021010 软扩刀头
Soft Reamer Head,Φ13mm



92222202-92003001 子钉孔定位器
Locator For Sub Screw Hole, 95mm

92220001-92003101 空箱
Box

