



**\* \* KTO 滑动加压伽玛钉安装器械包小号加强型(加长版) 92239904**

**KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)**



H



\*\*\*\*\*

91123201-91020001 螺纹导针  
Threaded Guide Pin,  $\Phi 3.2 \times 400\text{mm}$




\*\*

91124001-91010001 导针  
Guide Pin,  $\Phi 4.0 \times 300\text{mm}$



92230101-92061001 导针测深器  
Depth Gauge For Guide Pin



92230201-92001001 显影量尺  
Radiographic Ruler



92100702-92023001 开孔器  
Awl



92230401-92042001 保护套筒  
Protection Sleeve,  $\Phi 15\text{mm}$

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object. 20220629



**\* \* KTO 滑动加压伽玛钉安装器械包小号加强型(加长版) 92239904**

**KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)**




92230501-92081001 导针套管  
Guide Pin Sleeve ,Φ15mm



92230601-92065001 近端空心钻  
Proximal Cannulated Drill Bit,Φ15mm



92100101-92066001 导针把持器  
Guide Pin Gripper



92100201-92061001 球头导针  
Ball-tipped Head Guide Pin,Φ2.5×800



92100301-92004001 复位棒(导针引入器)  
Reduction Tool



92230701-92069001 手把  
Handle



92230804-92048001 近端导杆  
Proximal Aiming Arm



92230901-92080001 连接螺杆  
Connecting Rod,SW7.0



92231001-92001001 连接螺杆扳手  
Connecting Rod WrenchSW7.0



92231104-92006001 远端导杆  
Distal Aiming Arm



92231201-92080001 导杆锁轮  
Locking Bolt For Aiming Arm



92231301-92065001 近端套管  
Proximal Sleeve



92231401-92042001 近端皮质套管  
Proximal Cortical Sleeve



92231501-92080001 近端套管针  
Proximal Trocar



92231601-92061001 支持加压螺母  
Buttress Compression Nut



\* \* KTO 滑动加压伽玛钉安装器械包小号加强型(加长版) 92239904

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)



92231701-92061001 子钉开口钻  
Entwy Drill For Sub Screw Φ7.1mm



92231801-92002001 子钉钻头  
Drill Bit For Sub Screw, Φ6.5mm



92231901-92003001 子钉孔定位器  
Locator For Sub Screw Hole, 65mm



92232001-92065001 母钉皮质骨空心钻  
Cortical Bone Cannulated Bit For Parent Screw, Φ9.6mm



92232101-92046001 母钉限位空心钻  
Cannulated Drill Bit For Parent Screw Φ9.0mm



92232201-92001001 限位器  
Limiter



92232301-92046001 母钉空心丝锥  
Cannulated Tap For Parent Screw, Φ9mm



92232401-92021001 母钉持杆  
Holding Rod For Parent Screw, M5



92232501-92042001 母钉中空持杆  
Cannulated Holder Parent Screw, M5



92232601-92046001 母钉上钉器  
Parent Screw Screwdriver



92232701-92081001 母钉加压螺母  
Compression Nut For Parent Screw



92232801-92002001 母钉导针套管  
Guide Pin Sleeve For Parent Screw, Φ3.2mm



92232901-92084001 子钉上钉器  
Sub Screw Screwdriver, SW4.0/M2.5



91114001-91040001 骨钻  
Bone Drill Bit, Φ4.0×300mm



91150801-91030001 限位器  
Limiter, Φ4.0/SW3.5



**\* \* KTO 滑动加压伽玛钉安装器械包小号加强型(加长版) 92239904**

**KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)**



**\*\***

92233001-92002001 远端锁钉套管  
Distal Locking Screw Sleeve,Φ8.2mm



**\*\***

92233101-92002001 钻套  
Drill Sleeve,Φ4.2mm



92233201-92080001 套管针  
Trocarnail,Φ4mm



92233301-92002001 测深器  
Depth Gauge



92233401-92001001 锁钉自持杆  
Self Holding Rod For Locking Screw,SW3.5



92233501-92001001 T型六角扳手  
Hex T Wrench ,SW3.5



92233601-92001001 六角起子(快装式)  
Hex Screwdriver(Quick Coupling),,SW4.0



92233701-92005001 锁钉扳手  
Locking Screw Wrench,SW4.0/M2.5



92233801-92001001 连接螺栓/尾帽扳手  
Wrench For Connecting Bolt/End Cap,SW7.0



92233901-92027001 软杆尾帽万向扳手  
Hexible End Cap Wrench,SW7.0



92100401-92003001 快装手柄  
Quick Coupling Handle



92100501-92062001 转换接头  
Adapter



92234001-92062001 打入取出连接杆  
Inserter For Main Nail,M10



92234101-92062001 主钉取出杆  
Extractor Rod For Main Nail,M10



92100601-92027001 滑锤  
Mallet




**\* \* KTO 滑动加压伽玛钉安装器械包小号加强型(加长版) 92239904**

**KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)**




92101401-92090001 组合扳手  
Combination Wrench, SW11



92235401-92020001 清洁导针  
Cleaning Guide Pin,  $\Phi 3.0\text{mm}$   
曾用编码92213401



**\*\***  
92100801-92025101 软扩刀杆  
Soft Reamer Shaft



92100902-92005001 软扩刀盒  
Soft Reamer Head Box



91109801-91021001 软扩刀头  
Soft Reamer Head,  $\Phi 8.5\text{mm}$



91109802-91021002 软扩刀头  
Soft Reamer Head,  $\Phi 9.0\text{mm}$



91109803-91021003 软扩刀头  
Soft Reamer Head,  $\Phi 9.5\text{mm}$



91109804 -91021004 软扩刀头  
Soft Reamer Head,  $\Phi 10\text{mm}$



91109805-91021005 软扩刀头  
Soft Reamer Head,  $\Phi 10.5\text{mm}$



91109806-91021006 软扩刀头  
Soft Reamer Head,  $\Phi 11\text{mm}$



91109807-91021007 软扩刀头  
Soft Reamer Head,  $\Phi 11.5\text{mm}$



91109808-91021008 软扩刀头  
Soft Reamer Head,  $\Phi 12\text{mm}$



91115201-91040001 定位杆钻头  
Located Rod Drill Bit,  $\Phi 5.2 \times 300\text{mm}$



92234201-92001001 定位杆平头钻  
Flat Drill Bit For Located Rod,  $\Phi 5.2\text{mm}$



92234301-92021001 定位杆  
Located Rod



**\* \* KTO 滑动加压伽玛钉安装器械包小号加强型(加长版) 92239904**

**KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)**



92234401-91080001 定位杆套管针  
Trocar For Located Rod



92234501-92002001 定位杆套管  
Located Rod Sleeve



92234601-92002001 定位杆钻套  
Drill Sleeve For Located Rod, Φ5.2mm



92234701-92001001 定位卡块  
Located Plug



\*\*\*  
92234801-92080001 远端连接锁轮  
Distal Connecting Nut



92234901-92025001 股骨近端瞄准架  
Proximal Aiming Jig



92235001-92040101 长钉远端连接导杆  
Distal Connecting Guide Rod For Long Screw



92235101-92048001 长钉远端导杆  
Distal Guide Rod For Long Screw



92235201-92046001 长钉远端定位架左  
Distal Located Bracket For Long Screw(Left)



92235301-92046001 长钉远端定位架右  
Distal Located Bracket For Long Screw(Right)



92231902-92003001 子钉孔定位器  
Locator For Sub Screw Hole, 95mm



92235601-92048001 近端瞄准架  
Proximal Position Arm

92230004-92063104 空箱  
Box

