





KTO 肱骨带锁髓内钉安装器械包 92519901

KTO Humeral Intramedullary Nail Instrument Set



92510201-92065001 导针把持器
Guide Pin Holder



92510301-92001001 显影量尺
Developing Measure Gauge



92100705-92003001 中空开孔器
Hollow Awl



92510401-92002001 保护套筒
Protective Sleeve, $\Phi 9.5$



92510501-92064001 近端空心钻
Proximal Hollow Drill, $\Phi 9.5$



92100403-92062001 T型快装手柄
T Quick Coupling Handle



92100502-92042001 转换接头
Exchanging Joint



91122506-91010001 螺纹导针
Threaded Guide Pin, $\Phi 2.5 \times 200$



92100206-92060001 球头导针
Ball Head Guide Pin, $\Phi 2.0 \times 660\text{mm}$

Our company have final right to modify all the specification in this catalogue, If there is any inconformitybetween contant, model, color, please refer to the object.20220629



KTO 肱骨带锁髓内钉安装器械包 92519901

KTO Humeral Intramedullary Nail Instrument Set



92510801-92001001 埋头锥
Countersink Tap



92100302-92063001 复位棒
Reduction Rod



92510101-92004001 钉盒
Screw Box



91114002-91080001 骨钻
Bone Drill Bit, $\Phi 4.0 \times 120\text{mm}$



92511001-92063001 软扩
Soft Expansion, $\Phi 7.0\text{mm}$



92511201-92063002 软扩
Soft Expansion, $\Phi 8.5\text{mm}$



92510701-92001001 埋头钻
Countersink Drill Bit



92511101-92063003 软扩
Soft Expansion, $\Phi 8.0\text{mm}$



92511301-92063004 软扩
Soft Expansion, $\Phi 9.0\text{mm}$




92511401-92068001 手把
Handle



92511501-92067001 导杆
Guide Rod



92511601-92005001 近端钉导杆
Proximal Screw Guide Rod



92511701-92046001 远端定位架A
Distal Located Bracket



92511801-92046002 远端定位架B
Distal Located Bracket



92511901-92070001 连接锁轮
Connecting Locking Nut, M6/SW5



92512001-92021001 定位卡快
Located Plug



92512101-92042001 连接螺栓
Connecting Bolt



92512201-92001001 加压螺栓
Compression Bolt



KTO 肱骨带锁髓内钉安装器械包 92519901

KTO Humeral Intramedullary Nail Instrument Set



92512301-92001001 辅助定位杆
Auxiliary Located Rod



92101403-92001001 组合扳手
Combination Wrench, SW10



92512501-92021001 临时定位杆
Temporary Located Rod, $\Phi 3.0\text{mm}$



91150403-91060001 限位器扳手
Limiter Wrench, SW3.0



91150404-91060001 锁轮扳手
Locking Nut Wrench, SW5.0



92512601-92021001 定位杆
Located Rod, $\Phi 3.6/\Phi 7.1\text{mm}$



92512701-92002001 外套管
Outer Sleeve, $\Phi 7.1/\Phi 9\text{mm}$



92512801-92080001 套管针
Trocar, $\Phi 7.1\text{mm}$



92512901-92002001 定位杆钻套
Drill Sleeve For Located Rod, $\Phi 3.6\text{mm}$



91113601-91040001 定位杆骨钻
Bone Drill For Located Rod, $\Phi 3.6 \times 180\text{mm}$



92513101-92001001 平头钻
Flat Drill Bit, $\Phi 3.6\text{mm}$



92513201-92001001 梅花起子持杆
Torx Screwdriver Holding Rod, T3.4



92513301-92081001 梅花起子
Torx Screwdriver, T3.4



92513401-92001001 梅花起子 (快装式)
Torx Screwdriver (Quick Coupling), T3.4



91122003-92600001 导针
Guide Pin, $\Phi 2.0 \times 200\text{mm}$



KTO 肱骨带锁髓内钉安装器械包 92519901

KTO Humeral Intramedullary Nail Instrument Set



91122004-91010001 螺纹导针
Threaded Guide Pin, $\Phi 2.0 \times 200\text{mm}$



92513501-92021001 导针测深器
Depth Gauge For Guide Pin, $\Phi 2.0\text{mm}$



92513601-92064001 空心钻
Hollow Drill Bit



92513701-92064001 叶片插入器
Blade Inserter



92513801-92042001 叶片连接杆
Blade Connecting Rod, $\Phi 2.0/\text{M5}$



**

92513901-92002001 锁钉套管
Locking Sleeve, $\Phi 7.1/\Phi 13$



92514001-92061001 导针钻套
Drill Sleeve For Guide Pin, $\Phi 2.0/\Phi 7.1$



92514101-92060001 套管针
Trocac, $\Phi 2.0\text{mm}$



92100602-92046001 滑锤
Slide Hammer



91123003-91010001 导针
Guide Pin, $\Phi 3.0 \times 220$



91113003-91040001 骨钻
Bone Drill, $\Phi 3.0 \times 220\text{mm}$



91150815-92030001 限位器
Limiter, $\Phi 3.0/\text{SW}3.0$



92514401-92002001 钻套
Drill Sleeve, $\Phi 3.0/\Phi 7.1\text{mm}$



92514601-92023001 滑锤导向杆
Slide Hammer For Guide Rod



92514501-92021001 测深器
Depth Gauge

92510001-92000101 空箱
Box