




*** * PEEK 颈椎融合器安装器械包 93659901**

PEEK Thoracolumbar fusion cage instrument set

A



93650101-93065001 上钉器
Screwdriver



93650201-93021001 撑开螺钉
Distractor



93650301-93042001 推进器
Thruster



93650401-93044001 持取器
Holder



93650501-93042001 填充器
Applier



93651601-93063001 植骨台
Bone graft trolley



93650601-93023001 滑锤
Slide hammer



93650701-93003001 开孔器 (弯)
Awl (bend)



93650801-93003001 开孔器 (直)
Awl (straight)

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object.20220616



* * PEEK 颈椎融合器安装器械包 93659901

PEEK Thoracolumbar fusion cage instrument set

A



93650901-93021001 拉钩 (大)
Drag hook (big)



93650902-93021002 拉钩 (小)
Drag hook (little)



93651001-93065001 刮匙
Curet, 3mm



93651101-93065001 终板刮刀
End plate scraper, 4x7mm



93640801-93025001 植骨漏斗
Bone graft funnel
曾用编码 93651201



93640701-93082001 植骨器
Implant bone aretifacts
曾用编码 93651301



93651401-93000101 颈椎牵开器
Cervical retractor



93651506-93021001 试模
Trial Tewplate, 4mm



93651501-93021002 试模
Trial Tewplate, 5mm



93651502-93021003 试模
Trial Tewplate, 6mm



93651503-93021004 试模
Trial Tewplate, 7mm



93651504-93021005 试模
Trial Tewplate, 8mm



93651505-93021006 试模
Test mode, 9mm



93651402-93001101 撑开器 (选配)
Cervical retractor

93650001-93048001 空箱
Box