







*** * 股骨倒打带锁髓内钉安装器械包 92249901**

Femoral Intramedullary Nail Instrument Set



92240101-92001001 显影量尺
Developing Measure Gauge



91123203-91020001 螺纹导针
Threaded Guide Pin, $\Phi 3.2 \times 300\text{mm}$



92240201-92042001 保护套筒
Protective Sleeve



92240301-92065001 近端空心钻
Proximal Hollow Drill Bit, $\Phi 14\text{mm}$



92100401-92003001 快装手柄
Quick Coupling Handle



92100501-92062001 转换接头
Exchanging Joint



92100710-92003001 中空开孔器
Hollow Awl



92100101-92066001 导针把持器
Guide Pin Holder




92100304-92004001 复位棒 (导针引入器)
Reduction Rod (Guide Pin Guider), $\Phi 2.5 \times 320\text{mm}$

Our company have final right to modify all the specification in this catalogue, If there is any inconformity/between contant, model, color, please refer to the object.20220629



*** * 股骨倒打带锁髓内钉安装器械包 92249901**

Femoral Intramedullary Nail Instrument Set



92100201-92061001 球头导针
Ball Head Guide Pin, $\Phi 2.5 \times 800$



91109801-91021001 软扩刀头
Soft Reamer Head, $\Phi 8.5\text{mm}$



91109802-91021002 软扩刀头
Soft Reamer Head, $\Phi 9.0\text{mm}$



91109803-91021003 软扩刀头
SSoft Reamer Head, $\Phi 9.5\text{mm}$



91109805-91021005 软扩刀头
Soft Reamer Head, $\Phi 10.5\text{mm}$



91109807-91021007 软扩刀头
Soft Reamer Head, $\Phi 11.5\text{mm}$



91109804-91021004 软扩刀头
Soft Reamer Head, $\Phi 10\text{mm}$



91109806-91021006 软扩刀头
Soft Reamer Head, $\Phi 11\text{mm}$



91109808-91021008 软扩刀头
Soft Reamer Head, $\Phi 12\text{mm}$



91109809-91021009 软扩刀头
Soft Reamer Head $\Phi 12.5\text{mm}$



91109810-910212010 软扩刀头
Soft Reamer Head, $\Phi 13\text{mm}$



92100803-92063001 软扩刀杆
Soft Expansion Tool Bar



92240401-92069001 手把
Handle



92240501-92080001 连接螺杆
Connecting Rod, SW6.5



92240601-92021001 连接螺杆扳手
Connecting Rod Wrench, SW6.5



92240701-92062001 延长连接螺栓
Extended Connecting Bolt



92101405-92090001 开口扳手
Opening Wrench, SW11



92240801-92062001 主钉打拔杆
Inserter-Extractor Rod For Main Nail , M10



*** * 股骨倒打带锁髓内钉安装器械包 92249901**

Femoral Intramedullary Nail Instrument Set



92100601-92027001 滑锤
Inserter Extractor



92240901-92027001 近端导杆
Proximal Guide Rod



92241001-92060001 手把连接锁轮
Connecting Nut For Handle



91150404-91060001 锁轮扳手
Locking Nut Wrench, SW5.0



92241101-92048001 远端导杆
Distal Guide Rod



92241201-92080001 远端连接锁轮
Distal Locking Nut, M8



92241301-92046001 定位架
Located Bracket



92241401-92002001 外套管/定位杆套管
Outer Sleeve



92241501-92080001 软组织分离器
Soft Tissue Separator



**

92241601-92002001 远端钻套/定位杆钻套
Distal Drill Sleeve/Drill Sleeve For Located Rod, Φ5.2mm



92241701-92001001 T型平头钻
T Flat Drill Bit, Φ5.2mm



9115202-91040001 定位杆骨钻
Bone Drill For Located Rod, Φ5.2×250mm



92241801-92021001 定位杆
Located Rod



92241901-92001001 定位卡块
Located Rod Plug, Φ9mm、Φ10mm



92242001-92001001 定位卡块
Located Rod Plug, Φ11mm、Φ12mm



*** * 股骨倒打带锁髓内钉安装器械包 92249901**

Femoral Intramedullary Nail Instrument Set



92242101-92002001 近端钻套
Proximal Drill Sleeve, $\Phi 4.0\text{mm}$



91114001-91040001 近端骨钻
Proximal Bone Drill Bit, $\Phi 4.0 \times 300\text{mm}$



91150801-91030001 限位器
Limiter, $\Phi 4.0/\text{SW}3.5$



92242201-92021001 直角定位杆
Right-Angled Located Rod, $\Phi 4.0\text{mm}$



92242301-92002001 测深器
Depth Gauge



92242401-92001001 T型锁钉扳手
T Locking Screw Wrench, SW3.5



92242501-92001001 锁钉扳手
Locking Screw Wrench, SW3.5



91150803-91030001 限位器
Depth Gauge, $\Phi 5.2/\text{SW}3.5$



92242601-92081001 远端锁定丝锥
Distal Locking Tap



92242701-92081001 尾帽球头扳手
Ball Head Wrench For End Cap, SW3.5



92242801-92001001 尾帽持杆
Holding Rod For End Cap, SW3.5



92242901-92061001 主钉拔出连接杆
Inserter-Extractor Connecting Rod For Main Nail



91115201-91040001 远端骨钻
Bone Drill For Located Rod, $\Phi 5.2 \times 300\text{mm}$



91124001-91010001 导针
Guide Pin, $\Phi 4.0 \times 300\text{mm}$



92100901-92005001 软扩刀盒
Soft Reamer Head Box

92240001-92002101 空箱
Box