





* * 钛缆固定系统安装器械包 91729902

Titanium Cable Fixation System Instrument Set



91720101-91042001 直引针
Straight Guide Pin



91720201-91023001 弯引针-1
Curve Guide Pin



91720202-91023002 弯引针-2
Curve Guide Pin



91720301-91005001 前压缆锁
Cable Compression Lock(Front), $\Phi 1.3/\Phi 1.8$



91720401-91044001 后压缆锁
Cable Compression Lock(Back), $\Phi 1.3/\Phi 1.8$



91720501-91063001 紧缆钳
Cable Tighten Plier



91720602-91023101 剪缆钳
Cable Cut Plier



91720701-91006101 骨针剪
Skeletal Pins Scissor



91720801-91004101 锁紧扣压紧钳
Compression Plier For Fasten Buckle

Our company have final right to modify all the specification in this catalogue, If there is any inconformitybetween contant, model, color, please refer to the object.20211130



* * 钛缆固定系统安装器械包 91729902

Titanium Cable Fixation System Instrument Set



91720901-91044001 骨针折弯器
Skeletal Pins Bender



91721002-91063001 引缆钩
Cable Guiding Hook



91720302-91005001 前压缆锁(选配)
Cable Compression Lock(Front),Φ1.3



91720303-91005001 前压缆锁 (选配)
Cable Compression Lock(Front),Φ1.8



91720402-91044001 后压缆锁 (选配)
Cable Compression Lock(Back),Φ1.3



91720403-91044002 后压缆锁 (选配)
Cable Compression Lock(Back),Φ1.8

91720002-91027002 空箱
Box



G