


# Anterior Cervical Plate-II System

CE 0197 ISO 13485 FDA



**KANGLI ORTHOPAEDICS  
INSTRUMENT**


## Anterior Cervical Plate-II System



Lower cutting design, the thickness of plate is 1.8mm. Special window design, to be observed easily.



Provide the condition of 0°-12° in the coronal plane between the screw and plate.



Provide the condition of 0°-10° in the vertical plane between the screw and plate.

# Anterior Cervical Plate-II System

**Anterior Cervical Plate-II**



Product Name	Material	Holes	Length	Product No.
Anterior Cervical Plate-II	Pure Titanium	4	22	T6520422
			25	T6520425
			27	T6520427
			30	T6520430
		6	32	T6520632
			35	T6520635
			38	T6520638
			41	T6520641
			44	T6520644
			47	T6520647
		8	50	T6520650
			53	T6520853
			56	T6520856
			59	T6520859
			62	T6520862
			65	T6520865
			68	T6520868
			71	T6520871

**Anterior Cervical Plate-II**



**Fixed Screw  
T654**



Product Name	Material	Diameter (mm)	Length (mm)	Product No.
Fixed Screw	Titanium Alloy	4.0	13	T6544013
			14	T6544014
			15	T6544015
			16	T6544016
			17	T6544017
			18	T6544018

**Rescue Screw  
T655**



Product Name	Material	Diameter (mm)	Length (mm)	Product No.
Rescue Screw	Titanium Alloy	4.5	14	T6554514
			15	T6554515
			16	T6554516
			17	T6554517