

Femoral Intramedullary Nail System

CE 0197 ISO 13485 FDA



KangLi

KANGLI ORTHOPAEDICS
INSTRUMENT

Femoral Intramedullary Nail System



Femoral Intramedullary Nail-DTD



Femoral Reverse Intramedullary Nail Aiming Device

**Femoral Intramedullary Nail
T8194、T8110、T8111**



| Product Name | Material | Diameter | Length | Product No. |
|-----------------------------|----------------|----------|---------------------|-----------------------|
| Femoral Intramedullary Nail | Titanium Alloy | Φ9.4 | 320、340、360、380 | T8194320- T8194380 |
| | | Φ10 | 320、340、360、380、400 | T8110320- T8110400 |
| | | Φ11 | — — 360、380、400 | T8111360- T8111400 |

**Femoral Reverse Intramedullary Nail
T8290、T8210、T8211**



| Product Name | Material | Diameter | Length | Product No. |
|-------------------------------------|----------------|----------|-----------------------------|-----------------------|
| Femoral Reverse Intramedullary Nail | Titanium Alloy | Φ9 | — 260、280、300 | T8290260- T8290300 |
| | | Φ10 | 240、260、280、300、320、340、360 | T8210240- T8210360 |
| | | Φ11 | 300、320、340、360 | T8211300- T8211360 |

**Φ4.5 Full Threaded Locking Screw
T852**



| Product Name | Material | Diameter | Length | Product No. |
|-----------------------------|----------------|----------|----------------------------|-----------------------|
| Full Threaded Locking Screw | Titanium Alloy | Φ4.5 | 30 | T8524030 |
| | | | 32 | T8524032 |
| | | | 35 | T8524035 |
| | | | 38 | T8524038 |
| | | | 40、45、50、55、60、65、70、75、80 | T8524040- T8524080 |

**Φ6.5 Full Threaded Locking Screw
T853**



| Product Name | Material | Diameter | Length | Product No. |
|-----------------------------|----------------|----------|-------------------------|-----------------------|
| Full Threaded Locking Screw | Titanium Alloy | Φ6.5 | 60、65、70、75、80、85、90、95 | T8536560- T8536595 |