


# JY-Locking Compression Plate Series

CE 0197  
ISO 13485 FDA



KangLi

KANGLI ORTHOPAEDICS  
INSTRUMENT

## JY Locking Plate Series

### Indications

- Open and closed wedged Osteotomies of the medial proximal tibia.
- Medial or lateral Unicondylar osteoarthritis, malalignment of proximal tibia.
- Idiopathic or traumatic varus or valgus of proximal tibia.



### JY HTO Tibial Osteotomy Standard Locking Plate 232140



Product Name	Material	Holes	Length	Product No.
JY HTO Tibial Osteotomy Standard Locking Plate	Pure Titanium	4	115	23214040

Note: used for  $\Phi$ 5.0 locking screw(Anodized)、HA4.5 screw (Sample: 4 holes)

### JY HTO Tibial Osteotomy Small Locking Plate 232141



Product Name	Material	Holes	Length	Product No.
JY HTO Tibial Osteotomy Small Locking Plate	Pure Titanium	4	112	23214140

Note: used for  $\Phi$ 5.0 locking screw(Anodized)、HA4.5螺钉 (Sample: 4 holes)

## JY Locking Plate Series

**JY S Clavicle Locking Plate(Low Stress)(Anodized)  
23083**



Product Name	Material	Holes	Length	Product No.
JY S Clavicle Locking Plate (Low Stress)(Anodized)	Pure Titanium	7L	91	23083071
		7R	91	23083072
		8L	102	23083081
		8R	102	23083082
		9L	113	23083091
		9R	113	23083092
		10L	123	23083101
		10R	123	23083102
		12L	145	23083121
		12R	145	23083122

Note: used for  $\Phi$ 3.5 locking screw(Anodized)、HA3.5 screw(Sample: 9 holes)

**JY Clavicle Stem Locking Plate(Left,Right)  
23071、 23077**



Product Name	Material	Holes	Length	Product No.
JY Clavicle Stem Locking Plate (Left,Right)	Pure Titanium	short 7 L	90	23077071
		short 7 R	90	23077072
		7L	109	23071071
		7R	109	23071072
		8L	122	23071081
		8R	122	23071082
		9L	136	23071091
		9R	136	23071092
		11L	162	23071111
		11R	162	23071112

Note: used for  $\Phi$ 3.5 locking screw(Anodized)、HA3.5 screw(Sample: 8 holes)

**JY Distal Clavicle Locking Plate(Left,Right)  
23084**



Product Name	Material	Holes	Length	Product No.
Locking Plate (Left,Right)	Pure Titanium	5L	89	23084051
		5R	89	23084052
		6L	100	23084061
		6R	100	23084062
		8L	125	23084081
		8R	125	23084082
		10L	148	23084101
		10R	148	23084102

Note: used for  $\Phi$ 2.7、 $\Phi$ 3.5 locking screw(Anodized)、HA3.5 screw(Sample: 6 holes)

**JY Clavicle Hook Locking Plate(Left,Right)  
23075**



**JY Clavicle Hook Locking Plate(Left,Right)  
23079**



Product Name	Material	Height	Holes	Length	Product No.
JY Clavicle Hook Locking Plate(Left,Right)	Pure Titanium	14	8L	54	23075081
			8R	54	23075082
			9L	66	23075091
			9R	66	23075092
			11L	92	23075111
			11R	92	23075112
JY Clavicle Hook Locking Plate(Left,Right)	Pure Titanium	17	8L	54	23079081
			8R	54	23079082
			9L	66	23079091
			9R	66	23079092
			11L	92	23079111
			11R	92	23079112

Note: used for  $\Phi$ 2.7、 $\Phi$ 3.5 locking screw(Anodized)、HA3.5screw(Sample:11 holes)

## JY Locking Plate Series

**JY Humerus Tubercle Locking Plate(Left,Right)  
23073**



Product Name	Material	Holes	Length	Product No.
JY Humerus Tubercle Locking Plate(Left,Right)	Pure Titanium	5L	42	23073051
		5R	42	23073052
Note: used for $\Phi 2.7$ 、 $\Phi 3.5$ locking screw(Anodized)、HA2.7、HB4.0 screw (Sample: 5 holes)				

**JY Proximal Humerus Locking Plate  
23051**



Product Name	Material	Holes	Length	Product No.
JY Proximal Humerus Locking Plate	Pure Titanium	3	91	23051030
		5	115	23051050
Note: used for $\Phi 3.5$ locking screw(Anodized)、HA3.5、HB4.0 screw (Sample: 5 holes)				

**JY Proximal Humerus Locking Plate(Long)  
23052**



Product Name	Material	Holes	Length	Product No.
JY Proximal Humerus Locking Plate(Long)	Pure Titanium	5	144	23052050
		7	180	23052070
		9	216	23052090
		11	252	23052110
		13	288	23052130
Note: used for $\Phi 3.5$ locking screw(Anodized)、HA3.5、HB4.0 screw (Sample: 5 holes)				

**JY Distal Lateral Humerus Locking Plate(Left,Right)  
23081**



Product Name	Material	Holes	Length	Product No.
JY Distal Lateral Humerus Locking Plate(Left,Right)	Pure Titanium	5L	92	23081051
		5R	92	23081052
		7L	118	23081071
		7R	118	23081072
		9L	144	23081091
		9R	144	23081092
		11L	170	23081111
		11R	170	23081112
Note: used for $\Phi 3.5$ locking screw(Anodized)、HA3.5、HB4.0screw (Sample: 11 holes)				

**JY Distal humerus Medial Locking Plate  
23085**



Product Name	Material	Holes	Length	Product No.
JY Distal humerus Medial Locking Plate	Pure Titanium	4L	82	23085041
		4R	82	23085042
		5L	94	23085051
		5R	94	23085052
		7L	118	23085071
		7R	118	23085072
		9L	142	23085091
		9R	142	23085092
		11L	166	23085111
		11R	166	23085112
		Note: used for $\Phi 3.5$ locking screw(Anodized)、 $\Phi 2.7$ locking screw(Anodized)、HB4.0screw (Sample: 9 holes)		

## JY Locking Plate Series

**JY Distal Humerus Subcondyle Locking Plate(Left,Right)  
23074**



Product Name	Material	Holes	Length	Product No.
JY Distal Humerus Subcondyle Locking Plate(Left,Right)	Pure Titanium	12L	142	23074121
		12R	142	23074122
		15L	188	23074151
		15R	188	23074152
Note: used for $\Phi$ 3.5 locking screw(Anodized)、HA3.5 screw (Sample:15 holes)				

**JY Proximal Ulna Locking Plate (Left,Right)  
23081**



Product Name	Material	Holes	Length	Product No.
JY Proximal Ulna Locking Plate (Left,Right)	Pure Titanium	4L	112	23081041
		4R	112	23081042
		6L	140	23081061
		6R	140	23081062
		8L	168	23081081
		8R	168	23081082
Note: used for $\Phi$ 3.5 locking screw(Anodized)、HA3.5 screw(Sample: 6 holes)				

**JY 2.7 Distal Radius Buttress Locking Plate(Mini)(Left,Right)  
23053**



Product Name	Material	Holes	Length	Product No.
JY 2.7 Distal Radius Buttress Locking Plate(Mini)(Left,Right)	Pure Titanium	4L	55	23053041
		4R	55	23053042
		5L	63	23053051
		5R	63	23053052
		7L	79	23053071
		7R	79	23053072
		9L	95	23053091
		9R	95	23053092
		11L	111	23053111
		11R	111	23053112
		13L	127	23053131
		13R	127	23053132
Note: used for $\Phi$ 2.7 locking screw(Anodized)、HA2.7 screw (Sample: 4 holes)				

**JY 2.7 Distal Radius Buttress Locking Plate(Large)(Left,Right)  
23054**



Product Name	Material	Holes	Length	Product No.
JY 2.7 Distal Radius Buttress Locking Plate(Large)(Left,Right)	Pure Titanium	4L	55	23054041
		4R	55	23054042
		5L	63	23054051
		5R	63	23054052
		7L	79	23054071
		7R	79	23054072
		9L	95	23054091
		9R	95	23054092
Note: used for $\Phi$ 2.7 locking screw(Anodized)、HA2.7 screw (Sample: 5 holes)				

**JY Distal Back Radius Locking Plate(Left,Right)  
23055**



Product Name	Material	Holes	Product No.
JY Distal Back Radius Locking Plate (Left,Right)	Pure Titanium	14L	23055141
		14R	23055142
Note: used for $\Phi$ 2.7 locking screw(Anodized) (Sample: 14 holes)			

## JY Locking Plate Series

**JY Distal Ulna Locking Plate  
23080**



Product Name	Material	Holes	Length	Product No.
JY Distal Ulna Locking Plate	Pure Titanium	8	45	23080080

Note: used for  $\Phi$ 2.0 locking screw(Anodized)、HA2.0 screw (Sample: 8 holes)

**JY 2.7 Radius Head Locking Plate(Left,Right)  
23056**



Product Name	Material	Holes	Length	Product No.
JY 2.7 Radius Head Locking Plate(Left,Right)	Pure Titanium	3L	40	23056031
		3R	40	23056032
		4L	48	23056041
		4R	48	23056042
		5L	56	23056051
5R	56	23056052		

Note: used for  $\Phi$ 2.7 locking screw(Anodized)(Sample: 4 holes)

**JY 2.7 Distal Fibula Reconstruction Locking Plate(Left,Right)  
23201**



Product Name	Material	Holes	Length	Product No.
JY 2.7 Distal Fibula Reconstruction Locking Plate (Left,Right)	Pure Titanium	4L	68	23201041
		4R	68	23201042
		6L	88	23201061
		6R	88	23201062
		8L	108	23201081
		8R	108	23201082
		10L	128	23201101
		10R	128	23201102
		12L	148	23201121
		12R	148	23201122

Note: used for  $\Phi$ 2.7 locking screw(Anodized)、HA2.7 screw (Sample: 6 holes)

**JY 2.7 Distal Lateral Fibula Locking Plate- II (Left,Right)  
23226**



Product Name	Material	Holes	Length	Product No.
JY 2.7 Distal Lateral Fibula Locking Plate- II (Left,Right)	Pure Titanium	3L	86	23226031
		3R	86	23226032
		5L	112	23226051
		5R	112	23226052
		7L	138	23226071
		7R	138	23226072

Note: used for  $\Phi$ 2.7  $\Phi$ 3.5 locking screw(Anodized)、HA2.7、HA3.5 screw (Sample: 3 holes)

**JY Combo Minimally Invasive Calcaneal Locking Plate(Left,Right)  
23209**



Product Name	Material	Holes	Product No.
JY Combo Minimally Invasive Calcaneal Locking Plate (Left,Right)	Pure Titanium	9L	23209091
		9R	23209092

Note: used for  $\Phi$ 3.5 locking screw(Anodized)

## JY Locking Plate Series

### JY Combo Minimally Invasive Calcaneal Locking Plate(Max)(Left,Right)

23232



Product Name	Material	Holes	Product No.
JY Combo Minimally Invasive Calcaneal Locking Plate (Max)(Left,Right)	Pure Titanium	8L	23232081
		8R	23232082

Note: used for  $\Phi$ 3.5 locking screw(Anodized)

### JY Combo Minimally Invasive Calcaneal Locking Plate(Mini)(Left,Right)

23212



Product Name	Material	Holes	Product No.
JY Combo Minimally Invasive Calcaneal Locking Plate (Mini)(Left,Right)	Pure Titanium	8L	23212081
		8R	23212082

Note: used for  $\Phi$ 3.5 locking screw(Anodized)

### JY Posterior Tuberosity Minimally Invasive Calcaneal Locking Plate(Left,Right)

23208



Product Name	Material	Holes	Product No.
JY Posterior Tuberosity Minimally Invasive Calcaneal Locking Plate(Left,Right)	Pure Titanium	6L	23208061
		6R	23208062

Note: used for  $\Phi$ 3.5 locking screw(Anodized)

### JY Posterior Tuberosity Minimally Invasive Calcaneal Locking Plate (Mini)(Left,Right)

23211



Product Name	Material	Holes	Product No.
JY Posterior Tuberosity Minimally Invasive Calcaneal Locking Plate (Mini)(Left,Right)	Pure Titanium	5L	23211051
		5R	23211052

Note: used for  $\Phi$ 3.5 locking screw(Anodized)

### JY Anterior Process Minimally Invasive Calcaneal Locking Plate(Left,Right)

23210



Product Name	Material	Holes	Product No.
JY Anterior Process Minimally Invasive Calcaneal Locking Plate (Left,Right)	Pure Titanium	5L	23210051
		5R	23210052

Note: used for  $\Phi$ 3.5 locking screw(Anodized)

## JY Locking Plate Series


**JY Anterior Process Minimally Invasive Calcaneal Locking Plate (Mini)(Left,Right)  
23213**



Product Name	Material	Holes	Product No.
JY Anterior Process Minimally Invasive Calcaneal Locking Plate (Mini)(Left,Right)	Pure Titanium	5L	23213051
		5R	23213052

Note: used for  $\Phi$ 3.5 locking screw(Anodized)

**JY Calcaneus Anatomical Locking Plate(Mini)(Left,Right)  
23227**



Product Name	Material	Holes	Length	Product No.
JY Calcaneus Anatomical Locking Plate(Mini)(Left,Right)	Pure Titanium	11L	58	23227111
		11R	58	23227112

Note: used for  $\Phi$ 3.5 locking screw(Anodized)


**JY Calcaneus Anatomical Locking Plate(Mid)(Left,Right)  
23228**



Product Name	Material	Holes	Length	Product No.
JY Calcaneus Anatomical Locking Plate(Mid)(Left,Right)	Pure Titanium	11L	64	23228111
		11R	64	23228112

Note: used for  $\Phi$ 3.5 locking screw(Anodized)

**JY Calcaneus Anatomical Locking Plate(Max)(Left,Right)  
23230**



Product Name	Material	Holes	Length	Product No.
JY Calcaneus Anatomical Locking Plate(Max)(Left,Right)	Pure Titanium	11L	71	23230111
		11R	71	23230112

Note: used for  $\Phi$ 3.5 locking screw(Anodized)

**JY 3.5 Proximal Medial Tibial Locking Plate(Left,Right)  
23206**



Product Name	Material	Holes	Length	Product No.
JY 3.5 Proximal Medial Tibial Locking Plate(Left,Right)	Pure Titanium	5L	106	23206051
		5R	106	23206052
		7L	134	23206071
		7R	134	23206072
		9L	162	23206091
		9R	162	23206092
		11L	190	23206111
		11R	190	23206112
		13L	218	23206131
		13R	218	23206132

Note: used for  $\Phi$ 3.5 locking screw(Anodized) , HA3.5screw (Sample: 7 holes)



## JY Locking Plate Series

**JY 5.0 Proximal Medial Tibial Locking Plate(Left,Right)**  
**23207**



Product Name	Material	Holes	Length	Product No.
JY 5.0 Proximal Medial Tibial Locking Plate(Left,Right)	Pure Titanium	5L	125	23207051
		5R	125	23207052
		7L	161	23207071
		7R	161	23207072
		9L	197	23207091
		9R	197	23207092
		11L	233	23207111
		11R	233	23207112
		13L	269	23207131
		13R	269	23207132
		15L	305	23207151
		15R	305	23207152

Note: used for  $\Phi$ 5.0 locking screw(Anodized)、HA4.5screw (Sample: 9 holes)

**JY KSL Proximal Lateral Tibial Locking Plate(Left,Right)**  
**23205**



Product Name	Material	Holes	Length	Product No.
JY KSL Proximal Lateral Tibial Locking Plate(Left,Right)	Pure Titanium	5L	98	23205051
		5R	98	23205052
		7L	130	23205071
		7R	130	23205072
		9L	162	23205091
		9R	162	23205092
		11L	194	23205111
		11R	194	23205112
		13L	226	23205131
		13R	226	23205132

Note: used for  $\Phi$ 3.5 locking screw(Anodized)、HA3.5 screw (Sample: 9 holes)

**JY Proximal Lateral Tibial Locking Plate(Left,Right)**  
**23231**



Product Name	Material	Holes	Length	Product No.
JY Proximal Lateral Tibial Locking Plate(Left,Right)	Pure Titanium	4L	148	23231041
		4R	148	23231042
		6L	184	23231061
		6R	184	23231062
		8L	220	23231081
		8R	220	23231082
		10L	256	23231101
		10R	256	23231102
		12L	292	23231121
		12R	292	23231122
		14L	328	23231141
		14R	328	23231142

Note: used for  $\Phi$ 5.0 locking screw(Anodized)、HA4.5 screw (Sample: 8 holes)

**JY Distal Tibial Medial Malleolus Locking Plate**  
**23231**



Product Name	Material	Holes	Length	Product No.
JY Distal Tibial Medial Malleolus Locking Plate	Pure Titanium	4	88	23231040
		6	111	23231060

Note: used for  $\Phi$ 3.5 locking screw(Anodized)、HA3.5 screw (Sample: 6 holes)

## JY Locking Plate Series

**JY Distal Medial Tibial Locking Plate(Left,Right)  
23203**



Product Name	Material	Holes	Length	Product No.
JY Distal Medial Locking Plate(Left,Right)	Pure Titanium	5L	148	23203051
		5R	148	23203052
		7L	180	23203071
		7R	180	23203072
		9L	212	23203091
		9R	212	23203092
		11L	244	23203111
		11R	244	23203112
		13L	276	23203131
13R	276	23203132		

Note: used for  $\Phi$ 3.5 locking screw(Anodized)、HA3.5 screw (Sample: 7 holes)

**JY Distal Anterolateral Tibial Locking Plate(Left,Right)  
23204**



Product Name	Material	Holes	Length	Product No.
JY Distal Anterolateral Tibial Locking Plate(Left,Right)	Pure Titanium	5L	123	23204051
		5R	123	23204052
		7L	155	23204071
		7R	155	23204072
		9L	187	23204091
		9R	187	23204092
		11L	219	23204111
		11R	219	23204112
		13L	251	23204131
		13R	251	23204132

Note: used for  $\Phi$ 3.5 locking screw(Anodized)、HA3.5 screw (Sample: 7 holes)

**JY Distal Tibial Posterior Malleolus Locking Plate  
23234**



Product Name	Material	Holes	Length	Product No.
JY Distal Tibial Posterior Malleolus Locking Plate	Pure Titanium	4	66	23234040
		6	86	23234060

Note: used for  $\Phi$ 3.5 locking screw、HA3.5 screw (Sample: 6 holes)

**JY Pelvic Reconstruction Locking Compression Plate  
(Curved) (Anodized)  
23235**



Product Name	Material	Holes	Length	Product No.
VA Pelvic Reconstruction Locking Compression Plate(Curved) (Anodized)	Pure Titanium	5	62	23235050
		6	74	23235060
		8	98	23235080
		10	122	23235100
		12	146	23235120
		14	170	23235140

Note: used for  $\Phi$ 3.5 locking screw(Anodized) (Sample: 10 holes)

**JY Pelvic Reconstruction Locking Compression Plate  
(Straight) (Anodized)  
23014**



Product Name	Material	Holes	Length	Product No.
JY pelvic Reconstruction Locking Compression Plate (Straight) (Anodized)	Pure Titanium	6	76	23014060
		8	102	23014080
		10	128	23014100
		12	154	23014120
		14	180	23014140
		16	206	23014160
		18	232	23014180
		20	258	23014200

Note: used for  $\Phi$ 3.5 locking screw (Anodized)、HA3.5 screw (Sample: 8 holes)

## JY Locking Plate Series

**JY Femoral Prosthesis Revision Locking Plate (Straight)**  
**23013**



Product Name	Material	Holes	Length	Product No.
JY Femoral Prosthesis Revision Locking Plate (Straight)	Pure Titanium	10	257	23013100
		12	301	23013120

Note: used for  $\Phi$ 5.0 locking screw(Anodized) (Sample: 10 holes)

**JY Femoral Prosthesis Revision Locking Plate**  
**23233**



Product Name	Material	Holes	Length	Product No.
JY Femoral Prosthesis Revision Locking Plate	Pure Titanium	5	125	23233050
		14L	271	23233141
		14R	271	23233142
		16L	297	23233161
		16R	297	23233162

Note: used for  $\Phi$ 5.0 locking screw(Anodized) (Sample: 14 holes)

**JY KUH Proximal Femur Locking Plate**  
**23239**



Product Name	Material	Holes	Length	Product No.
JY KUH Proximal Femur Locking Plate	Pure Titanium	3L	121	23239031
		3R	121	23239032
		4L	141	23239041
		4R	141	23239042
		5L	161	23239051
		5R	161	23239052
		7L	201	23239071
		7R	201	23239072
		9L	241	23239091
		9R	241	23239092
		11L	281	23239111
		11R	281	23239112
		14L	340	23239141
14R	340	23239142		

Note: used for  $\Phi$ 6.5locking screw in head (Anodized),other  $\Phi$ 5.0 locking screw (Anodized) and HA4.5 screw (Sample: 11 holes)

**JY Distal Lateral Femoral Locking Plate(Left,Right)**  
**23202**



Product Name	Material	Holes	Length	Product No.
JY Distal Lateral Femoral Locking Plate(Left,Right)	Pure Titanium	5L	178	23202051
		5R	178	23202052
		7L	222	23202071
		7R	222	23202072
		9L	266	23202091
		9R	266	23202092
		11L	310	23202111
		11R	310	23202112
		13L	353	23202131
		13R	353	23202132

Note: used for  $\Phi$ 5.0 locking screw(Anodized)、HA4.5 screw (Sample: 11 holes)

## JY Locking Plate Series

**JY Femoral Locking Plate**  
23012、23004



Product Name	Material	Holes	Length	Product No.
JY Femoral Locking Plate	Pure Titanium	8short	194	23012080
		8	213	23004080
		9	235	23004090
		10	257	23004100
		12	301	23004120
		14	345	23004140

Note: used for  $\Phi$ 5.0 locking screw(Anodized)、HA4.5 screw (Sample: 9 holes)

**JY Tibial Metaphyseal Locking Plate**  
23216



Product Name	Material	Holes	Length	Product No.
JY Tibial Metaphyseal Locking Plate	Pure Titanium	10	158	23216010
		12	194	23216012
		14	230	23216014
		16	266	23216016

Note: used for  $\Phi$ 3.5 locking screw(Anodized)、HA3.5 screw (Sample: 12 holes)

**JY Tibial Locking Plate**  
23003



Product Name	Material	Holes	Length	Product No.
JY Tibial Locking Plate	Pure Titanium	7	159	23003070
		8	177	23003080
		9	195	23003090
		10	213	23003100
		12	249	23003120
		14	285	23003140

Note: used for  $\Phi$ 5.0 locking screw(Anodized)、HA4.5 screw (Sample: 8 holes)

**JY Humeral Metaphyseal Locking Plate**  
23070



Product Name	Material	Holes	Length	Product No.
JY Humeral Metaphyseal Locking Plate	Pure Titanium	7	112	23070070
		9	144	23070090
		11	176	23070110
		13	208	23070130

Note: used for  $\Phi$ 3.5 locking screw(Anodized)、HA3.5 screw (Sample: 9 holes)

**JY Humeral Locking Plate**  
23002



Product Name	Material	Holes	Length	Product No.
JY Humeral Locking Plate	Pure Titanium	7	131	23002070
		8	147	23002080
		9	163	23002090
		10	179	23002100
		12	211	23002120
		14	243	23002140

Note: used for  $\Phi$ 3.5 locking screw(Anodized)、HA3.5 screw (Sample: 8 holes)

**JY Ulna Radius Locking Plate**  
23005



Product Name	Material	Holes	Length	Product No.
JY Ulna Radius Locking Plate	Pure Titanium	6	106	23005060
		7	122	23005070
		8	138	23005080
		9	154	23005090
		10	170	23005100

Note: used for  $\Phi$ 3.5 locking screw(Anodized)、HA3.5 screw (Sample: 6 holes)

## JY Locking Plate Series

**JY Reconstruction Locking Plate  
23001**



Product Name	Material	Holes	Length	Product No.
JY Reconstruction Locking Plate	Pure Titanium	5	87	23001050
		6	103	23001060
		7	119	23001070
		8	135	23001080
		9	151	23001090
		10	167	23001100
		12	199	23001120
		14	231	23001140

Note: used for  $\Phi$ 3.5 locking screw(Anodized)、HA3.5 screw (Sample: 7 holes)

**JY Scapula Acromion Locking Plate(Left,Right)  
23077**



Product Name	Material	Holes	Product No.
JY Scapula Acromion Locking Plate(Left,Right)	Pure Titanium	9L	23077091
		9R	23077092

Note: used for  $\Phi$ 2.7 locking screw(Anodized) (Sample: 9 holes)

**JY Scapula Lateral Margin Locking Plate(Left,Right)  
23078**



Product Name	Material	Holes	Product No.
JY Scapula Lateral Margin Locking Plate(Left,Right)	Pure Titanium	9L	23078091
		9R	23078092

Note: used for  $\Phi$ 2.7 locking screw(Anodized) (Sample: 9 holes)

**JY Scapula Glenoid Locking Plate(Left,Right)  
23076**



Product Name	Material	Holes	Product No.
JY Scapula Glenoid Locking Plate(Left,Right)	Pure Titanium	4L	23076041
		4R	23076042

Note: used for  $\Phi$ 2.7 locking screw(Anodized) (Sample: 4 holes)


**JY Scapula Medial Margin Locking Plate(Left,Right)  
23082**



Product Name	Material	Holes	Product No.
JY Scapula Medial Margin Locking Plate(Left,Right)	Pure Titanium	10L	23082101
		10R	23082102

Note: used for  $\Phi$ 2.7 locking screw(Anodized) (Sample: 10 holes)

## JY Locking Plate Series



**LCC Locking composite pressure connecting plate**  
23236



Product Name	Material	Holes	Product No.
LCC Locking composite pressure connecting plate	Pure Titanium	3	23236030
		4	23236040

**LCC Dynamic Rod**  
33507



Product Name	Material	Diameter	Length	Product No.
LCC Dynamic Rod	Titanium Alloy	10	75-130 The specification is 5mm	33507065- 33507120
LCC Anti-rotation screw			6.5	75-130 The specification is 5mm

**LCC Anti-rotation screw**  
33506



Product Name	Material	Diameter	Length	Product No.
Torx Locking Screw (Anodized)	Titanium Alloy	5.0	14-50 The specification is 2mm	T4935014- T4935050
			55-100 The specification is 5mm	T4935055- T4935100

**Torx Locking Screw (Anodized)**  
T493



## JY Locking Plate Series

### Sterile JY rib locking fixation system




### Character Of The System

- No need to remove periosteum, ensure adequate blood supply.
- Lockong plate provides stable fixation.
- Pre-plastion plate,reduce surgery time.
- Low profile design.
- Depth limit drill,Prevent drilling too deep.

# JY Locking Plate Series

## Sterile JY rib locking fixation system



### JY rib locking fixation system (Sterile) 25615



Product Name	Material	Product specifications	Product No.	Part used
Sterile JY rib locking fixation system	Pure Titanium	8	25615080	Universal
		15L ( 3L )	25615151	Left 3rd rib
		15R ( 3R )	25615152	Right 3rd rib
		16L(4-5L)	25615161	Left 4th/5th rib
		16R(4-5R)	25615162	Right 4th/5th rib
		17L(6-7L)	25615171	Left 6th/7th rib
		17R(6-7R)	25615172	Right 6th/7th rib
		18L(8-9L)	25615181	Left 8th/9th rib
		18R(8-9R)	25615182	Right 8th/9th rib

### Sterile crossed screw (Anodized) (Sterile) 35008



Product Name	Material	Diameter	Length	Product No.
Sterile crossed screw (Anodized)	Titanium Alloy	3.0	6	35008006
			8	35008008
			10	35008010
			12	35008012
			14	35008014



## Screw Series (Anodized)

**VA Locking Screw(Anodized)**  
**33157**



Product Name	Material	Diameter	Length	Product No.
Universal Locking Screw(Anodized)	Titanium Alloy	1.5	6、7、8、10、12、14、16、18、20、22	33157006-33157022

**VA Locking Screw(Anodized)**  
**33158**



Product Name	Material	Diameter	Length	Product No.
Universal Locking Screw(Anodized)	Titanium Alloy	2.0	6、7、8、10、12、14、16、18、20、22	33158006-33158022

**VA Locking Screw(Anodized)**  
**T480**



Product Name	Material	Diameter	Length	Product No.
Universal Locking Screw(Anodized)	Titanium Alloy	2.7	10、12、14、16、18、20、22、24、26、28、	T4802710-T4802728
			30、35、40	T4802730-T4802740

**Universal Locking Screw(Anodized)**  
**33151**



Product Name	Material	Diameter	Length	Product No.
Universal Locking Screw(Anodized)	Titanium Alloy	3.5	10-50 The specification is 2mm	33151010-33151050
			55-80 The specification is 5mm	33151055-33151080

**Triangle Locking Screw(Anodized)**  
**33009**



Product Name	Material	Diameter	Length	Product No.
Triangle Locking Screw (Anodized)	Titanium Alloy	2.0	6	33009006
			7	33009007
			8	33009008
			10-22 The specification is 2mm	33009010-33009022

**Torx Locking Screw(Anodized)**  
**T491**



Product Name	Material	Diameter	Length	Product No.
Torx Locking Screw (Anodized)	Titanium Alloy	2.7	6-60 The specification is 2mm	T4912706-T4912760

## Screw Series (Anodized)

### Torx Locking Screw (Anodized) T492



Product Name	Material	Diameter	Length	Product No.
Torx Locking Screw (Anodized)	Titanium Alloy	3.5	10-50 The specification is 2mm	T4923510- T4923550
			55-80 The specification is 5mm	T4923555- T4923580

### Torx Locking Screw (Anodized) T493



Product Name	Material	Diameter	Length	Product No.
Torx Locking Screw (Anodized)	Titanium Alloy	5.0	14-50 The specification is 2mm	T4935014- T4935050
			55-100 The specification is 5mm	T4935055- T4935100

### Metal Bone Fracture Screw (Anodized) T52235



Product Name	Material	Diameter	Length	Product No.
Metal Bone Fracture Screw (Anodized)	Titanium Alloy	HA3.5	12-50 The specification is 2mm	T5223512- T5223550

### Metal Bone Fracture Screw (Anodized) T52245



Product Name	Material	Diameter	Length	Product No.
Metal Bone Fracture Screw (Anodized)	Titanium Alloy	HA4.5	20-50 The specification is 2mm	T5224520- T5224550
			55、60	T5224555 T5224560

### Filling Screw (Anodized) 33005、33006



Product Name	Material	Specification	Product No.
Filling Screw (Anodized)	Titanium Alloy	3.5*1.5mm	33005104
		5.0*1.5mm	33006105