

# JZWA Spinal Fixation System (Minimally Invasive)

CE 0197 ISO 13485 FDA




KangLi

KANGLI ORTHOPAEDICS  
INSTRUMENT

## KAS Spinal Fixation Series

### JZWA 5.5 Spinal Fixation System(Minimally Invasive)(4370)

The pre-bending rod is more suitable for the normal physiological curvature of human body.



Double thread and quadruple thread design  
Faster screw insertion  
Consistent screw insertion speed  
Pedicle cortex fixation, increase screw strength

Hollow Type



JZWA Minimally Invasive  
U-rotating Pedicle  
Screw allows the screw body to rotate along the longitudinal axis of the screw

## KAS Spinal Fixation Series

**JZWA Minimally Invasive Multiaxial Pedicle Screw  
43703**



Product Name	Material	Diameter	Length	Product No.
JZWA Minimally Invasive Multiaxial Pedicle Screw	Titanium Alloy	5.5	35	43703235
			40	43703240
			45	43703245
		6.0	35	43703335
			40	43703340
			45	43703345
		6.5	50	43703350
			35	43703435
			40	43703440
			45	43703445
		7.0	50	43703450
			35	43703535
			40	43703540
			45	43703545
7.5	50	43703550		
	35	43703635		
			40	43703640

**JZWA Minimally Invasive Monoaxial Pedicle Screw  
43704**




Product Name	Material	Diameter	Length	Product No.
JZWA Minimally Invasive Monoaxial Pedicle Screw	Titanium Alloy	5.5	35	43704235
			40	43704240
			45	43704245
		6.0	35	43704335
			40	43704340
			45	43704345
		6.5	50	43704350
			35	43704435
			40	43704440
			45	43704445
		7.0	50	43704450
			35	43704535
			40	43704540
			45	43704545
7.5	50	43704550		
	35	43704635		
			40	43704640

**JZWA Arc Shaped Minimally Invasive Rod  
43703C**



Product Name	Material	Total Length	Effective Length	Product No.
JZWA Arc Shaped Minimally Invasive Rod	Titanium Alloy	45	35	43703045
		50	40	43703050
		55	45	43703055
		60	50	43703060
		65	55	43703065
		70	60	43703070
		75	65	43703075
		80	70	43703080
		85	75	43703085
		90	80	43703090
		95	85	43703095
		100	90	43703100
		105	95	43703105
		110	100	43703110
115	105	43703115		
125	115	43703125		
135	125	43703135		
145	135	43703145		
155	145	43703155		

**JZWA Straight Minimally Invasive Rod  
43706**



Product Name	Material	Total Length	Effective Length	Product No.
JZWA Straight Minimally Invasive Rod	Titanium Alloy	45	35	43706045
		50	40	43706050
		55	45	43706055
		60	50	43706060
		65	55	43706065
		70	60	43706070
		75	65	43706075
		80	70	43706080
		85	75	43706085
		90	80	43706090
		95	85	43706095
		100	90	43706100
		105	95	43706105
		110	100	43706110
115	105	43706115		
125	115	43706125		
135	125	43706135		
145	135	43706145		
155	145	43706155		