


KSS-I 6.0 Spinal Fixation System (Torx)

CE 0197 ISO 13485 FDA




KANGLI ORTHOPAEDICS
INSTRUMENT

KSS Spinal Fixation Series

Kss- I 6.0 Spinal Fixation System (Torx)



The top of the screw adopt the design of the Trapezoid type, to lower the external tension minimally, it is the best that the ability of anti-prolapsus and locking.




Up-to-date. design, that is easily operated for the Bayonet type and stretch hatch.



Special design of the torx nut can avoid slipping during the surgery.



The design for the cancellous bone, that is most suitable for the centrum. The increment diameter of the bottom can be less the risk of break.



KSS Spinal Fixation Series


**KSS 6.0 Monoaxial Pedicle Screw
T581**




**KSS 6.0 Monoaxial Reduction Pedicle Screw
T588**



**KSS 6.0 Multiaxial Reduction Pedicle Screw (Hexagon Type)
T587**



**KSS 6.0 Multiaxial Reduction Pedicle Screw (Torx Type)
T583**



**Rod 6.0
T6161**



**EF 6.0 Crosslinks
T619**



Product Name	Material	Diameter	Length	Product No.
KSS 6.0 Monoaxial Pedicle Screw	Titanium Alloy	4.5	30	T5814530
			35	T5814535
		5.0	30	T5815030
			35	T5815035
		5.5	35	T5815535
			40	T5815540
		6.0	35	T5816035
			40	T5816040
			45	T5816045
		6.5	35	T5816535
			40	T5816540
			45	T5816545
		7.0	50	T5816550
			35	T5817035
		KSS 6.0 Reduction Monoaxial Pedicle Screw	Titanium Alloy	4.5
30	T5885030			
5.0	35			T5885035
	35			T5885535
5.5	40			T5885540
	45			T5885545
6.0	35			T5886035
	40			T5886040
	45			T5886045
6.5	50			T5886050
	35			T5886535
	40			T5886540
7.0	45			T5886545
	50			T5886550
KSS 6.0 KSS 6.0 Multiaxial Reduction Pedicle Screw (Hexagon Type)	Titanium Alloy			4.5
		35	T5874535	
		5.0	30	T5875030
			35	T5875035
		5.5	40	T5875040
			30	T5875530
		6.0	35	T5875535
			40	T5875540
			45	T5875545
		6.5	30	T5876030
			35	T5876035
			40	T5876040
		7.0	45	T5876045
			50	T5876050
		KSS 6.0 KSS 6.0 Multiaxial Reduction Pedicle Screw (Torx Type)	Titanium Alloy	4.5
35	T5876535			
5.0	40			T5876540
	35			T5876545
5.5	50			T5876550
	30			T5877030
6.0	40			T5877040
	45			T5877045
	50			T5877050
6.5	30			T5834530
	35			T5834535
	30			T5835030
5.5	35			T5835035
	30			T5835530
6.0	35			T5835535
	40	T5835540		
	45	T5835545		
6.5	30	T5836030		
	35	T5836035		
	40	T5836040		
7.0	45	T5836045		
	30	T5836530		
	35	T5836535		
7.5	40	T5836540		
	45	T5836545		
8.0	50	T5836550		
	6.0	50-180 (The specification is 10mm)	T6161050- T6161180	
6.0	200-400 (The specification is 20mm)	T6161200- T6161400		
	EF 6.0 Crosslinks	Titanium Alloy	50	T6190050
60			T6190060	
70			T6190070	
80			T6190080	
90			T6190090	