

KTO-Sliding Compression Gamma Intramedullary Nail System(Asian Version)

CE 0197 ISO 13485 FDA



KangLi

KANGLI ORTHOPAEDICS
INSTRUMENT

KTO-Sliding Compression Gamma Intramedullary Nail System(Asian Version)

Character Of The System



- Proximal trapezoid design, better stability.
- 5° lateral offset, minimally invasive greater trochanter approach.
- Proximal diameter is small and more suitable for Asians.
- 127° is more suitable for Asians

- Integrated interlocking Screw: anti-rotation, anti-cutting, compression.

- Holding type intraoperative anti-off screw
- Screw speed of double thread is faster

Distal locking screw plan




Plan 1



Plan 2

KTO-Sliding Compression Gamma Intramedullary Nail System(Asian Version)

Character Of The System



KTO-Sliding Compression Gamma Intramedullary Nail System(Asian Version)


Character Of The System



Extended

Femoral model

KTO-Sliding Compression Gamma Intramedullary Nail System(Asian Version)



KTO-Sliding Compression Gamma Intramedullary Nail System(Asian Version)

KTO®-Sliding Compression Gamma Nail Small(Extended edition)
63037、63038、63039



Product Name	Material	Diameter	Length	Product No.
KTO®-Sliding Compression Gamma Nail	Titanium Alloy	Φ 9.4	300L	63037301
			300R	63037302
			320L	63037321
			320R	63037322
			340L	63037341
			340R	63037342
			360L	63037361
			360R	63037362
			380L	63037381
			380R	63037382
			400L	63037401
			400R	63037402
		Φ 10	320L	63038321
			320R	63038322
			340L	63038341
			340R	63038342
			360L	63038361
			360R	63038362
			380L	63038381
			380R	63038382
			400L	63038401
			400R	63038402
			420L	63038421
			420R	63038422
		Φ 11	320L	63039321
			320R	63039322
			340L	63039341
			340R	63039342
			360L	63039361
			360R	63039362
			380L	63039381
			380R	63039382
			400L	63039401
			400R	63039402
			420L	63039421
			420R	63039422

KTO-Sliding Compression Gamma Intramedullary Nail System(Asian Version)

**KTO®-Sliding Compression Gamma Nail Small
63035**



Product Name	Material	Diameter	Length	Product No.
KTO®-Sliding Compression Gamma Nail Small	Titanium Alloy	Φ 9.4	180、200、220	63035118- 63035122
		Φ 10	180、200、220	63035218- 63035222
		Φ 11	180、200、220	63035318- 63035322
		Φ 12	180、200、220	63035418- 63035422

**KTO®-Sliding Compression Lag Screw(Parent Screw) Small
63041**



Product Name	Material	Diameter	Length	Product No.
KTO®-Sliding Compression Lag Screw (Parent Screw)	Titanium Alloy	Φ 9.0	70、75、80、85、 90、95、100、105、 110、115、120	63041070- 63041120

**KTO®-Sliding Compression Lag Screw(Sub Screw) Small
63042**



Product Name	Material	Diameter	Length	Product No.
KTO®-Sliding Compression Lag Screw (Sub Screw)	Titanium Alloy	Φ 6.0	65、70、75、80、 85、90、95、100、 105、110、115	63042065- 63042115

**Φ4.5 Holding Type Full Threaded Locking Screw
63009**




Product Name	Material	Diameter	Length	Product No.
Φ4.5 Holding Type Thread Locking Screw	Titanium Alloy	Φ 4.5	30	63009030
			32	63009032
			35	63009035
			38	63009038
			40、45、50、55、 60、65、70、75、 80、85	63009040- 63009085

**Locking End Cap
63043**



Product Name	Material	Specification	Product No.
Locking End Cap	Titanium Alloy	0	63043000
		5	63043005
		10	63043010
		15	63043015

**Locking End Cap
63306**



Product Name	Material	Length	Product No.
End Cap	Titanium Alloy	0	63306100
		5	63306105
		10	63306110
		15	63306115