

KTO PFNA Gamma Intramedullary Nail System

CE 0197 ISO 13485 FDA




KangLi

KANGLI ORTHOPAEDICS
INSTRUMENT

KTO PFNA Gamma Intramedullary Nail System



KTO[®]-PFNA Gamma Intramedullary Nail System



KTO[®]-PFNA Gamma Intramedullary Nail System (Extended Edition)

KTO PFNA Gamma Intramedullary Nail System

KTO® -PFNA Mail Nail
T8794、T8710、T8711、T8712



Product Name	Material	Diameter	Length	Product No.
KTO®-PFNA Mail Nail	Titanium Alloy	Φ9.4	180、200、220、240	T8794180-T8794240
		Φ10	180、200、220、240	T8710180-T8710240
		Φ11	180、200、220、240	T8711180-T8711240
		Φ12	180、200、220、240	T8712180-T8712240

KTO® -PFNA Mail Nail (Extended Edition)
T8794、T8710、T8711



Product Name	Material	Diameter	Length	Product No.
KTO®-PFNA Mail Nail (Extended Edition)	Titanium Alloy	Φ9.4	320、340、360、380、400	T8794320R/L-T8794400R/L
		Φ10	320、340、360、380、400、420	T8710320R/L-T8710420R/L
		Φ11	320、340、360、380、400	T8711320R/L-T8711400R/L

KTO® -PFNA Blade Nail
T8541



Product Name	Material	Diameter	Length	Product No.
KTO®-PFNA Blade Nail	Titanium Alloy	Φ10.5	70、75、80、85、90、95、100、105、110、115、120	T8541070-T8541120


Φ4.5 Full Threaded Locking Screw
T852



Product Name	Material	Diameter	Length	Product No.
Full Threaded Locking Screw	Titanium Alloy	Φ4.5	30	T8524030
			32	T8524032
			35	T8524035
			38	T8524038
			40、45、50、55、60、65、70、75、80	T8524040-T8524080

Note: Choose a SW3.5 Hollow Hexagonal Screwdriver

Locking End Cap
T8571



Product Name	Material	Specification	Product No.
Locking End Cap	Titanium Alloy	0	T8571000
		5	T8571005
		10	T8571010
		15	T8571015

Normal End Cap
T858



Product Name	Material	Specification	Product No.
Normal End Cap	Titanium Alloy	0	T8580000
		5	T8580005
		10	T8580010
		15	T8580015