

康力·脊柱·创伤产品目录

Kangli Spinal Trauma Products Catalogue

2025 Edition

INTRODUCTION

康力医疗（以下简称康力）成立于2002年，是专注骨科植入物及手术器械研发制造的国家高新技术企业。公司投资资金3.6亿元，注册资金1.3亿元，并拥有4万平方米智能化生产基地。拥有FDA、ISO 13485、CE等国际权威认证，公司获评“江苏省示范智能车间”、“江苏省企业研究生工作站单位”、“江苏省工程技术研究中心”等荣誉，产品质量达到全球主流市场标准。

公司以创新驱动发展，配备德国EOS激光烧结平台、医学影像处理系统等先进设备，与上海交通大学、南京航空航天大学建立产学研合作，累计获授权专利100余项。针对中国骨科临床需求，已形成创伤、脊柱、关节、颅颌面、足踝、胸外及小儿骨科七大产品线，覆盖骨科全领域解决方案，并布局神经外科、胸外科器械领域。

在中国骨科器械国产替代加速的背景下，康力构建了覆盖全国25个省级行政区的营销网络，产品同步出口欧美、东南亚及非洲市场。依托智能化生产体系和成本控制能力，公司可为医疗机构提供高性价比的骨科器械选择，同时支持定制化临床需求开发。

作为民族医疗企业代表，康力始终践行“诚信、务实、服务创造价值”的理念，依托本土化服务优势，已为数百家医院提供骨科手术技术支持及3D打印服务。我们诚邀各级医疗机构、专家莅临考察，共同推进骨科器械技术创新与临床应用，助力中国骨科医疗事业发展。

KangliOrthopaedics (Kangli), founded in 2002, is a leading high-tech innovator specializing in orthopedic implants and surgical solutions. With \$51.4 million in total investment and \$18.6 million in registered capital, our sprawling 40,000-square-meter manufacturing facility operates under stringent FDA, ISO13485, and CE certifications. We integrate German EOS laser sintering systems with advanced medical imaging technologies, accumulated 100+ patents and established strategic R&D partnerships with institutions such as Shanghai Jiao Tong University.

Kangli delivers comprehensive solutions across trauma, spine, craniomaxillofacial, foot/ankle, pediatric orthopedics, and thoracic surgery, including specialized xenograft bone implants and minimally invasive balloon systems. Vertically integrated smart manufacturing enables the cost-effective production of standardized and customized devices to meet diverse clinical needs.

Through an extensive network of distributors across China, Europe, the Americas, Southeast Asia, and Africa, Kangli combines localized technical support with responsive OEM/ODM services and clinical co-development programs. Trusted by over 500+ hospitals worldwide, we balance precision engineering with competitive accessibility, ensuring reliable solutions for enhanced surgical outcomes. Our distributor-centric model prioritizes agile logistics, regulatory compliance, and collaborative growth opportunities in evolving global markets.





INTRODUCTION



德玛吉加工中心群组
DMG Machine Tool Center



纵切车削中心群组
Turning and Milling Center



电子动静态万能材料试验机
INSTRON: Dynamic Simulation Experiment Machine




化学实验室
Chemistry Lab



材料表面检测中心
Material Surface Testing Center



EOS M290植入物3D打印
EOS M290 Selective Laser Melting (SLS) Equipment



尼康工业CT测试仪
Nikon industry CT testing device



火花直读光谱仪
氧氮氢气体分析仪
Direct-reading Spectrometer
Oxygen, Nitrogen & Hydrogen Analyzer

IDEA OF KANGLI

Mission


关爱生命、呵护健康！
Caring for life and health!

Vision

创民族一流企业！
Create a national first-class enterprise!

Values

诚信、务实、服务创造价值！
Integrity, pragmatism, and service create value.



CE认证
CE Certification



祝贺康力公司零483通过美国FDA审核
 Congratulations To KangLi Passed(USA)FDA
 Inspection with ZERO 483

CONTENTS

目 录



A

脊柱内固定系列
Spinal Fixation Series

01-23



B

颅颌面固定系列
Craniomaxillofacial Fixation Series

24-34



C

髓内钉系列
Intramedullary Nail Series

35-64



D

骨针、钛缆、髌骨、空心钉
Skeletal Pins, Titanium Cable, Patella, Cannulated Screw 65-75



E

微型锁定指掌骨系列
Micro Phalangeal Metacarpal Locking Series

76-82



F

JY、WX 锁定钢板螺钉系列
JY、WX Locking Plate Screw Series

83-105



G

(本色) 锁定钢板螺钉系列
(Base Color) Locking Plate Series

106-120


H


肋骨钢板系列
Rib Plate Series

121-124


I

普通钢板螺钉系列
Common Plate Screw Series

125-131


J

外固定支架系列
External Fixation Series

132-138

A

B

C

D

E

F

G

H

I

J


A



脊柱内固定系列 Spinal Fixation Series

脊柱（右侧视图）

Vertebral (Spinal) column (seen from right side)



JZWA 5.5 脊柱微创内固定系统

JZWA 5.5 Spinal Fixation System (Minimally Invasive)

A



视频在线播放
Video online play

预弯棒更符合人体正常生理曲度。

Pre-bending rod conforms to the normal physiological curvature of the human body.




双线和四线螺纹

- 螺钉拧入更快
- 一致的螺钉拧入速度
- 椎弓根皮质固定, 增加螺钉强度

Double thread and quadruple thread design

- Faster screw insertion speed
- Consistent screw insertion speed
- Pedicle cortex fixation, increase screw strength



JZWA 微创 U 型旋转椎弓根钉允许钉身沿螺钉纵轴转动。


JZWA minimally invasive U-rotating pedicle screw allows the screw body to rotate along the screw's longitudinal axis.

JZWA 脊柱微创内固定系统安装器械包 93519901

JZWA Minimally Invasive Spinal Fixation System Instrument Set




PDF文件
PDF File



JZWA 脊柱微创内植入物箱 706013

JZWA Minimally Invasive Spinal Implant Box



JZWA 5.5 脊柱微创内固定系统

JZWA 5.5 Spinal Fixation System (Minimally Invasive)

JZWA 微创U型万向椎弓根钉 43703
JZWA Minimally Invasive U Multiaxial Pedicle Screw




产品名称 Product Name	材质 Material	直径 Diameter	长度 Length	产品编号 Product No.
JZWA 微创U型万向椎弓根钉 JZWA Minimally Invasive U Multiaxial Pedicle Screw	钛合金 Titanium Alloy	5.5	35	43703235
			40	43703240
			45	43703245
		6.0	35	43703335
			40	43703340
			45	43703345
		6.5	50	43703350
			35	43703435
			40	43703440
		7.0	45	43703445
			50	43703450
			35	43703535
		7.5	40	43703540
			45	43703545
			50	43703550
			35	43703635
			40	43703640

JZWA 微创U型旋转椎弓根钉 43704
JZWA Minimally Invasive U Rotating Pedicle Screw



产品名称 Product Name	材质 Material	直径 Diameter	长度 Length	产品编号 Product No.
JZWA 微创U型旋转椎弓根钉 JZWA Minimally Invasive U Rotating Pedicle Screw	钛合金 Titanium Alloy	5.5	35	43704235
			40	43704240
			45	43704245
		6.0	35	43704335
			40	43704340
			45	43704345
		6.5	50	43704350
			35	43704435
			40	43704440
		7.0	45	43704445
			50	43704450
			35	43704535
		7.5	40	43704540
			45	43704545
			50	43704550
			35	43704635
			40	43704640

JZWA 5.5 弧形微创棒 43703C
JZWA 5.5 Arc Shaped Minimally Invasive Rod



产品名称 Product Name	材质 Material	总长 Total Length	有效长度(mm) Effective Length	产品编号 Product No.
JZWA 5.5 弧形微创棒 JZWA 5.5 Arc Shaped Minimally Invasive Rod	钛合金 Titanium Alloy	45	35	43703045
		50	40	43703050
		55	45	43703055
		60	50	43703060
		65	55	43703065
		70	60	43703070
		75	65	43703075
		80	70	43703080
		85	75	43703085
		90	80	43703090
		95	85	43703095
		100	90	43703100
		105	95	43703105
		110	100	43703110
		115	105	43703115
		125	115	43703125
		135	125	43703135
		145	135	43703145
		155	145	43703155

JZWA 5.5 直形微创棒 43706
JZWA 5.5 Straight Minimally Invasive Rod



产品名称 Product Name	材质 Material	总长 Total Length	有效长度(mm) Effective Length	产品编号 Product No.
JZWA 5.5 直形微创棒 JZWA 5.5 Straight Minimally Invasive Rod	钛合金 Titanium Alloy	45	35	43706045
		50	40	43706050
		55	45	43706055
		60	50	43706060
		65	55	43706065
		70	60	43706070
		75	65	43706075
		80	70	43706080
		85	75	43706085
		90	80	43706090
		95	85	43706095
		100	90	43706100
		105	95	43706105
		110	100	43706110
		115	105	43706115
		125	115	43706125
		135	125	43706135
		145	135	43706145
		155	145	43706155

JZPA 5.5 脊柱内固定系统

JZPA 5.5 Spinal Fixation System

A



视频在线播放
Video online play



双线和四线螺纹

- 螺钉拧入更快
- 一致的螺钉拧入速度
- 椎弓根皮质固定，增加螺钉强度

Double thread and quadruple thread design

- Faster screw insertion speed
- Consistent screw insertion speed
- Pedicle cortex fixation, increase screw strength



JZPA 5.5 U 型旋转折断椎弓根钉允许钉身沿螺钉纵轴转动

JZPA 5.5 monoaxial reduction pedicle screw (rotating) allows the screw body to rotate along the longitudinal axis of the screw

JZPA 5.5/6.0 脊柱内固定系统安装器械包 93359901

JZPA 5.5/6.0 Spinal Fixation System Instrument Set



PDF文件
PDF File



钉棒内固定盒 (空箱、钉架) 706003

Screw Rod Implant Box (Empty box, screw frame)



JZPA 5.5 脊柱内固定系统

JZPA 5.5 Spinal Fixation System

JZPA 5.5 U型万向折断椎弓根钉 43010

JZPA 5.5 U Multiaxial Break-Off Pedicle Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
JZPA 5.5 U型万向折断椎弓根钉 JZPA 5.5 U Multiaxial Break-Off Pedicle Screw	钛合金 Titanium Alloy	4.5	30	43010130
			35	43010135
		5.0	30	43010230
			35	43010235
		5.5	35	43010335
			40	43010340
		6.0	30	43010430
			35	43010435
			40	43010440
			45	43010445
		6.5	50	43010450
			30	43010630
			35	43010635
			40	43010640
		7.0	45	43010645
			50	43010650
			35	43010735
			40	43010740
			45	43010745
			50	43010750

JZPA 5.5 U型旋转折断椎弓根钉 43011

JZPA 5.5 U Rotating Break-Off Pedicle Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
JZPA 5.5 U型旋转折断椎弓根钉 JZPA 5.5 U Rotating Break-Off Pedicle Screw	钛合金 Titanium Alloy	4.5	30	43011130
			35	43011135
		5.0	30	43011230
			35	43011235
		5.5	35	43011335
			40	43011340
		6.0	35	43011435
			40	43011440
			45	43011445
			50	43011450
		6.5	35	43011535
			40	43011540
			45	43011545
			50	43011550
		7.0	35	43011635
			40	43011640
			45	43011645
			50	43011650

JZPA 5.5 连接棒 43020

JZPA 5.5 Rod



JZPA EF 5.5 横连器 43030

JZPA EF 5.5 Crosslinks



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
JZPA 5.5 连接棒 JZPA 5.5 Rod	钛合金 Titanium Alloy	5.5	50	43020050
			60	43020060
			70	43020070
			80	43020080
			90	43020090
			100	43020100
			110	43020110
			120	43020120
			140	43020140
			160	43020160
			180	43020180
			200	43020200
			300	43020300
			50	43030050
60	43030060			
70	43030070			
80	43030080			
90	43030090			
JZPA EF 5.5 横连器 JZPA EF 5.5 Crosslinks		5.5		

JZPA 6.0 脊柱内固定系统

JZPA 6.0 Spinal Fixation System



视频在线播放
Video online play




双线和四线螺纹

- 螺钉拧入更快
- 一致的螺钉拧入速度
- 椎弓根皮质固定, 增加螺钉强度

Double thread and quadruple thread design

- Faster screw insertion speed
- Consistent screw insertion speed
- Pedicle cortex fixation, increase screw strength



此浅孔为 JZPA 6.0 系统标志孔
This Shallow Hole Is The JZPA 6.0 System Mark Hole

JZPA 6.0 U 型旋转折断椎弓根钉允许钉身沿螺钉纵轴转动


JZPA 6.0 monoaxial reduction pedicle screw (rotating) allows the screw body to rotate along the longitudinal axis of the screw

JZPA 5.5/6.0 脊柱内固定系统安装器械包 93359901

JZPA 5.5/6.0 Spinal Fixation System Instrument Set



PDF文件
PDF File



钉棒内固定盒 (空箱、钉架) 706003

Screw Rod Inplant Box (Empty box, screw frame)



JZPA 6.0 脊柱内固定系统

JZPA 6.0 Spinal Fixation System



A

JZPA 6.0 U型万向折断椎弓根钉 43112
JZPA 6.0 U Multiaxial Break-Off Pedicle Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
JZPA 6.0 U型万向折断椎弓根钉 JZPA 6.0 U Multiaxial Break-Off Pedicle Screw	钛合金 Titanium Alloy	4.5	30	43112130
			35	43112135
		5.0	30	43112230
			35	43112235
		5.5	35	43112335
			40	43112340
		6.0	35	43112435
			40	43112440
		6.5	45	43112445
			50	43112450
		7.0	35	43112535
			40	43112540
		7.5	45	43112545
			50	43112550
		8.0	35	43112635
			40	43112640
		8.5	45	43112645
			50	43112650

JZPA 6.0 U型旋转折断椎弓根钉 43113
JZPA 6.0 U Rotating Break-Off Pedicle Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
JZPA 6.0 U型旋转折断椎弓根钉 JZPA 6.0 U Rotating Break-Off Pedicle Screw	钛合金 Titanium Alloy	4.5	30	43113130
			35	43113135
		5.0	30	43113230
			35	43113235
		5.5	35	43113335
			40	43113340
		6.0	35	43113435
			40	43113440
		6.5	45	43113445
			50	43113450
		7.0	40	43113540
			45	43113545
		7.5	50	43113550
			45	43113645
		8.0	50	43113650
			50	43113650

JZPA 6.0 连接棒 43121
JZPA 6.0 Rod



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
JZPA 6.0 连接棒 JZPA 6.0 Rod	钛合金 Titanium Alloy	6.0	50	43121050
			60	43121060
			70	43121070
			80	43121080
			90	43121090
			100	43121100
			110	43121110
			120	43121120
			140	43121140
			160	43121160
			180	43121180
			200	43121200
			300	43121300
			50	43122050
JZPA EF 6.0 横连器 JZPA EF 6.0 Crosslinks	钛合金 Titanium Alloy	6.0	60	43122060
			70	43122070
			80	43122080
			90	43122090

JZPA EF 6.0 横连器 43122
JZPA EF 6.0 Crosslinks



DTL 5.5 双螺纹低切迹脊柱系统



DTL 5.5 Double Thread Low-profile Spinal System

A

DTL 5.5 U型折断万向椎弓根钉 4515
DTL 5.5 U Multiaxial Break-off Pedicle Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
DTL 5.5 U型折断万向椎弓根钉 DTL 5.5 U Multiaxial Break-off Pedicle Screw	钛合金 Titanium Alloy	4.0	25	45154025
			30	45154030
		4.5	30	45154530
			35	45154535
		5.0	30	45155030
			35	45155035
		5.5	35	45155535
			40	45155540
		6.0	30	45156030
			35	45156035
			40	45156040
			45	45156045
		6.5	50	45156050
			35	45156535
			40	45156540
		7.0	45	45156545
			50	45156550
		7.5	45	45157045
			50	45157050
		7.5	45	45157545
			50	45157550

DTL 5.5 U型折断旋转椎弓根钉 4516
DTL 5.5 U Rotating Break-off Pedicle Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
DTL 5.5 U型折断旋转椎弓根钉 DTL 5.5 U Rotating Break-off Pedicle Screw	钛合金 Titanium Alloy	4.0	25	45164025
			30	45164030
		4.5	30	45164530
			35	45164535
		5.0	30	45165030
			35	45165035
		5.5	35	45165535
			40	45165540
		6.0	35	45166035
			40	45166040
			45	45166045
			50	45166050
		6.5	35	45166535
			40	45166540
			45	45166545
		7.0	50	45166540
			45	45167045
		7.5	50	45167050
			45	45167545
			50	45167550

JZPA 5.5 连接棒 43020
JZPA 5.5 Rod



JZPA EF 5.5 横连器 43030
JZPA EF 5.5 Crosslinks




产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
JZPA 5.5 连接棒 JZPA 5.5 Rod	钛合金 Titanium Alloy	5.5	50	43020050
			60	43020060
			70	43020070
			80	43020080
			90	43020090
			100	43020100
			110	43020110
			120	43020120
			140	43020140
			160	43020160
JZPA EF 5.5 横连器 JZPA EF 5.5 Crosslinks	钛合金 Titanium Alloy	5.5	180	43020180
			200	43020200
			300	43020300
			50	43030050
			60	43030060
			70	43030070
			80	43030080
			90	43030090

JZPA 5.5/6.0 脊柱内固定系统安装器械包 93359901
JZPA 5.5/6.0 Spinal Fixation System Instrument Set



PDF文件
PDF File



DTL 6.0 双螺纹低切迹脊柱系统

DTL 6.0 Double Thread Low-profile Spinal System

DTL 6.0 U型折断万向椎弓根钉 4333
DTL 6.0 U Multiaxial Break-Off Pedicle Screw



DTL 6.0 U型折断旋转椎弓根钉 4334
DTL 6.0 U Rotating Break-off Pedicle Screw



DTL 6.0 U型短尾万向椎弓根钉 4331
DTL 6.0 U Multiaxial Pedicle Screw (Short Tail)



DTL 6.0 U型短尾旋转椎弓根钉 4332
DTL 6.0 U Rotating Pedicle Screw (Short Tail)




产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
DTL 6.0 U型折断万向椎弓根钉 DTL 6.0 U Multiaxial Break-Off Pedicle Screw	钛合金 Titanium Alloy	4.0	25	43334025
			30	43334030
		4.5	30	43334530
			35	43334535
		5.0	30	43335030
			35	43335035
		5.5	35	43335535
			40	43335540
		6.0	30	43336030
			35	43336035
			40	43336040
			45	43336045
			50	43336050
		6.5	35	43336535
			40	43336540
			45	43336545
DTL 6.0 U型折断旋转椎弓根钉 DTL 6.0 U Rotating Break-off Pedicle Screw	钛合金 Titanium Alloy	4.0	25	43344025
			30	43344030
		4.5	30	43344530
			35	43344535
		5.0	30	43345030
			35	43345035
		5.5	35	43345535
			40	43345540
		6.0	35	43346035
			40	43346040
			45	43346045
			50	43346050
		6.5	35	43346535
			40	43346540
			45	43346545
			50	43346550
DTL 6.0 U型短尾万向椎弓根钉 DTL 6.0 U Multiaxial Pedicle Screw (Short Tail)	钛合金 Titanium Alloy	4.0	25	43314025
			30	43314030
		4.5	30	43314530
			35	43314535
		5.0	30	43315030
			35	43315035
		5.5	35	43315535
			40	43315540
		6.0	30	43316030
			35	43316035
			40	43316040
			45	43316045
			50	43316050
		6.5	35	43316535
			40	43316540
			45	43316545
DTL 6.0 U型短尾旋转椎弓根钉 DTL 6.0 U Rotating Pedicle Screw (Short Tail)	钛合金 Titanium Alloy	4.0	25	43324025
			30	43324030
		4.5	30	43324530
			35	43324535
		5.0	30	43325030
			35	43325035
		5.5	35	43325535
			40	43325540
		6.0	35	43326035
			40	43326040
			45	43326045
			50	43326050
		6.5	35	43326535
			40	43326540
			45	43326545
			50	43326550

JZKA 6.0 脊柱内固定系统

JZKA 6.0 Spinal Fixation System

A



JZKA 6.0 脊柱内固定系统

JZKA 6.0 Spinal Fixation System

JZKA 6.0 U型万向折断椎弓根钉 (骨质加强型) 43206
JZKA 6.0 U Multiaxial Break-Off Pedicle Screw
(Bone enhancement type)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
JZKA 6.0 U型万向折断椎弓根钉 (骨质加强型) JZKA 6.0 U Multiaxial Break-Off Pedicle Screw (Bone enhancement type)	钛合金 Titanium Alloy	5.5	35	43206035
			40	43206040
		6.0	35	43206135
			40	43206140
			45	43206145
		6.5	50	43206150
			35	43206235
			40	43206240
		7.0	45	43206245
			50	43206250
			35	43206335
			40	43206340
			45	43206345
			50	43206350

JZPA 6.0 连接棒 43121
JZPA 6.0 Rod



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
JZPA 6.0 连接棒 JZPA 6.0 Rod	钛合金 Titanium Alloy	6.0	50	43121050
			60	43121060
			70	43121070
			80	43121080
			90	43121090
			100	43121100
			110	43121110
			120	43121120
			140	43121140
			160	43121160
			180	43121180
			200	43121200
			300	43121300

JZPA EF 6.0 横连器 43122
JZPA EF 6.0 Crosslinks




产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
JZPA EF 6.0 横连器 JZPA EF 6.0 Crosslinks	钛合金 Titanium Alloy	6.0	50	43122050
			60	43122060
			70	43122070
			80	43122080
			90	43122090

JZPA 5.5/6.0 脊柱内固定系统安装器械包 93359901
JZPA 5.5/6.0 Spinal Fixation System Instrument Set



PDF文件
PDF File



钉棒内固定盒 (空箱、钉架) 706003

Screw Rod Implant Box (Empty box, screw frame)



KSS-I 6.0 脊柱内固定系统（梅花型）

KSS-I 6.0 Spinal Fixation System (Torx)

A




视频在线播放
Video online play




螺钉头部采用了梯形螺纹设计，将螺钉头的外张力降到最低，抗脱出能力最强，锁紧能力最强。

The Top Of The Screw Adopt The Design Of The Trapezoid Type, To Lower The External Tension Minimally, It Is The Best That The Ability Of Anti-prolapsus And Locking.



最新设计，上开口弹性固定操作极为方便，固定牢固。

Up-to-date. Design, That Is Easily Operated For The Bayonet Type And Stretch Hatch.




独特的梅花顶丝设计可有效防止顶丝手术中打滑现象。

Special Design Of The Torx Nut Can Avoid Slipping During The Surgery.



松质骨螺纹设计，最适合椎体。增大的底径减少折断风险。

The Design For The Cancellous Bone, That Is Most Suitable For The Centrum. The Increment Diameter Of The Bottom Can Be Less The Risk Of Break.




JZPA 5.5/6.0 脊柱内固定系统安装器械包 93359901

JZPA 5.5/6.0 Spinal Fixation System Instrument Set





PDF文件
PDF File



钉棒内固定盒（空箱、钉架）706003

Screw Rod Implant Box (Empty box, screw frame)



KSS-I 6.0 脊柱内固定系统（梅花型）

KSS-I 6.0 Spinal Fixation System (Torx)

A

KSS 6.0 U 型椎弓根钉 T581

KSS 6.0 U Pedicle Screw



KSS 6.0 U 型可折断椎弓根钉 T588

KSS 6.0 U Break-Off Pedicle Screw



KSS 6.0 U 型万向可折断椎弓根钉（六角万向头）T587

KSS 6.0 U Multiaxial Break-Off Pedicle Screw (Hexagon Type)



KSS 6.0 U 型万向可折断椎弓根钉（梅花万向头）T583

KSS 6.0 U Multiaxial Break-Off Pedicle Screw (Torx Type)



连接棒 (6.0) T6161

Rod 6.0



EF 6.0 横连器 T6190

EF 6.0 Crosslinks



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KSS 6.0 U型 椎弓根钉 KSS 6.0 U Pedicle Screw	钛合金 Titanium Alloy	4.5	30	T5814530
			35	T5814535
		5.0	30	T5815030
			35	T5815035
		5.5	35	T5815535
			40	T5815540
		6.0	35	T5816035
			40	T5816040
			45	T5816045
			35	T5816535
		6.5	40	T5816540
			45	T5816545
			50	T5816550
		7.0	35	T5817035
			40	T5817040
KSS 6.0 U型 可折断椎弓根钉 KSS 6.0 U Break-Off Pedicle Screw	钛合金 Titanium Alloy	4.5	35	T5884535
			30	T5885030
		5.0	35	T5885035
			35	T5885535
		5.5	40	T5885540
			45	T5885545
		6.0	35	T5886035
			40	T5886040
			45	T5886045
			50	T5886050
		6.5	35	T5886535
			40	T5886540
			45	T5886545
		7.0	50	T5886550
			35	T5887035
KSS 6.0 U型万向 可折断椎弓根钉 (六角万向头) KSS 6.0 KSS 6.0 U Multiaxial Break-Off Pedicle Screw (Hexagon Type)	钛合金 Titanium Alloy	4.5	30	T5874530
			35	T5874535
		5.0	30	T5875030
			35	T5875035
		5.5	40	T5875040
			30	T5875530
		6.0	35	T5875535
			40	T5875540
		6.5	45	T5875545
			30	T5876030
			35	T5876035
		7.0	40	T5876040
			45	T5876045
KSS 6.0 U型万向 可折断椎弓根钉 (梅花万向头) KSS 6.0 U Multiaxial Break-Off Pedicle Screw (Torx Type)	钛合金 Titanium Alloy	4.5	30	T5834530
			35	T5834535
		5.0	30	T5835030
			35	T5835035
		5.5	30	T5835530
			35	T5835535
		6.0	40	T5835540
			45	T5835545
			30	T5836030
			35	T5836035
		6.5	40	T5836040
			45	T5836045
			30	T5836530
连接棒 Rod	钛合金 Titanium Alloy	4.5	30	T5834530
			35	T5834535
		5.0	30	T5835030
			35	T5835035
		5.5	30	T5835530
			35	T5835535
		6.0	40	T5835540
			45	T5835545
			30	T5836030
			35	T5836035
		6.5	40	T5836040
			45	T5836045
			30	T5836530
EF 6.0 横连器 EF 6.0 Crosslinks	钛合金 Titanium Alloy	6.0	50-180 (每10mm为一档) The specification is 10mm)	T6161050- T6161180
			200-400 (每20mm为一档) The specification is 20mm)	T6161200- T6161400
EF 6.0 横连器 EF 6.0 Crosslinks	钛合金 Titanium Alloy	6.0	50	T6190050
			60	T6190060
			70	T6190070
			80	T6190080
			90	T6190090

KSS-II 5.5 脊柱内固定系统（后路固定模式）

KSS-II 5.5 Spinal Fixation System (Posterior Fixation)

A



KSS 5.5 U 型椎弓根钉 T621
KSS 5.5 U Pedicle Screw



KSS 5.5 U 型可折断椎弓根钉 T622
KSS 5.5 U Break-Off Pedicle Screw



KSS 5.5 U 型万向可折断椎弓根钉 T623
KSS 5.5 U Multiaxial Break-Off Pedicle Screw



连接棒(5.5) T6241
Rod 5.5




EF 5.5 横连器 T6290
EF 5.5 Crosslinks



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KSS 5.5 U型 椎弓根钉 KSS 5.5 U Pedicle Screw	钛合金 Titanium Alloy	4.0	25	T6214025
			30	T6214030
		4.5	30	T6214530
			35	T6214535
		5.0	30	T6215030
			35	T6215035
			40	T6215040
		5.5	30	T6215530
			35	T6215535
			40	T6215540
		6.0	35	T6216035
			40	T6216040
			45	T6216045
			50	T6216050
		6.5	35	T6216535
			40	T6216540
			45	T6216545
			50	T6216550
		7.0	55	T6216555
			35	T6217035
			40	T6217040
KSS 5.5 U型 可折断椎弓根钉 KSS 5.5 U Break-Off Pedicle Screw	钛合金 Titanium Alloy	4.0	25	T6224025
			30	T6224030
		4.5	30	T6224530
			35	T6224535
		5.0	30	T6225030
			35	T6225035
			35	T6225535
		5.5	40	T6225540
			45	T6225545
			35	T6226035
		6.0	40	T6226040
			45	T6226045
			50	T6226050
			35	T6226535
		6.5	40	T6226540
			45	T6226545
			50	T6226550
			35	T6227035
		7.0	40	T6227040
			35	T6227035
KSS 5.5 U型万向 可折断椎弓根钉 KSS 5.5 U Multiaxial Break-Off Pedicle Screw	钛合金 Titanium Alloy	4.0	25	T6234025
			30	T6234030
		4.5	30	T6234530
			35	T6234535
		5.0	30	T6235030
			35	T6235035
			40	T6235040
		5.5	30	T6235530
			35	T6235535
			40	T6235540
			30	T6236030
		6.0	35	T6236035
			40	T6236040
			45	T6236045
			35	T6236535
		6.5	40	T6236540
			45	T6236545
			50	T6236550
			50	T6236550
连接棒 Rod	钛合金 Titanium Alloy	5.5	60、70、80、90、100、110、120	T6241060-T6241120
			140、160、180、200	T6241140-T6241200
EF 5.5 横连器 EF 5.5 Crosslinks	钛合金 Titanium Alloy	5.5	50	T6290050
			60	T6290060
			70	T6290070
			80	T6290080

KSS-II 5.5（前路固定模式）

KSS-II 5.5 (Anterior Fixation)



KSS 5.5 U 型椎弓根钉 T621
KSS 5.5 U Pedicle Screw



连接棒(5.5) T6241
Rod 5.5



双棒垫片 T627
Double Rods Shim



单棒垫片 T628
Single Rod Shim



前路横连器(5.5) T626
Anterior Crosslink 5.5



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KSS 5.5 U型 椎弓根钉 KSS 5.5 U Pedicle Screw	钛合金 Titanium Alloy	4.0	25	T6214025
			30	T6214030
		4.5	30	T6214530
			35	T6214535
		5.0	30	T6215030
			35	T6215035
			40	T6215040
		5.5	30	T6215530
			35	T6215535
			40	T6215540
		6.0	35	T6216035
			40	T6216040
			45	T6216045
			50	T6216050
		6.5	35	T6216535
			40	T6216540
			45	T6216545
			50	T6216550
连接棒 Rod	钛合金 Titanium Alloy	5.5	60、70、80、90 100、110、120	T6241060- T6241120
			140、160、180、 200	T6241140- T6241200
双棒垫片 Double Rods Shim	钛合金 Titanium Alloy	头向（小/中/大） Cephalad (Small/Mid/Large)		T6270001
				T6270002
				T6270003
		尾向（小/中/大） Caudad (Small/Mid/Large)		T6270011
				T6270012
单棒垫片 Single Rod Shim				T6270013
前路横连器 Anterior Crosslink	钛合金 Titanium Alloy	5.5	13	T6260013
			14	T6260014
			15	T6260015
			16	T6260016
			17	T6260017
			18	T6260018

KSS-II 5.5mm 脊柱固定器安装器械包 93339901
KSS-II 5.5 Spinal Fixation Instrument Case



PDF 文件
PDF File



JZP 6.0 脊柱内固定系统

JZP 6.0 Spinal Fixation System

A

JZP 6.0 U型万向折断椎弓根钉 42225
JZP 6.0 U Multiaxial Break-Off Pedicle Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
JZP 6.0 U型万向折断椎弓根钉 JZP 6.0 U Multiaxial Break-Off Pedicle Screw	钛合金 Titanium Alloy	4.5	30	42225130
			35	42225135
		5.0	30	42225230
			35	42225235
		5.5	35	42225335
			40	42225340
		6.0	30	42225430
			35	42225435
			40	42225440
			45	42225445
			50	42225450
		6.5	35	42225535
			40	42225540
			45	42225545
			50	42225550
		7.0	45	42225645
			50	42225650

JZP 6.0 U型旋转折断椎弓根钉 42226
JZP 6.0 U Rotating Break-Off Pedicle Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
JZP 6.0 U型旋转折断椎弓根钉 JZP 6.0 U Rotating Break-Off Pedicle Screw	钛合金 Titanium Alloy	4.5	30	42226130
			35	42226135
		5.0	30	42226230
			35	42226235
		5.5	35	42226335
			40	42226340
		6.0	40	42226440
			45	42226445
			35	42226535
		6.5	40	42226540
			45	42226545
			50	42226550
		7.0	45	42226645
			50	42226650

连接棒 (6.0) T6161
Rod 6.0



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
连接棒 Rod	钛合金 Titanium Alloy	6.0	50-180 (每10mm为一档) The specification is 10mm	T6161050- T6161180
			200-400 (每20mm为一档) The specification is 20mm	T6161200- T6161400

EF 6.0 横连器 T6190
EF 6.0 Crosslinks




产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
EF 6.0 横连器 EF 6.0 Crosslinks	钛合金 Titanium Alloy	6.0	50	T6190050
			60	T6190060
			70	T6190070
			80	T6190080
			90	T6190090

JZPA 5.5/6.0 脊柱内固定系统安装器械包 93359901
JZPA 5.5/6.0 Spinal Fixation System Instrument Set



PDF文件
PDF File




KLP Posterior Cervical Fixation System

KLP 椎板成型钢板 (防滑脱型) 52213
KLP Laminoplasty Plate (Anti-slipping)



产品名称 Product Name	材质 Material	高度 Height	跨度(mm) Span	产品编号 Product No.
KLP 椎板成型钢板 (防滑脱型) KLP Laminoplasty Plate (Anti-slipping)	纯钛 Pure Titanium	4	8	52213108
			10	52213110
			12	52213112
			14	52213114
		6	16	52213116
			8	52213208
			10	52213210
			12	52213212
			14	52213214
			16	52213216

KLP 椎板成型钢板 (防滑脱T型) 52214
KLP Laminoplasty Plate (Anti-slipping T)



产品名称 Product Name	材质 Material	高度 Height	跨度(mm) Span	产品编号 Product No.
KLP 椎板成型钢板 (防滑脱T型) KLP Laminoplasty Plate (Anti-slipping T)	纯钛 Pure Titanium	4	8	52214108
			10	52214110
			12	52214112
			14	52214114
		6	8	52214208
			10	52214210
			12	52214212
			14	52214214

KLP 椎板成型钢板 (单弯) T6584
KLP Laminoplasty Plate (Single Bending)



产品名称 Product Name	材质 Material	孔数 Holes	跨度(mm) Span	产品编号 Product No.
KLP 椎板成型钢板 (单弯) KLP Laminoplasty Plate (Single Bending)	纯钛 Pure Titanium	5	8	T6584508
			10	T6584510
			12	T6584512
			14	T6584514

KLP 椎板成型钢板 (双弯) T6583
KLP Laminoplasty Plate (Double Bending)



产品名称 Product Name	材质 Material	孔数 Holes	跨度(mm) Span	产品编号 Product No.
KLP 椎板成型钢板 (双弯) KLP Laminoplasty Plate (Double Bending)	纯钛 Pure Titanium	5	8	T6583508
			10	T6583510
			12	T6583512
			14	T6583514

KLP 颈椎螺钉 52305
KLP Cervical Screw



钉帽形式 Head of Nali Type	材质 Material	螺纹直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
十字型 Cross Type	钛合金 Titanium Alloy	2.0	5	52305005
			6	52305006
			8	52305008
			10	52305010
			12	52305012

KLP 颈椎修正螺钉 52306
KLP Cervical Correction Screw




钉帽形式 Head of Nali Type	材质 Material	螺纹直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
十字型 Cross Type	钛合金 Titanium Alloy	2.5	6	52306006
			8	52306008
			10	52306010
			12	52306012

颈椎后路椎板成形系统器械包 93229901
Posterior Cervical Laminar Shaping System Instrument Set



PDF 文件
PDF File






KCO 3.2 钉棒系统

KCO 3.2 Screw Rod System

A



视频在线播放
Video online play



KCO 3.2mm 钉棒器械包 93219901
KCO 3.2 Screw Rod Instrument Set



PDF 文件
PDF File



KCO 3.2 钉棒系统

KCO 3.2 Screw Rod System

KCO U型万向椎弓根钉 T641
KCO Multiaxial Pedicle Screw




产品名称 Product Name	材质 Material	直径×长度 Diameter×Length	产品编号 Product No.	直径×长度 Diameter×Length	产品编号 Product No.
KCO U型万向 椎弓根钉 KCO U-Multi- Axial Pedicle Screw	钛合金 Titanium Alloy	φ3.5×10	T64113510	φ4.0×10	T64114010
		φ3.5×12	T64113512	φ4.0×12	T64114012
		φ3.5×14	T64113514	φ4.0×14	T64114014
		φ3.5×16	T64113516	φ4.0×16	T64114016
		φ3.5×18	T64113518	φ4.0×18	T64114018
		φ3.5×20	T64113520	φ4.0×20	T64114020
		φ3.5×22	T64113522	φ4.0×22	T64114022
		φ3.5×24	T64113524	φ4.0×24	T64114024
		φ3.5×26	T64113526	φ4.0×26	T64114026
		φ3.5×28	T64113528	φ4.0×28	T64114028
		φ3.5×30	T64113530	φ4.0×30	T64114030
		φ3.5×32	T64113532	φ4.0×32	T64114032
			φ4.0×34	T64114034	

KCO 连接棒 T644
KCO Rod



KCO 枕骨连接棒 T643
KCO Occipital Bone Rod



EF 3.2 横连器 T647
EF 3.2 Crosslinks



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KCO 连接棒 KCO Rod	钛合金 Titanium Alloy	3.2	50、60、70、	T64432050-70
			80、90、100	T64432080-100
			200	T64432200
KCO 枕骨连接棒 KCO Occipital Bone Rod			200	T64332200
EF 3.2 横连器 EF 3.2 Crosslinks			40	T6470040
			50	T6470050
			60	T6470060

KCO 椎板钩 T642
KCO Laminar Hook



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
KCO 椎板钩 KCO Laminar Hook	钛合金 Titanium Alloy	4.5	T6423245
		6.0	T6423260

KCO 多米诺固定块 T645
KCO Domino Bolts



KCO 侧方固定块 T646
KCO Lateral Bolts



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
KCO 多米诺固定块 KCO Domino Bolts	钛合金 Titanium Alloy	φ3.2×5.5	T6453255
		φ3.2×6.0	T6453260
KCO 侧方固定块 KCO Lateral Bolts		φ3.2×3.2	T6463232

KCO-S 枕骨板 T6481
KCO-S Occipital Plate



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
KCO-S 枕骨板 KCO-S Occipital Plate	钛合金 Titanium Alloy	6	T6481006

KCO-R 枕骨板 T6481C
KCO-R Occipital Plate



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
KCO-R 枕骨板 KCO-R Occipital Plate	钛合金 Titanium Alloy	5	T6481005

KCO 枕骨螺钉 T649
KCO Occipital Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KCO 枕骨螺钉 KCO Occipital Screw	钛合金 Titanium Alloy	4.0	6	T6494006
			8	T6494008
			10	T6494010
			12	T6494012

Intervertebral Fusion Cage Series

A



适应症

Indications

适用于（C2-C7）节段，颈椎椎间盘退变性疾病。

系统特点

Character Of The System


中线零切迹，降低吞咽困难及食管瘘的发生。两枚螺钉固定，相对传统钉板系统较少的螺钉，可变角度，灵活置钉，一步锁定技术，操作简单，可变角度置钉，头尾端成角27-44°，内外成角15-29°，安全限深板设计，保证产品前缘和椎体缘齐平，避免产品过深进入椎间隙，避免置钉时内植入物位移，能更好的维持颈椎生理曲度，PEEK材料与钛合金结合，提供最佳的固定效果与应变分享，不透射线钽金属标记物，可在成像过程中作为后方结构的标记，PEEK材质融合器，表面的齿状结构提供初始稳定性，锁定螺钉结构，提供安全的锁定机制。融合器弧形外形设计，解剖贴合，利于融合。避免对相邻节段的干扰，预防邻近阶段骨化。

零切迹融合器器械包 93249901

Zero-profile Fusion Cage Instrument Set



PDF 文件
PDF File



Intervertebral Fusion Cage Series

A

颈椎零切迹 椎间融合器（弧形） 58020
Cervical Zero Notch Interbody Fusion Cage (Curve)



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
颈椎零切迹 椎间融合器（弧形） Cervical Zero Notch Interbody Fusion Cage (Curve)	PEEK	5	58020050
		6	58020060
		7	58020070
		8	58020080
		9	58020090
		10	58020100

4.0 自钻自攻型螺钉 58204
4.0 Self Drilling And Tapping Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
4.0 自钻自攻型螺钉 4.0 Self Drilling And Tapping Screw	钛合金 Titanium Alloy	4.0	12	58204012
			14	58204014
			16	58204016

4.0 自攻型螺钉 58214
4.0 Self Tapping Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
4.0 自攻型螺钉 4.0 Self Tapping Screw	钛合金 Titanium Alloy	4.0	12	58214012
			14	58214014
			16	58214016
			18	58214018

颈椎零切迹 椎间融合器（弧形小号） 58021
Cervical Zero Notch Interbody Fusion Cage (Curve Small)



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
颈椎零切迹 椎间融合器（弧形小号） Cervical Zero Notch Interbody Fusion Cage (Curve Small)	PEEK	5	58021050
		6	58021060
		7	58021070
		8	58021080
		9	58021090
		10	58021100

3.5 自钻自攻型螺钉 58203
3.5 Self Drilling And Tapping Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
3.5 自钻自攻型螺钉 3.5 Self Drilling And Tapping Screw	钛合金 Titanium Alloy	3.5	12	58203512
			14	58203514
			16	58203516

3.5 自攻型螺钉 58213
3.5 self tapping screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
3.5 自攻型螺钉 3.5 Self tapping Screw	钛合金 Titanium Alloy	3.5	12	58213512
			14	58213514
			16	58213516

椎间椎体融合器系列

Intervertebral Fusion Cage Series



A

PEEK 无菌 颈椎椎间融合器 57501
PEEK Cervical Fusion Cage (Sterile)



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
PEEK 无菌 颈椎椎间融合器 PEEK Cervical Fusion Cage (Sterile)	PEEK	4	57501040
		5	57501050
		6	57501060
		7	57501070
		8	57501080
		9	57501090

PEEK 无菌 腰椎椎间融合器 57401
PEEK Lumbar Fusion Cage (Sterile)





产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
PEEK 无菌 腰椎椎间融合器 PEEK Lumbar Fusion Cage (Sterile)	PEEK	8×22	57401022
		8×26	57401026
		9×22	57401422
		9×26	57401426
		10×22	57401122
		10×26	57401126
		11×22	57401522
		11×26	57401526
		12×22	57401222
		12×26	57401226
		13×22	57401622
		13×26	57401626
		14×22	57401322
		14×26	57401326

PEEK 颈椎融合器安装器械包 93659901
PEEK Cervical Fusion Cage Instrument Set




PDF 文件
PDF File



Intervertebral Fusion Cage Series

颈椎椎间融合器 T682
Cervical Fusion Cage



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
颈椎椎间融合器 Cervical Fusion Cage	钛合金 Titanium Alloy	15×4.5、5.5、6.5、7.5、8.5	T6821545-T6821585

腰椎椎间融合器 T681
Lumbar Fusion Cage



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
腰椎椎间融合器 Lumbar Fusion Cage	钛合金 Titanium Alloy	8×20、22、26、32、36	T6810820-T6810836
		10×20、22、26、32、36	T6811020-T6811036
		12×20、22、26、32、36	T6811220-T6811236
		14×20、22、26、32、36	T6811420-T6811436

钛笼 T671
Titanium Cage



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
钛笼 Titanium Cage	纯钛 Pure Titanium	10	60	T6711060
		12	60	T6711260
		14	60	T6711460
		16	60	T6711660
		19	60	T6711960
		24	60	T6712460

钛笼盖板 T672
Cover Of Titanium Cage





产品名称 Product Name	材质 Material	直径(mm) Diameter	产品编号 Product No.
钛笼盖板 Cover Of Titanium Cage	纯钛 Pure Titanium	10	T6720010
		12	T6720012
		14	T6720014
		16	T6720016
		19	T6720019
		24	T6720024

颈椎椎间融合器安装器械包 93619902
Cervical Fusion Device Instrument Set




PDF 文件
PDF File



腰椎椎间融合器安装器械包 93629902
Lumbar Fusion Device Instrument Set






PDF 文件
PDF File



钛笼系统安装器械包 93639901
Titanium Mesh Cage Instrument Case



PDF 文件
PDF File




Anterior Cervical Plate-II System




PDF 文件
PDF File


A



低切迹设计，板厚仅1.8mm。
独特的开窗设计，更利于观察
植骨的情况。
Lower Cutting Design, The
Thickness Of Plate Is 1.8mm.
Special Window Design, To Be
Observed Easily.



螺钉与钢板间在横断面上可提
供0°-12°范围内可调节的角度
固定。
Provide The Condition Of
0°-12°in The Coronal Plane
Between The Screw And
Plate.



螺钉与钢板间在矢状面上
可提供0°-10°范围内可调
节的角度固定。
Provide The Condition
Of 0°-10°in The Vertical
Plane Between The
Screw And Plate.



颈椎前路钢板器械包 93239901
Anterior Cervical Plate Instrument Set

颈椎前路钢板 II 型 T652
Anterior Cervical Plate-II



固定螺钉 T654
Fixed Screw




修正螺钉 T655
Correction screw



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
颈椎前路钢板 II 型 Anterior Cervical Plate-II	纯钛 Pure Titanium	4	22	T6520422
			25	T6520425
			27	T6520427
			30	T6520430
		6	32	T6520632
			35	T6520635
			38	T6520638
			41	T6520641
			44	T6520644
			47	T6520647
			50	T6520650
		8	53	T6520853
			56	T6520856
			59	T6520859
			62	T6520862
			65	T6520865
			68	T6520868
			71	T6520871

产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
固定螺钉 Fixed Screw	钛合金 Titanium Alloy	4.0	13	T6544013
			14	T6544014
			15	T6544015
			16	T6544016
			17	T6544017
			18	T6544018
修正螺钉 Correction screw	钛合金 Titanium Alloy	4.5	14	T6554514
			15	T6554515
			16	T6554516
			17	T6554517



颅颌面固定系列

Craniomaxillofacial Fixation Series

B

螺钉可持功能，临床手术更方便
Screw With Holding Function, More Convenient For Surgery

Φ1.5mm


自攻型刃口
Self tapping blade

十字螺钉防打滑
Cross Screw Antislipping



MCF中芯十字起子杆 91130206

MCF Filler Cross Screwdriver Rod, 40mm



MCF Multifunction Craniomaxillofacial Fixation System

MCF 直型桥板 T705A、T707A
MCF Straight Plate




产品名称 Product Name	材质 Material	厚度 Thickness(mm)	孔数×长度 Bar Length(mm)	产品编号 Product No.
MCF 直型桥板 MCF straight Plate	纯钛 Pure Titanium	0.6薄型	2×15	T705A0200
			4×20	T705A0400
			6×29	T705A0600
			8×38	T705A0800
		0.8薄型	2×16	T707A0200
			4×21	T707A0400
			6×30	T707A0600
			8×39	T707A0800

MCF 直型板 T705A、T707A
MCF Straight Plate




产品名称 Product Name	材质 Material	厚度 Thickness(mm)	孔数×长度 Bar Length(mm)	产品编号 Product No.
MCF 直型板 MCF straight Plate	纯钛 Pure Titanium	0.6薄型	12×53	T705A1200
		0.8薄型	12×54	T707A1200

MCF 斜L型板 (左右) T705B、T707B
MCF Oblique L Plate (Left,Right)




产品名称 Product Name	材质 Material	厚度 Thickness(mm)	孔数×长度 Bar Length(mm)	产品编号 Product No.
MCF 斜L型板 (左右) MCF Oblique L Plate (Left,Right)	纯钛 Pure Titanium	0.6薄型	4×17左	T705B0401
			4×17右	T705B0402
		0.8薄型	4×17左	T707B0401
			4×17右	T707B0402

MCF 斜L型板 (左右) T705B、T707B
MCF Oblique L Plate (Left,Right)



产品名称 Product Name	材质 Material	厚度 Thickness(mm)	孔数×长度 Bar Length(mm)	产品编号 Product No.
MCF 斜L型板 (左右) MCF Oblique L Plate (Left,Right)	纯钛 Pure Titanium	0.6薄型	5×21左	T705B0501
			5×21右	T705B0502
		0.8薄型	5×22左	T707B0501
			5×22右	T707B0502

MCF 斜L型板 (左右) T705B、T707B
MCF Oblique L Plate (Left,Right)



产品名称 Product Name	材质 Material	厚度 Thickness(mm)	孔数×长度 Bar Length(mm)	产品编号 Product No.
MCF 斜L型板 (左右) MCF Oblique L Plate (Left,Right)	纯钛 Pure Titanium	0.6薄型	7×26左	T705B0701
			7×26右	T705B0702
		0.8薄型	7×27左	T707B0701
			7×27右	T707B0702

MCF Multifunction Craniomaxillofacial Fixation System

B

MCF T型板 T705C、T707C
MCF T Plate




产品名称 Product Name	材质 Material	厚度 Thickness(mm)	孔数×长度 Bar Length(mm)	产品编号 Product No.
MCF T型板 MCF T Plate	纯钛 Pure Titanium	0.6薄型	7×25	T705C0700
		0.8薄型	7×26	T707C0700

MCF Y型板 T705D、T707D
MCF Y Plate




产品名称 Product Name	材质 Material	厚度 Thickness(mm)	孔数×长度 Bar Length(mm)	产品编号 Product No.
MCF Y型板 MCF Y Plate	纯钛 Pure Titanium	0.6薄型	5×20	T705D0500
		0.8薄型	5×21	T707D0500

MCF 弧型板 T705E、T707E
MCF Arc-Shaped Plate



产品名称 Product Name	材质 Material	厚度 Thickness(mm)	孔数×长度 Bar Length(mm)	产品编号 Product No.
MCF 弧形板 MCF Arc-Shaped Plate	纯钛 Pure Titanium	0.6薄型	6×30	T705E0600
			8×39	T705E0800
		0.8薄型	6×30	T707E0600
			8×39	T707E0800

MCF 双Y型板 T705F、T707F
MCF Double Y Plate



产品名称 Product Name	材质 Material	厚度 Thickness(mm)	孔数×长度 Bar Length(mm)	产品编号 Product No.
MCF 双Y型板 MCF Double Y Plate	纯钛 Pure Titanium	0.6薄型	6×20	T705F0600
		0.8薄型	6×21	T707F0600


MCF X型板 T705G、T707G
MCF X plate



产品名称 Product Name	材质 Material	厚度 Thickness(mm)	孔数×长度 Bar Length(mm)	产品编号 Product No.
MCF X型板 MCF X plate	纯钛 Pure Titanium	0.6薄型	4×10	T705G0400
		0.8薄型	4×11	T707G0400

MCF Multifunction Craniomaxillofacial Fixation System

MCF解剖型L型板 正颌 T707H
MCF Anatomic L Plate Orthognathic



MCF 斜L型板 正颌 T707J
MCF Oblique L Plate Orthognathic



MCF 双弯颌板 正颌 T710C
MCF Double Curved Jaw Plate Orthognathic



MCF 直矢状劈开板 正颌 T710B
MCF Straight Sagittal Split Plate Orthognathic



MCF 弯曲矢状劈开板 正颌 T710A
MCF Curved Sagittal Split Plate Orthognathic



MCF 眶底解剖重建板 T705H
MCF Orbital Floor Anatomical Reconstruction Plate



产品名称 Product Name	材质 Material	厚度 Thickness(mm)	孔数 Bar Length(mm)	桥长 Long bridge	产品编号 Product No.
MCF 解剖型L型板 正颌 MCF Anatomic L Plate Orthognathic	纯钛 Pure Titanium	0.8薄型	6孔左	4	T707H6041
			6孔右	4	T707H6042
			6孔左	8	T707H6081
			6孔右	8	T707H6082
			6孔左	12	T707H6121
MCF 斜L型板 正颌 MCF Oblique L Plate Orthognathic	纯钛 Pure Titanium	0.8薄型	4孔左	2	T707J4021
			4孔右	2	T707J4022
			4孔左	6	T707J4061
			4孔右	6	T707J4062
			4孔左	11	T707J4111
			4孔右	11	T707J4112
			6孔左	4	T707J6041
			6孔右	4	T707J6042
			6孔左	8	T707J6081
			6孔右	8	T707J6082
			6孔左	13	T707J6131
			6孔右	13	T707J6132
MCF 双弯颌板 正颌 MCF Double Curved Jaw Plate Orthognathic	纯钛 Pure Titanium	0.8薄型	5孔	4	T710C0504
				6	T710C0506
				8	T710C0508
				10	T710C0510

产品名称 Product Name	材质 Material	厚度 Thickness(mm)	孔数 Bar Length(mm)	桥长 Long bridge	产品编号 Product No.
MCF 直矢状劈开板 正颌 MCF Straight Sagittal Split Plate Orthognathic	纯钛 Pure Titanium	1.0标准型	4孔	6	T710B0406
				8	T710B0408
				10	T710B0410
				12	T710B0412
MCF 弯曲矢状 劈开板 正颌 MCF Curved Sagittal Split Plate Orthognathic	纯钛 Pure Titanium	1.0标准型	6孔	6	T710A0660
				8	T710A0860
				10	T710A1060
				12	T710A1260

产品名称 Product Name	材质 Material	厚度 Thickness(mm)	规格 Bar Length(mm)	型号 Material	产品编号 Product No.
MCF 眶底解剖重建板 MCF Orbital Floor Anatomical Reconstruction Plate	纯钛 Pure Titanium	0.6薄型	大号左	42x37	T705H001L
			大号右	42x37	T705H001R
			小号左	42x36	T705H002L
			小号右	42x36	T705H002R

MCF Multifunction Craniomaxillofacial Fixation System

MCF 0.6 薄型 颌面板系列

MCF 0.6 Thin Maxillofacial Plate Series



配套螺钉

Matching screws

MCF $\Phi 1.5$ 自钻自攻螺钉

MCF $\Phi 1.5$ Self Drilling And Tapping Screw



MCF $\Phi 1.5$ 自攻螺钉

MCF $\Phi 1.5$ Self Tapping Screw



MCF $\Phi 1.7$ 修正螺钉

MCF $\Phi 1.7$ Self Tapping Correction Screw



MCF 0.8 薄型 颌面板系列

MCF 0.8 Thin Maxillofacial Plate Series



配套螺钉

Matching screws

MCF $\Phi 2.0$ 自钻自攻螺钉

MCF $\Phi 2.0$ Self Drilling And Tapping Screw




MCF $\Phi 2.0$ 自攻螺钉

MCF $\Phi 2.0$ Self Tapping Screw



MCF 0.8 薄型 正颌颌面板系列

MCF 0.8 Standard Orthognathic Maxillofacial Plate Series



MCF 眶底解剖重建板

Mcf Orbital Floor Anatomical Reconstruction Plate



MCF 1.0 标准型 正颌颌面板系列

MCF 1.0 Standard Orthognathic Maxillofacial Plate Series



配套螺钉

Matching screws

MCF $\Phi 2.0$ 自钻自攻螺钉

MCF $\Phi 2.0$ Self Drilling And Tapping Screw



MCF $\Phi 2.0$ 自攻螺钉

MCF $\Phi 2.0$ Self Tapping Screw



MCF Multifunction Craniomaxillofacial Fixation System

MCF $\Phi 1.5$ 自钻自攻螺钉 32325
MCF $\Phi 1.5$ Self Drilling And Tapping Screw



产品名称 Product Name	材质 Material	直径 Bar Length(mm)	长度 Thickness(mm)	产品编号 Product No.
MCF $\Phi 1.5$ 自钻自攻螺钉 MCF $\Phi 1.5$ Self Drilling And Tapping Screw	钛合金 Titanium Alloy	$\Phi 1.5$	4	32325004
			5	32325005
			6	32325006
			7	32325007

MCF $\Phi 1.5$ 自攻螺钉 32321
MCF $\Phi 1.5$ Self Tapping Screw



产品名称 Product Name	材质 Material	直径 Bar Length(mm)	长度 Thickness(mm)	产品编号 Product No.
MCF $\Phi 1.5$ 自攻螺钉 MCF $\Phi 1.5$ Self Tapping Screw	钛合金 Titanium Alloy	$\Phi 1.5$	4	32321004
			5	32321005
			6	32321006
			7	32321007

MCF $\Phi 1.7$ 自攻修正螺钉 32322
MCF $\Phi 1.7$ Self Tapping Correction Screw



产品名称 Product Name	材质 Material	直径 Bar Length(mm)	长度 Thickness(mm)	产品编号 Product No.
MCF $\Phi 1.7$ 修正螺钉 MCF $\Phi 1.7$ Self Tapping Correction Screw	钛合金 Titanium Alloy	$\Phi 1.7$	5	32322005
			6	32322006
			7	32322007
			8	32322008

MCF $\Phi 2.0$ 自钻自攻螺钉 32326
MCF $\Phi 2.0$ Self Drilling And Tapping Screw



产品名称 Product Name	材质 Material	直径 Bar Length(mm)	长度 Thickness(mm)	产品编号 Product No.
MCF $\Phi 2.0$ 自钻自攻螺钉 MCF $\Phi 2.0$ Self Drilling And Tapping Screw	钛合金 Titanium Alloy	$\Phi 2.0$	5	32326005
			6	32326006
			7	32326007
			8	32326008

MCF $\Phi 2.0$ 自攻螺钉 32323
MCF $\Phi 2.0$ Self Tapping Screw





产品名称 Product Name	材质 Material	直径 Bar Length(mm)	长度 Thickness(mm)	产品编号 Product No.
MCF $\Phi 2.0$ 自攻螺钉 MCF $\Phi 2.0$ Self Tapping Screw	钛合金 Titanium Alloy	$\Phi 2.0$	5	32323005
			6	32323006
			7	32323007
			8	32323008

MCF 多功能颅颌面固定系统 91279901
MCF Multifunction Craniomaxillofacial Fixation System



PDF 文件
PDF File





MCF 直型板（无菌方案）
MCF straight Plate (Sterile Solution)



无菌包装
Sterile Packaging




无菌包装产品数量图示
Sterile Packaging Products Quantity Graphic



配套工具 Supporting Tools

91130221 MCF 一体式起子 十字 138mm
MCF All In One Screwdriver Cross 138mm



MCF 直型桥板 T705A
MCF Straight Plate



0.6mm

产品名称 Product Name	材质 Material	厚度 Thickness(mm)	孔数×长度 Bar Length(mm)	产品编号 Product No.
MCF 直型桥板 MCF Straight Plate	纯钛 Pure Titanium	0.6	2×15	T705A0200

MCF $\Phi 1.5$ 自钻自攻螺钉 32325
MCF $\Phi 1.5$ Self Drilling And Tapping Screw




产品名称 Product Name	材质 Material	直径 Bar Length(mm)	长度 Thickness(mm)	产品编号 Product No.
MCF $\Phi 1.5$ 自钻自攻螺钉 MCF $\Phi 1.5$ Self Drilling And Tapping Screw	钛合金 Titanium Alloy	$\Phi 1.5$	4	32325004
			5	32325005
			6	32325006
			7	32325007

Titanium Mesh System

数字化三维钛网 (定制)
Digital Three Dimension Titanium Mesh (Customization)




患者CT.MRI.DSA数据
CT.MRI.DSA Data Of Patient



计算机颅骨成像重建
Peconstruction Image



计算机钛网定制设计
Titanium Net Custom-designed




数字化三维钛网成型
Multi-faceted Shaped Of Restoration Of
Data Applications Titanium Mesh



视频在线播放
Video online play

钛网 T708
Titanium Mesh



MCFΦ2.0 自钻自攻螺钉 32326
MCF Φ2.0 Self Drilling And Tapping Screw






产品名称 Product Name	材质 Material	厚度 (mm) Thickness	规格 Specification	产品编号 Product No.
钛网 Titanium Mesh	纯钛 Pure Titanium	0.8	90×90	T708W0090
			80×100	T708W0100
			100×119	T708W0119
			112×132	T708W0132
			145×150	T708W0150
			142×158	T708W0158

产品名称 Product Name	材质 Material	直径 Bar Length(mm)	长度 Thickness(mm)	产品编号 Product No.
MCF Φ2.0 自钻自攻螺钉 MCF Φ2.0 Self Drilling And Tapping Screw	钛合金 Titanium Alloy	Φ2.0	5	32326005
			6	32326006
			7	32326007
			8	32326008


Maxillofacial Bone Plate Series

B



Maxillofacial Bone Plate Series

类别指示图:
Category directive



T710N0411



T710N0822



T710N1235



T710N1646

B

MCFΦ2.0 自钻自攻螺钉 32326
MCF Φ2.0 Self Drilling And Tapping Screw



产品名称 Product Name	材质 Material	直径 Bar Length(mm)	长度 Thickness(mm)	产品编号 Product No.
MCF Φ2.0 自钻自攻螺钉 MCF Φ2.0 Self Drilling And Tapping Screw	钛合金 Titanium Alloy	Φ2.0	5	32326005
			6	32326006
			7	32326007
			8	32326008

配套工具 Accessory tools

91130221 MCF 一体式起子 十字 138mm
Mcf All In One Screwdriver Cross 138mm



产品特点: 自攻自钻, 螺钉可持取 Product Features: Self-Tapping And Drilling, Screw Can Be Retained

金属接骨螺钉 (十字型) T5022
Metal Bone Fracture Screw (Cross Type)



产品名称 Product Name	材质 Material	直径 Bar Length(mm)	长度 Thickness(mm)	产品编号 Product No.
金属接骨螺钉 (十字型) Metal Bone Fracture Screw (Cross Type)	钛合金 Titanium Alloy	Φ2.0	5	T5022105
			6	T5022106
			7	T5022107
			8	T5022108

配套工具 Accessory tools

91111501 限位钻头
Limit Drill Bit, Φ1.5×8mm

91130201 十字起子
Crosshead Screwdriver




产品特点: 非自攻自钻, 螺钉不可持取 Product Features: Non Self-Tapping And Drilling, Screw Can Not Be Retained



髓内钉系列
Intramedullary Nail Series

弹性髓内针（无菌包装）

Elastic Intramedullary Pins (Sterile Packaging)



金属髓内针 T894
Metal Intramedullary Pins

产品名称 Product Name	包装形式 Packaging Form	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
金属髓内针 Metal Intramedullary Pins	无菌包装 Sterile Packaging	钛合金 Titanium Alloy	2.0	320	T89420320-S
			2.5	400	T89425400-S
			3.0	400	T89430400-S
			3.5	400	T89435400-S
			4.0	400	T89440400-S

尾帽 T895
End Cap




产品名称 Product Name	包装形式 Packaging Form	材质 Material	规格 Specification	产品编号 Product No.
尾帽 End Cap	无菌包装 Sterile Packaging	钛合金 Titanium Alloy	小 (Mini)	T8950000-S
			大 (Large)	T8950001-S

备注：Φ2.0、Φ2.5金属髓内针选用小号尾帽，Φ3.0、3.5、4.0金属髓内针选用大号尾帽。
Note: Φ2.0, Φ2.5 metal intramedullary pins use small end caps, Φ3.0、3.5、4.0 metal intramedullary pins use large end caps.

金属弹性髓内针安装器械包 92619901
Metal Elastic Intramedullary Pins Instrument Set



PDF 文件
PDF File



QML 多维尺、桡骨髓内钉系统（无菌包装）

QML Multi-Dimensional Ulna & Radius Intramedullary Nail System (Sterile Packaging)



QML 尺桡骨髓内钉安装器械包 92629901

QML Ulna&Radius Intramedullary Nail System Instrument Set



PDF 文件
PDF File



QML 多维尺、桡骨髓内钉系统（无菌包装）



QML Multi-Dimensional Ulna & Radius Intramedullary Nail System(Sterile Packaging)

QML 多维尺骨带锁髓内钉 65059、65060
QML Multi-Dimensional Ulna Intramedullary Nail



产品名称 Product Name	产品形式 Packaging From	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
QML 多维尺骨带锁髓内钉 QML Multi-Dimensional Ulna Intramedullary Nail	无菌包装 Sterile Packaging	钛合金 Titanium Alloy	3.0	210左	65059211
				210右	65059212
				220左	65059221
				220右	65059222
				230左	65059231
				230右	65059232
				240左	65059241
				240右	65059242
				250左	65059251
				250右	65059252
			3.5	260左	65059261
				260右	65059262
				270左	65059271
				270右	65059272
				210左	65060211
				210右	65060212
				220左	65060221
				220右	65060222
				230左	65060231
				230右	65060232

QML 多维桡骨带锁髓内钉 65061、65062
QML Multi-Dimensional Radius Intramedullary Nail



产品名称 Product Name	产品形式 Packaging From	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
QML 多维桡骨带锁髓内钉 QML Multi-Dimensional Radius Intramedullary Nail	无菌包装 Sterile Packaging	钛合金 Titanium Alloy	3.0	190左	65062191
				190右	65062192
				200左	65062201
				200右	65062202
				210左	65062211
				210右	65062212
				220左	65062221
				220右	65062222
				230左	65062231
				230右	65062232
			3.5	190左	65061191
				190右	65061192
				200左	65061201
				200右	65061202
				210左	65061211
				210右	65061212
				220左	65061221
				220右	65061222
				230左	65061231
				230右	65061232

QML 尺、桡骨尾帽 65064
QML Ulna & Radius End Cap



产品名称 Product Name	产品形式 Packaging From	材质 Material	长度(mm) Length	产品编号 Product No.
QML 尺、桡骨尾帽 QML Ulna & Radius End Cap	无菌包装 Sterile Packaging	钛合金 Titanium Alloy	0	65064100

QMLΦ3.0 全螺纹锁钉 65063
QML Φ3.0 Full Thread Locking Screw




产品名称 Product Name	产品形式 Packaging From	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
QML Φ3.0 全螺纹锁钉 Φ3.0 Full Thread Locking Screw	无菌包装 Sterile Packaging	钛合金 Titanium Alloy	3.0	10~60 每2mm为一档 The specification is 2mm)	65063010 ~65063060

QML 多维肱骨带锁髓内钉系统

QML Multi-Dimensional Humeral Intramedullary Nail System



视频在线播放
Video online play



钉中钉设计, 改善骨质疏松
骨骼的固定。
Screw In Screw Design,
Improve Osteoporosis And
Bone Fixation.



短型
Short

上行螺钉提供内侧稳定支撑
The Uplink Screw Provides
Innerstability Support



加长型
Extended




多平面交锁
Multiplanar
locking

QML 多维肱骨带锁髓内钉安装器械包 92529901

QML Multi-Dimensional Humeral Intramedullary Nail System Instrument Set



PDF 文件
PDF File



QML 多维肱骨带锁髓内钉系统

QML Multi-Dimensional Humeral Intramedullary Nail System

QML 多维肱骨带锁髓内钉 (短型) 63044、63045
QML Multi-Dimensional Humeral Intramedullary Nail (Short)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
QML 多维肱骨带锁髓内钉 (短型) QML Multi-Dimensional Humeral Intramedullary Nail (Short)	钛合金 Titanium Alloy	8.5	160左	63044161
			160右	63044162
		9.5	160左	63045161
			160右	63045162

QML 多维肱骨带锁髓内钉 (加长型) 63049、63050
QML Multi-Dimensional Humeral Intramedullary Nail (Extended)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
QML 多维肱骨带锁髓内钉 (加长型) Multi-Dimensional Humeral Intramedullary Nail (Extended)	钛合金 Titanium Alloy	7.0	200左	63049201
			200右	63049202
			220左	63049221
			220右	63049222
			240左	63049241
			240右	63049242
			260左	63049261
			260右	63049262
			280左	63049281
			280右	63049282
		8.0	300左	63049301
			300右	63049302
			200左	63050201
			200右	63050202
			220左	63050221
			220右	63050222
			240左	63050241
			240右	63050242
			260左	63050261
			260右	63050262
			280左	63050281
			280右	63050282
			300左	63050301
			300右	63050302

QML Φ 4.5 全螺纹锁钉 63046
QML Φ 4.5 Full Thread Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
QML Φ 4.5 全螺纹锁钉 QML Φ 4.5 Full Thread Locking Screw	钛合金 Titanium Alloy	4.5	20-60 每2mm为一档 The specification is 2mm)	63046020 -63046060

QML Φ 3.5 全螺纹锁钉 63047
QML Φ 3.5 Full Thread Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
QML Φ 3.5 全螺纹锁钉 QML Φ 3.5 Full Thread Locking Screw	钛合金 Titanium Alloy	3.5	20-46 每2mm为一档 The specification is 2mm)	63047020 -63047046

QML 锁定尾帽 63048
QML Locking End Cap



产品名称 Product Name	材质 Material	长度(mm) Length	产品编号 Product No.
QML 锁定尾帽 QML Locking End Cap	钛合金 Titanium Alloy	0	63048100
		0 (加长型)	63048200
		5	63048105
		5 (加长型)	63048215
		10	63048110

Φ 3.5 梅花型锁定螺钉(阳极氧化) T492
Torx Locking Screw (Anodized)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型锁定螺钉 (阳极氧化) Torx Locking Screw (Anodized)	钛合金 Titanium Alloy	3.5	10-50 (每2mm为一档) The specification is 2mm	T4923510- T4923550
			55-80 (每5mm为一档) The specification is 5mm	T4923555- T4923580

QML 多维胫骨髌上入路髓内钉系统

QML Multi-Dimensional Tibial Suprapatellar approach Intramedullary Nail System



视频在线播放
Video online play



短型
Short

- 可持式 $\Phi 5.0$ 全螺纹锁钉
- 松质骨螺纹提供更好的把持力
- Holding Type Full Thread Locking Screw
- Can Cellous Bone Thread Provide Better Holding Force



加长型
Extended



- 可持式 $\Phi 4.5$ 双线全螺纹锁钉。
- 双线螺纹，拧钉更快。
- $\Phi 4.5$ Holding Type Double Thread Locking Screw
- Screw Speed Of Double Thread Is Faster

QML 多维胫骨髌上入路髓内钉安装器械包 92339901

QML Multi-Dimensional Tibial Supratellar approach Intramedullart Nail System Instrument Set



PDF 文件
PDF File



QML 多维胫骨髌上入路髓内钉系统

QML Multi-Dimensional Tibial Suprapatellar approach Intramedullary Nail System

QML 多维胫骨髌上入路髓内钉 63201、63202、63203
QML Multi-Dimensional Tibial Suprapatellar approach Intramedullary Nail



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
QML 多维胫骨髌上入路髓内钉 QML Multi-Dimensional Tibial Suprapatellar approach Intramedullary Nail	钛合金 Titanium Alloy	9.0	160	63201160
			180	63201180
			200	63201200
		10	160	63202160
			180	63202180
			200	63202200
		11	160	63203160
			180	63203180
			200	63203200

QML 多维胫骨髌上入路髓内钉 (加长型) 63204、63205、63206
QML Multi-Dimensional Tibial Suprapatellar approach Intramedullary Nail (Extended)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
QML 多维胫骨髌上入路髓内钉 (加长型) QML Multi-Dimensional Tibial Suprapatellar approach Intramedullary Nail (Extended)	钛合金 Titanium Alloy	8.4	255	63204255
			270	63204270
			285	63204285
			300	63204300
			315	63204315
			330	63204330
			345	63204345
			360	63204360
			375	63204375
		9.0	255	63205255
			270	63205270
			285	63205285
			300	63205300
			315	63205315
			330	63205330
			345	63205345
			360	63205360
			375	63205375
		10	255	63206255
			270	63206270
			285	63206285
			300	63206300
			315	63206315
			330	63206330
			345	63206345
			360	63206360
			375	63206375

QML Φ5.0全螺纹锁钉 63207
QML Φ5.0 Full Thread Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
QML 5.0全螺纹锁钉 QML Φ5.0 Full Thread Locking Screw	钛合金 Titanium Alloy	5.0	25-85 每5mm为一档 (in 5mm increments)	63207025 -63207085

Φ4.5可持式双线全螺纹锁钉 63009
Φ4.5 Holding Type Double Thread Locking Screw




产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
Φ4.5 可持式双线程 全螺纹锁钉 Φ4.5 Holding Type Double Thread Locking Screw	钛合金 Titanium Alloy	4.5	30	63009030
			32	63009032
			35	63009035
			38	63009038
			40、45、50、55、 60、65、70、75、 80、85	63009040- 63009085
备注：选用SW4.0中空六角起子 Note: Choose a SW4.0 Hollow Hexagonal Screwdriver				

QML 锁定尾帽 63208
QML Locking End Cap



产品名称 Product Name	材质 Material	长度(mm) Length	产品编号 Product No.
QML 锁定尾帽 QML Locking End Cap	钛合金 Titanium Alloy	0	63208100

QML 尾帽 63209
QML End Cap



产品名称 Product Name	材质 Material	长度(mm) Length	产品编号 Product No.
QML 尾帽 End Cap	钛合金 Titanium Alloy	5	63209005
		10	63209010
		15	63209015

QML Multi-Dimensional Femoral Reverse Intramedullary Nail System




视频在线播放
Video online play






QML 多维股骨倒置髓内钉安装器械包 92279901

QML Multi-Dimensional Femoral Reverse Intramedullart Nail System Instrument Set



PDF 文件
PDF File



QML 多维股骨倒置带锁髓内钉系统

QML Multi-Dimensional Femoral Reverse Intramedullary Nail System

QML 多维股骨倒置带锁髓内钉 63301、63302、63303
QML Multi-Dimensional Femoral Reverse
Intramedullary Nail



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
QML 多维股骨倒置带锁髓内钉 QML Multi-Dimensional Femoral Reverse Intramedullary Nail	钛合金 Titanium Alloy	9.4	180	63301180
			220	63301220
			260	63301260
		10	180	63302180
			220	63302220
			260	63302260
		11	180	63303180
			220	63303220
			260	63303260

QML 多维股骨倒置带锁髓内钉 63301、63302、63303
QML Multi-Dimensional Femoral Reverse
Intramedullary Nail



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
QML 多维股骨倒置带锁髓内钉 QML Multi-Dimensional Femoral Reverse Intramedullary Nail	钛合金 Titanium Alloy	9.4	300	63301300
			320	63301320
			340	63301340
			360	63301360
			380	63301380
			400	63301400
		10	300	63302300
			320	63302320
			340	63302340
			360	63302360
			380	63302380
			400	63302400
		11	300	63303300
			320	63303320
			340	63303340
			360	63303360
			380	63303380
			400	63303400

QML Φ5.0 全螺纹锁钉 63207
QML Φ5.0 Full Thread Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
QML 5.0 全螺纹锁钉 QML Φ5.0 Full Thread Locking Screw	钛合金 Titanium Alloy	5.0	30-85 每5mm为一档 in 5mm increments)	63207030 -63207085

Φ4.5 可持式双线全螺纹锁钉 63009
Φ4.5 Holding Type Double Thread Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
Φ4.5 可持式双线全螺纹锁钉 Φ4.5 Holding Type Double Thread Locking Screw	钛合金 Titanium Alloy	4.5	30	63009030
			32	63009032
			35	63009035
			38	63009038
			40、45、50、55、60、65、70、75、80、85	63009040-63009085


备注：选用SW4.0中空六角起子
Note: Choose a SW4.0 Hollow Hexagonal Screwdriver

QML 锁定尾帽 63304
QML Locking End Cap



产品名称 Product Name	材质 Material	长度(mm) Length	产品编号 Product No.
QML 锁定尾帽 QML Locking End Cap	钛合金 Titanium Alloy	0	63304100

QML 尾帽 63305
QML End Cap



产品名称 Product Name	材质 Material	长度(mm) Length	产品编号 Product No.
QML 尾帽 QML End Cap	钛合金 Titanium Alloy	0	63305100
		5	63305105
		10	63305110
		15	63305115

KTO®-滑动加压伽玛钉系统 (亚洲版)



KTO®-Sliding Compression Gamma Intramedullary Nail System (Asian Version)

系统特点

Character Of The System



视频在线播放
Video online play

C

远端锁钉方案

Distal Locking Screw Plan



方案一
Plan 1



方案二
Plan 2

- 联合交锁钉：抗旋转、防切出、可加压。
- Integrated Interlocking Screw: Anti-rotation, anti-cutting, Compression.



- 近端直径小，更适合亚洲人。
- 5°外偏角，大粗隆顶点微创入路。
- 近端梯型设计，稳定性更好。
- 127°颈干角更适合亚洲人。
- Proximal Trapezoid Design, better Stability.
- 5° lateral Offset, minimally Invasive Greater Trochanter Approach.
- Proximal Diameter Is Small And More Suitable For Asians.
- 127° is More Suitable For Asians

- 可持术中防脱钉
- 双线螺纹拧钉速度更快。
- Holding Type Intraoperative Anti-off Screw
- Screw Speed Of Double Thread Is Faster




KTO 滑动加压伽玛钉安装器械包(小号) 92239901

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Small)



PDF 文件
PDF File

- 一个手术切口，完成两处锁钉（仅180mm主钉）。
- One Surgery Incision, two Locking Screws Completed (Only 180mm Main Nail)



KTO®-滑动加压伽玛钉系统 (亚洲版)

KTO®-Sliding Compression Gamma Intramedullary Nail System (Asian Version)

系统特点

Character Of The System



加长型
Extended

股骨钉模式
Femoral model


- 可持式 $\Phi 4.5$ 双线全螺纹锁钉
- 多平面交锁
- Holding Type Full Thread Locking Screw
- Multi Plane Cross Locjing

KTO 滑动加压伽玛钉安装器械包小号(加长版) 92239902

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended version)



PDF 文件



KTO®-滑动加压伽玛钉系统（亚洲版）

KTO®-Sliding Compression Gamma Intramedullary Nail System (Asian Version)

KTO®-滑动加压伽玛钉小号 63035

KTO®-Sliding Compression Gamma Nail Small



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KTO®-滑动加压伽玛钉 小号 KTO®-Sliding Compression Gamma Nail Small	钛合金 Titanium Alloy	9.4	180、200、220	63035118- 63035122
		10	180、200、220	63035218- 63035222
		11	180、200、220	63035318- 63035322
		12	180、200、220	63035418- 63035422

KTO®-滑动加压拉力钉（母钉）小号 63041

KTO®-Sliding Compression Lag Screw (Parent Screw) Small



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KTO®-滑动加压拉力钉 (母钉)小号 KTO®-Sliding Compression Lag Screw (Parent Screw) Small	钛合金 Titanium Alloy	9.0	70、75、80、85、 90、95、100、105、 110、115、120	63041070- 63041120

KTO®-滑动加压拉力钉（子钉）小号 63042

KTO®-Sliding Compression Lag Screw (Sub Screw) Small



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KTO®-滑动加压拉力钉 (子钉) KTO®-Sliding Compression Lag Screw (Sub Screw)	钛合金 Titanium Alloy	6.0	65、70、75、80、 85、90、95、100、 105、110、115、	63042065- 63042115

Φ4.5 可持式双线全螺纹锁钉 63009

Φ4.5 Holding Type Double Thread Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
Φ4.5 可持式双线 全螺纹锁钉 Φ4.5 Holding Type Double Thread Locking Screw	钛合金 Titanium Alloy	4.5	30	63009030
			32	63009032
			35	63009035
			38	63009038
			40、45、50、55、 60、65、70、75、 80、85	63009040- 63009085

备注：选用SW4.0中空六角起子
Note: Choose a SW4.0 Hollow Hexagonal Screwdriver

锁定尾帽 小号 63043


Locking End Cap Small



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
锁定尾帽 Locking End Cap Small	钛合金 Titanium Alloy	0	63043000
		5	63043005
		10	63043010
		15	63043015

普通尾帽 63306

Normal End Cap



产品名称 Product Name	材质 Material	长度(mm) Length	产品编号 Product No.
普通尾帽 Normal End Cap	钛合金 Titanium Alloy	0	63306100
		5	63306105
		10	63306110
		15	63306115

KTO®-滑动加压伽玛钉系统（亚洲版）

KTO®-Sliding Compression Gamma Intramedullary Nail System (Asian Version)

KTO®-滑动加压伽玛钉小号(加长型) 63037 、63038、63039

KTO®-Sliding Compression Gamma Nail Small(Extended edition)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KTO®-滑动加压伽玛钉 小号(加长型) KTO®-Sliding Compression Gamma Nail Small (Extended edition)	钛合金 Titanium Alloy	9.4	300左	63037301
			300右	63037302
			320左	63037321
			320右	63037322
			340左	63037341
			340右	63037342
			360左	63037361
			360右	63037362
			380左	63037381
			380右	63037382
			400左	63037401
			400右	63037402
		10	320左	63038321
			320右	63038322
			340左	63038341
			340右	63038342
			360左	63038361
			360右	63038362
			380左	63038381
			380右	63038382
			400左	63038401
			400右	63038402
			420左	63038421
			420右	63038422
			440左	63038441
			440右	63038442
		11	320左	63039321
			320右	63039322
			340左	63039341
			340右	63039342
			360左	63039361
			360右	63039362
			380左	63039381
			380右	63039382
			400左	63039401
			400右	63039402
			420左	63039421
			420右	63039422

C

KTO®-滑动加压伽玛玛钉系统（加强型）

KTO®-Sliding Compression Gamma Intramedullary Nail System
(Enhanced Type)



KTO®-滑动加压伽玛钉系统（加强型）



KTO®-Sliding Compression Gamma Intramedullary Nail System
(Enhanced Type)

KTO®-滑动加压伽玛钉小号（加强型） 6306
KTO®-Sliding Compression Gamma Nail Small
(Enhanced Edition)



KTO®-滑动加压伽玛钉小号(加长加强型)
63065、63066、63067
KTO®-Sliding Compression Gamma Nail Small
(Extended Enhanced Edition)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KTO®-滑动加压伽玛钉 小号（加强型） KTO®-Sliding Compression Gamma Nail Small (Enhanced Edition)	钛合金 Titanium Alloy	9.4	180、200、220、240	63060180- 63060240
		10	180、200、220、240	63061180- 63061240
		11	180、200、220、240	63062180- 63062240
		12	180、200、220、240	63063180- 63063240
		13	180、200、220、240	63064180- 63064240
产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KTO®-滑动加压伽玛钉 小号(加长加强型) KTO®-Sliding Compression Gamma Nail Small (Extended Enhanced Edition)	钛合金 Titanium Alloy	9.4	300左	63065301
			300右	63065302
			320左	63065321
			320右	63065322
			340左	63065341
			340右	63065342
			360左	63065361
			360右	63065362
			380左	63065381
			380右	63065382
			400左	63065401
			400右	63065402
		10	320左	63066321
			320右	63066322
			340左	63066341
			340右	63066342
			360左	63066361
			360右	63066362
			380左	63066381
			380右	63066382
			400左	63066401
			400右	63066402
			420左	63066421
			420右	63066422
		11	320左	63067321
			320右	63067322
			340左	63067341
			340右	63067342
			360左	63067361
			360右	63067362
			380左	63067381
			380右	63067382
			400左	63067401
			400右	63067402
			420左	63067421
			420右	63067422



KTO®-滑动加压伽玛钉系统

KTO®-Sliding Compression Gamma Intramedullary Nail System



视频在线播放
Video online play

C




联合交锁钉组合:


最佳的稳定性, 抗旋转, 防切出,
拧入加压钉时产生显著的加压效果。

Integrated Interlocking Screws
Combination:

Superior Stability, Anti-rotation,
Anti-cutting, Significant Compression
Effects Can Be Produced When Screwing
Into The Compression Screw.



远端钉孔可选择动态或静态交锁
Distal Locking Slot Can Choose
Dynamic Or Static Locking.



5° 外偏角提供大粗隆顶点微创入路

5° lateral Offset For Minimally Invasive
Trochanteric Tip Entry

15.8 mm



较小的近端设计, 有利于保留正常骨质和肌腱
等软组织。

近端梯形的横断面犹如关节假体柄设计, 能增
强稳定性和力学优势。

Smaller Proximal Design Can Help To
Preserve Normal Bone And Tendon
And Other Soft Tissue
Proximal Trapezoidal Cross-section Just
Like Joint Prosthesis Stem, Which Can
Enhance Stability And Mechanical Superiority.

可持, 术中防脱钉。
双线螺纹, 拧钉更快。

Holding Type Intraoperative Anti-off Screw
Screw Speed Of Double Thread Is Faster

KTO®-Sliding Compression Gamma Intramedullary Nail System

KTO®-滑动加压伽玛钉 T831

KTO®-Sliding Compression Gamma Nail



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KTO®-滑动加压伽玛钉 KTO®-Sliding Compression Gamma Nail	钛合金 Titanium Alloy	9.4	180、200、220	T83194180- T83194220
		10	180、200、220	T83110180- T83110220
		11	180、200、220	T83111180- T83111220
		12	180、200	T83112180- T83112200

KTO®-滑动加压伽玛钉（母钉） T832

KTO®-Sliding Compression Gamma Nail (Parent Screw)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KTO®-滑动加压伽玛钉 (母钉) KTO®-Sliding Compression Gamma Nail (Parent Screw)	钛合金 Titanium Alloy	10.5	70、75、80、85、 90、95、100、105、 110、115、120	T83200070- T83200120

KTO®滑动加压拉力钉（子钉） T833

KTO®Sliding Compression Lag Screw(Sub Screw)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KTO®-滑动加压拉力钉 (子钉) KTO®-Sliding Compression Lag Screw (Sub Screw)	钛合金 Titanium Alloy	6.5	65、70、75、80、 85、90、95、100、 105、110、115	T83300065- T83300115

Φ4.5 可持式双线全螺纹锁钉 63009

Φ4.5 Holding Type Double Thread Locking Screw




产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
Φ4.5 可持式双线 全螺纹锁钉 Φ4.5 Holding Type Double Thread Locking Screw	钛合金 Titanium Alloy	4.5	30	63009030
			32	63009032
			35	63009035
			38	63009038
			40、45、50、55、 60、65、70、75、 80、85	63009040- 63009085

备注：选用SW4.0中空六角起子
Note: Choose a SW4.0 Hollow Hexagonal Screwdriver

锁定尾帽 T8340

Locking End Cap



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
锁定尾帽 Locking End Cap	钛合金 Titanium Alloy	0	T8340000
		5	T8340005
		10	T8340010
		15	T8340015

KTO 滑动加压伽玛钉安装器械包 92229901

KTO Gamma Intramedullary Nail Instrument Set

PDF 文件
PDF File

KTO®Femoral Intramedullary Nail System




视频在线播放
Video online play

C




股骨钉模式
Femoral Model






重建钉模式
Reconstruction Model

KTO 股骨带锁髓内钉安装器械包 92259901
KTO Femoral Intramedullary Nail Instrument Set



PDF 文件
PDF File



KTO®Femoral Intramedullary Nail System

KTO®股骨带锁髓内钉 T820
KTO®Femoral Intramedullary Nail



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KTO®股骨 带锁髓内钉 KTO®Femoral Intramedullary Nail	钛合金 Titanium Alloy	Φ9.4	320、340、360、 380、400	T82094320R/L- T82094400R/L
		Φ10	320、340、360、 380、400、420	T82010320R/L- T82010420R/L
		Φ11	320、340、360、 380、400	T82011320R/L- T82011400R/L

防旋钉 T856
Anti-Rotation Screw




产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
防旋钉 Anti-Rotation Screw	钛合金 Titanium Alloy	Φ6.5	60、65、70、75、 80、85、90、95、 100、105、110、 115	T8566460- T85664115

Φ4.5 全螺纹锁钉 T852
Φ4.5 Full Thread Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
全螺纹锁钉 Φ4.5 Full Thread Locking Screw	钛合金 Titanium Alloy	Φ4.5	30	T8524030
			32	T8524032
			35	T8524035
			38	T8524038
			40、45、50、55、 60、65、70、75、 80	T8524040- T8524080

尾帽 T8201
End Cap




产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
尾帽 End Cap	钛合金 Titanium Alloy	0	T8201001
		5	T8201002
		10	T8201003
		15	T8201004
		20	T8201005


KTO®Tibial Intramedullary Nail System






视频在线播放
Video online play



KTO 胫骨带锁髓内钉安装器械包 92319901
KTO Tibial Intramedullary Nail Instrument Set



PDF 文件
PDF File



KTO®Tibial Intramedullary Nail System

KTO®胫骨带锁髓内钉(空心)T8684、T8690、T8610
KTO®Tibial Intramedullary Nail(Cannulated)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KTO®胫骨 带锁髓内钉(空心) KTO®Tibial Intramedullary Nail (Cannulated)	钛合金 Titanium Alloy	Φ8.4	255、270、285、 300、315、330、 345、360、375	T8684255- T8684375
		Φ9.0	255、270、285、 300、315、330、 345、360、375	T8690255- T8690375
		Φ10	255、270、285、 300、315、330、 345、360、375	T8610255- T8610375

Φ4.5 全螺纹锁钉 T852
Φ4.5 Full Thread Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
全螺纹锁钉 Full Thread Locking Screw	钛合金 Titanium Alloy	Φ4.5	30	T8524030
			32	T8524032
			35	T8524035
			38	T8524038
			40、45、50、55、 60、65、70、75、 80	T8524040- T8524080

Φ5.0 全螺纹锁钉 T851
Φ5.0 Full Thread Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
全螺纹锁钉 Full Thread Locking Screw	钛合金 Titanium Alloy	Φ5.0	40、45、50、55、 60、65、70、75、 80	T8515040- T8515080

普通尾帽 T8600001
Normal End Cap



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
普通尾帽 Normal End Cap	钛合金 Titanium Alloy	0	T8600001

锁定尾帽 T8600002
Locking End Cap



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
锁定尾帽 Locking End Cap	钛合金 Titanium Alloy	0	T8600002

C


KTO®PFNA Gamma Intramedullary Nail System



视频在线播放
Video online play



KTO-PFNA伽玛钉系统
KTO-PFNA Gamma Intramedullary Nail System






KTO-PFNA伽玛钉系统(加长型)
KTO-PFNA Gamma Intramedullary Nail System (Extended Edition)

KTO-PFNA伽玛带锁髓内钉安装器械包 92219901
KTO-PFNA Gamma Intramedullary Nail Instrument Set



PDF 文件
PDF File



KTO®PFNA Gamma Intramedullary Nail System

KTO®-PFNA主钉 T8794、T8710、T8711、T8712
KTO®-PFNA Main Nail



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KTO®-PFNA主钉 KTO®PFNA Main Nail	钛合金 Titanium Alloy	Φ9.4	180、200、 220、240	T8794180- T8794240
		Φ10	180、200、 220、240	T8710180- T8710240
		Φ11	180、200、 220、240	T8711180- T8711240
		Φ12	180、200、 220、240	T8712180- T8712240

KTO®-PFNA主钉 (加长版) T8794、T8710、T8711
KTO®-PFNA Main Nail (Extended Edition)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KTO®-PFNA主钉 (加长版) KTO®-PFNA Main Nail (Extended Edition)	钛合金 Titanium Alloy	Φ9.4	320、340、360、 380、400	T8794320R/L- T8794400R/L
		Φ10	320、340、360、 380、400、420	T8710320R/L- T8710420R/L
		Φ11	320、340、360、 380、400	T8711320R/L- T8711400R/L

KTO®-PFNA刀片钉 T8541
KTO®-PFNA Blade Nail




产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KTO®-PFNA刀片钉 KTO®-PFNA Blade Nail	钛合金 Titanium Alloy	Φ10.5	70、75、80、85、 90、95、100、 105、110、115、 120	T8541070- T8541120

Φ4.5 全螺纹锁钉 T852
Φ4.5 Full Thread Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
全螺纹锁钉 Full Thread Locking Screw	钛合金 Titanium Alloy	Φ4.5	30	T8524030
			32	T8524032
			35	T8524035
			38	T8524038
			40、45、50、55、 60、65、70、75、 80	T8524040- T8524080

锁定尾帽 T8571
Locking End Cap




产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
锁定尾帽 Locking End Cap	钛合金 Titanium Alloy	0	T8571000
		5	T8571005
		10	T8571010
		15	T8571015


普通尾帽 T858
Normal End Cap



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
普通尾帽 Normal End Cap	钛合金 Titanium Alloy	0	T8580000
		5	T8580005
		10	T8580010
		15	T8580015




KTO®-PFNA伽玛加强型系统
KTO®-PFNA Gamma Intramedullary Nail System





KTO®-PFNA伽玛钉系统(加长加强型)
KTO®-PFNA Gamma Intramedullary Nail System (Extended Edition)

KTO®- PFNA安装器械包加强型 短款 92219906

KTO®- PFNA Intramedullary Nail Instrument Set Enhanced Type (Short)



PDF 文件
PDF File



KTO®PFNA伽玛钉系统（加强型）

KTO® PFNA Gamma IntramedullaryNail System (Enhanced Edition)

KTO®-PFNA主钉（加强型） 63020
KTO®-PFNA Main Nail (Enhanced Edition)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KTO®PFNA主钉 (加强型) KTO®-PFNA Main Nail (Enhanced Edition)	钛合金 Titanium Alloy	Φ9.4	180、200、 220、240	63020118- 63020124
		Φ10	180、200、 220、240	63020218- 63020224
		Φ11	180、200、 220、240	63020318- 63020324
		Φ12	180、200、 220、240	63020418- 63020424

KTO®-PFNA加强型主钉（加长加强型） 63020
KTO®-PFNA Enhanced Main Nail
(Extended Enhanced Edition)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KTO®-PFNA主钉 (加长加强型) KTO®-PFNA Enhanced Main Nail (Extended Enhanced Edition)	钛合金 Titanium Alloy	Φ9.4	320、340、360、 380、400	63020132- 63020140
		Φ10	320、340、360、 380、400、420	63020232- 63020242
		Φ11	320、340、360、 380、400	63020332- 63020340

KTO®-PFNA刀片钉 T8541
KTO®-PFNA Blade Nail




产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KTO®-PFNA刀片钉 KTO®-PFNA Blade Nail	钛合金 Titanium Alloy	Φ10.5	70、75、80、85、 90、95、100、 105、110、115、 120	T8541070- T8541120

Φ4.5 全螺纹锁钉 T852
Φ4.5 Full Thread Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
全螺纹锁钉 Full Thread Locking Screw	钛合金 Titanium Alloy	Φ4.5	30	T8524030
			32	T8524032
			35	T8524035
			38	T8524038
			40、45、50、55、 60、65、70、75、 80	T8524040- T8524080

锁定尾帽 T8571
Locking End Cap



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
锁定尾帽 Locking End Cap	钛合金 Titanium Alloy	0	T8571000
		5	T8571005
		10	T8571010
		15	T8571015

普通尾帽 T8580
Normal End Cap



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
普通尾帽 Normal End Cap	钛合金 Titanium Alloy	0	T8580000
		5	T8580005
		10	T8580010
		15	T8580015

PF Small Hole Compression Gamma Nail System

PF 小孔伽玛钉 DT8410、DT8411
PF Small Hole Gamma Nail



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
PF 小孔伽玛钉 PF Small Hole Gamma Nail	钛合金 Titanium Alloy	10	180	DT8410180
			200	DT8410200
		11	180	DT8411180
			200	DT8411200

PF 小孔伽玛型拉力锁钉 DT8550
PF Small Hole Gamma Lag Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
PF 小孔伽玛型拉力锁钉 PF Small Hole Gamma Lag Screw	钛合金 Titanium Alloy	Φ9.5	70、75、80、85、 90、95、100、 105、110、115	DT8550070- DT8550115

防旋钉 T856
Anti-Rotation Screw




产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
防旋钉 Anti-Rotation Screw	钛合金 Titanium Alloy	Φ6.5	60、65、70、75、 80、85、90、95、 100、105、110、 115	T8566460- T85664115

Φ4.5 全螺纹锁钉 T852
Φ4.5 Full Thread Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
全螺纹锁钉 Full Thread Locking Screw	钛合金 Titanium Alloy	Φ4.5	30	T8524030
			32	T8524032
			35	T8524035
			38	T8524038
			40、45、50、55、 60、65、70、75、 80	T8524040- T8524080

锁定尾帽 T8571
Locking End Cap



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
锁定尾帽 Locking End Cap	钛合金 Titanium Alloy	0	T8571000
		5	T8571005
		10	T8571010
		15	T8571015

普通尾帽 T8580
Normal End Cap



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
普通尾帽 Normal End Cap	钛合金 Titanium Alloy	0	T8580000
		5	T8580005
		10	T8580010
		15	T8580015

KTO®Humerus Intramedullary Nail System

KTO®肱骨带锁髓内钉 T8870、T8875、T8880
KTO®Humerus Intramedullary Nail



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KTO®肱骨带锁髓内钉 KTO®Humerus Intramedullary Nail	钛合金 Titanium Alloy	Φ7.0	180、200、220、240、260、280	T8870180-T8870280
		Φ7.5	180、200、220、240、260、280	T8875180-T8875280
		Φ8.0	180、200、220、240、260、280	T8880180-T8880280

Φ3.5 全螺纹锁钉 T8920
Φ3.5 Full Thread Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
全螺纹锁钉 Full Thread Locking Screw	钛合金 Titanium Alloy	Φ3.5	20、22、24、26、28、30、32、34、36、38、40、42、44、46	T8920020-T8920046

主钉尾帽 T8800
Main End Cap




产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
主钉尾帽 Main End Cap	钛合金 Titanium Alloy	0	T8800001
		5	T8800005
		10	T8800010
		15	T8800015

KTO 肱骨带锁髓内钉安装器械包 92519901
KTO Humeral Intramedullary Nail Instrument Set



PDF 文件
PDF File



Tibial Intramedullary Nail System

胫骨带锁髓内钉瞄准装置

Tibial Intramedullary Nail- Aiming device



胫骨带锁髓内钉 T8384、T8390、T8310

Tibial Intramedullary Nail



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
胫骨带锁髓内钉 Tibial Intramedullary Nail	钛合金 Titanium Alloy	Φ8.4	260、280、300、 320	T8384260- T8384320
		Φ9.0	— 280、300、 320、340	T8390280- T8390340
		Φ10	— — 300、 320、340	T8310300- T8310340

Φ4.5 全螺纹锁钉 T852

Φ4.5 Full Thread Locking Screw





产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
全螺纹锁钉 Full Thread Locking Screw	钛合金 Titanium Alloy	Φ4.5	30	T8524030
			32	T8524032
			35	T8524035
			38	T8524038
			40、45、50、55、 60、65、70、75、 80	T8524040- T8524080

胫骨带锁髓内钉安装器械包 92329901

Tibial Intramedullary Nail Instrument Set



PDF 文件
PDF File



Femoral Intramedullary Nail System

股骨带锁髓内钉 T8194、T8110、T8111
Femoral Intramedullary Nail



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
股骨带锁髓内钉 Femoral Intramedullary Nail	钛合金 Titanium Alloy	Φ9.4	320、340、360、380	T8194320-T8194380
		Φ10	320、340、360、380、400	T8110320-T8110400
		Φ11	— — 360、380、400	T8111360-T8111400

股骨倒置带锁髓内钉 T8290、T8210、T8211
Femoral Reverse Intramedullary Nail



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
股骨倒置带锁髓内钉 Femoral Reverse Intramedullary Nail	钛合金 Titanium Alloy	Φ9	— 260、280、300	T8290260-T8290300
		Φ10	240、260、280、300、320、340、360	T8210240-T8210360
		Φ11	300、320、340、360	T8211300-T8211360

Φ4.5 全螺纹锁钉 T852
Φ4.5 Full Thread Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
全螺纹锁钉 Full Thread Locking Screw	钛合金 Titanium Alloy	Φ4.5	30	T8524030
			32	T8524032
			35	T8524035
			38	T8524038
			40、45、50、55、60、65、70、75、80	T8524040-T8524080

Φ6.5 全螺纹锁钉 T853
Φ6.5 Full Thread Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
全螺纹锁钉 Full Thread Locking Screw	钛合金 Titanium Alloy	Φ6.5	60、65、70、75、80、85、90、95	T8536560-T8536595

股骨带锁髓内钉安装器械包 92269901
Femoral Intramedullary Nail Instrument Set



PDF 文件
PDF File



股骨带锁髓内钉瞄准装置
Femoral Intramedullary Nail-DTD

股骨倒置带锁髓内钉安装器械包 92249901
Femoral Reverse Intramedullary Nail System Instrument Set




PDF 文件
PDF File



股骨倒置带锁髓内钉瞄准装置
Femoral Reverse Intramedullary Nail Aiming Device

D



骨针、钛缆、髓
Skeletal Pins, Tit
Patella, Cannu

固定示例

鹰嘴骨折示例



Kirschner Wire (Sterile Packaging)



D

系统特点
Character Of The System

- 国家三类注册证，植入更安全。
- 钛合金材质，生物相容性好。
- 锋利三角尖，植入更方便。
- National Class Iii Registration Certificate, more Safe To Implant.
- Titanium Alloy Material, good Biocompatibility.
- Sharp Triangle Tip, more Convenient To Implant.

金属骨针（双头尖） T896
Metal Skeletal Pins(Double End Tips)




产品名称 Product Name	包装形式 Packaging Form	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product NO.
金属骨针（双头尖） Metal Skeletal Pins (Double End Tips)	无菌包装 Sterile Packaging	钛合金 Titanium Alloy	0.8	150	T89608150-S
			1.0	150	T89610150-S
			1.2	150	T89612150-S
			1.5	150	T89615150-S
			1.8	150	T89618150-S
			2.0	150	T89618150-S
			2.5	150	T89618150-S

金属骨针 T897
Metal Skeletal Pins



产品名称 Product Name	包装形式 Packaging Form	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product NO.
金属骨针 Metal Skeletal Pins	无菌包装 Sterile Packaging	钛合金 Titanium Alloy	2.0	250	T89720250-S
			2.5	250	T89725250-S
			3.0	250	T89730250-S
			3.5	250	T89735250-S
			4.0	250	T89740250-S

无菌包装
Sterile Packaging



Titanium Cable Fixation System (Sterile Packaging)



系统特点

Character Of The System

- 由多根钛合金丝捻制而成，柔韧性好，具有出众的静态强度和更高的疲劳强度。
Made Of Multiple Twisted Titanilum Alloy Wire,good Flexibility,outstanding Static Sterength And Higher Fatigye Strength.
- 钛缆具有很强的抗应变能力，可以理想的拉紧和固定。
Titanium Cable Has Strong Strain Resistance Capability,ideal For Tensioning And Fixing.
- 锁扣机制锁定，无需扭转钛缆。
Fasten Buckle Mechanism Locking,no Need Of Twisting Wire Cable.
- 通过收紧器可精确控制张力。
Control Tension Precisely With Tightener.




钛缆固定系统安装器械包 91729902

Titanium Cable Fixation System Instrument Set



PDF 文件
PDF File



e Packaging)

产品：普通无针缆索（Φ1.3）1根

穿针4支

环形锁紧扣（枕形锁紧扣）

配套工具：钻头

多功能收紧器（或缆索双向收紧器+锁紧扣夹紧钳）

缆索剪

趾骨骨折示例



产品名称 Product Name	包装形式 Packaging Form	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
金属缆锁 Metal Titanium Cable	无菌包装 Sterile Packaging	钛合金 Titanium Alloy	1.3	520	T7713520-S
			1.8	520	T7718520-S

产品名称 Product Name	包装形式 Packaging Form	材质 Material	直径(mm) Diameter	产品编号 Product No.
锁紧扣 Fasten Buckle	无菌包装 Sterile Packaging	钛合金 Titanium Alloy	1.4	T7721408-S
			1.9	T7721908-S

产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
缆索填塞帽 Cable filling cap	钛合金 Titanium Alloy	3.5	T4863501
		5.0	T4864500

产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
锁紧接骨螺钉（缆索） 填塞型 Locking Screw(Cable) Filling Type	钛合金 Titanium Alloy	5.0	33072015

产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 股骨假体翻修锁定钢板(直型) JY Femoral Prosthesis Revision Locking Plate (Straight)	纯钛 Pure Titanium	10孔	257	23013100
		12孔	301	23013120
备注：选用Φ5.0锁定螺钉(阳极氧化)（图例为10孔） Note: Used For Φ5.0 Locking Screw(Anodized) (Sample: 10 Holes)				

产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 股骨假体翻修锁定钢板 JY Femoral Prosthesis Revision Locking Plate	纯钛 Pure Titanium	5孔	125	23233050
		14孔左	271	23233141
		14孔右	271	23233142
		16孔左	297	23233161
		16孔右	297	23233162
备注：选用Φ5.0锁定螺钉(阳极氧化)（图例为14孔） Note: Used For Φ5.0 Locking Screw(Anodized) (Sample: 14 Holes)				


髌骨固定系列

Patella Fixation Series

WX 2.7 髌骨锁定板
VA 2.7 Patella Locking Plate



WX 2.7 髌骨锁定板 23572、23571、23570
VA 2.7 Patella Locking Plate



Φ2.7 万向锁定螺钉(阳极氧化) T48027
Φ2.7 VA Locking Screw(Anodized)



91150119 锁定式直型手柄(微型) (D型)
Locking type straight handle(Micro)(Type D)



91132301 快装起子 (D型快装)
Quick Coupling Screwdriver(Type D Quick Coupling), T2.3



钩形张力带锁定接骨板

顶紧式锁紧扣 (可根据需要选择)

配套工具: 缆索单向收紧器 (加附件头)

缆索剪

大、中、小号引线器

内六角扳手 (SW2.5)

产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
WX 2.7 髌骨锁定板 VA 2.7 Patella Locking Plate	纯钛 Pure Titanium	大号 Large	23572160
		中号 Middle	23571160
		小号 Small	23570160

产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
万向锁定螺钉 (阳极氧化) VA Locking Screw(Anodized)	钛合金 Titanium Alloy	2.7	10、12、14、16、 18、20、22、24、 26、28、	T4802710- T4802728
			30	T4802730
			35	T4802735
			40	T4802740

接骨螺钉1支

穿针2支（不使用接骨螺钉情况下可3支）

环形锁紧扣（枕形锁紧扣）

配套工具：钻头

多功能收紧器（或缆索双向收紧器+锁紧扣夹紧钳）

缆索剪

产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
张力带别针 Tension Band Pins	钛合金 Titanium Alloy	2.0	60	T5220060
			70	T5220070
			80	T5220080
			100	T5220100
			120	T5220120
		2.5	60	T5225060
			70	T5225070
			80	T5225080
			100	T5225100
			120	T5225120



鹰嘴骨折



髌骨骨折



远端桡骨骨折

金属缆索内固定系统/06



91222001 髌骨钳
Patella Forceps



91133402 梅花起子
Torx Screwdriver, T3.4



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
髌骨爪 Patella Hook	纯钛 Pure Titanium	中号 Mid	T3000002
		小号 Small	T3000003
		特小号 Smaller	T3000004
		特小号II型 Smaller-II	T3000005
髌骨爪(无菌包装) Patella Hook (Sterile Packaging)	纯钛 Pure Titanium	中号 Mid	14301011
		小号 Small	14301012
		特小号 Smaller	14301013
		特小号II型 Smaller-II	14301014

产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
髌骨钳 Patella Forceps	/	/	91222001

产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
梅花起子 Torx Screwdriver	/	T3.4	91133402

Headless Cannulated Compression Screw System

D



- 中空设计，提高经皮植入的准确性。
- 变螺距螺纹和锥度设计，螺距有大变小，头部进入骨质的速度比尾部大，从而螺钉拧入时，对骨折面进行加压。
- 全螺纹设计，在受到周期性负载时能更好地保持加压功能。
- 无头设计，在关节区域使用中，减少碰撞以及对软组织的刺激。
- Hollow Design, Improve The Accuracy Of The Percutaneous Implant.
- Variable Pitch Thread And Taper Design Pitch From Big To Small, The Speed Of The Head Into The Bone Quicker Than The Speed Of The Tail, When The Screw Is Implanted, The Fracture Surface Is Well Pressured.
- Full Threaded Design, Better Keep The Pressure Function When Receiving Periodic Load.
- Headless Design, To Reduce Collisions And Stimulation To The Soft Tissue In The Joint Area In Use.

Headless Cannulated Compression Screw System

Φ2.5 中空无头加压螺钉 T57025

Φ2.5 Headless Cannulated Compression Screw



Φ3.5 中空无头加压螺钉 T57035

Φ3.5 Headless Cannulated Compression Screw



Φ4.5 中空无头加压螺钉 T57045

Φ4.5 Headless Cannulated Compression Screw



Φ6.5 中空无头加压螺钉 33505

Φ6.5 Headless Cannulated Compression Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
中空无头加压螺钉 Headless Cannulated Compression Screw	钛合金 Titanium Alloy	Φ2.5	8、10、12、 14、16、18、 20、22、24、 26、28、30	T5702508- T5702530
		Φ3.5	20、22、24、 26、28、30、 32、34、36、 38、40	T5703520- T5703540
		Φ4.5	30、32、34、 36、38、40、 42、44、46、 48、50、52、 54、56、58、 60	T5704530- T5704560
		Φ6.5	40、45、50、 55、60、65、 70、75、80、 85、90、95、 100、105、110、 115、120	33505040- 33505120

Φ2.5 无头加压空心钉安装器械包 91619902

Φ2.5 Headless Compression Cannulated Screw Instrument Set



PDF 文件
PDF File



托盘/Tray 91611101

Φ3.5 无头加压空心钉安装器械包 91629902

Φ3.5 Headless Compression Cannulated Screw Instrument Set



PDF 文件
PDF File



托盘/Tray 91621101

Φ4.5 无头加压空心钉安装器械包 91639902

Φ4.5 Headless Compression Cannulated Screw Instrument Set



PDF 文件
PDF File



托盘/Tray 91631101

Φ6.5 无头加压空心钉安装器械包 91649901

Φ6.5 Headless Compression Cannulated Screw Instrument Set



PDF 文件
PDF File



托盘/Tray 91640501

Counter Sink Cannulated Compression Screw System

Φ2.5 中空双头加压螺钉阳极氧化 33507
Φ2.5 Double-headed Cannulated Compression Screw



Φ3.0 中空双头加压螺钉 T5503
Φ3.0 Double-headed Cannulated Compression Screw



Φ4.0 中空双头加压螺钉 T5504
Φ4.0 Double-headed Cannulated Compression Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
中空双头 加压螺钉 Double-headed Cannulated Compression Screw	钛合金 Titanium Alloy	Φ2.5	10	33507010
			12	33507012
			14	33507014
			16	33507016
			18	33507018
			20	33507020
			22	33507022
			24	33507024
			26	33507026
			28	33507028
中空双头 加压螺钉 Double-headed Cannulated Compression Screw	钛合金 Titanium Alloy	Φ3.0	30	33507030
			14	T5503014
			16	T5503016
			18	T5503018
			20	T5503020
			22	T5503022
			24	T5503024
			26	T5503026
			28	T5503028
			30	T5503030
中空双头 加压螺钉 Double-headed Cannulated Compression Screw	钛合金 Titanium Alloy	Φ4.0	32	T5503032
			34	T5503034
			36	T5503036
			38	T5503038
			40	T5503040
			26	T5504026
			28	T5504028
			30	T5504030
			32	T5504032
			34	T5504034
中空双头 加压螺钉 Double-headed Cannulated Compression Screw	钛合金 Titanium Alloy	Φ4.0	36	T5504036
			38	T5504038
			40	T5504040
			42	T5504042
			44	T5504044
			46	T5504046
			48	T5504048
			50	T5504050

Φ2.5 双头加压空心钉安装器械包 91519901
Φ2.5 Double-End Compression Cannulated Screw Instrument Set



PDF 文件
PDF File



托盘/Tray 91510601

Φ3.0 双头加压空心钉安装器械包 91529901
Φ3.0 Double-End Compression Cannulated Screw Instrument Set



PDF 文件
PDF File



托盘/Tray 91520601

Φ4.0 双头加压空心钉安装器械包 91539901
Φ4.0 Double-End Compression Cannulated Screw Instrument Set



PDF 文件
PDF File



托盘/Tray 91530601

Cannulated Compression Screw Series

中空加压螺钉 (Φ2.0) 33500
Cannulated Compression Screw (Φ2.0)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
中空加压螺钉 Cannulated Compression Screw	钛合金 Titanium Alloy	Φ2.0	6、8、10、12、 14、16、18、 20、22	33500006- 33500022

中空加压螺钉 (Φ3.5) T551
Cannulated Compression Screw (Φ3.5)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
中空加压螺钉 Cannulated Compression Screw	钛合金 Titanium Alloy	Φ3.5	20、22、24、 26、28、30、 32、34、36、 38、40	T5513520- T5513540

中空加压螺钉 (Φ4.0) T552
Cannulated Compression Screw (Φ4.0)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
中空加压螺钉 Cannulated Compression Screw	钛合金 Titanium Alloy	Φ4.0	26、28、30、 32、34、36、 38、40、42、 44、46、48、 50	T5524026- T5524050

中空加压螺钉 (Φ4.5) T553
Cannulated Compression Screw (Φ4.5)




产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
中空加压螺钉 Cannulated Compression Screw	钛合金 Titanium Alloy	Φ4.5	24、26、28、 30、32、34、 36、38、40、 42、44、46、 48、50、52、 54、56、58、 60	T5534524- T5534560

Φ3.5 空心钉安装器械包 91449902
Φ3.5 Cannulated Screw Instrument Set



PDF 文件
PDF File



Φ4.0 空心钉安装器械包 91439902
Φ4.0 Cannulated Screw Instrument Set




PDF 文件
PDF File



Φ4.5 空心钉安装器械包 91429902
Φ4.5 Cannulated Screw Instrument Set



PDF 文件
PDF File



中空加压螺钉系列

Cannulated Compression Screw Series

股骨颈空心加压螺钉 (Φ6.5) T562

Femoral Neck Cannulated Compression Screw (Φ6.5)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
股骨颈空心 加压螺钉 Femoral Neck Cannulated Compression Screw	钛合金 Titanium Alloy	Φ6.5	60、65、70、 75、80、85、 90、95、100、 105	T56265060- T56265105
备注: 适用于股骨颈偏小患者 Note: Apply To The Patient Of Smaller Femur				

股骨颈空心加压螺钉 (Φ7.3) T563

Femoral Neck Cannulated Compression Screw (Φ7.3)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
股骨颈空心 加压螺钉 Femoral Neck Cannulated Compression Screw	钛合金 Titanium Alloy	Φ7.3	60、65、70、 75、80、85、 90、95、100、 105、110、115	T56373060- T56373115

双头空心加压螺钉 (Φ6.5) T561

Double-headed Cannulated Compression Screw (Φ6.5)




产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
双头空心加压螺钉 Double-headed Cannulated Compression Screw	钛合金 Titanium Alloy	Φ6.5	65、70、75、 80、85、90、 95、100、105、 110、115	T56165065- T56165115

Φ7.3 空心钉安装器械包 91419901

Φ7.3 Cannulated Screw Instrument Set



PDF 文件
PDF File



自钻自攻螺钉 T5012

Self-tapping & Drilling Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
自钻自攻螺钉 Self-tapping & Drilling Screw	钛合金 Titanium Alloy	Φ2.0	11、12、13、 14、15、16、 18、20	T5012011- T5012020
		Φ2.5	11、12、13、 14、15、16、 18、20、22	T5012511- T5012522

Φ2.0 / Φ2.5 自钻自攻螺钉安装器械包 91249901



PDF 文件
PDF File





E

微型锁定指掌骨系列
Micro Phalangeal Metacarpal
Locking Series

TML 1.5 微型锁定指掌骨系列

TML 1.5 Micro Phalangeal Metacarpal Locking Series

TML 快装起子杆
TML Quick Coupling Screw
Screwdriver Rod



螺钉可持功能
Screw holding

梅花螺钉防打滑 T1.3
Torx screw anti-slipping T1.3



TML 1.5 微型锁定钉
TML 1.5 Micro Locking Screw

自攻型刃口
Self tapping Cutting edge
 $\Phi 1.5\text{mm}$

孔距4.5mm
Hole spacing 4.5mm



板宽4.0mm
plate width 4.0mm



板厚1.0mm
plate thickness 1.0mm

TML 1.5 微型锁定指掌骨安装器械包 91859901

TML 1.5 Micro phalangeal metacarpal locking series instrument Set



PDF 文件
PDF File




91850701



91850801




91850901



91851001



91851101



91851201

TML 1.5 微型锁定指掌骨系列

TML 1.5 Micro Phalangeal Metacarpal Locking Series

TML 1.5 微型直型锁定板 235811
TML 1.5 Micro Straight Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度 Holes	产品编号 Product No.
TML 1.5 微型直型锁定板 TML 1.5 Micro Straight Locking Plate	纯钛 Pure Titanium	5孔	22	23581105
		7孔	31	23581107
		9孔	40	23581109
		11孔	49	23581111
备注: TML 1.5锁定螺钉(阳极氧化) (图例为9孔) Note: Used For TML Locking Screw (Anodized) 、(Sample: 9 Holes)				

TML 1.5 微型T型锁定 头部2孔 235812
TML 1.5 Micro T Locking Plate (2 holes in the head)



产品名称 Product Name	材质 Material	孔数 Holes	长度 Holes	产品编号 Product No.
TML 1.5 微型T型锁定钢板 头部2孔 TML Micro T Locking Plate (2 holes in the head)	纯钛 Pure Titanium	5孔	19	23581205
		6孔	24	23581206
		7孔	28	23581207
		8孔	33	23581208
		9孔	37	23581209
备注：选用TML 1.5锁定螺钉(阳极氧化)（图例为6孔） Note: Used For TML 1.5 Locking Screw (Anodized) (Sample: 6 Holes)				

TML 1.5 微型髌部锁定板 235813
TML 1.5 Micro Condylar Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度 Holes	产品编号 Product No.
TML 1.5 微型髌部锁定板 TML 1.5 Micro Condylar Locking Plate	纯钛 Pure Titanium	5孔	19	23581305
		6孔	24	23581306
		7孔	28	23581307
		8孔	33	23581308
		9孔	37	23581309
备注：选用TML 1.5锁定螺钉(阳极氧化)（图例为8孔） Note: Used For TML 1.5 Locking Screw (Anodized)(Sample: 8 Holes)				

TML 1.5 微型L型锁定板 235814
TML Micro L Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度 Holes	产品编号 Product No.
TML 1.5 微型L型锁定板 TML 1.5 Micro L Locking Plate	纯钛 Pure Titanium	5孔左	18	23581451
		5孔右	18	23581452
		6孔左	23	23581461
		6孔右	23	23581462
		7孔左	27	23581471
		7孔右	27	23581472
		8孔左	32	23581481
		8孔右	32	23581482
备注：选用TML 1.5 锁定螺钉(阳极氧化) (图例为7孔) Note: Used For TML 1.5 Locking Screw (Anodized) (Sample: 7 Holes)				

TML 1.5 微型锁定钉 (阳极氧化) 33020
TML 1.5 Micro Locking Screw (Anodized)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
TML 1.5 微型锁定钉 (阳极氧化) TML 1.5 Micro Locking Screw(Anodized)	钛合金 Titanium Alloy	1.5	6	33020006
			7	33020007
			8	33020008
			10	33020010
			12	33020012
			14	33020014


TML 1.5 微型螺钉 (阳极氧化) T5201
TML 1.5 Micro Screw (Anodized)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
TML 1.5 微型螺钉 (阳极氧化) TML 1.5 Micro Screw (Anodized)	钛合金 Titanium Alloy	1.5	5	T5201505
			6	T5201506
			7	T5201507
			8	T5201508
			10	T5201510
			12	T5201512
			14	T5201514

TML 2.0 微型锁定指掌骨系列

TML 2.0 Micro Phalangeal Metacarpal Locking Series





TML 2.0 微型锁定指掌骨安装器械包 91829901

TML 2.0 Micro phalangeal metacarpal locking series instrument Set



PDF 文件
PDF File



91821301



91821401




91821501



91821601



91270302



91821801

TML 2.0 微型锁定指掌骨系列

TML 2.0 Micro Phalangeal Metacarpal Locking Series

TML 2.0 微型直型锁定板 235801
TML 2.0 Micro Straight Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度 Holes	产品编号 Product No.
TML 2.0 微型直型锁定板 TML 2.0 Micro Straight Locking Plate	纯钛 Pure Titanium	5孔	25	23580105
		7孔	35	23580107
		9孔	45	23580109
		11孔	55	23580111
备注：TML 2.0锁定螺钉(阳极氧化)（图例为9孔） Note: Used For TML 2.0 Locking Screw (Anodized) 、(Sample: 9 Holes)				

TML 2.0 微型T型锁定 头部2孔 235802
TML Micro T Locking Plate (2 holes in the head)



产品名称 Product Name	材质 Material	孔数 Holes	长度 Holes	产品编号 Product No
TML 2.0 微型T型锁定钢板 头部2孔 TML Micro T Locking Plate (2 holes in the head)	纯钛 Pure Titanium	5孔	22	23580205
		6孔	27	23580206
		7孔	32	23580207
		8孔	37	23580208
		9孔	42	23580209
备注：选用TML 2.0锁定螺钉(阳极氧化)（图例为6孔） Note: Used For TML 2.0 Locking Screw (Anodized) (Sample: 6 Holes)				

TML 2.0 微型髌部锁定板 235803
TML Micro Condylar Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度 Holes	产品编号 Product No
TML 2.0 微型髌部锁定板 TML Micro Condylar Locking Plate	纯钛 Pure Titanium	5孔	21	23580305
		6孔	26	23580306
		7孔	31	23580307
		8孔	36	23580308
		9孔	41	23580309
备注：选用TML 2.0锁定螺钉(阳极氧化)（图例为8孔） Note: Used For TML 2.0 Locking Screw (Anodized)(Sample: 8 Holes)				

TML 2.0 微型L型锁定板 235804
TML Micro L Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度 Holes	产品编号 Product No.
TML 2.0 微型L型锁定板 TML Micro L Locking Plate	纯钛 Pure Titanium	5孔左	21	23580451
		5孔右	21	23580452
		6孔左	26	23580461
		6孔右	26	23580462
		7孔左	31	23580471
		7孔右	31	23580472
		8孔左	36	23580481
		8孔右	36	23580482
备注： 选用TML 2.0锁定螺钉(阳极氧化)（图例为7孔） Note: Used For TML2.0 Locking Screw (Anodized) (Sample: 7 Holes)				

TML 2.0 微型锁定钉 (阳极氧化) 33010
TML 2.0 Micro Locking Screw (Anodized)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
TML 2.0 微型锁定钉 (阳极氧化) TML 2.0 Micro Locking Screw	钛合金 Titanium Alloy	2.0	6、8、10、12、 14、16、18、20	33010006- 33010020

TML 2.0 微型螺钉 (阳极氧化) T5202
TML 2.0 Micro Screw (Anodized)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
TML 2.0 微型螺钉 (阳极氧化) TML 2.0 Micro Screw (Anodized)	钛合金 Titanium Alloy	2.0	6、8、10、12、 14、16、18、20	T5202006- T5202020

KML 2.0 Locking Plate Series

2.0 T型锁定钢板 头部2孔 T41201
2.0 T-Type Locking Plate(Two Holes In The Head)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.0 T型锁定钢板 头部2孔 2.0 T-Type Locking Plate (Two Holes In The Head)	纯钛 Pure Titanium	5/6/7/8/9	T4120105-T4120109
备选: 选用Φ2.0锁定螺钉 (三角型), HA2.0螺钉 (图例为7孔) Note:used For Φ2.0 Locking Screw (Triangle),HA2.0 Screw (Sample:7holes)			

2.0 T型锁定钢板 头部3孔 T41202
2.0 T-Type Locking Plate(Three Holes In The Head)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.0 T型锁定钢板 头部3孔 2.0 T-Type Locking Plate (Three Holes In The Head)	纯钛 Pure Titanium	6/7/8/9	T4120206-T4120209
备选: 选用Φ2.0锁定螺钉 (三角型), HA2.0螺钉 (图例为8孔) Note:used For Φ2.0 Locking Screw (Triangle),HA2.0 Screw (Sample:8holes)			

2.0 踝部锁定钢板 T41204
2.0 Ankle Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.0 踝部锁定钢板 2.0 Ankle Locking Plate	纯钛 Pure Titanium	5/6/7/8/9	T4120405-T4120409
备选: 选用Φ2.0锁定螺钉 (三角型), HA2.0螺钉 (图例为9孔) Note:used For Φ2.0 Locking Screw (Triangle),HA2.0 Screw (Sample:9holes)			

2.0 L型锁定钢板 (左右) T41206
2.0 L-Type Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.0 L型锁定钢板 (左右) 2.0 L-Type Locking Plate (Left,Right)	纯钛 Pure Titanium	5/6/7/8	T4120605-T4120608R/L
备选: 选用Φ2.0锁定螺钉 (三角型), HA2.0螺钉 (图例为8孔) Note:used For Φ2.0 Locking Screw (Triangle),HA2.0 Screw (Sample:8holes)			

2.0 直型锁定钢板 T41200
2.0 Straight Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.0 直型锁定钢板 2.0 Straight Locking Plate	纯钛 Pure Titanium	4/5/6/7/8/9/10	T4120004-T4120010
备选: 选用Φ2.0锁定螺钉 (三角型) (图例为6孔) Note:used For Φ2.0 Locking Screw (Triangle) (Sample:6holes)			

三角形锁定螺钉 T480
Triangle Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
三角形锁定螺钉 Triangle Locking Screw	钛合金 Titanium Alloy	2.0	6、8、10、12、 14、16、18、20、 22	T4802006- T4802022

金属接骨螺钉 (三角形) T541
Metal Bone Fracture Screw (Triangle)



钉帽形式 Head of Nail Type	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
三角形 Triangle	钛合金 Titanium Alloy	2.0	6、8、10、12、 14、16、18、20、 22	T5412006- T5412022

KML 锁定钢板2.5系列

KML 2.5 Locking Plate Series

2.5 T型锁定钢板 头部2孔 T41251
2.5 T-Type Locking Plate (Two Holes In The Head)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.5 T型锁定钢板 头部2孔 2.5 T-Type Locking Plate (Two Holes In The Head)	纯钛 Pure Titanium	5/6/7/8/9	T4125105-T4125109
备选：选用Φ2.5锁定螺钉，HA2.5螺钉（图例为7孔） Note:used For Φ2.5 Locking Screw,HA2.5 Screw (Sample:7holes)			

2.5 T型锁定钢板 头部3孔 T41252
2.5 T-Type Locking Plate (Three Holes In The Head)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.5 T型锁定钢板 头部3孔 2.5 T-Type Locking Plate (Three Holes In The Head)	纯钛 Pure Titanium	6/7/8/9/10	T4125206-T4125210
备选：选用Φ2.5锁定螺钉，HA2.5螺钉（图例为8孔） Note:used For Φ2.5 Locking Screw,HA2.5 Screw(Sample:8holes)			

2.5 踝部锁定钢板 T41254
2.5 Ankle Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.5 踝部锁定钢板 2.5 Ankle Locking Plate	纯钛 Pure Titanium	5/6/7/8/9	T4125405-T4125409
备选：选用Φ2.5锁定螺钉，HA2.5螺钉（图例为8孔） Note:used For Φ2.5 Locking Screw,HA2.5 Screw (Sample:8holes)			

2.5 L型锁定钢板（左右） T41256
2.5 L-Type Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.5 L型锁定钢板（左右） 2.5 L-Type Locking Plate (Left,Right)	纯钛 Pure Titanium	5/6/7/8/9	T4125605-T4125609R/L
备选：选用Φ2.5锁定螺钉，HA2.5螺钉（图例为8孔） Note:used For Φ2.5 Locking Screw,HA2.5 Screw (Sample:8holes)			

2.5 直型锁定钢板 T41250
2.5 Straight Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.5 直型锁定钢板 2.5 Straight Locking Plate	纯钛 Pure Titanium	4/5/6/7/8/9/10	T4125004-T4125010
备选：选用Φ2.5锁定螺钉（图例为6孔） Note:used For Φ2.5 Locking Screw (Sample:6holes)			

梅花型锁定螺钉 T481/T482
Torx Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型锁定螺钉 Torx Locking Screw	钛合金 Titanium Alloy	2.5	6、8、10、12、 14、16、18、20、 22、24、26、28、 30	T4812506- T4812530



JY、WX 锁定钢板螺钉系列
JY、WX Locking Plate Screw Series

JY. WX Locking Plate Series

WX 2.7 桡骨远端锁定钢板 VA 2.7 Distal Radius Locking Plate



视频在线播放
Video online play

专用螺钉可用于桡骨茎突
固定和支撑月骨关节面和
下尺桡关节

Special Screw Is Used
For The Fixation Of Radial
Styloid Process And Support
Of Lunate Articular Surface
And Distal Ulna Radial Articular.

克氏针孔用于接骨板的临时固定
High Kirschner Wire Hole Is Used
For The Temporary Fixation Of Bone Plate.

解剖型设计
高抛光的平滑边缘可靠近掌侧嵴，
减少软组织刺激



High Polished Design High
Polished And Smooth Edge
Can Be Near To Volar Crest,
Reduce Soft Tissue Irritation.

双柱固定，增加稳定性
High Double Post Fixation,
Enhance Stability.

加长结合孔用于接骨板位置调整
Lengthened Combination
Hole Is Used For The Bone
Plate Position Adjustment.



2.1 万向套筒
2.1 VA Sleeve



头部万向结构提供30°旋转，
有效固定粉碎骨折块。
Universal Head Structure Provides
30° Rotation, Fix The Comminuted
Fracture Fragments Effectively.

JY. WX Locking Plate Series

WX 2.7 桡骨远端锁定钢板(小号)(左右) T4326
VA 2.7 Distal Radius Locking Plate(Mini)(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
WX 2.7 桡骨远端锁定 钢板(小号)(左右) VA 2.7 Distal Radius Locking Plate(Mini)(Left,Right)	纯钛 Pure Titanium	3	55	T4326003R/L
		4	67	T4326004R/L
		5	76	T4326005R/L
备注: 选用Φ2.7万向锁定螺钉(阳极氧化)、Φ2.7锁定螺钉(阳极氧化)、HA2.7螺钉 (图例为4孔) Note: Used For Φ2.7 Variable Angle Locking Screw(Anodized)、 Φ2.7 Locking Screw(Anodized)、HA2.7 Screw (Sample: 4 Holes)				

WX 2.7 桡骨远端锁定钢板(大号)(左右) T4327
VA 2.7 Distal Radius Locking Plate(Large)(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
WX 2.7 桡骨远端锁定 钢板(大号)(左右) VA 2.7 Distal Radius Locking Plate(Large)(Left,Right)	纯钛 Pure Titanium	3	56	T4327003R/L
		4	68	T4327004R/L
		5	77	T4327005R/L
备注: 选用Φ2.7万向锁定螺钉(阳极氧化)、Φ2.7锁定螺钉(阳极氧化)、HA2.7螺钉 (图例为4孔 Note: Used For Φ2.7 Variable Angle Locking Screw(Anodized)、 Φ2.7 Locking Screw(Anodized)、HA2.7 Screw (Sample: 4 Holes)				

F

无菌 WX 2.7 桡骨远端锁定导向板 25070
VA 2.7 Distal Radius Locking Guide Plate (Sterile)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
无菌 WX 2.7 桡骨远端锁定导向板 VA 2.7 Distal Radius Locking Guide Plate (Sterile)	纯钛 Pure Titanium	4孔左	57	25070041
		4孔右	57	25070042
		5孔左	65	25070051
		5孔右	65	25070052
		7孔左	81	25070071
		7孔右	81	25070072
		13孔左	129	25070131
		13孔右	129	25070132
		17孔左	160	25070171
		17孔右	160	25070172
备注：选用Φ2.7万向锁定螺钉(阳极氧化)、Φ2.7锁定螺钉(阳极氧化)、HA2.7螺钉（图例为4孔 Note: Used For Φ2.7 Variable Angle Locking Screw(Anodized)、 Φ2.7 Locking Screw(Anodized)、HA2.7 Screw (Sample:4 Holes)				

JY 2.7 桡骨远端支持锁定钢板小号 (左右) 23053
JY 2.7 Distal Radius Buttress Locking Plate (Mini)
(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 2.7 桡骨远端支持锁定 钢板小号 (左右) JY 2.7 Distal Radius Buttress Locking Plate(Mini)(Left,Right)	纯钛 Pure Titanium	4孔左	55	23053041
		4孔右	55	23053042
		5孔左	63	23053051
		5孔右	63	23053052
		7孔左	79	23053071
		7孔右	79	23053072
		9孔左	95	23053091
		9孔右	95	23053092
		11孔左	111	23053111
		11孔右	111	23053112
		13孔左	127	23053131
		13孔右	127	23053132
备注：选用Φ2.7锁定螺钉(阳极氧化)、HA2.7螺钉 (图例为4孔) Note: Used For Φ2.7 Locking Screw(Anodized)、HA2.7 Screw (Sample: 4 Holes)				

JY. WX Locking Plate Series

JY 2.7 桡骨远端支持锁定钢板大号 (左右) 23054
JY 2.7 Distal Radius Buttress Locking Plate (Large)
(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 2.7 桡骨远端支持锁定 钢板大号 (左右) JY 2.7 Distal Radius Buttress Locking Plate(Large)(Left,Right)	纯钛 Pure Titanium	4孔左	55	23054041
		4孔右	55	23054042
		5孔左	63	23054051
		5孔右	63	23054052
		7孔左	79	23054071
		7孔右	79	23054072
		9孔左	95	23054091
		9孔右	95	23054092
备注：选用Φ2.7锁定螺钉(阳极氧化)、HA2.7螺钉 (图例为5孔) Note: Used For Φ2.7 Locking Screw(Anodized)、HA2.7 Screw (Sample: 5 Holes)				

JY 桡骨远端背侧锁定钢板 (左右) 23055
JY Distal Back Radius Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
JY 桡骨远端背侧锁定钢板（左右） JY Distal Back Radius Locking Plate (Left,Right)	纯钛 Pure Titanium	14孔左	23055141
		14孔右	23055142
备注：选用Φ2.7锁定螺钉(阳极氧化) Note: Used For Φ2.7 Locking Screw(Anodized)			

JY 2.7 桡骨头锁定钢板 (左右) 23056
JY 2.7 Radius Head Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 2.7 桡骨头锁定钢板 (左右) JY 2.7 Radius Head Locking Plate(Left,Right)	纯钛 Pure Titanium	3孔左	40	23056031
		3孔右	40	23056032
		4孔左	48	23056041
		4孔右	48	23056042
		5孔左	56	23056051
		5孔右	56	23056052
备注: 选用Φ2.7锁定螺钉(阳极氧化) (图例为4孔) Note: Used For Φ2.7 Locking Screw(Anodized)(Sample: 4 Holes)				

WX 桡骨远端中间柱锁定钢板 (左右) 23467
VA Distal Radius Intermediate Column Locking Plate
(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
WX 桡骨远端中间柱锁定钢板 (左右) VA Distal Radius Intermediate Column Locking Plate (Left,Right)	纯钛 Pure Titanium	6孔左	50	23467061
		6孔右	50	23467062
备注: 选用Φ2.7万向锁定螺钉(阳极氧化)、HA2.7螺钉 Note: Used For Φ2.7 Variable Angle Locking Screw(Anodized)、HA2.7 Screw				

WX 桡骨远端桡侧柱锁定钢板 23468
VA Distal Radial Column Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
WX 桡骨远端桡侧柱锁定钢板 VA Distal Radial Column Locking Plate	纯钛 Pure Titanium	5孔	50	23468050
		6孔	55	23468060
备注：选用Φ2.7万向锁定螺钉(阳极氧化)、HA2.7螺钉(图例为6孔) Note: Used For Φ2.7 Variable Angle Locking Screw(Anodized)、HA2.7 Screw (Sample: 6 Holes)				

JY. WX Locking Plate Series

JY 尺骨远端锁定钢板 23080

JY Distal Ulna Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 尺骨远端锁定钢板 JY Distal Ulna Locking Plate	纯钛 Pure Titanium	8孔	45	23080080
备注：选用Φ2.0锁定螺钉(阳极氧化)、HA2.0螺钉 Note: Used For Φ2.0 Locking Screw(Anodized)、HA2.0 Screw				

JY 尺骨近端锁定钢板（左右） 23081

JY Proximal Ulna Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 尺骨近端锁定钢板（左右） JY Proximal Ulna Locking Plate （Left,Right）	纯钛 Pure Titanium	4孔左	112	23081041
		4孔右	112	23081042
		6孔左	140	23081061
		6孔右	140	23081062
		8孔左	168	23081081
		8孔右	168	23081082
备注：选用Φ3.5锁定螺钉(阳极氧化)，HA3.5螺钉(阳极氧化)（图例为6孔） Note: Used For Φ3.5 Locking Screw(Anodized)，HA3.5 Screw(Anodized)(Sample: 6 Holes)				

WX 尺骨近端锁定加压钩板 23474

VA Proximal Ulna Locking Compression Hook Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
WX 尺骨近端锁定加压钩板 VA Proximal Ulna Locking Compression Hook Plate	纯钛 Pure Titanium	5孔	64	23474050
备注：选用Φ3.5万向锁定螺钉(阳极氧化)、Φ3.5锁定螺钉(阳极氧化)、HA3.5螺钉(阳极氧化) Note: Used For Φ3.5 Variable Angle Locking Screw(Anodized)、Φ3.5 Locking Screw (Anodized)、HA3.5 Screw(Anodized)				

JY 肱骨大结节锁定钢板（左右） 23073

JY Humerus Tubercle Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 肱骨大结节锁定钢板（左右） JY Humerus Tubercle Locking Plate(Left,Right)	纯钛 Pure Titanium	5孔左	42	23073051
		5孔右	42	23073052
备注：选用Φ2.7、Φ3.5锁定螺钉（阳极氧化）、HA2.7螺钉（阳极氧化）、HB4.0螺钉 Note: Used For Φ2.7、Φ3.5 Locking Screw(Anodized)、HA2.7 Screw(Anodized)、 HB4.0 Screw				

JY 肱骨近端锁定钢板 23051

JY Proximal Humerus Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 肱骨近端锁定钢板 JY Proximal Humerus Locking Plate	纯钛 Pure Titanium	3孔	91	23051030
		5孔	115	23051050
备注：选用Φ3.5锁定螺钉(阳极氧化)、HA3.5螺钉（阳极氧化）、HB4.0螺钉（图例为5孔） Note: Used For Φ3.5 Locking Screw(Anodized)、HA3.5 Screw(Anodized)、HB4.0 Screw (Sample: 5 Holes)				

JY. WX Locking Plate Series

JY 肱骨近端锁定钢板 (长) 23052

JY Proximal Humerus Locking Plate (Long)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 肱骨近端锁定钢板（长） JY Proximal Humerus Locking Plate(Long)	纯钛 Pure Titanium	5孔	144	23052050
		7孔	180	23052070
		9孔	216	23052090
		11孔	252	23052110
		13孔	288	23052130
备注：选用Φ3.5锁定螺钉（阳极氧化）、HA3.5螺钉（阳极氧化）、HB4.0螺钉（图例为5孔） Note: Used For Φ3.5 Locking Screw(Anodized)、HA3.5 Screw(Anodized)、HB4.0 Screw (Sample: 5 Holes)				

JY 肱骨远端外侧锁定钢板 (左右) 23081

JY Distal Lateral Humerus Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 肱骨远端外侧锁定钢板（左右） JY Distal Lateral Humerus Locking Plate (Left,Right)	纯钛 Pure Titanium	5孔左	92	23081051
		5孔右	92	23081052
		7孔左	118	23081071
		7孔右	118	23081072
		9孔左	144	23081091
		9孔右	144	23081092
		11孔左	170	23081111
		11孔右	170	23081112
备注：选用Φ3.5锁定螺钉（阳极氧化）、HA3.5螺钉（阳极氧化）、HB4.0螺钉（图例为11孔） Note: Used For Φ3.5 Locking Screw(Anodized)、HA3.5 Screw(Anodized)、HR4.0screw (Sample: 11 Holes)				

JY 肱骨远端内侧锁定钢板 (左右) 23085

JY Distal Humerus Medial Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 锁定肱骨远端内侧锁定钢板 (左右) JY Distal Humerus Medial Locking Plate (Left,Right)	纯钛 Pure Titanium	4孔左	82	23085041
		4孔右	82	23085042
		5孔左	94	23085051
		5孔右	94	23085052
		7孔左	118	23085071
		7孔右	118	23085072
		9孔左	142	23085091
		9孔右	142	23085092
		11孔左	166	23085111
		11孔右	166	23085112
备注：选用Φ3.5锁定螺钉（阳极氧化）、Φ2.7锁定螺钉(阳极氧化)、HB4.0螺钉（图例为9孔） Note: Used For Φ3.5 Locking Screw(Anodized)、Φ2.7 Locking Screw(Anodized)、HB4.0screw (Sample: 9 Holes)				

JY 肱骨远端亚髁锁定钢板 (左右) 23074

JY Distal Humerus Subcondyle Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 肱骨远端亚髁锁定钢板 (左右) JY Distal Humerus Subcondyle Locking Plate(Left,Right)	纯钛 Pure Titanium	12孔左	142	23074121
		12孔右	142	23074122
		15孔左	188	23074151
		15孔右	188	23074152
备注: 选用Φ3.5锁定螺钉 (阳极氧化)、HA3.5螺钉 (阳极氧化) (图例为15孔) Note: Used For Φ3.5 Locking Screw(Anodized)、HA3.5 Screw(Anodized) (Sample:15 Holes)				

JY、WX 锁定钢板系列

JY. WX Locking Plate Series

JY 锁骨S型锁定低应力板 23083
JY S Clavicle Locking Plate (Anodized)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 锁骨S型锁定低应力板 JY S Clavicle Locking Plate (Anodized)	纯钛 Pure Titanium	7孔左	91	23083071
		7孔右	91	23083072
		8孔左	102	23083081
		8孔右	102	23083082
		9孔左	113	23083091
		9孔右	113	23083092
		10孔左	123	23083101
		10孔右	123	23083102
		12孔左	145	23083121
		12孔右	145	23083122
备注：选用Φ3.5锁定螺钉（阳极氧化）、HA3.5螺钉（阳极氧化）（图例为9孔） Note: Used For Φ3.5 Locking Screw(Anodized)、HA3.5 Screw (Anodized)(Sample: 9 Holes)				

JY 锁骨干锁定钢板（左右） 23071、23077
JY Clavicle Stem Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 锁骨干锁定钢板（左右） JY Clavicle Stem Locking Plate (Left,Right)	纯钛 Pure Titanium	短7孔左	90	23077071
		短7孔右	90	23077072
		7孔左	109	23071071
		7孔右	109	23071072
		8孔左	122	23071081
		8孔右	122	23071082
		9孔左	136	23071091
		9孔右	136	23071092
		11孔左	162	23071111
		11孔右	162	23071112
备注：选用Φ3.5锁定螺钉（阳极氧化）、HA3.5螺钉（阳极氧化）（图例为8孔） Note: Used For Φ3.5 Locking Screw(Anodized)、HA3.5 Screw (Anodized)(Sample: 8 Holes)				

JY 锁骨干端锁定低应力板（左右） 23084
JY Distal Clavicle Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 锁骨远端锁定低应力板 (左右) JY Distal Clavicle Locking Plate (Left,Right)	纯钛 Pure Titanium	5孔左	89	23084051
		5孔右	89	23084052
		6孔左	100	23084061
		6孔右	100	23084062
		8孔左	125	23084081
		8孔右	125	23084082
		10孔左	148	23084101
		10孔右	148	23084102
备注：选用Φ2.7、Φ3.5锁定螺钉（阳极氧化）、HA3.5螺钉（阳极氧化）（图例为6孔） Note: Used For Φ2.7、Φ3.5 Locking Screw(Anodized)、HA3.5 Screw (Anodized) (Sample: 6 Holes)				

JY 肩锁钩锁定钢板（左右） 23075
JY Clavicle Hook Locking Plate (Left,Right)



JY 肩锁钩锁定钢板（左右） 23079
JY Clavicle Hook Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	高度(mm) Height	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 肩锁钩锁定钢板 (左右) JY Clavicle Hook Locking Plate (Left,Right)	纯钛 Pure Titanium	14	8孔左	54	23075081
			8孔右	54	23075082
			9孔左	66	23075091
			9孔右	66	23075092
			11孔左	92	23075111
			11孔右	92	23075112
JY 肩锁钩锁定钢板 (左右) JY Clavicle Hook Locking Plate (Left,Right)	纯钛 Pure Titanium	17	8孔左	54	23079081
			8孔右	54	23079082
			9孔左	66	23079091
			9孔右	66	23079092
			11孔左	92	23079111
			11孔右	92	23079112
备注：选用Φ2.7、Φ3.5锁定螺钉（阳极氧化）、HA3.5螺钉（阳极氧化）（图例为11孔） Note: Used For Φ2.7、Φ3.5 Locking Screw(Anodized)、HA3.5screw (Anodized) (Sample:11 Holes)					

JY、WX 锁定钢板系列

JY. WX Locking Plate Series

JY 肩胛骨肩峰锁定钢板（左右） 23077
JY Scapula Acromion Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
JY 肩胛骨肩峰锁定钢板（左右） JY Scapula Acromion Locking Plate(Left,Right)	纯钛 Pure Titanium	9孔左	23077091
		9孔右	23077092
备注：选用Φ2.7锁定螺钉(阳极氧化) Note: Used For Φ2.7 Locking Screw (Anodized)			

JY 肩胛骨外侧缘锁定钢板（左右） 23078
JY Scapula Lateral Margin Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
JY 肩胛骨外侧缘锁定钢板（左右） JY Scapula Lateral Margin Locking Plate(Left,Right)	纯钛 Pure Titanium	9孔左	23078091
		9孔右	23078092
备注：选用Φ2.7锁定螺钉(阳极氧化) Note: Used For Φ2.7 Locking Screw (Anodized)			

JY 肩胛骨关节盂锁定钢板（左右） 23076
JY Scapula Glenoid Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
JY 肩胛骨关节盂锁定钢板（左右） JY Scapula Glenoid Locking Plate(Left,Right)	纯钛 Pure Titanium	4孔左	23076041
		4孔右	23076042
备注：选用Φ2.7锁定螺钉(阳极氧化) Note: Used For Φ2.7 Locking Screw (Anodized)			


JY 肩胛骨内侧缘锁定钢板（左右） 23082
JY Scapula Medial Margin Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
JY 肩胛骨内侧缘锁定钢板（左右） JY Scapula Medial Margin Locking Plate(Left,Right)	纯钛 Pure Titanium	10孔左	23082101
		10孔右	23082102
备注：选用Φ2.7锁定螺钉(阳极氧化) Note: Used For Φ2.7 Locking Screw (Anodized)			



视频在线播放
Video online play



F

LCC锁定复合加压连接板 23236

LCC Locking Composite Pressure Connecting Plate



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
LCC锁定复合加压连接板 LCC Locking Composite Pressure Connecting Plate	纯钛 Pure Titanium	3孔	23236030
		4孔	23236040

LCC 动力棒 33508

LCC Dynamic Rod



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
LCC 动力棒 LCC Dynamic Rod	钛合金 Titanium Alloy	10	75-130 (每5mm为一档) The specification is 5mm	33508065-33508120
LCC 抗旋螺钉 LCC Anti-rotation screw		6.5	75-130 (每5mm为一档) The specification is 5mm	33506060-33506115

LCC 抗旋螺钉 33506

LCC Anti-rotation screw



Φ5.0 梅花型锁定螺钉(阳极氧化) T493

Torx Locking Screw (Anodized)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型锁定螺钉 (阳极氧化) Torx Locking Screw (Anodized)	钛合金 Titanium Alloy	5.0	14-50 (每2mm为一档) The specification is 2mm	T4935014-T4935050
			55-100 (每5mm为一档) The specification is 5mm	T4935055-T4935100

LCC 锁定复合加压系统器械包 91749901

LCC Locking compound pressurization system Instrument Set (FNS)



PDF 文件
PDF File



JY、WX 锁定钢板系列

JY. WX Locking Plate Series



适应症

Indications

- 胫骨近端内侧开放式闭合楔型截骨术；
 - 内侧或外侧单髁膝关节病，胫骨近端对线不良；
 - 胫骨近端特发性或创伤性内翻或外翻畸形。
- Open and closed wedged Osteotomies of the medial proximal tibia.
 - Medial or lateral Unicompartmental osteoarthritis, malalignment of proximal tibia.
 - Idiopathic or traumatic varus or valgus of proximal tibia.



F

JY HTO 胫骨截骨标准型锁定钢板 232140

JY HTO Tibial Osteotomy Standard Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY HTO 胫骨截骨 标准型锁定钢板 JY HTO Tibial Osteotomy Standard Locking Plate	纯钛 Pure Titanium	4孔	115	23214040
备注：选用Φ5.0锁定螺钉(阳极氧化)、HA4.5螺钉（图例为4孔） Note: used for Φ5.0 locking screw(Anodized)、HA4.5 screw (Sample: 4 holes)				

JY HTO 胫骨截骨小型锁定钢板 232141

JY HTO Tibial Osteotomy Small Locking Plate





产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY HTO 胫骨截骨 小型锁定钢板 JY HTO Tibial Osteotomy Small Locking Plate	纯钛 Pure Titanium	4孔	112	23214140
备注：选用Φ5.0锁定螺钉(阳极氧化)、HA4.5螺钉（图例为4孔） Note: used for Φ5.0 locking screw(Anodized)、HA4.5 screw (Sample: 4 holes)				

HTO锁定胫骨截骨安装器械包 91379901

HTO Tibial Osteotomy Locking Plate Instrument Set



PDF 文件
PDF File




JY、WX 锁定钢板系列



JY. WX Locking Plate Series

JY 3.5 胫骨近端内侧锁定钢板（左右） 23206
JY 3.5 Proximal Medial Tibial Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 3.5 胫骨近端内侧锁定钢板（左右） JY 3.5 Proximal Medial Tibial Locking Plate (Left,Right)	纯钛 Pure Titanium	5孔左	106	23206051
		5孔右	106	23206052
		7孔左	134	23206071
		7孔右	134	23206072
		9孔左	162	23206091
		9孔右	162	23206092
		11孔左	190	23206111
		11孔右	190	23206112
		13孔左	218	23206131
		13孔右	218	23206132
备注：选用Φ3.5锁定螺钉（阳极氧化）、HA3.5螺钉（阳极氧化）（图例为7孔） Note: Used For Φ3.5 Locking Screw(Anodized)、HA3.5screw(Anodized) (Sample: 7 Holes)				

JY 5.0 胫骨近端内侧锁定钢板（左右） 23207
JY 5.0 Proximal Medial Tibial Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 5.0 胫骨近端内侧锁定 钢板（左右） JY 5.0 Proximal Medial Tibial Locking Plate (Left,Right)	纯钛 Pure Titanium	5孔左	125	23207051
		5孔右	125	23207052
		7孔左	161	23207071
		7孔右	161	23207072
		9孔左	197	23207091
		9孔右	197	23207092
		11孔左	233	23207111
		11孔右	233	23207112
		13孔左	269	23207131
		13孔右	269	23207132
		15孔左	305	23207151
		15孔右	305	23207152
备注：选用Φ5.0锁定螺钉（阳极氧化）、HA4.5螺钉（阳极氧化）（图例为9孔） Note: Used For Φ5.0 Locking Screw(Anodized)、HA4.5screw(Anodized) (Sample: 9 Holes)				

JY KSL 胫骨近端外侧锁定钢板（左右） 23205
JY KSL Proximal Lateral Tibial Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY KSL 胫骨近端外侧锁定钢板 JY KSL Proximal Lateral Tibial Locking Plate (Left,Right)	纯钛 Pure Titanium	5孔左	98	23205051
		5孔右	98	23205052
		7孔左	130	23205071
		7孔右	130	23205072
		9孔左	162	23205091
		9孔右	162	23205092
		11孔左	194	23205111
		11孔右	194	23205112
		13孔左	226	23205131
		13孔右	226	23205132
备注：选用Φ3.5锁定螺钉（阳极氧化）、HA3.5螺钉（阳极氧化）（图例为9孔） Note: Used For Φ3.5 Locking Screw(Anodized)、HA3.5 Screw(Anodized) (Sample: 9 Holes)				

JY 锁定胫骨平台前内侧（左右） 23409
JY Anteromedial tibial plateau locking plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 锁定胫骨平台前内侧（左右） JY Anteromedial tibial plateau locking plate (Left,Right)	纯钛 Pure Titanium	7孔左	39	23409071
		7孔右	39	23409072
		8孔左	48	23409081
		8孔右	48	23409082
备注：选用Φ3.5锁定螺钉（阳极氧化）（图例为8孔） Note: Used For Φ3.5 Locking Screw(Anodized)(Sample: 8 Holes)				

JY、WX 锁定钢板系列




JY. WX Locking Plate Series

JY 胫骨近端外侧锁定钢板（左右） 23231
JY Proximal Lateral Tibial Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 胫骨近端外侧锁定钢板 (左右) JY Proximal Lateral Tibial Locking Plate (Left,Right)	纯钛 Pure Titanium	4孔左	148	23231041
		4孔右	148	23231042
		6孔左	184	23231061
		6孔右	184	23231062
		8孔左	220	23231081
		8孔右	220	23231082
		10孔左	256	23231101
		10孔右	256	23231102
		12孔左	292	23231121
		12孔右	292	23231122
		14孔左	328	23231141
		14孔右	328	23231142
备注：选用Φ5.0锁定螺钉（阳极氧化）、HA4.5螺钉（阳极氧化）（图例为8孔） Note: Used For Φ5.0 Locking Screw(Anodized)、HA4.5 Screw(Anodized) (Sample: 8 Holes)				

JY 锁定胫骨近端外侧II型（左右） 23408
JY Proximal Lateral Tibial Locking Plate-II (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Le-ngth	产品编号 Product No.
JY 锁定胫骨近端外侧II型 (左右) JY Proximal Lateral Tibial Locking Plate-II (Left,Right)	纯钛 Pure Titanium	5孔左	85	23408051
		5孔右	85	23408052
		7孔左	111	23408071
		7孔右	111	23408072
		9孔左	137	23408091
		9孔右	137	23408092
		11孔左	163	23408111
		11孔右	163	23408112
		13孔左	189	23408131
		13孔右	189	23408132
		15孔左	215	23408151
		15孔右	215	23408152
		16孔左	228	23408161
		16孔右	228	23408162
备注：选用Φ3.5锁定螺钉（阳极氧化）HA3.5螺钉（阳极氧化）（图例为13孔） Note: Used For Φ3.5 Locking Screw(Anodized)、HA3.5 Screw(Anodized) (Sample: 13 Holes)				

JY 胫骨远端内踝锁定钢板 23231
JY Distal Tibial Medial Malleolus Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 胫骨远端内踝锁定钢板 JY Distal Tibial Medial Malleolus Locking Plate	纯钛 Pure Titanium	4孔	88	23231040
		6孔	111	23231060
备注：选用Φ3.5锁定螺钉(阳极氧化)、HA3.5螺钉(阳极氧化)（图例为6孔） Note: Used For Φ3.5 Locking Screw(Anodized)、HA3.5 Screw(Anodized) (Sample: 6 Holes)				

JY 胫骨远端内侧锁定钢板（左右） 23203
JY Distal Medial Tibial Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 胫骨远端内侧锁定钢板 (左右) JY Distal Medial Locking Plate(Left,Right)	纯钛 Pure Titanium	5孔左	148	23203051
		5孔右	148	23203052
		7孔左	180	23203071
		7孔右	180	23203072
		9孔左	212	23203091
		9孔右	212	23203092
		11孔左	244	23203111
		11孔右	244	23203112
		13孔左	276	23203131
13孔右	276	23203132		
备注：选用Φ3.5锁定螺钉(阳极氧化)、HA3.5螺钉(阳极氧化)（图例为7孔） Note: Used For Φ3.5 Locking Screw(Anodized)、HA3.5 Screw (Anodized)(Sample: 7 Holes)				

F

JY、WX 锁定钢板系列

JY. WX Locking Plate Series

WX 胫骨远端前内侧锁定钢板 (左右) 23567
VA Distal Anteromedial Tibial Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
WX 胫骨远端前内侧锁定钢板、 (左右) VA Distal Anteromedial Tibial Locking Plate(Left,Right)	纯钛 Pure Titanium	4孔左	115	23567041
		4孔右	115	23567042
		6孔左	143	23567061
		6孔右	143	23567062
		8孔左	171	23567081
		8孔右	171	23567082
		10孔左	199	23567101
		10孔右	199	23567102
		12孔左	227	23567121
12孔右	227	23567122		
备注: 选用Φ3.5万向锁定螺钉(阳极氧化)、Φ3.5锁定螺钉(阳极氧化)、HA3.5螺钉(阳极氧化) (图例为8孔) Note: Used For Φ3.5 Variable Angle Locking Screw(Anodized)、 Φ3.5 Locking Screw(Anodized)、HA3.5 Screw(Anodized) (Sample: 8 Holes)				

WX 胫骨远端内侧锁定钢板 (左右) 23568
VA Distal Medial Tibial Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
WX 胫骨远端内侧锁定钢板（左右） VA Distal Medial Tibial Locking Plate(Left,Right)	纯钛 Pure Titanium	4孔左	115	23568041
		4孔右	115	23568042
		6孔左	143	23568061
		6孔右	143	23568062
		8孔左	171	23568081
		8孔右	171	23568082
		10孔左	199	23568101
		10孔右	199	23568102
		12孔左	227	23568121
		12孔右	227	23568122
备注：选用Φ3.5万向锁定螺钉(阳极氧化)、Φ3.5锁定螺钉(阳极氧化)、HA3.5螺钉(阳极氧化) (图例为8孔) Note: Used For Φ3.5 Variable Angle Locking Screw(Anodized)、 Φ3.5 Locking Screw(Anodized)、HA3.5 Screw(Anodized) (Sample: 8 Holes)				

JY 胫骨远端前外侧锁定钢板 (左右) 23204
JY Distal Anterolateral Tibial Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 胫骨远端前外侧锁定钢板 (左右) JY Distal Anterolateral Tibial Locking Plate(Left,Right)	纯钛 Pure Titanium	5孔左	123	23204051
		5孔右	123	23204052
		7孔左	155	23204071
		7孔右	155	23204072
		9孔左	187	23204091
		9孔右	187	23204092
		11孔左	219	23204111
		11孔右	219	23204112
		13孔左	251	23204131
		13孔右	251	23204132
备注：选用Φ3.5锁定螺钉(阳极氧化)、HA3.5螺钉(阳极氧化)（图例为7孔） Note: Used For Φ3.5 Locking Screw(Anodized)、HA3.5 Screw(Anodized) (Sample: 7 Holes)				

JY 胫骨远端后踝锁定板 23234
JY Distal Tibial Posterior Malleolus Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 胫骨远端后踝锁定板 JY Distal Tibial Posterior Malleolus Locking Plate	纯钛 Pure Titanium	4孔	66	23234040
		6孔	86	23234060
备注：选用Φ3.5锁定螺钉(阳极氧化)、HA3.5螺钉(阳极氧化)（图例为6孔） Note：Used For Φ3.5 Locking Screw(Anodized)、HA3.5 Screw(Anodized)(Sample: 6 Holes)				

JY. WX Locking Plate Series

JY 2.7 腓骨远端重建锁定钢板(左右) 23201

JY 2.7 Distal Fibula Reconstruction Locking Plate
(Left,Right)

产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 2.7 腓骨远端重建锁定钢板 (左右) JY 2.7 Distal Fibula Reconstruction Locking Plate (Left,Right)	纯钛 Pure Titanium	4孔左	68	23201041
		4孔右	68	23201042
		6孔左	88	23201061
		6孔右	88	23201062
		8孔左	108	23201081
		8孔右	108	23201082
		10孔左	128	23201101
		10孔右	128	23201102
		12孔左	148	23201121
		12孔右	148	23201122
备注：选用Φ2.7锁定螺钉(阳极氧化)、HA2.7螺钉（图例为6孔） Note: Used For Φ2.7 Locking Screw(Anodized)、HA2.7 Screw (Sample: 6 Holes)				

JY 2.7 腓骨远端外侧锁定钢板-II型(左右) 23226

JY 2.7 Distal Lateral Fibula Locking Plate-II(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 2.7 腓骨远端外侧锁定 钢板-II(左右) JY 2.7 Distal Lateral Fibula Locking Plate-II(Left,Right)	纯钛 Pure Titanium	3孔左	86	23226031
		3孔右	86	23226032
		5孔左	112	23226051
		5孔右	112	23226052
		7孔左	138	23226071
		7孔右	138	23226072
备注：选用Φ2.7、Φ3.5锁定螺钉(阳极氧化)、HA2.7、HA3.5螺钉(阳极氧化) (图例为3孔) Note: Used For Φ2.7/Φ3.5 Locking Screw(Anodized)、HA2.7、HA3.5 Screw(Anodized) (Sample: 3 Holes)				

JY 跟骨结合微创锁定钢板 (左右) 23209

JY Combo Minimally Invasive Calcaneal Locking
Plate(Left,Right)

产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
JY 跟骨结合微创锁定钢板(左右) JY Combo Minimally Invasive Calcaneal Locking Plate (Left,Right)	纯钛 Pure Titanium	9孔左	23209091
		9孔右	23209092
备注：选用Φ3.5锁定螺钉(阳极氧化) Note: Used For Φ3.5 Locking Screw(Anodized)			

JY 跟骨结合微创锁定钢板大号 (左右) 23232

JY Combo Minimally Invasive Calcaneal Locking
Plate(Max)(Left,Right)

产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
JY 跟骨结合微创锁定钢板大号(左右) JY Combo Minimally Invasive Calcaneal Locking Plate (Max)(Left,Right)	纯钛 Pure Titanium	8孔左	23232081
		8孔右	23232082
备注：选用Φ3.5锁定螺钉(阳极氧化) Note: Used For Φ3.5 Locking Screw(Anodized)			

JY 跟骨结合微创锁定钢板小号 (左右) 23212

JY Combo Minimally Invasive Calcaneal Locking
Plate(Mini)(Left,Right)

产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
JY 跟骨结合微创锁定钢板小号(左右) JY Combo Minimally Invasive Calcaneal Locking Plate (Mini)(Left,Right)	纯钛 Pure Titanium	8孔左	23212081
		8孔右	23212082
备注：选用Φ3.5锁定螺钉(阳极氧化) Note: Used For Φ3.5 Locking Screw(Anodized)			

JY、WX 锁定钢板系列

JY. WX Locking Plate Series

JY 跟骨后粗隆微创锁定钢板（左右） 23208
JY Posterior Tuberosity Minimally Invasive Calcaneal Locking Plate(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
JY 跟骨后粗隆微创锁定钢板(左右) JY Posterior Tuberosity Minimally Invasive Calcaneal Locking Plate(Left,Right)	纯钛 Pure Titanium	6孔左	23208061
		6孔右	23208062
备注：选用Φ3.5锁定螺钉(阳极氧化) Note: Used For Φ3.5 Locking Screw(Anodized)			

JY 跟骨后粗隆微创锁定钢板小号（左右） 23211
JY Posterior Tuberosity Minimally Invasive Calcaneal Locking Plate (Mini)(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
JY 跟骨后粗隆微创锁定钢板小号(左右) JY Posterior Tuberosity Minimally Invasive Calcaneal Locking Plate (Mini)(Left,Right)	纯钛 Pure Titanium	5孔左	23211051
		5孔右	23211052
备注：选用Φ3.5锁定螺钉(阳极氧化) Note: Used For Φ3.5 Locking Screw(Anodized)			

JY 跟骨前突微创锁定钢板（左右） 23210
JY Anterior Process Minimally Invasive Calcaneal Locking Plate(Left,Right)




产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
JY 跟骨前突微创锁定钢板(左右) JY Anterior Process Minimally Invasive Calcaneal Locking Plate (Left,Right)	纯钛 Pure Titanium	5孔左	23210051
		5孔右	23210052
备注：选用Φ3.5锁定螺钉(阳极氧化) Note: Used For Φ3.5 Locking Screw(Anodized)			

JY 跟骨前突微创锁定钢板小号（左右） 23213
JY Anterior Process Minimally Invasive Calcaneal Locking Plate (Mini)(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
JY 跟骨前突微创锁定钢板小号(左右) JY Anterior Process Minimally Invasive Calcaneal Locking Plate (Mini)(Left,Right)	纯钛 Pure Titanium	5孔左	23213051
		5孔右	23213052
备注：选用Φ3.5锁定螺钉(阳极氧化) Note: Used For Φ3.5 Locking Screw(Anodized)			

JY 跟骨解剖锁定钢板 小号（左右） 23227
JY Calcaneus Anatomical Locking Plate(Mini)(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 跟骨解剖锁定钢板 小号(左右) JY Calcaneus Anatomical Locking Plate(Mini)(Left,Right)	纯钛 Pure Titanium	11孔左	58	23227111
		11孔右	58	23227112
备注： 选用Φ3.5锁定螺钉(阳极氧化) Note: Used For Φ3.5 Locking Screw(Anodized)				


JY. WX Locking Plate Series

JY 跟骨解剖锁定钢板 中号 (左右) 23228
JY Calcaneus Anatomical Locking Plate(Mid)(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 跟骨解剖锁定钢板 中号(左右) JY Calcaneus Anatomical Locking Plate(Mid)(Left,Right)	纯钛 Pure Titanium	11孔左	64	23228111
		11孔右	64	23228112
备注：选用Φ3.5锁定螺钉(阳极氧化) Note: Used For Φ3.5 Locking Screw(Anodized)				

JY 跟骨解剖锁定钢板 大号 (左右) 23230
JY Calcaneus Anatomical Locking Plate(Max)(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 跟骨解剖锁定钢板 大号(左右) JY Calcaneus Anatomical Locking Plate(Max)(Left,Right)	纯钛 Pure Titanium	11孔左	71	23230111
		11孔右	71	23230112
备注：选用Φ3.5锁定螺钉(阳极氧化) Note: Used For Φ3.5 Locking Screw(Anodized)				

JY 重建骨盆锁定加压板弧形 23235
JY Pelvic Reconstruction Locking Compression Plate
(Curved) (Anodized)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 重建骨盆锁定加压板弧形 VA Pelvic Reconstruction Locking Compression Plate(Curved) (Anodized)	纯钛 Pure Titanium	5孔	62	23235050
		6孔	74	23235060
		8孔	98	23235080
		10孔	122	23235100
		12孔	146	23235120
		14孔	170	23235140
备注: 选用Φ3.5锁定螺钉（阳极氧化）（图例为10孔） Note: Used For Φ3.5 Locking Screw(Anodized) (Sample: 10 Holes)				

JY 重建骨盆锁定加压板直型 23014
JY Pelvic Reconstruction Locking Compression Plate
(Straight) (Anodized)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No
JY 重建骨盆锁定加压板直型 JY pelvic Reconstruction Locking Compression Plate (Straight) (Anodized)	纯钛 Pure Titanium	6孔	76	23014060
		8孔	102	23014080
		10孔	128	23014100
		12孔	154	23014120
		14孔	180	23014140
		16孔	206	23014160
		18孔	232	23014180
		20孔	258	23014200
备注：选用Φ3.5锁定螺钉（阳极氧化）（图例为8孔） Note: Used For Φ3.5 Locking Screw (Anodized) (Sample: 8 Holes)				

JY、WX 锁定钢板系列

JY. WX Locking Plate Series

JY KUH股骨近端锁定钢板 23239
JY KUH Proximal Femur Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY KUH股骨近端锁定钢板 JY KUH Proximal Femur Locking Plate	纯钛 Pure Titanium	3孔左	121	23239031
		3孔右	121	23239032
		4孔左	141	23239041
		4孔右	141	23239042
		5孔左	161	23239051
		5孔右	161	23239052
		7孔左	201	23239071
		7孔右	201	23239072
		9孔左	241	23239091
		9孔右	241	23239092
		11孔左	281	23239111
		11孔右	281	23239112
		14孔左	340	23239141
14孔右	340	23239142		
备注：头部选用Φ6.5中空锁定螺钉(阳极氧化)，其余选用Φ5.0锁定螺钉（阳极氧化）和HA4.5螺钉(阳极氧化)（图例为11孔） Note: Used For Φ6.5locking Screw In Head (Anodized), OtherΦ5.0 Locking Screw (Anodized) And HA4.5 Screw (Anodized)(Sample: 11 Holes)				

JY 股骨远端外侧锁定钢板（左右） 23202
JY Distal Lateral Femoral Locking Plate(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 股骨远端外侧锁定 钢板（左右） JY Distal Lateral Femoral Locking Plate (Left,Right)	纯钛 Pure Titanium	5孔左	178	23202051
		5孔右	178	23202052
		7孔左	222	23202071
		7孔右	222	23202072
		9孔左	266	23202091
		9孔右	266	23202092
		11孔左	310	23202111
		11孔右	310	23202112
		13孔左	353	23202131
13孔右	353	23202132		
备注：选用Φ5.0锁定螺钉（阳极氧化）、HA4.5螺钉（阳极氧化）（图例为11孔） Note: Used For Φ5.0 Locking Screw(Anodized)、HA4.5 Screw(Anodized) (Sample: 11 Holes)				

JY 锁定股骨远端内侧II型（左右） 23238
JY Distal Medial Femoral Locking Plate-II (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 锁定股骨远端内侧II型 (左右) JY Distal Medial Femoral Locking Plate-II (Left,Right)	纯钛 Pure Titanium	5孔左	79	23238105L
		5孔右	79	23238105R
		7孔左	111	23238107L
		7孔右	111	23238107R
备注：选用Φ3.5锁定螺钉（阳极氧化）、HA3.5螺钉（阳极氧化）（图例为7孔） Note: Used For Φ3.5 Locking Screw(Anodized)、HA3.5 Screw(Anodized) (Sample: 7 Holes)				

JY 股骨锁定钢板 23012、23004
JY Femoral Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 股骨锁定钢板 JY Femoral Locking Plate	纯钛 Pure Titanium	8孔短	194	23012080
		8孔	213	23004080
		9孔	235	23004090
		10孔	257	23004100
		12孔	301	23004120
		14孔	345	23004140
备注：选用Φ5.0锁定螺钉（阳极氧化）、HA4.5螺钉（阳极氧化）（图例为9孔） Note: Used For Φ5.0 Locking Screw(Anodized)、HA4.5 Screw(Anodized) (Sample: 9 Holes)				

JY、WX 锁定钢板系列

JY. WX Locking Plate Series

JY 重建锁定钢板 23001

JY Reconstruction Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 重建锁定钢板 JY Reconstruction Locking Plate	纯钛 Pure Titanium	5孔	87	23001050
		6孔	103	23001060
		7孔	119	23001070
		8孔	135	23001080
		9孔	151	23001090
		10孔	167	23001100
		12孔	199	23001120
		14孔	231	23001140
备注：选用Φ3.5锁定螺钉(阳极氧化)、HA3.5螺钉(阳极氧化)（图例为7孔） Note: Used For Φ3.5 Locking Screw(Anodized)、HA3.5 Screw(Anodized) (Sample: 7 Holes)				

JY 尺桡骨锁定钢板 23005

JY Ulna Radius Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 尺桡骨锁定钢板 JY Ulna Radius Locking Plate	纯钛 Pure Titanium	6孔	106	23005060
		7孔	122	23005070
		8孔	138	23005080
		9孔	154	23005090
		10孔	170	23005100
备注：选用Φ3.5锁定螺钉(阳极氧化)、HA3.5螺钉(阳极氧化)（图例为6孔） Note: Used For Φ3.5 Locking Screw(Anodized)、 HA3.5 Screw (Anodized)(Sample: 6 Holes)				

JY 肱骨干骺端锁定钢板 23070

JY Humeral Metaphyseal Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 肱骨干骺端锁定钢板 JY Humeral Metaphyseal Locking Plate	纯钛 Pure Titanium	7孔	112	23070070
		9孔	144	23070090
		11孔	176	23070110
		13孔	208	23070130
备注：选用Φ3.5锁定螺钉(阳极氧化)、HA3.5螺钉(阳极氧化)（图例为9孔） Note: Used For Φ3.5 Locking Screw(Anodized)、HA3.5 Screw(Anodized) (Sample: 9 Holes)				

JY 肱骨锁定钢板 23002

JY Humeral Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 肱骨锁定钢板 JY Humeral Locking Plate	纯钛 Pure Titanium	7孔	131	23002070
		8孔	147	23002080
		9孔	163	23002090
		10孔	179	23002100
		12孔	211	23002120
		14孔	243	23002140
备注：选用Φ3.5锁定螺钉(阳极氧化)、HA3.5螺钉(阳极氧化)（图例为8孔） Note: Used For Φ3.5 Locking Screw(Anodized)、HA3.5 Screw (Anodized)(Sample: 8 Holes)				

JY 胫骨干骺端锁定钢板 23216

JY Tibial Metaphyseal Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 胫骨干骺端锁定钢板 JY Tibial Metaphyseal Locking Plate	纯钛 Pure Titanium	10孔	158	23216010
		12孔	194	23216012
		14孔	230	23216014
		16孔	266	23216016
备注：选用Φ5.0、Φ3.5锁定螺钉(阳极氧化)、HA3.5、HA4.5螺钉(阳极氧化)（图例为12孔） Note: Used ForΦ5.0、Φ3.5 Locking Screw(Anodized)、HA3.5、HA4.5 Screw (Anodized) (Sample: 12 Holes)				

JY 胫骨锁定钢板 23003

JY Tibial Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
JY 胫骨锁定钢板 JY Tibial Locking Plate	纯钛 Pure Titanium	7孔	159	23003070
		8孔	177	23003080
		9孔	195	23003090
		10孔	213	23003100
		12孔	249	23003120
		14孔	285	23003140
备注：选用Φ5.0锁定螺钉(阳极氧化)、HA4.5螺钉(阳极氧化)（图例为8孔） Note: Used For Φ5.0 Locking Screw(Anodized)、HA4.5 Screw(Anodized) (Sample: 8 Holes)				

JY. WX Locking Plate Series



系统特点
Character Of The System

- 加压技术
- 万向锁定技术
- 解剖型和特定术式内植入物
- Compression Technology
- Variable Angled Locking Technology
- Anatomical And Special Surgical Implant

WX 第一跖跗关节锁定融合板 T47981
VA First TMT(Tarsometatarsal) Locking Fusion Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
WX 第一跖跗关节锁定融合板 VA First TMT(Tarsometatarsal) Locking Fusion Plate	纯钛 Pure Titanium	6	39	T4798106

备注：选用Φ 2.7万向锁定螺钉(阳极氧化)、HA2.7螺钉
Note: Used For Φ2.7 Variable Angle Locking Screw(Anodized)、HA2.7 Screw

WX 跖跗关节锁定融合板 T47982
VA TMT(Tarsometatarsal) Locking Fusion Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
WX 跖跗关节锁定融合板 VA TMT(Tarsometatarsal) Locking Fusion Plate	纯钛 Pure Titanium	4	44	T4798204

备注：选用Φ 2.7万向锁定螺钉(阳极氧化)、HA2.7螺钉
Note: Used For Φ2.7 Variable Angle Locking Screw(Anodized)、HA2.7 Screw

JY. WX Locking Plate Series

WX 开放楔形锁定板 T47983
VA Opening Wedge Variable Angle Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	间隔块宽度(mm) Spacer Block	产品编号 Product No.
WX 开放楔形锁定板 VA Opening Wedge Variable Angle Locking Plate	纯钛 Pure Titanium	4	3	T4798303
			4	T4798304
			5	T4798305
			6	T4798306
			7	T4798307
备注：选用Φ2.7万向锁定螺钉(阳极氧化) Note: Used For Φ2.7 Variable Angle Locking Screw(Anodized)				

WX 第一跖趾关节锁定融合板（左右） T47984
VA First MTP(Metatarsophalangeal) Locking
Fusion Plate(Left,Right)



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
WX 第一跖趾关节锁定融合板（左右） VA First MTP (Metatarsophalangeal)Locking Fusion Plate(Left,Right)	纯钛 Pure Titanium	5*6孔左	T4798406L
		5*6孔右	T4798406R
备注：选用Φ2.7万向锁定螺钉(阳极氧化)、HA2.7螺钉 Note: Used For Φ2.7 Variable Angle Locking Screw(Anodized)、HA2.7 Screw			

WX X型通用锁定融合板 T47985
VA X Type Universal Locking Fusion Plate




产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
WX X型通用锁定融合板 VA X Type Universal Locking Fusion Plate	纯钛 Pure Titanium	4	27	T4798504
备注：选用Φ2.7万向锁定螺钉(阳极氧化) Note: Used For Φ2.7 Variable Angle Locking Screw(Anodized)				

Φ2.7 万向锁定螺钉（阳极氧化） T48027
Φ2.7 VA Locking Screw(Anodized)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
Φ2.7 万向锁定螺钉 (阳极氧化) Φ2.7 VA Locking Screw(Anodized)	钛合金 Titanium Alloy	2.7	10、12、14、16、 18、20、22、24、 26、28、	T4802710- T4802728
			30	T4802730
			35	T4802735
			40	T4802740

JY. WX Locking Plate Series



JY 3.5 锁定连接板（阳极氧化）23016
JY 3.5 Locking Connecting Plate (Anodized)



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
JY 3.5 锁定连接板 (阳极氧化) JY 3.5 Locking Connecting Plate (Anodized)	纯钛 Pure Titanium	4x28	23016428
		4x31	23016431
		4x34	23016434
备注：选用Φ3.5锁定连接钉（阳极氧化）、Φ3.5锁定螺钉（阳极氧化） Note:use Φ3.5 Locking Connecting Screw (Anodized).Φ 3.5 Locking Screw (Anodized) .			

JY 5.0 锁定连接板（阳极氧化）23017
JY 5.0 Locking Connecting Plate (Anodized)



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
JY 5.0 锁定连接板 (阳极氧化) JY 5.0 Locking Connecting Plate (Anodized)	纯钛 Pure Titanium	4x35	23017435
		4x38	23017438
		4x41	23017441
备注：选用φ5.0锁定连接钉（阳极氧化）、φ5.0锁定螺钉（阳极氧化） Note:use φ5.0 Locking Connecting Screw (Anodized)φ 5.0 Locking Screw (Anodized).			

锁定连接钉（阳极氧化）33004
Locking Connecting Screw (Anodized)



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
3.5 锁定连接钉 (阳极氧化) 3.5 Locking Connecting Screw (Anodized)	钛合金 Titanium Alloy	Φ3.5x10	33004035
5.0 锁定连接钉（阳极氧化） 5.0 Locking Connecting Screw (Anodized)		Φ5.0x14	33004050

阳极氧化螺钉系列

Screw Series (Anodized)



Φ2.7 梅花型锁定螺钉（阳极氧化） T491
Torx Locking Screw (Anodized)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型锁定螺钉 (阳极氧化) Torx Locking Screw (Anodized)	钛合金 Titanium Alloy	2.7	6-60 (每2mm为一档) The specification is 2mm	T4912706- T4912760

Φ3.5 梅花型锁定螺钉（阳极氧化） T492
Torx Locking Screw (Anodized)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型锁定螺钉 (阳极氧化) Torx Locking Screw (Anodized)	钛合金 Titanium Alloy	3.5	10-50 (每2mm为一档) The specification is 2mm	T4923510- T4923550
			55-80 (每5mm为一档) The specification is 5mm	T4923555- T4923580

Φ5.0 梅花型锁定螺钉（阳极氧化） T493
Torx Locking Screw (Anodized)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型锁定螺钉 (阳极氧化) Torx Locking Screw (Anodized)	钛合金 Titanium Alloy	5.0	14-50 (每2mm为一档) The specification is 2mm	T4935014- T4935050
			55-100 (每5mm为一档) The specification is 5mm	T4935055- T4935100

填塞钉（阳极氧化） 33005、33006
Filling Screw (Anodized)



产品名称 Product Name	材质 Material	规格(mm) Specification	产品编号 Product No.
填塞钉(阳极氧化) Filling Screw (Anodized)	钛合金 Titanium Alloy	3.5*1.5mm	33005104
		5.0*1.5mm	33006105

F

Screw Series (Anodized)

Φ2.7 万向锁定螺钉(阳极氧化) T480
VA Locking Screw(Anodized)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
万向锁定螺钉阳极氧化 Universal Locking Screw(Anodized)	钛合金 Titanium Alloy	2.7	10、12、14、16、 18、20、22、24、 26、28、	T4802710- T4802728
			30、35、40	T4802730- T4802740

Φ3.5 万向锁定螺钉阳极氧化 33151
Universal Locking Screw(Anodized)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
万向锁定螺钉阳极氧化 Universal Locking Screw(Anodized)	钛合金 Titanium Alloy	3.5	10-50 (每2mm为一档) The specification is 2mm	33151010- 33151050
			55-80 (每5mm为一档) The specification is 5mm	33151055- 33151080

HA3.5 金属接骨螺钉 (阳极氧化) T52235
Metal Bone Fracture Screw (Anodized)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
金属接骨螺钉(阳极氧化) Metal Bone Fracture Screw(Anodized)	钛合金 Titanium Alloy	HA3.5	12-50 (每2mm为一档) The specification is 2mm	T5223512- T5223550

HA4.5 金属接骨螺钉 (阳极氧化) T52245
Metal Bone Fracture Screw (Anodized)



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
金属接骨螺钉(阳极氧化) Metal Bone Fracture Screw(Anodized)	钛合金 Titanium Alloy	HA4.5	20-50 (每2mm为一档) The specification is 2mm	T5224520- T5224550
			55、60	T5224555 T5224560



(本色) 锁定钢板螺钉系列
(Base Color) Locking Plate Series

Locking Plate Series

2.7 桡骨远端支持锁定钢板(大号)(左右) T4320

2.7 Distal Radius Buttress Locking Plate (Large)(Left,Right)



视频在线播放
Video online play

产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
2.7桡骨远端支持 锁定钢板（大号）（左右） 2.7Distal Radius Buttress Locking Plate (Large) (Left,Right)	纯钛 Pure Titanium	4	55	T4320004R/L
		5	63	T4320005R/L
		7	79	T4320007R/L
		9	95	T4320009R/L
备注：选用Φ2.7锁定螺钉，HA2.7螺钉（图例为5孔） Note:used For Φ3.5 Locking Screw, HA2.7 Screw (Sample:5 Holes)				

2.7 桡骨远端支持锁定钢板(小号)(左右) T4324

2.7 Distal Radius Buttress Locking Plate (Mini)(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
2.7桡骨远端支持 锁定钢板（小号）（左右） 2.7Distal Radius Buttress Locking Plate (Mini) (Left,Right)	纯钛 Pure Titanium	4	55	T4324004R/L
		5	63	T4324005R/L
		7	79	T4324007R/L
		9	95	T4324009R/L
备注：选用Φ2.7锁定螺钉，HA2.7螺钉（图例为4孔） Note:used For Φ2.7 Locking Screw, HA2.7 Screw (Sample:4 Holes)				

2.7 桡骨远端支持锁定钢板 (左右) T4321

2.7 Distal Radius Buttress Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.7桡骨远端背侧 锁定钢板 (左右) 2.7 Distal Back Radius Locking Plate (Left,Right)	纯钛 Pure Titanium	14	T4321014R/L
备注: 选用Φ2.7锁定螺钉 (图例为14孔) Note:used For Φ2.7 Locking Screw (Sample:14 Holes)			

小T型重建锁定钢板 T4330

Mini-T Reconstruction Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
小T型重建锁定钢板 Mini-T Reconstruction Locking Plate	纯钛 Pure Titanium	6	43	T4330006
		7	54	T4330007
		8	65	T4330008
		9	76	T4330009
备注：选用Φ3.5锁定螺钉（图例为7孔） Note:used For Φ3.5 Locking Screw (Sample:7 Holes)				

桡骨斜T型重建锁定钢板 T4340

Oblique T Radius Reconstruction Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
桡骨斜T型重建 锁定钢板（左右） Oblique T Radius Reconstruction Locking Plate (Left,Right)	纯钛 Pure Titanium	6	45	T4340006R/L
		7	56	T4340007R/L
		8	67	T4340008R/L
		9	78	T4340009R/L
备注：选用Φ3.5锁定螺钉（图例为7孔） Note:used For Φ3.5 Locking Screw (Sample:7 Holes)				

2.7 桡骨头锁定钢板 (左右) T4322

2.7 Radius Head Locking Plate (Left,Right)



视频在线播放
Video online play

产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
2.7桡骨头锁定钢板 (左右) 2.7 Radius Head Locking Plate (Left,Right)	纯钛 Pure Titanium	3	40	T43220003R/L
		4	48	T43220004R/L
		5	56	T43220005R/L
备注: 选用Φ2.7锁定螺钉 (图例为5孔) Note:used For Φ2.7 Locking Screw (Sample:5 Holes)				

Locking Plate Series

鹰嘴钩重建锁定钢板 (左右) T4310
Olecranon Hook Reconstruction Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
鹰嘴钩重建 锁定钢板 (左右) Olecranon Hook Reconstruction Locking Plate (Left,Right)	纯钛 Pure Titanium	4	87	T4310004R/L
		6	113	T4310006R/L
		8	139	T4310008R/L
备注：选用Φ3.5锁定螺钉 (图例为5孔) Note:used For Φ3.5 Locking Screw (Sample:5 Holes)				

尺骨近端锁定钢板 (左右) T4370
Proximal Ulna Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
尺骨近端锁定 钢板 (左右) Proximal Ulna Locking Plate (Left,Right)	纯钛 Pure Titanium	4	119	T4370004R/L
		6	145	T4370006R/L
		8	171	T4370008R/L
备注: 选用Φ3.5锁定螺钉 (图例为5孔) Note:used For Φ3.5 Locking Screw (Sample:5 Holes)				

锁骨远端锁定钢板 (左右) T4381
Distal Clavicle Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
锁骨远端锁定 钢板 (左右) Distal Clavicle Locking Plate (Left,Right)	纯钛 Pure Titanium	4	72	T4381004R/L
		5	84	T4381005R/L
		6	96	T4381006R/L
		7	108	T4381007R/L
		8	119	T4381008R/L
		10	141	T4381010R/L
备注：头部选用Φ2.7锁定螺钉，其余选用Φ3.5锁定螺钉（图例为4孔） Note:used For Φ2.7 Locking Screw In Head,other Φ3.5 Locking Screw (Sample: 4 Holes)				

3.0 锁骨S型锁定钢板 (左右) T4382
3.0 S Clavicle Locking Plate (Left,Right)




产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
3.0锁骨S型 锁定钢板 (左右) 3.0 S Clavicle Locking Plate (Left,Right)	纯钛 Pure Titanium	6	70	T4382006
		7	82	T4382007R/L
		8	94	T4382008R/L
		9	105	T4382008R/L
		10	116	T4382010R/L
		12	137	T4382012R/L
备注: 选用Φ3.5锁定螺钉 (图例为7孔) Note:used For Φ3.5 Locking Screw (Sample:7 Holes)				

锁定经皮锁骨钢板 (左右) T4383
Percutaneous Clavicle Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
锁定经皮锁骨 钢板（左右） Percutaneous Clavicle Locking Plate (Left,Right)	纯钛 Pure Titanium	7	82	T4383007R/L
		8	94	T4383008R/L
		9	105	T4383009R/L
		10	116	T4383010R/L
		12	137	T4383012R/L
备注：选用Φ3.5锁定螺钉（图例为9孔） Note:used For Φ3.5 Locking Screw (Sample:9 Holes)				

肩锁钩锁定钢板 (左右) T4390、T4391
Clavicle Hook Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	高度(mm) Height	孔数 Holes	长度(mm) Length	产品编号 Product No.
肩锁钩锁定 钢板 (左右) Clavicle Hook Locking Plate (Left,Right)	纯钛 Pure Titanium	14	6	57	T4390006R/L
			7	69	T4390007R/L
			8	81	T4390008R/L
			9	93	T4390009R/L
			10	105	T4390010R/L
		17	6	57	T4391006R/L
			7	69	T4391007R/L
			8	81	T4391008R/L
			9	93	T4391009R/L
			10	105	T4391010R/L
备注：选用Φ3.5锁定螺钉 (图例为7孔) Note:used For Φ3.5 Locking Screw (Sample:7 Holes)					

肱骨近端锁定钢板III型（左右） T4230
Proximal Humerus Locking Plate-III (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
肱骨近端锁定 钢板Ⅲ型（左右） Proximal Humerus Locking Plate-Ⅲ (Left,Right)	纯钛 Pure Titanium	4	103	T4230004R/L
		5	117	T4230005R/L
		6	131	T4230006R/L
		8	159	T4230008R/L
		10	187	T4230010R/L
		12	215	T4230012R/L
备注：选用Φ3.5锁定螺钉,HA3.5螺钉，HB4.0螺钉（图例为6孔） Note:used For Φ3.5 Locking Screw, HA3.5 Screw, HB4.0 Screw (Sample:6holes)				

肱骨近端锁定钢板II型 T4220
Proximal Humerus Locking Plate-II



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
肱骨近端锁定 钢板II型 Proximal Humerus Locking Plate-II	纯钛 Pure Titanium	3	96	T4220003
		4	108	T4220004
		5	120	T4220005
		6	132	T4220006
		8	156	T4220008
		10	180	T4220010
		12	204	T4220012
备注：选用Φ3.5锁定螺钉,HA3.5螺钉，HB4.0螺钉（图例为6孔） Note:used For Φ3.5 Locking Screw, HA3.5 Screw, HB4.0 Screw (Sample:6holes)				

肱骨远端外侧锁定钢板（左右） T4270
Distal Lateral Humerus Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
肱骨远端外侧 锁定钢板（左右） Distal Lateral Humerus Locking Plate (Left,Right)	纯钛 Pure Titanium	5	110	T4270005R/L
		7	124	T4270007R/L
		9	148	T4270009R/L
		11	172	T4270011R/L
备注：头部选用Φ2.7锁定螺钉， 其余选用Φ3.5锁定螺钉和HA3.5螺钉（图例为7孔） Note:used ForΦ2.7 Locking Screw In Head,other Φ3.5 Locking Screw And HA3.5 Screw(Sample: 7 Holes)				

肱骨远端内侧锁定钢板（左右） T4280
Distal Medial Humerus Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
肱骨远端内侧 锁定钢板（左右） Distal Medial Humerus Locking Plate (Left,Right)	纯钛 Pure Titanium	4	82	T4280004R/L
		5	94	T4280005R/L
		7	118	T4280007R/L
		9	142	T4280009R/L
备注：头部选用Φ2.7锁定螺钉，其余选用Φ3.5锁定螺钉和HA3.5螺钉（图例为7孔） Note:used ForΦ2.7 Locking Screw In Head,other Φ3.5 Locking Screw And HA3.5 Screw(Sample: 7 Holes)				

锁定钢板系列

Locking Plate Series

重建锁定钢板 T4040
Reconstruction Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
重建锁定钢板 Reconstruction Locking Plate	纯钛 Pure Titanium	5	63	T4040005
		6	76	T4040006
		7	89	T4040007
		8	102	T4040008
		9	115	T4040009
		10	128	T4040010
		12	154	T4040012
备注：选用Φ3.5锁定螺钉（图例为7孔） Note:used For Φ3.5 Locking Screw (Sample:7 Holes)				

1/3锁定钢板 T4100
One-third Tubular Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
1/3 锁定钢板 One-third Tubular Locking Plate	纯钛 Pure Titanium	5	64	T4100005
		6	80	T4100006
		7	92	T4100007
		8	104	T4100008
		10	128	T4100010
备注：选用Φ3.5锁定螺钉，HA3.5螺钉（图例为6孔） Note:used For Φ3.5 Locking Screw, HA3.5 Screw (Sample:6 Holes)				

尺桡骨锁定钢板 T4090
Ulna Radius Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
尺桡骨锁定钢板 Ulna Radius Locking Plate	纯钛 Pure Titanium	5	69	T4090005
		6	85	T4090006
		7	91	T4090007
		8	109	T4090008
		10	133	T4090010
		12	157	T4090012
备注：选用Φ3.5锁定螺钉，HA 3.5螺钉（图例为7孔） Note:used For Φ3.5 Locking Screw, HA3.5 Screw (Sample:7 Holes)				

肱骨锁定钢板 T4050
Humeral Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
肱骨锁定钢板 Humeral Locking Plate	纯钛 Pure Titanium	6	105	T4050006
		7	119	T4050007
		8	133	T4050009
		9	147	T4050009
		10	161	T4050010
		12	189	T4050012
		14	217	T4050014
备注：选用Φ3.5锁定螺钉，HA3.5螺钉（图例为7孔） Note:used For Φ3.5 Locking Screw, HA3.5 Screw (Sample:7 Holes)				

胫骨锁定钢板 T4071
Tibial Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
胫骨锁定钢板 Tibial Locking Plate	纯钛 Pure Titanium	6	127	T4071006
		7	143	T4071007
		8	159	T4071008
		10	191	T4071010
		12	223	T4071012
		14	255	T4071014
备注：选用Φ5.0锁定螺钉，HA4.5螺钉（图例为7孔） Note:used For Φ5.0 Locking Screw, HA4.5 Screw (Sample:7 Holes)				

股骨锁定钢板 T4083
Femoral Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
股骨锁定钢板 Femoral Locking Plate	纯钛 Pure Titanium	8	195	T4083008
		9	215	T4083009
		10	235	T4083010
		12	275	T4083012
		14	315	T4083014
		16	355	T4083016
备注：选用Φ5.0锁定螺钉，HA4.5螺钉（图例为10孔） Note:used For Φ5.0 Locking Screw, HA4.5 Screw (Sample:10 Holes)				

2.0 Locking Plate Series

2.0 T型锁定钢板 头部2孔 T41201
2.0 T-Type Locking Plate(Two Holes In The Head)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.0 T型锁定钢板 头部2孔 2.0 T-Type Locking Plate (Two Holes In The Head)	纯钛 Pure Titanium	5/6/7/8/9	T4120105-T4120109
备注: 选用Φ2.0锁定螺钉(三角形), HA2.0螺钉(三角形) (图例为7孔) Note:used For Φ2.0 Locking Screw (Triangle) , HB2.0screw (Triangle) (Sample:7 Holes)			

2.0 T型锁定钢板 头部3孔 T41202
2.0 T-Type Locking Plate(Three Holes In The Head)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.0 T型锁定钢板 头部3孔 2.0 T-Type Locking Plate (Three Holes In The Head)	纯钛 Pure Titanium	6/7/8/9	T4120206-T4120209
备注: 选用Φ2.0锁定螺钉(三角形), HA2.0螺钉(三角形) (图例为8孔) Note:used For Φ2.0 Locking Screw (Triangle) , HB2.0screw (Triangle) (Sample:8 Holes)			

2.0 踝部锁定钢板 T41204
2.0 Ankle Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.0 踝部锁定钢板 2.0 Ankle Locking Plate	纯钛 Pure Titanium	5/6/7/8/9	T4120405-T4120409
备注: 选用Φ2.0锁定螺钉(三角形), HA2.0螺钉(三角形) (图例为9孔) Note:used For Φ2.0 Locking Screw (Triangle) , HB2.0screw (Triangle) (Sample:9 Holes)			

2.0 L型锁定钢板 (左右) T41206
2.0 L-Type Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.0 L型锁定钢板 (左右) 2.0 L-Type Locking Plate (Left,Right)	纯钛 Pure Titanium	5/6/7/8	T4120605-T4120609
备注: 选用Φ2.0锁定螺钉(三角形), HA2.0螺钉(三角形) (图例为8孔) Note:used For Φ2.0 Locking Screw (Triangle) , HB2.0screw (Triangle) (Sample:8 Holes)			

2.0 直型锁定钢板 T41200
2.0 Straight Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.0 直型锁定钢板 2.0 Straight Locking Plate	纯钛 Pure Titanium	4/5/6/7/8/9/10	T4120004-T4120010
备注: 选用Φ2.0锁定螺钉(三角形) (图例为6孔) Note:used For Φ2.0 Locking Screw (Triangle) (Triangle) (Sample:6 Holes)			

三角形锁定螺钉 T480
Triangle Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
三角形锁定螺钉 Triangle Locking Screw	钛合金 Titanium Alloy	2.0	6、8、10、12、 14、16、18、20 22	T4802006- T4802022

金属接骨螺钉 (三角形) T541
Metal Bone Fracture Screw (Triangle)



钉帽形式 Head of Nail Type	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
三角形 Triangle	钛合金 Titanium Alloy	2.0	6、8、10、12、 14、16、18、20 22	T5412006- T5412022

锁定钢板2.5系列

2.5 Locking Plate Series

2.5 T型锁定钢板 头部2孔 T41251
2.5 T-Type Locking Plate (Two Holes In The Head)



2.5 T型锁定钢板 头部3孔 T41252
2.5 T-Type Locking Plate (Three Holes In The Head)



2.5 踝部锁定钢板 T41254
2.5 Ankle Locking Plate



2.5 L型锁定钢板 (左右) T41256
2.5 L-Type Locking Plate (Left,Right)



2.5 直型锁定钢板 T41250
2.5 Straight Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.5 T型锁定钢板 头部2孔 2.5 T-Type Locking Plate (Two Holes In The Head)	纯钛 Pure Titanium	5/6/7/8/9	T4125105-T4125109
备注: 选用Φ2.5锁定螺钉, HA2.5螺钉 (图例为7孔) Note:used For Φ2.5 Locking Screw, HA2.5 Screw (Sample:7 Holes)			

产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.5 T型锁定钢板 头部3孔 2.5 T-Type Locking Plate (Three Holes In The Head)	纯钛 Pure Titanium	6/7/8/9/10	T4125206-T4125210
备注: 选用Φ2.5锁定螺钉, HA2.5螺钉 (图例为8孔) Note:used For Φ2.5 Locking Screw, HA2.5 Screw (Sample:8 Holes)			

产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.5 踝部锁定钢板 2.5 Ankle Locking Plate	纯钛 Pure Titanium	5/6/7/8/9	T4125405-T4125409
备注: 选用Φ2.5锁定螺钉, HA2.5螺钉 (图例为9孔) Note:used For Φ2.5 Locking Screw, HA2.5 Screw (Sample:9 Holes)			

产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.5 L型锁定钢板 (左右) 2.5 L-Type Locking Plate (Left,Right)	纯钛 Pure Titanium	5/6/7/8/9	T4125605-T4125609R/L
备注: 选用Φ2.5锁定螺钉, HA2.5螺钉 (图例为8孔) Note:used For Φ2.5 Locking Screw, HA2.5 Screw (Sample:8 Holes)			

产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.5 直型锁定钢板 2.5 Straight Locking Plate	纯钛 Pure Titanium	4/5/6/7/8/9	T4125004-T4125010
备注: 选用Φ2.5锁定螺钉 (图例为6孔) Note:used For Φ2.5 Locking Screw (Sample:6 Holes)			

锁定钢板2.7系列

2.7 Locking Plate Series

2.7 H型锁定钢板 T41277
2.7 H-Type Locking Plate



2.7 X型锁定钢板 T41278
2.7 X-Type Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.7 H型锁定钢板 2.7 H-Type Locking Plate	纯钛 Pure Titanium	5	T4127705
		8	T4127708
备注: 选用Φ2.7锁定螺钉 Note:used For Φ2.7 Locking Screw			

产品名称 Product Name	材质 Material	宽度(mm) Width	长度(mm) Length	产品编号 Product No.
2.7 X型锁定钢板 2.7 X-Type Locking Plate	纯钛 Pure Titanium	8	16	T4127816
		12	18	T4127818
		14	24	T4127824
			30	T4127830
备注: 选用Φ2.7锁定螺钉 Note:used For Φ2.7 Locking Screw				

Locking Plate Series

KUS 股骨近端外侧锁定钢板（左右） T4424
KUS-Proximal Lateral Femoral Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
KUS 股骨近端外侧 锁定钢板(左右) KUS-Proximal Lateral Femoral Locking Plate (Left,Right)	纯钛 Pure Titanium	4	149	T4424004R/L
		5	169	T4424005R/L
		6	189	T4424006R/L
		8	229	T4424008R/L
		10	269	T4424010R/L
		12	309	T4424012R/L
备注：头部选用Φ6.5中空锁定螺钉，其余选用Φ5.0锁定螺钉和HA4.5螺钉(图为5孔) Note: Used For Φ6.5 Locking Screw In Head, Other Φ5.0 Locking Screw And HA4.5 Screw (Sample: 5 Holes)				

KUH 股骨近端锁定钢板（左右） T4430
KUH-Proximal Femoral Locking Plate (Left,Right)




产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
KUH 股骨近端 锁定钢板(左右) KUH-Proximal Femoral Locking Plate (Left,Right)	纯钛 Pure Titanium	4	128	T4430004R/L
		5	148	T4430005R/L
		6	168	T4430006R/L
		8	208	T4430008R/L
		10	248	T4430010R/L
		12	288	T4430012R/L
备注：头部选用Φ6.5中空锁定螺钉，其余选用Φ5.0锁定螺钉和HA4.5螺钉(图为5孔) Note: Used For Φ6.5 Locking Screw In Head, Other Φ5.0 Locking Screw And HA4.5 Screw (Sample: 5 Holes)				

KHS-I 股骨近端锁定钢板（左右） T4425
KHS-I Proximal Femoral Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
KHS-I 股骨近端锁定钢板(左右) KHS-I Proximal Femoral Locking Plate (Left,Right)	纯钛 Pure Titanium	4	T4425004R/L
备注：头部选用Φ7.0中空锁定螺钉，其余选用Φ5.0锁定螺钉 Note: Used For Φ7.0 Locking Screw In Head, Other Φ5.0 Locking Screw			

KHS-II 股骨近端加压锁定板（左右） T4426/T4427
KHS-II Proximal Femoral Locking Compression Plate (Left,Right)



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
KHS-II 股骨近端加压锁定板(左右) KHS-II Proximal Femoral Locking Compression Plate(Left,Right)	纯钛 Pure Titanium	大号3孔 Large 3 Holes	T4426003R/L
		小号3孔 Small 3 Holes	T4427003R/L

KHS-II 股骨颈固定螺钉 T4907
KHS-II Femoral Neck Fixing Screw



产品名称 Product Name	材质 Material	螺纹直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KHS-II 股骨颈固定螺钉 KHS-II Femoral Neck Fixing Screw	钛合金 Titanium Alloy	7.0	60-115 (每5mm为一档 The specification is 5mm)	T49070060 -T49070115

KHS-II 支头螺丝 T4910
KHS-II Nut



产品名称 Product Name	材质 Material	螺纹直径(mm) Diameter	产品编号 Product No.
KHS-II支头螺丝 KHS-II Nut	钛合金 Titanium Alloy	9.0	T4910000

锁定钢板系列

Locking Plate Series

骨盆重建锁定钢板 T4790
Pelvic Reconstruction Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
骨盆重建锁定钢板 Pelvic Reconstruction Locking Plate	纯钛 Pure Titanium	6	T4790006
		8	T4790008
		10	T4790010
		12	T4790012
		14	T4790014
备注：选用Φ3.5锁定螺钉（图例为12孔） Note:used For Φ3.5 Locking Screw (Sample:12holes)			

130°DHS 鹅头钉锁定钢板 T4500
130°DHS Compression Hip Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
130° DHS 鹅头钉 锁定钢板 130° DHS Compression Hip Locking Plate	纯钛 Pure Titanium	4	125	T4500004
		5	146	T4500005
		6	165	T4500006
		8	203	T4500008
备注：选用Φ5.0锁定螺钉，HA4.5螺钉，HB6.5螺钉(图例为5孔) Note:used For Φ5.0 Locking Screw,HA4.5 Screw,HB6.5 Screw (Sample:5holes)				

股骨远端外侧锁定钢板 (左右) T4411
Distal Lateral Femoral Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
股骨远端外侧 锁定钢板（左右） Distal Lateral Femoral Locking Plate (Left,Right)	纯钛 Pure Titanium	5	166	T4411005R/L
		7	206	T4411007R/L
		9	246	T4411009R/L
		11	285	T4411011R/L
		13	325	T4411013R/L
备注：选用Φ5.0锁定螺钉，HA4.5螺钉（图例为9孔） Note:used For Φ5.0 Locking Screw,HA4.5 Screw (Sample:9holes)				

胫骨近端外侧锁定钢板IV型 (左右) T4700
Proximal Lateral Tibial Locking Plate-IV(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
胫骨近端外侧锁定 钢板IV型（左右） Proximal Lateral Tibial Locking Plate-IV （Left,Right）	纯钛 Pure Titanium	4	113	T4700004R/L
		5	129	T4700005R/L
		7	161	T4700007R/L
		9	193	T4700009R/L
		11	225	T4700011R/L
		13	257	T4700013R/L
备注：选用Φ5.0锁定螺钉，HA4.5螺钉（图例为7孔） Note:used For Φ5.0 Locking Screw,HA4.5 Screw (Sample:7holes)				

KSL 胫骨近端外侧锁定钢板 (左右) T4711
KSL Proximal Lateral Tibial Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
KSL 胫骨近端外侧 锁定钢板 (左右) KSL Proximal Lateral Tibial Locking Plate (Left,Right)	纯钛 Pure Titanium	5	87	T4711005R/L
		7	116	T4711007R/L
		9	144	T4711009R/L
		11	172	T4711011R/L
		13	200	T4711013R/L
		15	228	T4711015R/L
备注：选用Φ3.5锁定螺钉，HA3.5螺钉（图例为9孔） Note:used For Φ3.5 Locking Screw,HA3.5 Screw (Sample:9holes)				

Locking Plate Series

3.5胫骨近端内侧锁定钢板 (左右) T4770
3.5 Proximal Medial Tibial Locking Plate
(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
3.5胫骨近端内 侧锁定钢板 (左右) 3.5 Proximal Medial Tibial Locking Plate (Left,Right)	纯钛 Pure Titanium	5	101	T4770005R/L
		7	129	T4770007R/L
		9	157	T4770009R/L
		11	185	T4770011R/L
备注: 选用Φ3.5锁定螺钉 (图例为7孔) Note:used For Φ3.5 Locking Screw (Sample:7holes)				

下肢T型重建锁定钢板 (左右) T4730
T Reconsyruction Locking Plate (Lower Limb) (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
下肢T型重建 锁定钢板（左右） T Reconstruction Locking Plate (Lower Limb) (Left,Right)	纯钛 Pure Titanium	7	91	T4730007R/L
		9	123	T4730009R/L
		11	155	T4730011R/L
备注：选用Φ5.0锁定螺钉,HA4.5螺钉（图例为9孔） Note:used For Φ5.0 Locking Screw, ha4.5 Screw (Sample:9holes)				

下肢L型重建锁定钢板 (左右) T4740
L Reconsyruction Locking Plate (Lower Limb)
(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
下肢L型重建 锁定钢板 (左右) L Reconstruction Locking Plate (Lower Limb) (Left,Right)	纯钛 Pure Titanium	7	101	T4740007R/L
		9	133	T4740009R/L
		11	165	T4740011R/L
备注: 选用Φ5.0锁定螺钉, HA4.5螺钉 (图例为7孔) Note:used For Φ5.0 Locking Screw,HA4.5 Screw (Sample:7holes)				

胫骨平台后柱锁定钢板 (左右) T4750
Tibial Posterior Column Platform Locking Plate
(Left,Right)




产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
胫骨平台后柱 锁定钢板（左右） Tibial Posterior Column Platform Locking Plate （Left,Right）	纯钛 Pure Titanium	3	65	T4750003R/L
		5	91	T4750005R/L
		7	117	T4750007R/L
备注：选用Φ3.5锁定螺钉（图例为7孔） Note:used For Φ3.5 Locking Screw (Sample:7holes)				

胫骨平台内柱锁定钢板 (左右) T4760
Tibial Interior Column Platform Locking Plate
(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
胫骨平台内柱 锁定钢板（左右） Tibial Interior Column Platform Locking Plate （Left,Right）	纯钛 Pure Titanium	4	61	T4760004R/L
		6	87	T4760006R/L
		8	113	T4760008R/L
备注：选用Φ3.5锁定螺钉（图例为8孔） Note:used For Φ3.5 Locking Screw (Sample:8holes)				

Locking connecting plate series



3.5 锁定连接板 23018
3.5 Locking Connecting Plate



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
3.5 锁定连接板 3.5 Locking Connecting Plate	纯钛 Pure Titanium	4x28	23018428
		4x31	23018431
		4x34	23018434
备注：选用Φ3.5锁定连接钉、Φ3.5锁定螺钉 Note:use Φ3.5 Locking Connecting Screw,Φ 3.5 Locking Screw .			

5.0 锁定连接板 23019
5.0 Locking Connecting Plate



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
5.0 锁定连接板 5.0 Locking Connecting Plate	纯钛 Pure Titanium	4x35	23019435
		4x38	23019438
		4x41	23019441
备注：选用Φ5.0锁定连接钉、Φ5.0锁定螺钉 Note:use Φ5.0 Locking Connecting Screw ,Φ 5.0 Locking Screw .			

锁定连接钉 33072
Locking Connecting Screw



产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
3.5 锁定连接钉 3.5 Locking Connecting Screw	钛合金 Titanium Alloy	Φ3.5x10	33072035
5.0 锁定连接钉 5.0 Locking Connecting Screw		Φ5.0x14	33072050

Locking Plate Series

胫骨远端前外侧锁定钢板（左右） T4630
Distal Anterolateral Tibial Locking Plate
(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
胫骨远端前外侧 锁定钢板（左右） Distal Anterolateral Tibial Locking Plate (Left,Right)	纯钛 Pure Titanium	5	104	T4630005R/L
		7	130	T4630007R/L
		9	156	T4630009R/L
		11	182	T4630011R/L
		13	208	T4630013R/L
备注：选用Φ3.5锁定螺钉，HJA3.5螺钉（图例为11孔） Note:used For Φ3.5 Locking Screw, HA3.5 Screw (Sample:11holes)				

胫骨远端外侧重建锁定钢板（左右） T4600
Distal Lateral Tibial Reconstuction Locking Plate
(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
胫骨远端外侧重建 锁定钢板（左右） Distal Lateral Tibial Reconstruction Locking Plate (Left,Right)	纯钛 Pure Titanium	5	125	T4600005R/L
		7	161	T4600007R/L
		9	197	T4600009R/L
		11	233	T4600011R/L
		13	269	T4600013R/L
备注：选用Φ5.0锁定螺钉（图例为7孔） Note:used For Φ5.0 Locking Screw (Sample:7holes)				

胫骨远端内侧锁定钢板（左右） T4620
Distal Medial Tibial Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
胫骨远端内侧 锁定钢板（左右） Distal Medial Tibial Locking Plate (Left,Right)	纯钛 Pure Titanium	5	139	T4620005R/L
		7	167	T4620007R/L
		9	195	T4620009R/L
		11	223	T4620011R/L
备注：选用Φ3.5锁定螺钉,HA3.5螺钉（图例为7孔） Note:Used For Φ3.5 Locking Screw, HA3.5 Screw (Sample:7holes)				

2.7 腓骨远端重建锁定钢板（左右） T4405
2.7 Distal Fibula Reconstruction Locking Plate
(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
2.7腓骨远端重建 锁定钢板（左右） 2.7Distal Fibula Reconstruction Locking Plate (Left,Right)	纯钛 Pure Titanium	4	68	T4405004R/L
		6	88	T4405006R/L
		8	108	T4405008R/L
		10	128	T4405010R/L
		12	148	T4405012R/L
备注：选用Φ2.7锁定螺钉,HA2.7螺钉（图例为6孔） Note:used For Φ2.7 Locking Screw, HA2.7 Screw (Sample:6holes)				

跟骨锁定钢板I型(左右) T4400
Calcaneus Locking Plate-I(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
跟骨锁定钢板 I 型 (左右) Calcaneus Locking Plate-I (Left,Right)	纯钛 Pure Titanium	8	49	T4400008R/L
		10	59	T4400010R/L

备注: 选用Φ3.5锁定螺钉 (图例为8孔)
Note:used For Φ3.5 Locking Screw (Sample:8 Holes)

跟骨锁定钢板II型(大号)(左右) T4402
Calcaneus Locking Plate-II(Max)(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
跟骨锁定钢板 II 型 (大号) (左右) Calcaneus Locking Plate-II (Max) (Left,Right)	纯钛 Pure Titanium	12	67	T4402012R/L

备注: 选用Φ3.5锁定螺钉
Note:used For Φ3.5 Locking Screw

跟骨锁定钢板II型(中号)(左右) T4401
Calcaneus Locking Plate-II(Mid)(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
跟骨锁定钢板 II 型 (中号) (左右) Calcaneus Locking Plate-II (Mid) (Left,Right)	纯钛 Pure Titanium	12	58	T4401012R/L

备注: 选用Φ3.5锁定螺钉
Note:used For Φ3.5 Locking Screw

跟骨锁定钢板II型(小号)(左右) T4401
Calcaneus Locking Plate-II(Mini)(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
跟骨锁定钢板 II 型 (小号) (左右) Calcaneus Locking Plate-II (Mini) (Left,Right)	纯钛 Pure Titanium	8	51	T4401008R/L

备注: 选用Φ3.5锁定螺钉
Note:used For Φ3.5 Locking Screw

跟骨锁定钢板II型(短尾)(左右) T4404
Calcaneus Locking Plate-II(short tail)(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
跟骨锁定钢板 II 型 (短尾) (左右) Calcaneus Locking Plate-II (short tail) (Left,Right)	纯钛 Pure Titanium	11	59	T4404011R/L
		13	69	T4404013R/L

备注: 选用Φ3.5锁定螺钉, HB4.0螺钉 (图例为13孔)
Note:used For Φ3.5 Locking Screw, HB4.0screw (Sample:13 Holes)

三角型锁定螺钉 T480
Triangle Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
三角型锁定螺钉 Triangle Locking Screw	钛合金 Titanium Alloy	2.0	6、8、10、12、 14、16、18、20 22	T4802006- T4802022

梅花型锁定螺钉 T481/T482
Torx Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型锁定螺钉 Torx Locking Screw	钛合金 Titanium Alloy	2.5	6、8、10、12、 14、16、18、20 22、24、26、28、 30	T4812506- T4812530
		2.7	6、8、10、12、 14、16、18、20 22、24、26、28	T4822706- T4822728
			30、35、40、45、 50、55、60	T4822730- T4822760

梅花型锁定螺钉 T483/T484
Torx Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型锁定螺钉 Torx Locking Screw	钛合金 Titanium Alloy	3.5	10、12、14、16、 18、20、22、24、 26、28、30、32、 34、36、38、40、 42、44、46、48	T4833510- T4833548
			50、55、60、65、 70、75、80	T4833550- T4833580
		5.0	10、12、14、16、 18、20、22、24、 26、28、30、32、 34、36、38、40、 42、44、46、48	T4848015- T4845048
			50、55、60、65、 70、75、80、85、 90、95、100	T4848050- T4845100

梅花型中空锁定螺钉 T487
Torx Hollow Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型中空 锁定螺钉 Torx Hollow Locking Screw	钛合金 Titanium Alloy	6.5	30、35、40、45、 50、55、60、65、 70、75、80、85、 90、95、100、 105、110	T48765030- T48765110

梅花型中空锁定螺钉 T488
Torx Hollow Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型中空 锁定螺钉 Torx Hollow Locking Screw	钛合金 Titanium Alloy	7.0	60、65、70、75、 80、85、90、95、 100、105、110、 115	T48870060- T48870115

Metallic Bone Fracture Screw Series

φ2.0 金属接骨螺钉（三角型） T54120
Metal Bone Fracture Screw (Triangle)



钉帽形式 Screw Head Type	材质 Material	螺纹直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
三角型 Triangle	钛合金 Titanium Alloy	2.0	6、8、10、12、 14、16、18、20、 22	T5412006- T5412022

φ2.5 金属接骨螺钉（梅花型） T54125
Metal Bone Fracture Screw (Torx Type)



钉帽形式 Screw Head Type	材质 Material	螺纹直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型 Torx Type	钛合金 Titanium Alloy	2.5	6、8、10、12、 14、16、18、20、 22、24、26	T5412506- T5412526

φ2.7 金属接骨螺钉（梅花型） T54127
Metal Bone Fracture Screw (Torx Type)



钉帽形式 Screw Head Type	材质 Material	螺纹直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型 Torx Type	钛合金 Titanium Alloy	2.7	6、8、10、12、 14、16、18、20、 22、24、26、28、 30	T5412706- T5412730

HA3.5 金属接骨螺钉（梅花型）（HA全牙） T54135
Metal Bone Fracture Screw (Torx Type))(HA Full Thread)



钉帽形式 Screw Head Type	材质 Material	螺纹直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型 Torx Type	钛合金 Titanium Alloy	3.5	12、14、16、18、 20、22、24、26、 28、30、32、34、 36、38、40、42、 44	T5413512- T5413544

HA4.5 金属接骨螺钉（梅花型）（HA全牙） T54345
Metal Bone Fracture Screw (Torx Type))(HA Full Thread)



钉帽形式 Screw Head Type	材质 Material	螺纹直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型 Torx Type	钛合金 Titanium Alloy	4.5	20、22、24、26、 28、30、32、34、 36、38、40、42、 44、46、48、50	T5434520- T5434550
			55、60	T5434555- T5434560

HB4.0 金属接骨螺钉（梅花型）（HB全牙） T545
Metal Bone Fracture Screw (Torx Type))(HB Full Thread)




钉帽形式 Screw Head Type	螺纹形式 Thread Type	材质 Material	螺纹直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型 Torx Type	HB全螺纹 HB Fully Threaded	钛合金 Titanium Alloy	4.0	12、14、16、18、 20、22、24、26、 28、30、32、34、 36、38、40、42、 44、46、48、50、 52、54、56、58、 60	T5454012- T54540605

HB6.5 金属接骨螺钉（梅花型）（HB全牙） T549
Metal Bone Fracture Screw (Torx Type))(HB Full Thread)



钉帽形式 Screw Head Type	螺纹形式 Thread Type	材质 Material	螺纹直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型 Torx Type	HB全螺纹 HB Fully Threaded	钛合金 Titanium Alloy	6.5	30、35、40、45、 50、55、60、65、 70、75、80、85	T5496530- T5496585




肋骨钢板系列 Rib Plate Series


纯钛肋骨板（无菌包装）



Pure Titanium Rib Plate (Sterile Packaging)



肋骨接骨板 T132
Rib Plate




产品名称 Product Name	包装形式 Packaging Form	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.	
肋骨接骨板 Rib Plate	洁净包装 Clean Packaging	纯钛 Pure Titanium	11	49	T1321149	
			15	34	T1321534	
			15	49	T1321545	
			15	59	T1321559	
			19	34	T1321934	
			19	49	T1321949	
	无菌包装 Sterile Packaging		19	59	T1321959	
			11	49	T1321149-S	
			15	34	T1321534-S	
			15	49	T1321545-S	
			15	59	T1321559-S	
			19	34	T1321934-S	
			19	49	T1321949-S	
			19	59	T1321959-S	

肋骨安装器械包 91719902
Rib Plate Instrument Set




PDF 文件
PDF File



JY Locking Plate Series (Sterile Packaging)

无菌 JY 锁定肋骨内固定系统

Sterile JY rib locking fixation system



系统特点

Character Of The System

- 无需去除骨膜，保证充足血供
- 锁定板提供稳定固定
- 预塑形板，减少手术时间
- 低切迹设计
- 限深钻头，防止钻孔过深
- No Need To Remove Periosteum, Ensure Adequate Blood Supply.
- Lockong Plate Provides Stable Fixtion.
- Pre-plastion Plate,reduce Surgery Time.
- Low Profile Design.
- Depth Limit Drill,prevent Drilling Too Deep.

JY 锁定钢板系列（无菌包装）

JY Locking Plate Series (Sterile Packaging)

无菌 JY 锁定肋骨内固定系统

Sterile JY rib locking fixation system




无菌十字型锁定螺钉（阳极氧化）35008

Sterile Crossed Screw (Anodized)



产品名称 Product Name	材质 Material	直径 Diameter	长度 Length	产品编号 Product No.
无菌十字型锁定螺钉(阳极氧化) Sterile Crossed Screw (Anodized)	钛合金 Titanium Alloy	3.0	6	35008006
			8	35008008
			10	35008010
			12	35008012
			14	35008014



普通钢板螺钉系列
Common Plate Screw Series

Metallic Bone Plate Series

指掌骨钢板系列
Phalangeal/Metacarpal Plate Series

微型钢板中号直型 T7510
Straight Micro Plate (Mid)



微型钢板中号T型 T7520
T Micro Plate (Mid)



微型钢板中号L型 (左右) T7530
L Micro Plate (Mid)(Left,Right)



微型钢板大号直型 T7610
Straight Micro Plate (Large)



微型钢板大号T型 T7620
T Micro Plate (Large)



微型钢板大号L型 (左右) T7630
L Micro Plate (Large)(Left,Right)



小T型重建钢板 T2010
Small T Reconstruction Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
微型钢板中号直型 Straight Micro Plate (Mid)	纯钛 Pure Titanium	4	23	T7510004
		5	29	T7510005
		6	35	T7510006
		7	41	T7510007
微型钢板中号T型 T Micro Plate (Mid)		4	18	T7520004
		5	24	T7520005
		6	30	T7520006
		7	36	T7520007
微型钢板中号L型 (左右) L Micro Plate (Mid)(Left,Right)		4	18	T7530004R/L
		5	24	T7530005R/L
		6	30	T7530006R/L

备注：选用HA2.0螺钉
Note:Used For HA2.0 Screw

产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
微型钢板大号直型 Straight Micro Plate (Large)	纯钛 Pure Titanium	4	31	T7610004
		5	39	T7610005
		6	47	T7610006
		8	63	T7610008
微型钢板大号T型 T Micro Plate (Large)		5	32	T7620005
		6	40	T7620006
		7	48	T7620007
		8	56	T7620008
微型钢板大号L型 (左右) L Micro Plate (Large)(Left,Right)		5	31	T7630005R/L
		6	39	T7630006R/L
		7	47	T7630007R/L

备注: 选用HA2.5螺钉
Note:Used For HA2.5 Screw

产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
小T型重建钢板 Small T Reconstruction Plate	纯钛 Pure Titanium	6	43	T2010006
		7	54	T2010007
		8	65	T2010008

备注: 头部选用HB4.0螺钉, 其余选用HA3.5螺钉 (图例为7孔)
Note: Used For HB4.0 Screw In Head, Other HA3.5 Screw (Sample: 7 Holes)

Metallic Bone Plate Series

骨盆钢板 T3820
Curved Pelvic Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
骨盆钢板 Curved Pelvic Plate	纯钛 Pure Titanium	5	60	T3820005
		6	72	T3820006
		8	96	T3820008
		10	120	T3820010
		12	144	T3820012
备注：选用HB4.0螺钉 (图例为10孔) Note: Used For HB4.0 Screw (Sample: 10 Holes)				

髌臼后壁钢板(左右) T3821
Posterior Wall Acetabular Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
髌臼后壁钢板(左右) Posterior Wall Acetabular Plate (Left,Right)	纯钛 Pure Titanium	9	T3821009R/L
备注：选用HB4.0螺钉 Note: Used For HB4.0 Screw			
专利产品 Patented Product			

弓状线耻骨钢板(左右) T3840
Pubis Lineae Actuate Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
弓状线耻骨钢板(左右) Pubis Lineae Actuate Plate (Left,Right)	纯钛 Pure Titanium	12	T3840012R/L
备注：选用HB4.0螺钉 Note: Used For HB4.0 Screw			
专利产品 Patented Product			

髌耻联合钢板(左右) T3850
Symphysis iliopectineal Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
髌耻联合钢板(左右) Symphysis iliopectineal Plate (Left,Right)	纯钛 Pure Titanium	13	T3850013R/L
备注：选用HB4.0螺钉 Note: Used For HB4.0 Screw			
专利产品 Patented Product			

Metallic Bone Plate Series

骶髂关节前路钢板 T3880
Anterior Articulation Sacroiliac Plate



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
骶髂关节前路钢板 Anterior Articulation Sacroiliac Plate	纯钛 Pure Titanium	4	T3880004
备注：选用HB4.0螺钉 Note: Used For HB4.0 Screw			

耻骨联合钢板 T3851
Symphysis Pubis Plate



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
耻骨联合钢板 Symphysis Pubis Plate	纯钛 Pure Titanium	8	T3851008
备注：选用HB4.0螺钉 Note: Used For HB4.0 Screw			

髋臼后壁钩钢板 T3821
Posterior Wall Hook Acetabular Plate



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
髋臼后壁钩钢板 Posterior Wall Hook Acetabular Plate	纯钛 Pure Titanium	3	T3821003
备注：选用HB4.0螺钉 Note: Used For HB4.0 Screw			
专利产品 Patented Product			

PM骨盆重建板 T1330
PM Pelvic Reconstruction Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
PM骨盆重建板 PM Pelvic Reconstruction Plate	纯钛 Pure Titanium	8	115	T1330008
		12	175	T1330012
备注：选用HB4.0螺钉（图例为8孔） Note: Used For HB4.0 Screw (Sample: 8 Holes)				

Metallic Bone Plate Series

限制型加压钢板(上肢) T1110
Restricted Compression Plate (Upper Limb)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
限制型加压钢板 (上肢) Restricted Compression Plate (Upper Limb)	纯钛 Pure Titanium	6	77	T1110006
		7	89	T1110007
		8	101	T1110008
		10	125	T1110010
		12	149	T1110012
备注: 选用HA3.5螺钉 (图例为6孔) Note: Used For HA3.5 Screw (Sample: 6 Holes)				

限制型加压钢板(胫骨) T1120
Restricted Compression Plate (Tibia)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
限制型加压钢板 (胫骨) Restricted Compression Plate (Tibia)	纯钛 Pure Titanium	7	130	T1120007
		8	146	T1120008
		10	178	T1120010
		12	210	T1120012
备注: 选用HA4.5螺钉 (图例为8孔) Note: Used For HA4.5 Screw (Sample: 8 Holes)				

LC-DCP 尺桡骨钢板 T1220
LC-DCP Ulna Radius Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
LC-DCP 尺桡骨钢板 LC-DCP Ulna Radius Plate	纯钛 Pure Titanium	5	62	T1220005
		6	74	T1220006
		7	86	T1220007
		8	98	T1220008
		10	122	T1220010
备注：选用HA3.5螺钉（图例为6孔） Note: Used For HA3.5 Screw (Sample: 6 Holes)				

1/3 圆钢板 T1410
One-third Tubular Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
1/3 圆钢板 One-third Tubular Plate	纯钛 Pure Titanium	5	61	T1410005
		6	73	T1410006
		7	85	T1410007
		8	97	T1410008
		10	121	T1410010
备注：选用HA3.5螺钉 (图例为6孔) Note: Used For HA3.5 Screw (Sample: 6 Holes)				

重建钢板 T1310
Reconstruction Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
重建钢板 Reconstruction Plate	纯钛 Pure Titanium	5	65	T1310005
		6	78	T1310006
		8	104	T1310008
		10	130	T1310010
		12	156	T1310012
备注：选用HA3.5、HB4.0螺钉（图例为6孔） Note: Used For HA3.5, HB4.0 Screw (Sample: 6 Holes)				

Metallic Bone Fracture Screw Series

φ1.5 金属接骨螺钉（三角型） T54115
Metal Bone Fracture Screw (Triangle)



钉帽形式 Screw Head Type	材质 Material	螺纹直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
三角型 Triangle	钛合金 Titanium Alloy	1.5	6、8、10、12、14	T5411506- T5411514

φ2.0 金属接骨螺钉（三角型） T54120
Metal Bone Fracture Screw (Triangle)



钉帽形式 Screw Head Type	材质 Material	螺纹直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
三角型 Triangle	钛合金 Titanium Alloy	2.0	6、8、10、12、14、16、18、20、22	T5412006- T5412022

φ2.5 金属接骨螺钉（梅花型） T54125
Metal Bone Fracture Screw (Torx Type)



钉帽形式 Screw Head Type	材质 Material	螺纹直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型 Torx Type	钛合金 Titanium Alloy	2.5	6、8、10、12、14、16、18、20、22、24、26	T5412506- T5412526

φ2.7 金属接骨螺钉（梅花型） T54127
Metal Bone Fracture Screw (Torx Type)



钉帽形式 Screw Head Type	材质 Material	螺纹直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型 Torx Type	钛合金 Titanium Alloy	2.7	6、8、10、12、14、16、18、20、22、24、26、28、30	T5412706- T5412730

HA3.5 金属接骨螺钉（梅花型）（HB全牙） T54135
Metal Bone Fracture Screw (Torx Type)(HA Full Thread)



钉帽形式 Screw Head Type	材质 Material	螺纹直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型 Torx Type	钛合金 Titanium Alloy	3.5	12、14、16、18、20、22、24、26、28、30、32、34、36、38、40、42、44	T5413512- T5413544

Metallic Bone Fracture Screw Series

HA4.5 金属接骨螺钉（梅花型）（HA全牙） T54345
Metal Bone Fracture Screw (Torx Type)(HA Full Thread)



钉帽形式 Screw Head Type	材质 Material	螺纹直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型 Torx Type	钛合金 Titanium Alloy	4.5	20、22、24、26、 28、30、32、34、 36、38、40、42、 44、46、48、50	T5434520- T5434550
			55、60	T5434555- T5434560

HB4.0 金属接骨螺钉（梅花型）（HB全牙） T545
Metal Bone Fracture Screw (Torx Type)(HB Full Thread)



钉帽形式 Screw Head Type	螺纹形式 Thread Type	材质 Material	螺纹直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型 Torx Type	HB全螺纹 HB Fully Threaded	钛合金 Titanium Alloy	4.0	12、14、16、18、 20、22、24、26、 28、30、32、34、 36、38、40、42、 44、46、48、50、 52、54、56、58、 60	T5454012- T54540605

HB6.5 金属接骨螺钉（梅花型）（HB全牙） T549
Metal Bone Fracture Screw (Torx Type)(HB Full Thread)



钉帽形式 Screw Head Type	螺纹形式 Thread Type	材质 Material	螺纹直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型 Torx Type	HB全螺纹 HB Fully Threaded	钛合金 Titanium Alloy	6.5	30、35、40、45、 50、55、60、65、 70、75、80、85	T5496530- T5496585

HB6.5 金属接骨螺钉（梅花型）（HB半牙） 548
Metal Bone Fracture Screw (Torx Type)(HB Half Thread)



钉帽形式 Screw Head Type	螺纹形式 Thread Type	材质 Material	螺纹直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型 Torx Type	HB半螺纹 HB Half Threaded	钛合金 Titanium Alloy	6.5	50、55、60、65、 70、75、80、85、 90	T5486550- T5486590




外固定支架系列 External Fixation Series

SSK 外固定支架系统（上肢综合包）



SSK External Fixation System(Upper limb integrated instrument set)




图号 Figure No.	产品编号 Product No.	产品名称 Product Name	规格型号 Specification	数量 Qty
1	89601	针杆固定夹 Pin-Rod fixation Clamp	5*4	12
2	89602	杆杆固定夹 Pin-Rod fixation Clamp	5*5	12
3	89603	多孔固定夹 Multi-holes fixation clamp	/	3
4	89604	支柱 Support	5*0°	2
5	89605	支柱 Support	5*30°	4
6	89610	碳纤维棒 Carbon Fiber Rod	5*250	4
6	89611	碳纤维棒 Carbon Fiber Rod	5*200	4
6	89612	碳纤维棒 Carbon Fiber Rod	5*180	4
6	89613	碳纤维棒 Carbon Fiber Rod	5*150	4
6	89614	碳纤维棒 Carbon Fiber Rod	5*100	4
7	89606	骨牵引针 Bone Traction Pin	4*130	8

图号 Figure No.	产品编号 Product No.	产品名称 Product Name	规格型号 Specification	数量 Qty
7	89607	骨牵引针 Bone Traction Pin	4*120	8
7	89608	骨牵引针 Bone Traction Pin	4*100	8
7	89609	骨牵引针 Bone Traction Pin	4*80	8
8	89622	L型扳手 L Wrench	SW2.5	1
9	89617	定位器 Locator	5*4*3	1
10	89621	弓形钻架 Arch Drill Frame	4	1
11	89620	T型扳手 L Wrench	上肢 Upper Limb	1
12	89619	上钉器 Screwdriver	4	1
13	89623	钻头 Drill Bit	3.0*150	2
14	89616	稳定复位扳手 Stable Reduction Wrench	上肢 Upper Limb	1
15	89615	肘关节活动器 Elbow joint movement device	5	1
16	89618	拇指轮 Speedy Screwdriver	/	1

SSK 外固定支架系统（下肢综合包）



SSK External Fixation System(Lower limb integrated instrument set)



图号 Figure No.	产品编号 Product No.	产品名称 Product Name	规格型号 Specification	数量 Qty
1	89701	针杆固定夹 Pin-Rod fixation Clamp	8*5	12
2	89702	杆杆固定夹 Pin-Rod fixation Clamp	8*8	12
3	89703	五孔固定夹 Five-holes fixation clamp	/	3
4	89704	支柱 Support	8*0°	2
5	89705	支柱 Support	8*30°	4
6	89710	碳纤维棒 Carbon Fiber Rod	8*150	4
6	89711	碳纤维棒 Carbon Fiber Rod	8*200	4
6	89712	碳纤维棒 Carbon Fiber Rod	8*250	4
6	89713	碳纤维棒 Carbon Fiber Rod	8*300	4
6	89714	碳纤维棒 Carbon Fiber Rod	8*350	4
6	89715	碳纤维棒 Carbon Fiber Rod	8*400	4
7	89706	骨牵引针 Bone Traction Pin	5*200	8
7	89707	骨牵引针 Bone Traction Pin	5*150	8

图号 Figure No.	产品编号 Product No.	产品名称 Product Name	规格型号 Specification	数量 Qty
7	89708	骨牵引针 Bone Traction Pin	5*130	8
7	89709	骨牵引针 Bone Traction Pin	5*120	8
8	89717	定位器 Locator	6*5*4	1
9	89720	弓形钻架 Arch Drill Frame	5	1
10	89716	稳定复位扳手 Stable Reduction Wrench	下肢 Lower Limb	1
11	89721	半环连接杆 Semi-ring Connecting Rod	8*160	1
11	89722	半环连接杆 Semi-ring Connecting Rod	8*180	1
11	89723	半环连接杆 Semi-ring Connecting Rod	8*200	1
12	89724	扭力扳手 Torque wrench	/	1
13	89726	上钉器扳手 Screwdriver Wrench	下肢多功能5 Lower limb multifunction5	1
14	91114001	钻头 Drill Bit	4*300	2
15	89618	拇指轮 Speedy Screwdriver	/	1

J

SSK External Fixation System - Application

胫骨近端骨折
Proximal Tibial
External Fixator



股骨干骨折
Femur Shaft External
Fixator



半环形 (标准型)
Half Ring (Standard)
External Fixator



股骨近端骨折
Proximal Femoral
External Fixator



膝关节固定
Knee Joint External
Fixator



骨盆固定
Pelvic External Fixator



踝关节固定 (A型)
Ankle External Fixator
(A Type)



桡骨关节端固定
Radius Joint External
Fixator



前臂架
Forearm External
Fixator



肘关节骨折
Elbow External Fixator



肱骨近端骨折
Proximal Humeral
External Fixator



指骨骨折
Mini Fragment Fixator
(Finger and palm)



SSk External Fixation System - Emergency Sterile Solution



系统特点

System features

- 无菌包装使用方便，减少临床消毒时间。
 - 三维固定模块，镶嵌式设计，易于安装。
 - 可根据术中要求自由装配。
 - 碳纤维棒材质连接杆具有强度高、耐腐蚀、质量轻等特点。
 - 适用于多种创伤和骨病治疗。
 - 骨针采用自攻自钻锥形螺纹设计，易于拧入。
 - 整体设计合理，术后减少患者负担。
- Sterile Packaging Is Easy To Use And Reduces Clinical Disinfection Time.
 - Three Dimensional Fixing Module, inlaid Design, Easy To Assemble.
 - Freely Assemble According To Surgical Requirement.
 - Carbon Fiber Rod Material Connection Rod Has High Strength\corrosion Resistance\light Weight And Other Characteristics.
 - Suitable For A Variety Of Trauma And Bone Disease Treatment.
 - Self-tapping & Drilling Taper Thread Design, Easy To Insert. The Overall Design Is Reasonable And The Postoperative Burden Of Patients Is Reduced.


SSK External Fixation System – Emergency Sterile Solution

肱骨固定应用方案
Humeral Fixation Application Plan




肱骨建议配置 Humeral Recommended configuration			
产品名称 Product Name	规格 Specification	配置数量 Configured Qty	产品编号 Product No.
针杆固定夹 Pin-To-Rod Clamp	Φ5×8	4个	84001001
骨牵引针 Bone Traction Pin	杆5钉4×120	4根	84040120-S
碳纤维棒 Carbon Fiber Rod	Φ8×200	1根	84004200

胫骨固定应用方案
Tibial Fixation Application Plan



胫骨建议配置 Tibial Recommended configuration			
产品名称 Product Name	规格 Specification	配置数量 Configured Qty	产品编号 Product No.
针杆固定夹 Pin-To-Rod Clamp	Φ5×8	4个	84001001
骨牵引针 Bone Traction Pin	Φ5×150	4根	84050150-S
碳纤维棒 Carbon Fiber Rod	Φ8×250	1根	84004250

股骨固定应用方案
Femoral Fixation Application Plan



股骨建议配置 Femoral Recommended configuration			
产品名称 Product Name	规格 Specification	配置数量 Configured Qty	产品编号 Product No.
针杆固定夹 Pin-To-Rod Clamp	Φ5×8	4个	84001001
骨牵引针 Bone Traction Pin	Φ5×180	4根	84050180-S
碳纤维棒 Carbon Fiber Rod	Φ8×300	1根	84004300

SSK External Fixation System - Emergency Sterile Solution

针杆固定夹
Pin-To-Rod Clamp



产品名称 Product Name	包装形式 Packaging Form	材质 Material	直径(mm) Diameter	产品编号 Product No.
针杆固定夹 Pin-To-Rod Clamp	无菌包装 Sterile Packaging	不锈钢 Stainless steel	Φ5×8	84001001

骨牵引针
Bone Traction Pin



产品名称 Product Name	包装形式 Packaging Form	材质 Material	直径(mm) Diameter	产品编号 Product No.
骨牵引针 Bone Traction Pin	无菌包装 Sterile Packaging	不锈钢 Stainless steel	杆5钉4×120	84040120-S
			Φ5×120	84050120-S
			Φ5×150	84050150-S
			Φ5×180	84050180-S

碳纤维棒
Carbon Fiber Rod



产品名称 Product Name	包装形式 Packaging Form	材质 Material	直径(mm) Diameter	产品编号 Product No.
碳纤维棒 Carbon Fiber Rod	无菌包装 Sterile Packaging	碳纤维 Carbon fiber	Φ8×150	84004150
			Φ8×200	84004200
			Φ8×250	84004250
			Φ8×300	84004300
			Φ8×350	84004350
			Φ8×400	84004400

上钉器
Screwdriver



产品名称 Product Name	材质 Material	直径(mm) Diameter	产品编号 Product No.
上钉器 Screwdriver	不锈钢 Stainless steel	Φ5	84100010

定位器
Locator



产品名称 Product Name	材质 Material	直径(mm) Diameter	产品编号 Product No.
定位器 Locator	不锈钢 Stainless steel	Φ6×5×4	84100020

钻头
Drill Bit



产品名称 Product Name	材质 Material	直径(mm) Diameter	产品编号 Product No.
钻头 Drill Bit	不锈钢 Stainless steel	Φ3.0×150	85100150
		Φ4.0×300	85100010

锁定螺钉系列

	Φ3.5锁定螺钉	Φ5.0锁定螺钉	Φ6.5锁定螺钉	Φ6.5空心锁定螺钉	Φ7.0空心锁定螺钉
钉帽					
钻头	Φ2.8×150	Φ4.3×280	Φ5.0×280	Φ5.0×250空心钻	Φ5.0×250空心钻
丝锥	/	/	/	/	/
起子	T3.4	T4.3	T5.3	T5.3	T5.3
导针	/	/	/	Φ3.0	Φ3.0

	Φ2.0锁定螺钉	Φ2.5锁定螺钉	Φ2.7锁定螺钉	TML 1.5 锁定螺钉	TML 2.0 锁定螺钉
钉帽					
钻头	Φ1.5×150	Φ2.0×100	Φ2.1×150	Φ1.0×85	Φ1.5×85
丝锥	/	/	/	/	/
起子	△Φ2.0	T2.3	T2.3	T1.3	T1.6
导针	/	/	/	/	/

金属接骨螺钉系列

	HA3.5皮质骨螺钉	HB4.0松质骨螺钉	HA4.5皮质骨螺钉	HB6.5松质骨螺钉
钉帽				
钻头	Φ2.5×125	Φ3.0×125	Φ3.2×150	Φ4.5×150
丝锥	HA3.5	HB4.0	HA4.5	HB6.5
起子	T3.4	T3.4	T4.3	T4.3
导针	/	/	/	/

	HA1.5皮质骨螺钉	HA2.0皮质骨螺钉	HA2.5皮质骨螺钉	HA2.7皮质骨螺钉
钉帽				
钻头	Φ1.1×80	Φ1.5×80	Φ1.8×80	Φ2.0×80
丝锥	HA1.5	HA2.0	HA2.5	HA2.7
起子	△2.0	△2.0	T2.3	T2.3
导针	/	/	/	/

空心螺钉系列

	Φ2.5中空无头加压螺钉	Φ3.5中空无头加压螺钉	Φ4.5中空无头加压螺钉
钉帽			
钻头	Φ1.7台阶空心钻	Φ2.1台阶空心钻	Φ2.0/Φ2.4×120台阶空心钻
丝锥	/	/	/
起子	△2.0	T2.3	T2.9
导针	Φ0.8×120	Φ1.0×120	Φ1.0×120

空心螺钉系列

	Φ2.5中空双头加压螺钉	Φ3.0中空双头加压螺钉	Φ4.0中空双头加压螺钉
钉帽			
钻头	Φ1.9空心钻	Φ2.0/Φ2.4×120台阶空心钻	Φ3.0×150空心钻
丝锥	/	/	Φ4.0
起子	T2.3	SW2.0	SW3.0
导针	Φ1.0×120	Φ1.0×120	Φ1.5×150

	Φ3.5中空加压螺钉	Φ4.0中空加压螺钉	Φ4.5中空加压螺钉
钉帽			
钻头	Φ2.5×150空心钻	Φ3.0×150空心钻	Φ3.2×150空心钻
丝锥	Φ3.5	Φ4.0	Φ4.5
起子	SW2.5	SW2.5	SW3.0
导针	Φ1.0×150	Φ1.2×150	Φ1.5×150

	Φ6.5双头中空加压螺钉	Φ6.5中空加压螺钉	Φ7.3中空加压螺钉
钉帽			
钻头	Φ5.0×250空心钻	Φ5.0×250空心钻	Φ5.0×250空心钻
丝锥	/	/	Φ7.3
起子	SW6.0套筒	SW4.0	SW4.0
导针	Φ2.5×250	Φ2.5×250	Φ2.5×250

阳极氧化螺钉系列

	Φ2.0锁定螺钉	Φ2.7锁定螺钉	Φ3.5锁定螺钉	Φ5.0锁定螺钉	Φ2.7万向锁定螺钉
钉帽					
钻头	Φ1.5×150	Φ2.1×150	Φ2.9×205	Φ4.3×280	Φ2.1×150
丝锥	/	/	/	/	/
起子	ΔΦ2.0	T2.3	T3.4	T4.3	T2.3
导针	/	/	/	/	/

阳极氧化螺钉系列

	Φ1.5 万向锁定螺钉	Φ2.0万向锁定螺钉	Φ3.5万向锁定螺钉	HA3.5皮质骨螺钉	HA4.5皮质骨螺钉
钉帽					
钻头	Φ1.1	Φ1.5	Φ2.8×150	Φ2.5×125	Φ3.2×150
丝锥	/	/	/	HA3.5	HA4.5
起子	ΔΦ2.0	ΔΦ2.0	T3.4	T3.4	T4.3
导针	/	/	/	/	/



康力微信公众号
WeChat Official Account



康力网站
Kangli Website



电子图本信息
Catalogue

 苏州市康力骨科器械有限公司
SUZHOU KANGLI ORTHOPAEDICS INSTRUMENT CO., LTD.

厂 址: 中国江苏张家港市南苑东路(塘桥镇)
邮 编: 215600
免费服务热线: 400-8282120 800-8282456
电 话: 0512-58955188 58560562 58358969
传 真: 0512-58561012
网 址: www.kangliMed.com、骨科器械.中国
邮 箱: ask@KangliMed.com

Address: Nanyuan East Road, Zhangjiagang,
Jiangsu, China 215600
Site: KangliMed.com
Email: ask@KangliMed.com



本目录所载之规格, 本公司保有变更之权利, 产品之内容、规格、颜色如有不符, 以实物为准。
The company reserves the right of final interpretation of this catalog.
In case of discrepancies, the actual product shall prevail.