

康力·器械产品目录

Kangli · Instrument Product Catalogue



20250530版

康力·器械产品目录

Kangli · Instrument Product Catalogue

康力·器械产品目录

A


脊柱系列安装器械包 Spinal Fixation Series



KCO 3.2mm 钉棒器械包 93219903·····	A 01
KCO 3.2mm 钉棒器械包 93219901·····	A 05
颈椎后路椎板成形系统器械包 93229901·····	A 08
颈椎前路钢板器械包 93239901·····	A 09
KSS-I6.0mm 脊柱固定器安装器械包 93319901·····	A 11
KSS-II 5.5mm 脊柱固定器安装器械包 93339901·····	A 15
KSS-IV 脊柱内固定系统安装器械包 93369901·····	A 19
JZPA 5.5/6.0 脊柱内固定系统安装器械包 93359901·····	A 23
JZPA 5.5 脊柱内固定系统安装器械包 93329901·····	A 29
JZWA 脊柱微创内固定系统安装器械包 93519901·····	A 33
零切迹融合器器械包 93249901·····	A 37
PEEK 胸腰椎融合器安装器械包 93649901·····	A 41
PEEK 颈椎融合器安装器械包 93659901·····	A 43
颈椎椎间融合器安装器械包 93619902·····	A 45
颈椎椎间融合器安装器械包 93619901·····	A 47
腰椎椎间融合器安装器械包 93629902·····	A 49
腰椎椎间融合器安装器械包 93629901·····	A 51
钛笼系统安装器械包 93639901·····	A 53
钛笼系统安装器械包 93639902·····	A 54

B

颌面系列安装器械包 Multifunction craniomaxillofacial fixation system



MCF 颌面钢板安装器械包 91279901·····	B 57
MCF 正颌面钢板安装器械包 91279901·····	B 59

C

微型板系列安装器械包

Multifunction craniomaxillofacial fixation system



TML 1.5 微型锁定指掌骨安装器械包 91859901·····	C 63
TML 2.0 微型锁定指掌骨安装器械包 91829901·····	C 65
2.0锁定螺钉系列器械包 91319901·····	C 67
Φ2.0锁定钢板安装器械包 PPSU 91319902·····	C 69
2.5锁定螺钉系列器械包 91389901·····	C 71
Φ2.5锁定钢板安装器械包 PPSU 91389902·····	C 73
2.7锁定螺钉系列器械包 91399901·····	C 75
微型钢板安装器械包 91219901·····	C 77
微型钢板安装器械包 91219902·····	C 79
微型钢板安装器械包 91219903·····	C 83

D

JY 锁定钢板螺钉系列

JY Locking Plate Screw Series




JY 2.7、3.5锁定钢板安装器械包 91819901·····	D 89
JY 2.7锁定钢板安装器械包 91349902·····	D 91
JY 2.7锁定钢板安装器械包 91349901·····	D 93
JY 3.5锁定钢板安装器械包 91359902·····	D 95
JY 3.5锁定钢板安装器械包 91359901·····	D 97
JY 5.0锁定钢板安装器械包 91369901·····	D 99
HTO锁定胫骨截骨安装器械包 91379901·····	D 101
LCC 锁定复合加压系统器械包 91749901·····	D 103

E

(原色) 锁定钢板螺钉系列

Kirschner Pin (Sterile Packaging)



上肢锁定安装器械包 91329901·····	E 107
上肢锁定接骨板安装器械包 (PPSU) 91329902·····	E 109
上肢锁定安装器械包 91329903·····	E 111
下肢锁定安装器械包 91339901·····	E 113
下肢锁定接骨板安装器械包 (PPSU) 91339902·····	E 115

F

普通钢板螺钉系列

Titanium Cable Fixation System(Sterile Packaging)



上肢骨折安装器械包 91229901·····	F 121
下肢骨折安装器械包 91239901·····	F 123
四肢骨折器械包 (重建) 91259901·····	F 125
鹅头钉安装器械包 91269901·····	F 129
断钉取出器安装器械包 91739901·····	F 131
断钉取出器安装器械包 91739902·····	F 135

G

空心钉、骨针、钛缆

Locking Connecting Plate Series



钛缆固定系统安装器械包 91729901·····	G 141
钛缆固定系统安装器械包 91729902·····	G 143
1.7 空心钉安装器械包 91459901·····	G 145
Φ7.3 空心钉安装器械包 91419901·····	G 147
Φ4.5 空心钉安装器械包 91429901·····	G 149
Φ4.5空心钉安装器械包 PPSU 91429902·····	G 150
Φ4.0 空心钉安装器械包 91439901·····	G 151
Φ4.0空心钉安装器械包 PPSU 91439902·····	G 152
Φ3.5 空心钉安装器械包 91449901·····	G 153
Φ3.5空心钉安装器械包 PPSU 91449902·····	G 154
Φ2.5 双头加压空心钉安装器械包 91519901·····	G 155
Φ3.0 双头加压空心钉安装器械包 91529901·····	G 156
Φ4.0 双头加压空心钉安装器械包 91539901·····	G 157
Φ4.0 双头加压空心钉安装器械包 91539902·····	G 158
Φ4.0 双头加压空心钉安装器械包 PPSU 91539903·····	G 159
Φ4.0 双头加压空心钉安装器械包 91539904·····	G 160
2.5 无头加压空心钉安装器械包 91619901·····	G 161
Φ2.5 无头加压空心钉安装器械包 91619902·····	G 162
3.5 无头加压空心钉安装器械包 91629901·····	G 163
Φ3.5 无头加压空心钉安装器械包 91629902·····	G 164
4.5 无头加压空心钉安装器械包 91639901·····	G 165
Φ4.5 无头加压空心钉安装器械包 91639902·····	G 166

Φ4.5 无头加压空心钉安装器械包 91639903·····	G 167
Φ6.5无头加压空心钉安装器械包 91649901·····	G 168
Φ2.0 / Φ2.5 自攻自钻螺钉安装器械包 91249901 ·····	G 169

H

髓内钉系列 Pure Titanium Rib Olate (Sterile Packaging)



QML 多维股骨倒打髓内钉安装器械包 92279901 ·····	D 173
KTO-PFNA安装器械包加强型 短款 92219906·····	D 177
KTO-PFNA安装器械包加强型 长款 92219905·····	D 181
KTO-PFNA伽玛带锁髓内钉安装器械包 92219901·····	D 187
KTO-PFNA伽玛带锁髓内钉安装器械包(加长版) 92219902·····	D 191
KTO-PFNA伽玛带锁髓内钉安装器械包(碳纤维) 92219904·····	D 197
KTO滑动加压伽玛钉安装器械包 92229901·····	D 201
KTO滑动加压伽玛钉安装器械包(加长版) 92229902·····	D 207
KTO滑动加压伽玛钉安装器械包加强型 (小号) 92239903·····	D 213
KTO 滑动加压伽玛钉安装器械包小号加强型 (加长版) 92239904·····	D 219
KTO滑动加压伽玛钉安装器械包小号加强型(加长)(加强型)(碳纤维款)92239905·····	D 225
KTO滑动加压伽玛钉安装器械包(小号) 92239901 ·····	D 231
KTO滑动加压伽玛钉安装器械包小号(加长版) 92239902·····	D 237
股骨倒打带锁髓内钉安装器械包 92249901 ·····	D 243
KTO 股骨带锁髓内钉安装器械包 92259901·····	D 247
股骨带锁髓内钉安装器械包 92269901·····	D 253
QML 多维胫骨髁上入路髓内钉安装器械包 92339901·····	D 257
KTO 胫骨带锁髓内钉安装器械包 92319901·····	D 263
胫骨带锁髓内钉安装器械包 92329901·····	D 267
QML 多维肱骨带锁髓内钉安装器械包 92529901·····	D 271
KTO 肱骨带锁髓内钉安装器械包 92519901·····	D 277
金属弹性髓内针安装器械包 92619901·····	D 281
尺桡骨髓内钉安装器械包系统 92629901·····	D 283

A

B

C

D

E

F

G

H

I

J

K

L

I

外固定支架系列
Emergency aseptic external fixation system



- 外固定支架器械包（上肢） 91759901.....I 289
- 外固定支架器械包（下肢） 91769901.....I 290

J

肋骨钢板系列
Pure Titanium Rib Olate (Sterile Packaging)



- 肋骨安装器械包 91719901.....J 293
- 肋骨安装器械包 91719902.....J 294

K

足踝钢板系列
Ankle Plate Series




- 足踝钢板安装器械包.....K 297

L

其它产品系列
other product lines



- 骨卡环骑缝钉安装器械包 97119902.....L 303




- 恒天腰椎椎间融合器安装器械包 95219901.....L 305
- BES 2.5螺钉安装器械包 95139901.....L 307
- BES 3.5螺钉安装器械包 95109901.....L 308
- BES 3.5空心螺钉安装器械包 95119901.....L 309

器械编码规则

器械编码总则311-358



脊柱系列安装器械包





* * KCO 3.2mm 钉棒器械包 93219903

KCO 3.2 Screw Rod Instrument Set

A



93210101 直探针
Straight Probe



93210202 开路锥
Awl



93210105 显影针 (柱形)
Developing Pin (Cylinder Shape)



93210306 显影针 (球形)
Developing Pin (Spherical Shape)



93210712 丝锥 (三角快装)
Drill Bit (Quick Coupling), HB3.5



93210713 丝锥 (三角快装)
Drill Bit (Quick Coupling), HB4.0



93212501 枕骨钻头 (三角快装)
Drill Bit (Quick Coupling), $\Phi 2.4$



93212601 枕骨丝锥 (三角快装)
Drill Bit (Quick Coupling), HB4.0




93211502 测深器
Depth Gauge



* * KCO 3.2mm 钉棒器械包 93219903

KCO3.2 Screw Rod Instrument Set

A



93211602 模棒
Trial Rod



93211902 抗旋转套筒扳手
Anti-Rotating Sleeve



91150126 直型快装手柄 (三角接口)
Quick Coupling Handle (Straight)



93210803 顶丝预锁紧起子
Curve Probe, T3.4



93211302 万向钉上钉器
Universal Screw Screwdriver



93210906 梅花起子 (快装三角) (短)
Universal Screw Screwdriver, T3.4



93210905 梅花起子 (快装三角) (长)
Universal Screw Screwdriver, T3.4



91132911 梅花起子 (快装三角)
Universal Screw Screwdriver, T2.9



93211702 持棒钳
Rod Holding Clamp




93211802 持钩钳
Hook Holding Clamp



93211405 KCO-DT 压棒器
Depth Gauge



93211403 压棒器
Trial Rod



93212102 剪棒钳
Rod Cutting Clamp



93212202 加压钳
Rod Bender (Right)



93212302 撑开钳
Rod Cutting Clamp



* * KCO 3.2mm 钉棒器械包 93219903

KCO 3.2 Screw Rod Instrument Set



93212401 枕骨导向器
Nut Holder, 6/8



93212402 枕骨导向器
Nut Holder, 10/12



93212702 弯棒 (板) 钳
Quick Coupling Handle (Straight)



93213002 钩开路器
Exchanging Joint



93213301 弯板钳
Tap Sleeve



93210505 直型开路器
Universal Screw Screwdriver



93210503 弯型开路器
Adjuster For Screw/Rod Pusher



93211701 持棒钳
Rod Holding Clamp



93211801 持钩钳
Hook Holding Clamp



93213501 3.5钉盒
Rod Bender , 14-16



93213601 3.5钉盒
Rod Bender , 18-34



93213701 4.0钉盒
Rod Cutting Clamp






93213801 横连钉盒
Nut Holder



* * KCO 3.2mm 钉棒器械包 93219901

KCO 3.2 Screw Rod Instrument Set

A



93210101 直探针
Straight Probe



93210102 弯探针
Curve Probe



93210201 开路锥
Awl



93210301 显影针 (柱形)
Developing Pin (Cylinder Shape)



93210302 显影针 (球形)
Developing Pin (Spherical Shape)



93210401 显影针植入器
Insert Device For Developing Pin



93210501 椎弓根探锥
Pedicle Prober



93210601 钻头 (三角快装式)
Drill Bit (Quick Coupling), $\Phi 2.4 \times 200$



93210701 丝攻 (三角快装式)
Tap (Quick Coupling), HB3.5



* * KCO 3.2mm 钉棒器械包 93219901

KCO3.2 Screw Rod Instrument Set



93210802 顶丝持杆
Nut Holder,T3.4



93210902 顶丝扳手 (三角快装式)
Nut Wrench (Quick Coupling),T3.4



91150105 快装手柄 (三角快装式)
Quick Coupling Handle (Straight)



91150201 转换接头 (三角快装式)
Exchanging Joint



93211201 丝攻套管
Tap Sleeve



93211301 万向钉上钉器
Universal Screw Screwdriver



93211401 钉头调整/推棒器
Adjuster For Screw/Rod Pusher



93211501 测深器
Depth Gauge



93211601 模棒
Trial Rod



93211701 持棒钳
Rod Holding Clamp




93211801 持钩钳
Hook Holding Clamp




93211901 抗旋转套筒
Anti-Rotataing Sleeve



93212001 弯棒器 (左)
Rod Bender (Left)



93212002 弯棒器 (右)
Rod Bender (Right)



93212101 剪棒钳
Rod Cutting Clamp



* * KCO 3.2mm 钉棒器械包 93219901

KCO 3.2 Screw Rod Instrument Set



93212201 加压钳
Compression Clamp



93212301 撑开钳
Distraction Clamp



93212401 枕骨导向器6/8
Occipital Guider




93212402 枕骨导向器10/12
Occipital Guider



93212501 枕骨钻头 (快装式)
Occipital Drill (Quick Coupling), Φ 2.4



93212601 枕骨丝锥 (快装式)
Occipital Tap (Quick Coupling),HB4.0



93212701 弯棒 (板) 钳
Rod (Plate) Bending Clamp



93213401 持钉套
Screw Holding Sleeve



91132902 梅花起子
Torx Screwdriver, T2.9



93213001 钩开路器
Hook Awl



93213101 带刻度导向器
Guider With Scale



93213201 枪式复位钳
Gun-Type Reduction Clamp




93213301 弯板钳
Plate Bending Clamp



91150115 锁定式直型手柄小号 (三角)
Locking type direct handle






91133405 三角形快装起子
Quick Coupling Screwdriver,T3.4 $\times$ 120mm



*** * 颈椎后路椎板成形系统器械包 93229901**

Posterior Cervical Laminar Shaping System Instrument Set

A



93220101 钉板盒
Screw&Plate Box

93220201 骨试模8/10
Bone Trial Implant

93220301 骨试模12/14
Bone Trial Implant



93220501 钻头
Drill Bit, $\Phi 1.7 \times 6\text{mm}$

93220503 钻头
Drill Bit, $\Phi 1.7 \times 10\text{mm}$

93220401 刮匙
Curette



93220502 钻头
Drill Bit, $\Phi 1.7 \times 8\text{mm}$

93220504 钻头
Drill Bit, $\Phi 1.7 \times 12\text{mm}$

91130201 十字起子
Crosshead Screwdriver



93220701 持板钳(镊)
Plate Holding Clamp (Forceps)


93220801 弯板剪板钳
Curved Plate Cutting Clamp

93220901 椎板撬钳
Lamina Pry Clamp



* * 颈椎前路钢板器械包 93239901

Anterior Cervical Plate Instrument Set



93230101 定位针
Located Pin



93230201 定位钉植入器
Insert Device For Located Screw



93230301 快装开孔器 (D型)
Awl, $\Phi 2.0$



93230401 快装钻头 (D型)
Drill Bit, $\Phi 2.5 \times 13\text{mm}$



93230402 快装钻头 (D型)
Drill Bit, $\Phi 2.5 \times 14\text{mm}$



93230403 快装钻头 (D型)
Drill Bit, $\Phi 2.5 \times 15\text{mm}$



93230404 快装钻头 (D型)
Drill Bit, $\Phi 2.5 \times 16\text{mm}$



93230501 快装丝锥 (D型)
Tap, $\Phi 4.0$



93230601 单孔导向器
Single Hole Guide



*** * 颈椎前路钢板器械包 93239901**

Anterior Cervical Plate Instrument Set



93230701 两孔导向器
Two Holes Guider



93650201 撑开钉
Distraction Screw



93230901 撑开钉起子
Distraction Screw Screwdriver,SW4.0



93231001 撑开器
Distraction Device



93231101 持钉套
Screw Holding Sleeve



93231201 快装内六角起子 (D型)
Inner Hexagonal Screwdriver,SW2.5



93231301 快装一字起子 (D型)
Straight Screwdriver



93231401 持板器
Plate Holder



93231501 弯板钳
Plate Bending Clamp



93231601 快装转换接头 (D型)
Exchanging Joint



93231701 快装手柄 (D型)
Quick Coupling Handle





93231801 钉板盒
Screw&Plate Box



* * KSS-I6.0mm 脊柱固定器安装器械包 93319901

KSS-I 6.0 Spinal Fixation Instrument Set



93310101 开口器
Awl



93310201 限位钻孔器
Limit Awl, $\Phi 3.0 \times 40$



93310301 直型开路器
Straight Awl, $\Phi 3.0$



93310302 直型开路器
Straight Awl, $\Phi 3.5$




93310401 显影针 (球形)
Developing Pin (Spherical Shape)




93310402 显影针 (柱形)
Developing Pin (Cylinder Shape)



93310501 显影针植入器
Insert Device For Developing Pin



93310601 直探针
Straight Probe



93310602 弯探针
Curve Probe



* * KSS-I6.0mm 脊柱固定器安装器械包 93319901

KSS-I 6.0 Spinal Fixation Instrument Set


A




93310701 丝锥
Tap,Φ5.0




93310801 固定钉上钉器 (短)
Screw Screwdrive (Short)




93310802 固定钉上钉器 (长)
Screw Screwdrive (Long)



93310901 多轴钉上钉器 (短)
Screwdriver For Polyaxial Pedicle Screw(Short),T3.4



93310902 多轴钉上钉器 (长)
Screwdriver For Polyaxial Pedicle Screw(Long),T3.4



93310903 多轴钉上钉器-梅花
Screwdriver For Polyaxial Pedicle Screw-torx,T3.4-M11



93311001 T型扳手 (短)
T Wrench(Short),T5.3



93311002 T型扳手 (长)
T Wrench (Long) ,T5.3



93311201 持锁紧螺塞器 (长)
Holding Set For Nut(Long),T5.3



93311202 持锁紧螺塞器 (短)
Holding Set For Nut(Short),T5.3



93311301 压棒器
Rod Pusher,Φ6.0



93311401 压棒钳
Rod Pushing Clamp,Φ2.5



93311501 模棒
Trial Rod,Φ6.0



93311601 持棒钳
Rod Holding Clamp



93311701 转棒钳
Rod Rotating Clamp



* * KSS-I6.0mm 脊柱固定器安装器械包 93319901

KSS-I 6.0 Spinal Fixation Instrument Set



93311801 转棒器
Rod Rotating Device, SW5.0



93311901 螺钉折断钳
Screw Cutter For Long Pedicle Screw



91133402 梅花起子
Torx Screwdriver, T3.4



93312001 抗旋转扳手 (直角)
Anti-Rotation Wrench (Right Angle)



93312101 抗旋转扳手
Anti-Rotation Wrench, M11



93312201 弯棒钳
Rod Bending Clamp



93312301 弯棒器
Rod Bender, Φ6.0



93312401 压缩钳
Compression Clamp





93312501 撑开钳
Distraction Clamp





* * KSS-II5.5mm 脊柱固定器安装器械包 93339901

KSS-II5.5 Spinal Fixation Instrument Set



93330101 开口器
Awl, $\Phi 4.5$



93330201 直型开路器
Straight Awl, $\Phi 3.0$



93330202 直型开路器
Straight Awl, $\Phi 3.5$



93310601 直探针
Straight Probe



93310602 弯探针
Curve Probe



93310401 显影针 (球形)
Developing Pin (Spherical Shape)

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object. 20220617



* * KSS-II5.5mm 脊柱固定器安装器械包 93339901

KSS-II5.5 Spinal Fixation Instrument Set



93310402 显影针 (柱形)
Developing Pin (Cylinder Shape)



93310502 显影针植入器
Insert Device For Developing Pin



93330501 丝锥
Tap,Φ3.5



93330502 丝锥
Tap,Φ4.5



93330503 丝锥
Tap,Φ5.5



93330504 丝锥
Tap,Φ6.5



93330601 固定钉上钉器 (短)
Screw Screwdriver (Short)



93330602 固定钉上钉器 (长)
Screw Screwdriver (Long)



93330701 多轴钉上钉器 (短)
Screwdriver For Polyaxial Pedicle Screw(Short),T3.4



93330702 多轴钉上钉器 (长)
Screwdriver For Polyaxial Pedicle Screw(Long),T3.4



93330801 模棒
Trial Rod,Φ5.5



93330901 持棒钳
Rod Holding Clamp



93331001 大力钳
Lock Wrench



93331101 转棒器
Rod Rotating Device,SW4.0



93331201 压棒器
Rod Pusher



* * KSS-II5.5mm 脊柱固定器安装器械包 93339901

KSS-II5.5 Spinal Fixation Instrument Set



93331301 体内弯棒器
Internal Rod Bender



93331401 可调节弯棒钳
Adjustable Rod Bending Clamp



93331501 持锁紧螺塞器 (短)
Holding Set For Nut(Short),T5.3/T3.4



93331502 持锁紧螺塞器 (长)
Holding Set For Nut(Long),T5.3/T3.4



93331601 直型快装手柄
Straight Handle With Quick Coupling



93331701 T型快装手柄
T Quick Coupling Handle



93331801 T型扳手 (短)
T Wrench (Short) ,T5.3



93331802 T型扳手 (长)
T Wrench (Long) ,T5.3



93331901 折断防松器
Cutter Locking Device



93332001 钉尾折断器
Screw Cutter For Long Pedicle Screw



93332101 蛙式复位钳
Frog-Type Reduction Clamp



93332201 撑开钳
Distraction Clamp



93332301 压缩钳
Compression Clamp



93332401 平行压缩钳
Parallel Compression Clamp



93332501 双棒垫片开路导向器
Double Rods Shim Open Guider



* * KSS-II5.5mm 脊柱固定器安装器械包 93339901

KSS-II5.5 Spinal Fixation Instrument Set



93332601 双棒垫片打入器
Double Rods Shim Impactor



93332701 导向器开口器 (快装式)
Guider Awl (Quick Coupling),Φ3.0



93332702 导向器开口器 (快装式)
Guider Awl (Quick Coupling),Φ4.5



93332801 单杆垫片打入器
Single Rod Shim Impactor



93332901 扩孔锥 (快装式)
Expanding Awl(Quick Coupling)



93333001 测深器
Depth Gauge




93333101 横连器起子
Crosslink Screwdriver,T3.4



93333201 横连规量13、14
Crosslink Gauge



93333202 横连规量15、16
Crosslink Gauge



93333203 横连规量17、18
Crosslink Gauge



93333301 持横连器钳
Clamp For Holding Crosslink



93333401 量规
Gauge



93333501 椎体撑开器
Centrum Distraction Device,Φ5.5

A








KSS-IV 脊柱内固定系统安装器械包 93369901

KSS-IV Spinal Fixation System Instrument Set

A



93320101 开路锥
Awl



93320201 弯型开路器
Curve Awl,3.5



93320202 直型开路器
Straight Awl,3.5



93320203 直型开路器
Straight Awl,3.0



93320204 直型开路器
Straight Awl,2.5



93320301 双头探针
Double-End Probe



93320401 直型探针(硬)
Straight Probe (Hard)



93320501 显影针植入器
Insert Device For Device For Developing Pin



93323101 显影针 (球形)
Developing Pin (Spherical Shape)

Our company have final right to modify all the specification in this catalogue, If there is any inconformitybetween contant, model, color, please refer to the object. 20211130



KSS-IV 脊柱内固定系统安装器械包 93369901

KSS-IV Spinal Fixation System Instrument Set

A



93323102 显影针 (柱形)
Developing Pin (Cylinder Shape)



93321601 模棒
Trial Rod, Φ5.5×400mm



93320601 双线丝锥
Tap, Φ4.5
曾用编码 93100205



93320602 双线丝锥
Tap, Φ5.0
曾用编码 93100201



93320603 双线丝锥
Tap, Φ5.5
曾用编码 93100202



93320604 双线丝锥
Tap, Φ6.0
曾用编码 93100203



93320605 双线丝锥
Tap, Φ6.5
曾用编码 93100204



93327101 万向钉调钉器
Universal Screw Adjuster



93320701 弯型持棒钳
Curve Rod Holding Clamp



93320801 对抗套筒
Counter Sleeve



93320901 长压棒器
Long Rod Pusher



93322001 锁紧起子
Nut Preloaded Wrench, T5.3



93331301 体内弯棒器(左右)
Internal Rod Bender (Right)



93321201 摇摆钳
Swing Clamp



93321301 弯棒钳
Rod Bending Clamp



KSS-IV 脊柱内固定系统安装器械包 93369901

KSS-IV Spinal Fixation System Instrument Set



93321401 长尾钉折断器
Screw Cutter For Long Pedicle Screw, 10mm



93321501 折断扶持器
Screw Cutter Support



93350101 模棒
Trial Rod, $\Phi 6.0 \times 400\text{mm}$



93321701 转棒钳
Rod Rotating Clamp




93321801 旋棒扳手 ($\Phi 5.5$ 棒用)
Rod Rotating Wrench, SW4.0




93311801 旋棒扳手 ($\Phi 6.0$ 棒用)
Rod Rotating Wrench, SW5.0



93321901 大力钳
Lock Wrench



93311201 顶丝预紧起子 (长)
Locking Screwdriver



93311202 顶丝预紧起子 (短)
Locking Screwdriver



93322101 万向钉退钉扳手
Universal Screw Removal Wrench



93322201 T型手柄(小)
T Handle (Small)



93322401 撑开钳
Distraction Clamp



93322501 加压钳
Compression Clamp



93322801 棘轮快换手柄
Quick Coupling Ratchet Handle



91133503 快装起子
Get the screwdriver, SW3.5



KSS-IV 脊柱内固定系统安装器械包 93369901

KSS-IV Spinal Fixation System Instrument Set



91150114 锁定式直型手柄大号 (三角)
Lock type straight handle large size (Triangular)



93323201 长尾万向钉上钉器(六角)
Long tail multiaxial pedicle screw screwdriver, SW3.5

曾用编码 93322601



93323301 短尾万向钉上钉器 (六角)
Short tail multiaxial pedicle screw screwdriver, SW3.5

曾用编码 93322602



93323401 短尾固定钉上钉器
Screw Screwdriver

曾用编码 93350201

A





* * JZPA 5.5/6.0 脊柱内固定系统安装器械包 93359901

JZPA 5.5/6.0 Spinal Fixation System Instrument Set



基础器械



93320101 开路锥
Awl



93320201 弯型开路器
Curve Awl,3.5



93320202 直型开路器
Straight Awl,3.5



93320203 直型开路器
Straight Awl,3.0



93320204 直型开路器
Straight Awl,2.5



93310601 直探针
Straight Probe



93310602 弯探针
Curve Probe



93320401 直型探针(硬)
Straight Probe (Hard)



93320501 显影针植入器
Insert Device For Device For Developing Pin

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object.20220617



* * JZPA 5.5/6.0 脊柱内固定系统安装器械包 93359901

JZPA 5.5/6.0 Spinal Fixation System Instrument Set

A



93323101 显影针 (球形)
Developing Pin (Spherical Shape)



93323102 显影针 (柱形)
Developing Pin (Cylinder Shape)



93320701 弯型持棒钳
Curve Rod Holding Clamp



93320801 对抗套筒
Counter Sleeve



93320901 长压棒器
Long Rod Pusher



93322001 锁紧起子
Nut Preloaded Wrench, T5.3



93321101 体内弯棒器(右)
Internal Rod Bender (Right)



93321102 体内弯棒器(左)
Internal Rod Bender (Left)



93321201 摇摆钳
Swing Clamp



93321301 弯棒钳
Rod Bending Clamp



93321401 长尾钉折断器
Screw Cutter For Long Pedicle Screw



93321501 折断扶持器
Screw Cutter Support



93321601 模棒
Trial Rod, Φ5.5×400mm



93350101 模棒
Trial Rod, Φ6.0×400mm



93321701 转棒钳
Rod Rotating Clamp



* * JZPA 5.5/6.0 脊柱内固定系统安装器械包 93359901

JZPA 5.5/6.0 Spinal Fixation System Instrument Set



93321801 旋棒扳手 (φ 5.5 棒用)
Rod Rotating Wrench, SW4.0




93311801 旋棒扳手 (φ 6.0 棒用)
Rod Rotating Wrench, SW5.0



93321901 大力钳
Lock Wrench



93311201 顶丝预紧起子 (长)
Locking Screwdriver (Long)



93311202 顶丝预紧起子 (短)
Locking Screwdriver (Short)



93322201 T型手柄(小)
T Handle (Small)



93322401 撑开钳
Distraction Clamp



93322501 加压钳
Compression Clamp



93322801 棘轮快换手柄 (选配)
Quick Coupling Ratchet Handle



91133503 快装起子
Get the screwdriver, SW3.5



93320902 DT 压棒器
DT Rod Pusher



91150114 锁定式直型手柄大号(三角)
Lock type straight handle large size (Triangular)
曾用编码 91150123

选配方案一 (适用于JZPA5.5/6.0系统) (双线)



93320601 双线丝锥
Tap, φ4.5
曾用编码 93100205



93320602 双线丝锥
Tap, φ5.0
曾用编码 93100201



93320603 双线丝锥
Tap, φ5.5
曾用编码 93100202



* * JZPA 5.5/6.0 脊柱内固定系统安装器械包 93359901

JZPA 5.5/6.0 Spinal Fixation System Instrument Set

A



93320604 双线丝锥
Tap,Φ6.0

曾用编码 93100203



93320605 双线丝锥
Tap,Φ6.5

曾用编码 93100204



93323201 长尾万向钉上钉器(六角)

Long tail multiaxial pedicle screw screwdriver, SW3.5

曾用编码 93322601



93323301 短尾万向钉上钉器(六角)

Short tail multiaxial pedicle screw screwdriver, SW3.5

曾用编码 93322602

选配方案二 (适用于KSS-I 6.0mm系统)



93100301 单线丝锥
Tap,Φ4.5



93100302 单线丝锥
Tap,Φ5.0



93100303 单线丝锥
Tap,Φ5.5



93100304 单线丝锥
Tap,Φ6.0



93100305 单线丝锥
Tap,Φ6.5



93323201 长尾万向钉上钉器(六角)

Long tail multiaxial pedicle screw screwdriver, SW3.5



93323301 短尾万向钉上钉器(六角)


Short tail multiaxial pedicle screw screwdriver, SW3.5



93323501 短尾固定钉上钉器

Screw Screwdriver

曾用编码 93350201



93323401 长尾固定钉上钉器

Screw Screwdriver



* * JZPA 5.5/6.0 脊柱内固定系统安装器械包 93359901

JZPA 5.5/6.0 Spinal Fixation System Instrument Set

T581



T582



T587



选配方案三 (适用于KSS-I 6.0mm系统)



93310301 单线丝锥
Tap,Φ4.5



93310302 单线丝锥
Tap,Φ5.0



93310303 单线丝锥
Tap,Φ5.5




93310304 单线丝锥
Tap,Φ6.0




93310305 单线丝锥
Tap,Φ6.5



93310901 短尾万向钉上钉器(梅花)
Short tail multiaxial pedicle screw screwdriver,T3.4



93310902 长尾万向钉上钉器(梅花)
Long tail multiaxial pedicle screw screwdriver,T3.4



93310903 多轴钉上钉器-梅花
Screwdriver For Polyaxial Pedicle Screw-torx,T3.4-M11



93353501 短尾固定钉上钉器
Screw Screwdrive



93353401 长尾固定钉上钉器
Screw Screwdrive

T581



T582



T583





* * JZPA 5.5/6.0 脊柱内固定系统安装器械包 93359901

JZPA 5.5/6.0 Spinal Fixation System Instrument Set

DT 专用器械



93320902 DT 压棒器 (选配)
DT Rod Pusher



93321502 DT 折断扶持器 (选配)
Screw Cutter Support




93320903 DT长尾钉折断器 (DT 11mm) (选配)
DT Screw Cutter For Long Pedicle Screw

A



* JZPA 5.5 脊柱内固定系统安装器械包 93329901

JZPA 5.5 Spinal Fixation System Instrument Set



93320101 开路锥
Awl



93320201 弯型开路器
Curve Awl,3.5



93320202 直型开路器
Straight Awl,3.5




93320203 直型开路器
Straight Awl,3.0



93320204 直型开路器
Straight Awl,2.5



93310601 直探针
Straight Probe



93310602 弯探针
Curve Probe



93320401 直型探针(硬)
Straight Probe (Hard)



93320501 显影针植入器
Insert Device For Device For Developing Pin

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object.20220617



* JZPA 5.5 脊柱内固定系统安装器械包 93329901


JZPA 5.5 Spinal Fixation System Instrument Set



93323101 显影针 (球形)
Straight




93323102 显影针 (柱形)
Straight



93320601 双线丝锥
Tap, $\Phi 4.5$



93320602 双线丝锥
Tap, $\Phi 5.0$



93320603 双线丝锥
Tap, $\Phi 5.5$



93320604 双线丝锥
Tap, $\Phi 6.0$



93320605 双线丝锥
Tap, $\Phi 6.5$



93320701 弯型持棒钳
Curve Rod Holding Clamp



93320801 对抗套筒
Counter Sleeve



93320901 长压棒器
Long Rod Pusher



93322001 锁紧起子
Nut Preloaded Wrench, T5.3



93321101 体内弯棒器(右)
Internal Rod Bender (Right)



93321102 体内弯棒器(左)
Internal Rod Bender (Left)



93321201 摇摆钳
Swing Clamp



93321301 弯棒钳
Rod Bending Clamp



* JZPA 5.5 脊柱内固定系统安装器械包 93329901

JZPA 5.5 Spinal Fixation System Instrument Set

A



93321401 长尾钉折断器
Screw Cutter For Long Pedicle Screw, 10mm



93321501 折断扶持器
Screw Cutter Support



93321601 模棒
Trial Rod, Φ5.5×400mm



93321602 模棒
Trial Rod, Φ5.5×250mm




93321701 转棒钳
Rod Rotating Clamp



93321801 旋棒扳手
Rod Rotating Wrench, SW4.0



93321901 大力钳
Lock Wrench



93311201 顶丝预紧起子 (长)
Locking Screwdriver



93311202 顶丝预紧起子 (短)
Locking Screwdriver



93322201 T型手柄(小)
T Handle (Small)



93322301 扭力手柄 (选配)
Torque Handle, 12NM



93322401 撑开钳
Distraction Clamp



93322501 加压钳
Compression Clamp



93323201 长尾万向钉上钉器(六角)
Long tail multi-axial pedicle screw screwdriver, SW3.5



93323301 短尾万向钉上钉器(六角)
Short tail multi-axial pedicle screw screwdriver, SW3.5



* JZPA 5.5 脊柱内固定系统安装器械包 93329901

JZPA 5.5 Spinal Fixation System Instrument Set



93322801 棘轮快换手柄(选配)
Quick Coupling Ratchet Handle



93322901 枪式复位钳(选配)
Gun-Type Reduction Clamp



93323001 延伸器(选配)
Extender




93320301 双头探针(选配)
Double-End Probe



93322701 万向钉调钉器(选配)
Universal Screw Adjuster




93322101 万向钉退钉扳手(选配)
Universal Screw Removal Wrench



* * JZWA脊柱微创固定系统安装器械包 93519901

JZWA Minimally Invasive Spinal Fixation System Instrument Set



93510101 螺纹导针
Threaded Guide Pin, $\Phi 1.5 \times 450\text{mm}$



93510201 预紧起子
Preloaded Screwdriver



93510301 六角起子
Hexagonal Screwdriver, SW2.5



93510401 顶丝扳手
Nut Wrench, T5.3



93510501 开口器
Awl



93510601 空心丝锥
Hollow Tap, $\Phi 5.0\text{mm}$




93510701 空心丝锥
Hollow Tap, $\Phi 6.0\text{mm}$



93510801 空心丝锥
Hollow Tap, $\Phi 7.0\text{mm}$



93510901 钉尾折断器
Screw Cutter



* * JZWA 脊柱微创内固定系统安装器械包 93519901

JZWA Minimally Invasive Spinal Fixation System Instrument Set

A



93511001 扩张套筒(大)
Expanding Sleeve(Large)



93511101 扩张套筒(中)
Expanding Sleeve(Medium)



93511201 扩张套筒(小)
Expanding Sleeve(Small)



93511301 对抗扳手
Counter Wrench



93511401 折断对抗器
Break-off Counter Device



93511501 置棒器
Rod Placing Device



93511601 测量器
Measurement Device



93511701 撑开钳
Distraction Clamp



93511801 加压钳
Compression Clamp



93511901 弯棒钳
Rod Bending Clamp



93512001 棘轮手柄
Ratchet Handle




93512101 扭力手柄
Torque Handle,12.5Nm




93512201 稳定套
Stable Sleeve



93512302 U型钉上钉器
Screwdriver For U Pedicle Screw,sw4.0



93512401 穿刺针 (选配)
Screwdriver For U Pedicle Screw



* * JZWA脊柱微创固定系统安装器械包 93519901

JZWA Minimally Invasive Spinal Fixation System Instrument Set



93512501 空心丝锥定制款 (选配)
Hollow Tap, $\Phi 4.5\text{mm}$



93512501 空心丝锥定制款 (选配)
Hollow Tap, $\Phi 5.0\text{mm}$



93512501 空心丝锥定制款 (选配)
Hollow Tap, $\Phi 6.0\text{mm}$





93512501 空心丝锥定制款 (选配)
Hollow Tap, $\Phi 7.0\text{mm}$





*** * 零切迹融合器器械包 93249901**

Zero - profile Fusion Cage Instrument Set



93240101 撑开钉
Distraction Screw



93240201 撑开钉起子
Distraction Screw Screwdriver



93240301 撑开器
Distraction Device



93240401 试模 (楔形)
Trial Template (Wedged), 5mm



93240402 试模 (楔形)
Trial Template (Wedged), 6mm



93240403 试模 (楔形)
Trial Template (Wedged), 7mm



93240404 试模 (楔形)
Trial Template (Wedged), 8mm



93240405 试模 (楔形)
Trial Template (Wedged), 9mm



93240406 试模 (楔形)
Trial Template (Wedged), 10mm

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object.20230928



* * 零切迹融合器器械包 93249901

Zero - profile Fusion Cage Instrument Set



93240501 试模 (弧形)
Trial Template (Curved), 5mm



93240502 试模 (弧形)
Trial Template (Curved), 6mm



93240503 试模 (弧形)
Trial Template (Curved), 7mm



93240504 试模 (弧形)
Trial Template (Curved), 8mm



93240505 试模 (弧形)
Trial Template (Curved), 9mm



93240506 试模 (弧形)
Trial Template (Curved), 10mm



93240601 植入物夹持器
Implant Holding Device



93240701 植骨仓
Bone Graft Case



93240801 松质骨压实器
Cancellous Bone Compactor



93240901 敲击锤
Knocking Hammer



93241001 持取打入钳
Screw Holding Blowing Plier



93241101 球形推进器
Spherical Pusher, Φ4.5mm



93241201 弧形推进器
Curved Pusher



93241301 限位开路椎
Limit Awl, Φ2.5mm



93241401 成角开角椎
Open Angle Awl, Φ2.5mm



* * 零切迹融合器器械包 93249901

Zero - profile Fusion Cage Instrument Set



93241501 导向器
Guider, $\Phi 2.8\text{mm}$



93241601 限位钻
Limit Drill, $\Phi 2.0 \times 12\text{mm}$



93241602 限位钻
Limit Drill, $\Phi 2.0 \times 14\text{mm}$



93241603 限位钻
Limit Drill, $\Phi 2.0 \times 16\text{mm}$



93241604 限位钻
Limit Drill, $\Phi 2.0 \times 18\text{mm}$



93241701 成角上钉器
Angled Screwdriver, SW2.5



93241801 上钉器把持柄
Screwdriver Holding Handle



93241901 六角起子
Hexagonal Screwdriver, SW2.5



93242001 持钉套筒
Screw Holding Sleeve



93242101 扭力快换
Torque Quick Coupling Joint, 1.2N.m



93242201 快换手柄
Quick Coupling Handle



93242301 PPSU 盒子
PPSU







* * PEEK 胸腰椎融合器安装器械包 93649901

PEEK Thoracolumbar fusion cage instrument set

A



93640101 试模
Trial Template, 7mm



93640103 试模
Trial Template, 9mm



93640105 试模
Trial Template, 11mm



93640102 试模
Trial Template, 8mm




93640104 试模
Trial Template, 10mm




93640106 试模
Trial Template, 12mm




93640107 试模
Trial Template, 13mm




93640108 试模
Trial Template, 14mm




93640201 铰刀
Reamer, 7mm



93640202 铰刀
Reamer, 8mm



93640203 铰刀
Reamer, 9mm




93640204 铰刀
Reamer, 10mm




*** * PEEK 胸腰椎融合器安装器械包 93649901**

PEEK Thoracolumbar fusion cage instrument set



93640205 铰刀
Reamer, 11mm



93640206 铰刀
Reamer, 12mm



93640301 脊柱植入物推送器
Spinal implant pusher



93640401 神经拉钩 (小)
Nerve hook(Small)




93640402 神经拉钩 (大)
Nerve hook(Large)



93640501 剥离器
Detacher



93640601 骨填充器
Bone Filler



93651601 植骨台
Bone graft device
曾用编码 93611001



93640701 植骨器
Bone graft device



93640801 植骨漏斗
Double chambers funnel



93640901 融合器把持器
Fusin cage holder



93641001 滑动锤
Slide hammer



93641101 方形刮刀 (直)
Square scraper(Curved)



93641201 水滴型刮刀弯 (右)
Drop-shaped scraper(Right)



93641202 水滴型刮刀弯 (左)
Drop-shaped scraper(Left)



93641301 带齿刮匙 (右)
Toothed curet(Right), 8mm





93641302 带齿刮匙 (左)
Toothed curet(Left), 8mm



* * PEEK 颈椎融合器安装器械包 93659901

PEEK Thoracolumbar fusion cage instrument set

A



93650101 上钉器
On the nail



93650201 撑开螺钉
Open screw



93650301 推进器
Thruster



93650401 持取器
A take device



93650501 填充器
Applier



93651601 植骨台
Bone graft taiwan



93650601 滑锤
Slide hammer



93650701 开孔器 (弯)
Hole saw (bend)



93650801 开孔器 (直)
Hole saw (straight)



* * PEEK 颈椎融合器安装器械包 93659901

PEEK Thoracolumbar fusion cage instrument set

A



93650901 拉钩 (大)
Drag hook (big)



93650902 拉钩 (小)
Drag hook (little)



93651001 刮匙
Curet, 3mm




93651101 终板刮刀
End plate scraper, 4x7mm



93640801 植骨漏斗
Bone graft funnel
曾用编码 93651201



93640701 植骨器
Implant bone aretifacts
曾用编码 93651301



93651401 颈椎牵开器
Cervical retractor



93651506 试模
Test mode, 4mm



93651501 试模
Test mode, 5mm



93651502 试模
Test mode, 6mm



93651503 试模
Test mode, 7mm



93651504 试模
Test mode, 8mm



93651505 试模
Test mode, 9mm



93651402 撑开器 (选配)
Cervical retractor








*** * 颈椎椎间融合器安装器械包 93619902**

Cervical Fusion Device Instrument Set

A



93650201 撑开钉
Distraction Screw



93610202 撑钉起子
Distraction Screw Screwdriver




93610302 开孔器
Awl



93610402 刮匙
Curette



93610502 定位器开孔锥
Awl For Located Device



93610602 刮刀
Scraper



93610701 平行定位器
Parallel Located Device



93610801 撑开器
Distraction Device



91150402 限位器扳手
Limiter Wrench,SW2.5

Our company have final right to modify all the specification in this catalogue, If there is any inconformitybetween contant, model, color, please refer to the object.20220617



*** * 颈椎椎间融合器安装器械包 93619902**

Cervical Fusion Device Instrument Set



93610901 试模
Trial Template, 12×15×4.5



93610902 试模
Trial Template, 12×15×5.5



93610903 试模
Trial Template, 12×15×6.5



93610904 试模
Trial Template, 12×15×7.5



93610905 试模
Trial Template, 12×15×8.5



93611001 融合器持取器
Fusion Holder



93611102 顶棒
Bucking Bar



93611201 漏斗
Funnel




93651601 植骨台
Bone Grafting Platform
曾用编码 93611302



93611401 植骨器
Bone Grafting Device



93611502 骨锤
Bone Hammer





93611601 颈椎拉钩
Cervical Pulling Hook



* * 颈椎椎间融合器安装器械包 93619901

Cervical Fusion Device Instrument Set

A



93650201 撑开钉
Distraction Screw



93610201 撑开钉起子
Distraction Screw Screwdriver



93610301 开孔器
Awl



93610401 刮匙
Curette



93610501 定位器开孔锥
Awl For Located Device



93610601 刮刀
Scraper



93610701 平行定位器
Parallel Located Device



93610801 撑开器
Distraction Device



91150402 限位器扳手
Limiter Wrench,SW2.5

Our company have final right to modify all the specification in this catalogue, If there is any inconformitybetween contant, model, color, please refer to the object.20220617



* * 颈椎椎间融合器安装器械包 93619901

Cervical Fusion Device Instrument Set



93610901 试模
Trial Template, 12×15×4.5



93610902 试模
Trial Template, 12×15×5.5



93610903 试模
Trial Template, 12×15×6.5



93610904 试模
Trial Template, 12×15×7.5



93610905 试模
Trial Template, 12×15×8.5



93611001 融合器持取器
Fusion Holder



93611101 顶棒
Bucking Bar



93611201 漏斗
Funnel




93651601 植骨台
Bone Grafting Platform



93611401 植骨器
Bone Grafting Device



93611501 骨锤
Bone Hammer






93611601 颈椎拉钩
Cervical Pulling Hook



* * 腰椎椎间融合器安装器械包 93629902

Lumbar Fusion Device Instrument Set

A



93620103 刮匙(小)
Curette (Small)



93620104 刮匙(大)
Curette (Large)



93620204 刮刀 (小)
Scraper (Small)



93620205 刮刀 (中)
Scraper (Mid)



93620206 刮刀 (大)
Scraper (Large)



93620302 拉式刮刀
Pull Type Scrapers



93620402 推式刮刀
Push Type Scraper



93620504 铰刀 (小)
Reamer (Small)

93620505 铰刀 (中)
Reamer (Mid)



93620610 试模
Trial Template, 6×24



93620611 试模
Trial Template, 8×22



93620612 试模
Trial Template, 8×26




93620613 试模
Trial Template, 10×22

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object. 20220617



*** * 腰椎椎间融合器安装器械包 93629902**


Lumbar Fusion Device Instrument Set




93620614 试模
Trial Template,10×26



93620615 试模
Trial Template,12×22



93620616 试模
Trial Template,12×26



93620617 试模
Trial Template,14×22




93620618 试模
Trial Template,14×26




93620702 融合器持取器
Fusion Holder



93640801 植骨漏斗
Bone Graft Funnel



93640701 植骨器
Bone Grafting Device



93620801 冲锤
Punch Hammer



93620901 剥离子
Detacher



93621001 神经根拉钩 (小)
Nerve Root Retractor(Small)



93621002 神经根拉钩 (大)
Nerve Root Retractor(Large)



93651601 植骨台
Bone Grafting Plate



93621201 冲板
Punch Plate



93621301 骨填充器
Bone Filler




A



* * 腰椎椎间融合器安装器械包 93629901

Lumbar Fusion Device Instrument Set

A



93620101 刮匙(小)
Curette (Small)



93620102 刮匙(大)
Curette (Large)



93620201 刮刀 (小)
Scraper (Small)



93620202 刮刀 (中)
Scraper (Mid)



93620203 刮刀 (大)
Scraper (Large)



93620301 拉式刮刀
Pull Type Scrapers



93620401 推式刮刀
Push Type Scraper



93620501 铰刀 (小)
Reamer (Small)



93620502 铰刀 (中)
Reamer (Mid)



93620601 试模
Trial Template,6×24



93620603 试模
Trial Template,8×26



93620503 铰刀 (大)
Reamer (Large)



93620602 试模
Trial Template,8×22



93620604 试模
Trial Template,10×22

Our company have final right to modify all the specification in this catalogue, If there is any inconformitybetween contant, model, color, please refer to the object.20220617



* * 腰椎椎间融合器安装器械包 93629901

Lumbar Fusion Device Instrument Set



93620605 试模
Trial Template,10×26



93620606 试模
Trial Template,12×22



93620607 试模
Trial Template,12×26



93620608 试模
Trial Template,14×22



93620609 试模
Trial Template,14×26




93620701 融合器持取器
Fusion Holder



93621401 植骨漏斗
Bone Graft Funnels



93621501 植骨器
Bone Grafting Device



93620801 冲锤
Punch Hammer



93620901 剥离子
Detacher



93621001 神经根拉钩 (小)
Nerve Root Retractor(Small)



93621002 神经根拉钩 (大)
Nerve Root Retractor(Large)



93651601 植骨台
Bone Grafting Plate



93621201 冲板
Punch Plate






93621301 骨填充器
Bone Filler



*** * 钛笼系统安装器械包 93639901**

Titaum Mesh Cage Instrument Set

A



93630501 植骨仓
Bone Graft Units



93630101 钛笼测量器
Titanium Mesh Gauge,10×12mm



93630102 钛笼测量器
Titanium Mesh Gauge,16×18mm



93630103 钛笼测量器
Titanium Mesh Gauge,22×25mm



93630601 钛笼修边挫
Titanium Mesh File



93630201 钛笼植骨器
Bone Grafting Device For Titanium Mesh,6mm



93630202 钛笼植骨器
Bone Grafting Device For Titanium Mesh,12mm



93630203 钛笼植骨器
Bone Grafting Device For Titanium Mesh,18mm



93630301 钛笼推入器 (直)
Titanium Mesh Pusher(Straight)



93630302 钛笼推入器 (弯)
Titanium Mesh Pusher(Cutved)



93630401 钛笼夹持钳
Titanium Mesh Holding Forceps



93630801 双关节钛笼剪
Double-Joint Titanium Mesh Scissor






93630701 钛笼剪
Titanium Mesh Scissor



* * 钛笼系统安装器械包 93639902

Titaium Mesh Cage Instrument Set

A



93630501 植骨仓
Bone Graft Units



93630101 钛笼测量器
Titanium Mesh Gauge,10×12mm



93630102 钛笼测量器
Titanium Mesh Gauge,16×18mm



93630103 钛笼测量器
Titanium Mesh Gauge,22×25mm



93630602 钛笼修边挫
Titanium Mesh File



93630204 钛笼植骨器
Bone Grafting Device For Titanium Mesh,6mm



93630205 钛笼植骨器
Bone Grafting Device For Titanium Mesh,12mm



93630206 钛笼植骨器
Bone Grafting Device For Titanium Mesh,18mm



93630301 钛笼推入器 (直)
Titanium Mesh Pusher(Straight)



93630302 钛笼推入器 (弯)
Titanium Mesh Pusher(Cutved)



93630401 钛笼夹持钳
Titanium Mesh Holding Forceps



93630801 双关节钛笼剪
Double-Joint Titanium Mesh Scissor



93630701 钛笼剪
Titanium Mesh Scissor


颌面、微型板系列 安装器械包





* * MCF 颌面钢板安装器械包 91279901

Lower Limb Fracture Instrument Set



91130206 MCF 中芯十字起子杆
MCF Filler Cross Screwdriver Rod,40mm



91130207 MCF 中芯十字起子杆
MCF Filler Cross Screwdriver Rod,65mm



91270101 MCF 开口器
MCF Awl,Φ1.4×5mm



91270102 MCF 开口器
MCF Awl,Φ1.4×8mm



91270103 MCF 开口器
MCF Awl, Φ1.4×10mm



91150124 MCF 快装手柄
MCF Quick Coupling Handle



91270201 MCF 弯板钳 (左)
MCF Plate Bending Plier(Left)



91270202 MCF 弯板钳 (右)
MCF Plate Bending Plier(Right)



91270301 MCF 螺钉盒
MCF Screw Box,Φ1.5螺钉



* * MCF 颌面钢板安装器械包 91279901

Lower Limb Fracture Instrument Set



91270302 MCF 螺钉盒
MCF Screw Box,Φ2.0螺钉



91270601 MCF 螺钉长度柱
MCF Screw Lenght Colum,4



91270602 MCF 螺钉长度柱
MCF Screw Lenght Colum,5



91270603 MCF 螺钉长度柱
MCF Screw Lenght Colum, 6



91270604 MCF 螺钉长度柱
MCF Screw Lenght Colum,7



91270605 MCF 螺钉长度柱
MCF Screw Lenght Colum,8



91270401 MCF 钢板圆盒
MCF Steel Plate Box, 0.6



91270402 MCF 钢板圆盒
MCF Steel Plate Box,0.8



91270702 MCF 螺钉架 1.7、2.0
MCF Screw frame



91270701 MCF 螺钉架 1.5
MCF Screw frame



91270502 MCF 钢板方盒
MCF Profiling Plate Box, PPSU 0.8



91270901 钢板圆盒硅胶垫
Silicone pad of screw box



91270505 MCF 钢板盒盖
MCF Profiling Plate Box, 0.6



91270506 MCF 钢板盒盖
MCF Profiling Plate Box, 0.6




91270507 MCF 钢板盒盖
MCF Profiling Plate Box, 0.8



91270508 MCF 钢板盒盖
MCF Profiling Plate Box, 0.8



B



* * MCF 正颌颌面钢板安装器械包 91289901

MCF Orthognathic Maxillofacial Plate Instrument Set

B



91130206 MCF 中芯十字起子杆
MCF Filler Cross Screwdriver Rod,40mm



91130207 MCF 中芯十字起子杆
MCF Filler Cross Screwdriver Rod,65mm



91270101 MCF 开口器
MCF Awl,Φ1.4×5mm



91270102 MCF 开口器
MCF Awl,Φ1.4×8mm



91270103 MCF 开口器
MCF Awl, Φ1.4×10mm



91150124 MCF 快装手柄
MCF Quick Coupling Handle



91270201 MCF 弯板钳 (左)
MCF Plate Bending Plier(Left)



91270202 MCF 弯板钳 (右)
MCF Plate Bending Plier(Right)



91270302 MCF 螺钉盒
MCF Screw Box,Φ2.0螺钉



* * MCF 正颌颌面钢板安装器械包 91289901

MCF Orthognathic Maxillofacial Plate Instrument Set



91280101 钢板塑形钳 (左)
Plate Shaping Plier(Left)



91280102 钢板塑形钳 (右)
Plate Shaping Plier(Right)




91280201 钢板折弯钳
Plate Bending Plier



91280301 钢板剪板钳
Plate Cutting Plier



91280401 钢板夹持镊
Plate Clamping Forceps




91270403 MCF 钢板圆盒
MCF Plate Round Box, 1.0



91270402 MCF 钢板圆盒
MCF Plate Round Box,0.8



91270503 MCF 钢板方盒
MCF Plate Square Box,PPSU 1.0



91270502 MCF 钢板方盒
MCF Plate Square Box, PPSU 0.8



91130221 MCF 一体式起子 十字 138mm
Mcf All In One Screwdriver Cross 138mm

微型板系列安装器械包



* * TML 1.5 微型锁定指掌骨安装器械包 91859901

TML 2.0 Micro phalangeal metacarpal locking series instrument Set



91820101 TML钢板折弯器
TML plate bender, 2×1mm



×2

91850101 TML 梅花起子
TML quick coupling handle, (D型)



×5

91111103 TML 快装钻头 (D型快装)
TML quick coupling drill bit, $\Phi 1.0 \times 85\text{mm}$



91850501 TML 1.5 埋头钻 (D型快装)
TML 1.5 countersink drill, HA1.5mm



×3

91850201 TML 1.5 锁定钻套
TML 1.5 locking drill sleeve, 30mm



×10

91850301 TML 1.5 mini 锁定钻套
TML 1.5 locking drill sleeve, 5mm



×1

91820401 TML 快装起子杆 (D型快装)
TML quick coupling screwdriver rod, T1.3×79mm



91161502 TML 1.5 丝锥
TML 1.5 tap, HA1.5mm



91850601 TML 1.5 测深器
TML 1.5 depth gauge, 30mm

* * TML 1.5 微型锁定指掌骨安装器械包 91859901


TML 2.0 Micro phalangeal metacarpal locking series instrument Set



91820601 TML 转接头 (D型快装)
TML exchanging joint



91850401 TML 1.5双头导向器
TML 1.5 double end guider, $\Phi 1.5 \times HA 1.5mm$




91850701
TML 微型直型锁定板盒 高11mm
TML micro straight locking plate box




91850801
TML 微型T型锁定板头部2孔盒 高11mm
TML micro locking plate head 2 holes box




91850901
TML 微型髌骨锁定板盒 高11mm
TML micro condylar locking plate box



91851001
TML 微型L型锁定板盒 高11mm
TML micro L locking plate box



91851101
TML 1.5 锁定钻套盒 高11mm
TML 1.5 locking drill box



91851201
MCF 1.5 螺钉盒 高24mm
MCF 1.5 screw box



91270601 MCF 螺钉长度柱
MCF Screw Length Column, 4



91270602 MCF 螺钉长度柱
MCF Screw Length Column, 5



91270603 MCF 螺钉长度柱
MCF Screw Length Column, 6



91270604 MCF 螺钉长度柱
MCF Screw Length Column, 7



91270605 MCF 螺钉长度柱
MCF Screw Length Column, 8



91270702 TML 螺钉架 1.7、1.5
MCF Screw frame



* * TML 2.0 微型锁定指掌骨安装器械包 91829901

TML 2.0 Micro phalangeal metacarpal locking series instrument Set



91820101 TML钢板折弯器
TML plate bender, 2×1mm



91150110 TML 快装手柄
TML quick coupling handle, (D字型)



91820301 TML 快装起子杆 (D型快装)
TML quick coupling screwdriver rod, T1.6×59mm



91820401 TML 快装起子杆 (D型快装)
TML quick coupling screwdriver rod, T1.6×79mm



91820501 TML 快装钻头 (D型快装)
TML quick coupling drill bit, Φ1.5×85mm



91820601 TML 转接头 (D型快装)
TML exchanging joint



91820701 TML 2.0 锁定钻套
TML 2.0 locking drill sleeve, 30mm



91820801 TML 2.0 mini 锁定钻套
TML 2.0 locking drill sleeve, 5mm



91820901 TML 2.0 双头导向器
TML 2.0 double end guider, Φ1.5×HA2.0mm

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object. 20220629




*** * TML 2.0 微型锁定指掌骨安装器械包 91829901**

TML 2.0 Micro phalangeal metacarpal locking series instrument Set




91821001 TML 2.0 埋头钻 (D型快装)
TML 2.0 countersink drill, HA2.0mm



91821101 TML 2.0 丝锥
TML 2.0 tap, HA2.0mm



91821201 TML 2.0 测深器
TML 2.0 depth gauge, 30mm




91821301 TML 微型直型锁定板盒
TML micro straght locking plate box, 高11mm



91821401 TML 微型T型锁定板头部2孔盒
TML micro locking plate head 2 hiles box, 高11mm



91821501 TML 微型髌部锁定板盒
TML micro condylar locking plate box, 高11mm



91821601 TML 微型L型锁定板盒
TML micro l locking plate box, 高11mm



91821801 TML 2.0 空盒
TML 2.0 empty box, 高11mm



91270302 MCF 2.0 螺钉盒
MCF 2.0 screw box, 高24mm



91270601 MCF 螺钉长度柱
MCF Screw Lenght Colum, 4



91270602 MCF 螺钉长度柱
MCF Screw Lenght Colum, 5



91270603 MCF 螺钉长度柱
MCF Screw Lenght Colum, 6



91270604 MCF 螺钉长度柱
MCF Screw Lenght Colum, 7



91270605 MCF 螺钉长度柱
MCF Screw Lenght Colum, 8



91270702 TML 螺钉架 1.7、2.0
MCF Screw frame



* * 2.0锁定螺钉系列器械包 91319901

2.0 Locking Screw Instrument Set



91310101 钻套
Drill Sleeve, Ø1.5



91210301 测深器
Depth Gauge, HA2.0



91162004 锁定丝锥
Tap, Ø2.0mm



91130101 三角起子
Triangle Screwdriver



91150903 持钉套
Screw Holding Sleeve





91111503 钻头
Drill Bit, Ø1.5×150





* $\Phi 2.0$ 锁定钢板安装器械包 PPSU 91319902

$\Phi 2.0$ Lock plate mountig Instrument Set



91121501 导针
Guide pin, $\Phi 1.5 \times 150\text{mm}$



91111504 骨钻 (三角快装式)
Bone drill(Quick coupling), $\Phi 1.5\text{mm}$



91310201 钻套
Drill bush, $\Phi 1.5\text{mm}$



91162003 丝锥 (三角快装式)
Tap(Quick coupling), $\Phi 2.0\text{mm}$



91310401 测深器
Bathometer



91310501 导向器
Guider, $\Phi 1.5\text{mm}$



91310601 骨撬
Bone Lever, $\Phi 6.5\text{mm}$



91310701 骨撬
Bone Lever, $\Phi 15\text{mm}$



91310801 剥离器 (圆)
detacher



* Φ2.0锁定钢板安装器械包 PPSU 91319902

Φ2.0Lock plate mountig Instrument Set



91310901 剥离器 (平)
detacher



91130105 三角起子
Triangle screwdriver, $\Delta 2.0\text{mm}$



91130108 三角起子 (三角快装式)
Triangle screwdriver(Quick coupling), $\Delta 2.0\text{mm}$



91150101 快装手柄 (三角快装)
Quick Coupling Handle



91151103 扭力手柄 (三角快装)
The torsion lever, 0.8N.m



91311001 螺钉把持套筒
Bone Lever, $\Delta 2.0\text{mm}$



91311101 尖头复位钳 (小号)
Tip reduction forceps(Small)



91121702 尖头复位钳
Tip reduction forceps



91121802 弯型复位钳
Bent reduction forceps




91311401 折弯钳
Bending pliers





* * 2.5锁定螺钉系列器械包 91389901


2.5 Locking Screw Instrument Set



91380101 钻套
Drill Sleeve, Ø2.0



91380201 测深器
Depth Gauge



91162504 锁定丝锥
Tap, Ø2.5mm



91132302 梅花起子
Torx Screwdriver, T2.3



91150904 持钉套
Screw Holding Sleeve




91112001 钻头
Drill Bit, Ø2.0×100





* $\Phi 2.5$ 锁定钢板安装器械包 PPSU 91389902

$\Phi 2.5$ Lock plate mountig Instrument Set



91121501 导针
Guide pin, $\Phi 1.5 \times 150\text{mm}$



91112001 骨钻
Bone drill, $\Phi 2.0 \times 100\text{mm}$



91150816 钻头限位器
Limiter For Drill Bit, $\Phi 2.0\text{mm}$



91380401 钻套
Drill bush, $\Phi 1.5\text{mm}$



91380501 钻套
Drill bush, $\Phi 2.0\text{mm}$



91622503 锁定丝锥 (三角快装式)
Tap (Quick coupling), $\Phi 2.5\text{mm}$



91310401 测深器
Bathometer



91380301 导向器
Guider, $\Phi 2.0\text{mm}/\Phi 2.5\text{mm}$



91310601 骨撬
Bone Lever, $\Phi 6.5\text{mm}$

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object.20220629



* $\Phi 2.5$ 锁定钢板安装器械包 PPSU 91389902

$\Phi 2.5$ Lock plate mounting Instrument Set



91310701 骨撬
Bone Lever, $\Phi 15\text{mm}$



91310801 剥离器 (圆)
detacher



91310901 剥离器 (平)
detacher



91132310 梅花起子
Quick Coupling Handle, T2.3



91132311 梅花起子 (三角快装式)
Quick Coupling Handle, T2.3



91151103 扭力手柄 (三角快装)
The torsion lever, 0.8N



91150101 快装手柄 (三角快装)
Quick Coupling Handle



91150402 限位器扳手
Quick Coupling Handle, SW2.5



91380701 螺钉把持套筒
Bone Lever, $\Delta 2.0\text{mm}$



91311101 尖头复位钳 (小号)
Tip reduction forceps (Small)



91221702 尖头复位钳
Tip reduction forceps




91221802 弯型复位钳
Bent reduction forceps



91311401 折弯钳
Bending pliers





* * 2.7锁定螺钉系列器械包 91399901

2.7 Locking Screw Instrument Set


C




91390101 钻套
Drill Sleeve,Φ2.1



91380201 测深器
Depth Gauge



91162703 锁定丝锥
Tap,Φ2.7mm



91132302 梅花起子
Torx Screwdriver,T2.3



91150904 持钉套
Screw Holding Sleeve





91112102 钻头
Drill Bit,Φ2.1×150



*** * 微型钢板安装器械包 91219901**

Micro Plate Instrument Set



91212101 钉盒
Screw Box ,HA2.0



91212102 钉盒
Screw Box ,HA2.5



91111501 限位钻头
Limit Drill Bit,Φ1.5×8mm



91111502 钻头
Drill Bit,Φ1.5×80mm



91111801 限位钻头
Limit Drill Bit,Φ1.8×8mm



91111802 钻头
Drill Bit,Φ1.8×80mm



91162001 丝锥
Tap,HA2.0



91162501 丝锥
Tap,HA2.5



91210301 测深器
Depth Gauge,HA2.0/HA2.5



91210101 双头导向器
Double-End Guider,HA2.0/Φ1.5



91210102 双头导向器
Double-End Guider,HA2.5/Φ1.8



91130101 三角起子
Triangle Screwdriver,△2.0



91150903 持钉套
Screw Holding Sleeve



91150904 持钉套
Screw Holding Sleeve, Φ2.3mm



91132302 梅花起子
Torx Screwdriver,T2.3



91130201 十字起子
Crosshead Screwdriver



* * 微型钢板安装器械包 91219901

Micro Plate Instrument Set



91150303 螺钉夹持钳
Screw Holding Clamp



91210601 埋头钻
Countersink Drill, HA2.0



91210701 埋头钻
Countersink Drill, HA2.5



91210801 双头剥离器
Double-End detacher, 3mm



91211001 双齿皮肤拉钩
Double Tooth Skin Pulling Hook



91211101 皮肤牵开器
Skin Retractor



91210802 双头剥离器
Double-End detacher, 5mm



91211201 锐钩
Sharp Hook



91211301 骨撬
Bone Lever



91211401 钢板折弯器
Plate Bender



91211501 弯板钳
Plate Bending Plier



91211601 弯型复位钳
Curved Reduction Plier



91211701 尖头复位钳
Reduction Plier With Points



91211801 持板钳
Plate Holding Plier



91211901 尖嘴钳
Sharp Nose Plier





91212001 剪板钳
Plate Cutting Plier



* * 微型钢板安装器械包 91219902

Micro Plate Instrument Set



91111102 钻头
Drill Bit, $\Phi 1.1 \times 80\text{mm}$



91111502 钻头
Drill Bit, $\Phi 1.5 \times 80\text{mm}$



91111802 钻头
Drill Bit, $\Phi 1.8 \times 80\text{mm}$



91210103 双头导向器
Double-End Guider, HA1.5



91210101 双头导向器
Double-End Guider, HA2.0




91210102 双头导向器
Double-End Guider, HA2.5



91210701 埋头钻
Countersink Drill, HA2.5



91210601 埋头钻
Countersink Drill, HA2.0




91210803 双头剥离器
Double-End detacher, 1.5mm

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object.20220629



* * 微型钢板安装器械包 91219902

Micro Plate Instrument Set



91210801 双头剥离器
Double-End detacher,3mm



91210802 双头剥离器
Double-End detacher,5mm



91211001 双齿皮肤拉钩
Double Tooth Skin Pulling Hook



91211101 皮肤牵开器
Skin Retractor



91150303 螺钉夹持钳
Screw Holding Plier



91211201 锐钩
Sharp Hook



91211301 骨撬
Bone Lever



91211401 钢板折弯器
Plate Bender



91211501 弯板钳
Plate Bending Plier



91211801 持板钳
Plate Holding Plier



91211601 弯型复位钳
Plate Holding Plier



91211701 尖头复位钳
Reduction Plier With Points



91211901 尖嘴钳
Sharp Nose Plier(Small)



91212001 剪板钳
Plate Cutting Plier



91162001 丝锥
Tap,HA2.0



* * 微型钢板安装器械包 91219902

Micro Plate Instrument Set



91162501 丝锥
Tap,HA2.5



91111101 限位钻头
Limit Drill Bit, $\Phi 1.1 \times 8\text{mm}$



91111501 限位钻头
Limit Drill Bit, $\Phi 1.5 \times 8\text{mm}$



91111801 限位钻头
Limit Drill Bit, $\Phi 1.8 \times 8\text{mm}$



91210302 测深器
Depth Gauge,HA1.5



91210301 测深器
Depth Gauge,HA2.0/HA2.5



91212103 钉盒
Screw Box ,HA1.5



91212101 钉盒
Screw Box ,HA2.0



91212102 钉盒
Screw Box ,HA2.5



91130101 三角起子
Triangle Screwdriver, $\Delta 2.0$



91150903 持钉套
Screw Holding Sleeve



91132302 梅花起子
Torx Screwdriver,T2.3



91150904 持钉套
Screw Holding Sleeve, $\Phi 1.3\text{mm}$




91130201 十字起子
Crosshead Screwdriver





* * 微型钢板安装器械包 91219903

Micro Plate Instrument Set



91111102 钻头
Drill Bit, $\Phi 1.1 \times 80\text{mm}$



91111502 钻头
Drill Bit, $\Phi 1.5 \times 80\text{mm}$



91111802 钻头
Drill Bit, $\Phi 1.8 \times 80\text{mm}$



91212501 双头导向器
Double-End Guider, HA $\Phi 1.5$



91212601 双头导向器
Double-End Guider, HA2.0



91212701 双头导向器
Double-End Guider, HA2.5

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between contant, model, color, please refer to the object.20220629



* * 微型钢板安装器械包 91219903

Micro Plate Instrument Set



91210302 测深器
Depth Gauge, HA1.5



91211601 弯型复位钳
Plate Holding Plier



91211401 钢板折弯器
Plate Bender



91161501 丝锥 (D型快装)
Tap, HA1.5



91162002 丝锥 (D型快装)
Tap, HA2.0



91162502 丝锥 (D型快装)
Tap, HA2.5



91130202 十字起子快装式 (D型快装)
Quick Coupling Crosshead Screwdriver



91130104 三角起子快装式 (D型快装)
Quick Coupling Triangle Screwdriver



91132305 梅花起子快装式 (D型快装)
Quick Coupling Torx Screwdriver, T2.3



91150910 持钉套
Screw Holding Sleeve



91150911 持钉套
Screw Holding Sleeve, Φ2.3mm



91150111 快装手柄 (D型快装)
Quick Coupling Handle



91242201 螺钉夹持钳
Screw Holding Clamp



91211102 骨撬
Bone Lever



91150103 T型快装手柄 (D型快装)
Quick Coupling Handle (Small T Type)



* * 微型钢板安装器械包 91219903

Micro Plate Instrument Set



91130105 三角起子
Triangle Screwdriver, Δ 2.0



91150903 持钉套
Screw Holding Sleeve



91150904 持钉套
Screw Holding Sleeve




91132306 梅花起子
Torx Screwdriver, T2.3



91130203 十字起子
Crosshead Screwdriver



91211001 双齿皮肤拉钩
Double Tooth Skin Pulling Hook



91210801 双头剥离器
Double-End detacher, 3mm



91210802 双头剥离器
Double-End detacher, 5mm



91211801 持板钳
Curved Reduction Plier



91212401 埋头钻 (D型快装)
Countersink Drill, HA2.5



91212402 埋头钻 (D型快装)
Countersink Drill, HA1.5/2.0



91212201 螺钉盒
Screw Box, HA1.5



91212202 螺钉盒
Screw Box, HA2.0



91212203 螺钉盒
Screw Box, HA2.5



91212301 钢板托盘
Plate Tray, HA1.5



* * 微型钢板安装器械包 91219903

Micro Plate Instrument Set



91212302 钢板托盘一
Plate Tray 1,HA2.0



91212303 钢板托盘二
Plate Tray 2,HA2.0



91212304 钢板托盘一
Plate Tray 1,HA2.5




91212305 钢板托盘二
Plate Tray 2,HA2.5



91212306 钢板托盘三
Plate Tray 3,HA2.5



JY锁定钢板螺钉系列





* * JY 2.7、3.5锁定钢板安装器械包 91819901

JY 2.7、3.5 Locking Plate Instrument Set



91810101 上肢锁定钉盒
Plate Bending Plier



91111901 钻头 (D型快装)
Drill Bit(Type D), Φ 1.9mm



91112101 钻头 (D型快装)
Drill Bit(Type D), Φ 2.1mm



91112701 钻头 (D型快装)
Drill Bit(Type D), Φ 2.7mm



91132301 快装起子 (D型快装)
Quick Coupling Screwdriver(Type D),T2.3



91150120 快装手柄 (D型快装)
quick mount handle(Type D)



91340301 钻套
Drill Sleeve, Φ 2.1mm



91340401 弯板钳(左)
Plate Bending Plier(Left)



91340402 弯板钳(右)
Plate Bending Plier(Right)



91150119 锁定式直型手柄(微型) (D型)
Locking type direct handle (Type D)



91340701 测深尺
Depth Gauge, Φ 2.7mm



91150405 L型扳手
L wrench,SW3.5



* * JY 2.7、3.5锁定钢板安装器械包 91819901

JY 2.7、3.5 Locking Plate Instrument Set



91340901 万向导向器
Universal Guider, Φ2.1mm



91351401 转换接头 (D型快装)
Exchanging Joint(Type D)



91112501 钻头 (D型快装)
Drill Bit(Type D),Φ2.5mm



91112803 钻头 (D型快装)
Drill Bit(Type D),Φ2.8mm



91113501 钻头 (D型快装)
Drill Bit(Type D),Φ3.5mm



91133401 快装起子 (D型快装)
Quick Coupling Screwdriver(Type D),T3.4



91150118 锁定式直型手柄(小号) (D型)
Locking type direct handle (Type D)



91150121 快装手柄 (D型快装)
Type D quick mount handle(Type D)



91162701 丝锥 (D型快装)
Tap,HA2.7mm



91163503 丝锥 (D型快装)
Tap,HA3.5mm



91164002 丝锥 (D型快装)
Tap,HB4.0mm



91340101 限位器
Limiter,Φ2.1mm



91350102 限位器
Limiter,Φ2.8mm



91340201 导向器
Guider,Φ1.9/Φ2.7mm



91350201 导向器
Guider,2.5/3.5mm



91350302 钻套
Drill Sleeve,Φ2.8mm



91350401 双头导钻
Double-End Guide Drill,Φ2.5mm



91350701 测深尺
Depth GaugeΦ3.5mm




91350901 万向导向器
Universal Guider,Φ2.9mm



91351001 弯板钳A
Plate Bending Plier




91351002 弯板钳B
Plate Bending Plier



* * JY 2.7锁定钢板安装器械包 91349902

JY 2.7 Locking Plate Instrument Set



91111901 钻头 (D型快装)
Drill Bit(Type D),Φ1.9mm



91112101 钻头 (D型快装)
Drill Bit(Type D),Φ2.1mm



91340101 限位器
Limiter,Φ2.1mm



91112701 钻头 (D型快装)
Drill Bit(Type D),Φ2.7mm



91340201 导向器
Guider,Φ1.9/Φ2.7mm



91340301 钻套
Drill Sleeve,Φ2.1mm



91162701 丝锥 (D型快装)
Tap,HA2.7mm



91340401 弯板钳(左)
Plate Bending Plier(Left)



91340402 弯板钳(右)
Plate Bending Plier(Right)

Our company have final right to modify all the specification in this catalogue, If there is any inconformitybetween contant, model, color, please refer to the object.20220706



*** * JY 2.7锁定钢板安装器械包 91349902**

JY 2.7 Locking Plate Instrument Set



91340701 测深尺
Depth Gauge, $\Phi 2.7\text{mm}$



91132301 快装起子 (D型 (D型快装))
Quick Coupling Screwdriver (Type D), T2.3



91340901 万向导向器
Universal Guider, $\Phi 2.1\text{mm}$



91351401 转换接头 (D型快装)
Exchanging Joint (Type D)



91341101 锁定钉盒
Locking Screw Box, $\Phi 2.7\text{mm}$



91341201 普通钉盒
Cortical & Cancellous Screw Box, HA2.7mm



91150405 L型扳手
L wrench, SW3.5



91150120 快装手柄 (D型快装)
quick mount handle (Type D)



91150119 锁定式直型手柄(微型) (D型)
Locking type direct handle (Type D)




91341301 加压导向器 (选配)
Compression Guider, $\Phi 1.9\text{mm}$





*** * JY 2.7锁定钢板安装器械包 91349901**

JY 2.7 Locking Plate Instrument Set



91111901 钻头(D型快装)
Drill Bit(Type D), Φ 1.9mm



91112101 钻头 (D型快装)
Drill Bit(Type D), Φ 2.1mm



91340101 限位器
Limiter, Φ 2.1mm



91112701 钻头 (D型快装)
Drill Bit(Type D), Φ 2.7mm



91340201 导向器
Guider, Φ 1.9/ Φ 2.7mm



91340301 钻套
Drill Sleeve, Φ 2.1mm



91162701 丝锥 (D型快装)
Tap,HA2.7mm



91340401 弯板钳(左)
Plate Bending Plier(Left)



91340402 弯板钳(右)
Plate Bending Plier(Right)



*** * JY 2.7锁定钢板安装器械包 91349901**

JY 2.7 Locking Plate Instrument Set



91340701 测深尺
Depth Gauge, Φ2.7mm



91132301 快装起子 (D型快装)
Quick Coupling Screwdriver(Type D),T2.3



91340901 万向导向器
Universal Guider, Φ2.1mm



91351401 转换接头 (D型快装)
Exchanging Joint(Type D)



91341101 锁定钉盒
Locking Screw Box,Φ2.7mm



91341201 普通钉盒
Cortical & Cancellous Screw Box,HA2.7mm



91150405 L型扳手
L wrench,SW3.5



91150120 快装手柄 (D型快装)
quick mount handle(Type D)




91150119 锁定式直型手柄(微型) (D型)
Locking type direct handle (Type D)




91341301 加压导向器 (选配)
Compression Guider, Φ1.9mm





* * JY 3.5锁定钢板安装器械包 91359902

JY 3.5 Locking Plate Instrument Set



91112501 钻头 (D型快装)
Drill Bit(Type D),Φ2.5mm



91112803 钻头 (D型快装)
Drill Bit(Type D),Φ2.8mm



91350102 限位器
Limiter,Φ2.8mm



91113501 钻头 (D型快装)
Drill Bit(Type D),Φ3.5mm



91350201 导向器
Guider,2.5/3.5mm



91350302 钻套
Drill Sleeve,Φ2.8mm



91350401 双头导钻
Double-End Guide Drill,Φ2.5mm



91163503 丝锥 (D型快装)
Tap,HA3.5mm



91164002 丝锥 (D型快装)
Tap,HB4.0mm



* * JY 3.5锁定钢板安装器械包 91359902

JY 3.5 Locking Plate Instrument Set



91350701 测深尺
Depth Gauge, Φ3.5mm



91133401 快装起子 (D型快装)
Quick Coupling Screwdriver (Type D), T3.4



91350901 万向导向器
Universal Guider, Φ2.9mm



91351401 转换接头 (D型快装)
Exchanging Joint (Type D)



91351001 弯板钳A
Plate Bending Plier



91351002 弯板钳B
Plate Bending Plier



91351201 锁定钉盒
Locking Screw Box, Φ3.5mm



91351301 普通钉盒
Cortical & Cancellous Screw Box, HA3.5/HB4.0mm



91150405 L型扳手
L wrench, SW3.5



91150121 快装手柄 (D型快装)
quick mount handle (Type D)



91150118 锁定式直型手柄(小号) (D型)
Locking type direct handle (Type D)




91351501 加压导向器 (选配)
Compression Guider, Φ2.5mm



* * JY 3.5锁定钢板安装器械包 91359901

JY 3.5 Locking Plate Instrument Set



91112501 钻头 (D型快装)
Drill Bit, $\Phi 2.5\text{mm}$ (Type D)



91112803 钻头 (D型快装)
Drill Bit, $\Phi 2.8\text{mm}$ (Type D)



91350102 限位器
Limiter, $\Phi 2.8\text{mm}$



91113501 钻头 (D型快装)
Drill Bit, $\Phi 3.5\text{mm}$ (Type D)



91350201 导向器
Guider, 2.5/3.5mm



91350302 钻套
Drill Sleeve, $\Phi 2.8\text{mm}$



91350401 双头导钻
Double-End Guide Drill, $\Phi 2.5\text{mm}$



91163503 丝锥 (D型快装)
Tap, HA3.5mm



91164002 丝锥 (D型快装)
Tap, HB4.0mm



*** * JY 3.5锁定钢板安装器械包 91359901**

JY 3.5 Locking Plate Instrument Set



91350701 测深尺
Depth Gauge, $\Phi 3.5\text{mm}$



91133401 快装起子 (D型快装)
Quick Coupling Screwdriver (Type D), T3.4



91350901 万向导向器
Universal Guider, $\Phi 2.9\text{mm}$



91351401 转换接头 (D型快装)
Exchanging Joint (Type D)



91351001 弯板钳A
Plate Bending Plier



91351002 弯板钳B
Plate Bending Plier



91351201 锁定钉盒
Locking Screw Box, $\Phi 3.5\text{mm}$



91351301 普通钉盒
Cortical & Cancellous Screw Box, HA3.5/HB4.0mm



91150405 L型扳手
L wrench, SW3.5



91150121 快装手柄 (D型快装)
quick mount handle (Type D)



91150118 锁定式直型手柄(小号) (D型)
Locking type direct handle (Type D)




91351501 加压导向器 (选配)
Compression Guider, $\Phi 2.5\text{mm}$



* * JY 5.0锁定钢板安装器械包 91369901

JY 5.0 Locking Plate Instrument Set



91114501 钻头 (D型快装)
Drill Bit(Type D),Φ4.5mm



91113201 钻头 (D型快装)
Drill Bit(Type D),Φ3.2mm



91114301 钻头 (D型快装)
Drill Bit(Type D),Φ4.3mm



91360101 限位器
Limiter,Φ4.3mm



91360201 限位器
Limiter,Φ5.2mm



91105201 空心钻头
Hollow Drill Bit,Φ5.2mm



91360401 钻套
Drill Sleeve,Φ4.3mm




91360301 导向器
Guider,3.2mm



91360501 两用导钻
Double-End Guide Drill,Φ3.2mm



91164502 丝锥 (D型快装)
Tap,HA4.5mm



91166503 丝锥 (D型快装)
Tap,HB6.5mm



91360601 测深尺
Depth Gauge,10-120mm, Φ5.0mm



91361501 钢板折弯器(左)
Plate Bender(Left)



91361601 钢板折弯器(右)
Plate Bender(Right)



* * JY 5.0锁定钢板安装器械包 91369901

JY 5.0 Locking Plate Instrument Set



91351401 转换接头 (D型快装)
Exchanging Joint



91145301 空心梅花起子
Hollow Torx Screwdriver,T5.3



91361001 钻套
Drill Sleeve,Φ5.2mm



91123001 螺纹导针
Thread Guide Pin,Φ3.0×230mm



91361201 导针导套
Guide Sleeve For Guide Pin,Φ3.0mm



91361301 测深尺
Depth Gauge,230mm



91134301 快装起子 (D型快装)
Quick Coupling Screwdriver(Type D),T4.3



91135301 快装起子 (D型快装)
Quick Coupling Screwdriver(Type D),T5.3



91361401 万向导向器
Universal Guider,Φ4.3mm



91150117 锁定式直型手柄 (大号) (D型)
Locking type direct handle (Type D)



91150122 快装手柄 (D型快装)
quick mount handle(Type D)



91150405 L型扳手
L wrench, SW3.5



91361701 锁定钉盒
Locking Screw Box,Φ5.0/6.5mm



91361801 普通钉盒
Cortical & Cancellous Screw Box,HA4.5/HB6.5mm



91361901 加压导向器 (选配)
Compression Guider, Φ3.2mm



* * HTO锁定胫骨截骨安装器械包 91379901

HTO Tibial Osteotomy Locking Plate Instrument Set



91122001 导针
Guide Pin, $\Phi 2.0 \times 280\text{mm}$



91122501 导针
Guide Pin, $\Phi 2.5 \times 280\text{mm}$



91113201 钻头 (D型快装)
Drill Bit, $\Phi 3.2\text{mm}$



91114304 钻头 (D型快装)
Drill Bit, $\Phi 4.3 \times 245\text{mm}$



91370101 限位器
Limiter, $\Phi 4.3\text{mm}$



91360401 钻套
Drill Sleeve, $\Phi 4.3\text{mm}$



91370201 导针导套
Guide Sleeve For Guide Pin, $\Phi 2.0\text{mm}$



91134301 快装起子 (D型快装)
Quick Coupling Screwdriver, T4.3



91133501 快装起子 (D型快装)
Quick Coupling Screwdriver, SW3.5



*** * HTO锁定胫骨截骨安装器械包 91379901**

HTO Tibial Osteotomy Locking Plate Instrument Set



91371601 导向器
Guider, $\Phi 3.2/\Phi 4.3\text{mm}$



91370301 快装手柄 (D型快装)
Quick Coupling Handle, T4.3



91370401 扭力手柄 (D型快装)
Torque Handle, 4.0Nm



91370501 撑开器
Distraction Device



91370601 骨凿
Osteotome, 10mm



91370701 骨凿
Osteotome, 15mm



91370801 骨凿
Osteotome, 20mm



91370901 骨凿
Osteotome, 25mm



91371001 导向块(大)
Guide Block (Large)




91371002 导向块(小)
Guide Block (Small)



91371101 软组织撑开器
Distraction Device For Soft Tissue



91371201 测深器
Depth Gauge



91371301 骨测量器
Bone Measurer



91371401 骨科定位片1
Orthopedics Locating Piece






91371501 骨科定位片2
Orthopedics Locating Piece



* * LCC 锁定复合加压系统器械包 91749901

LCC Locking compound pressurization systemInstrument Set (FNS)



91123201 螺纹导针
Threaded Guide Pin, $\phi 3.2 \times 400\text{mm}$



91740101 导向器
Guider



91740201 导针测深器
Depth Gauge Of Guide Pin, $\phi 3.2 \times 400\text{mm}$



91740301 校正器
Corrector, $\phi 3.2$



91740401 组合铰刀
Combined reamer



91740501 保护套筒
Protective Sleeve



91740601 套管测深器
Depth gauge for sleeve



91741601 梅花起子
Torx Screwdriver, T4.3



91740701 多功能杆
Multifunction rod





*** * LCC 锁定复合加压系统器械包 91749901**

LCC Locking compound pressurization system Instrument Set (FNS)



91740801 敲击套筒
Knocking sleeve



91740901 钻头
Drill Bit, $\phi 4.3\text{mm}$



91741001 限位器
Limiter



91741101 插入手柄
Insert handle



91741201 植入手柄
Implant handle



91741301 扭力接头
Torque joint, 4.0N.m



91741401 扭力快装手柄（下肢）
Torque fast loading handle(Lower limbs)




91741501 空心直型快装手柄
Cannulated straight fast loading handle



D




(原色) 锁定钢板螺钉系列





* * 上肢锁定安装器械包 91329901

Upper Limb Locking Plate Instrument Set



91121501 导针
Guide Pin, $\Phi 1.5 \times 150\text{mm}$



91112502 钻头
Drill Bit, $\Phi 2.5 \times 125\text{mm}$



91112801 钻头
Drill Bit, $\Phi 2.8 \times 150\text{mm}$



91150813 限位器
Limiter, $\Phi 2.8$



91112802 钻头(三角快装式)
Drill Bit (Quick Coupling), $\Phi 2.8 \times 150\text{mm}$



91163502 丝锥(三角快装式)
Tap (Quick Coupling), $\Phi 3.5\text{mm}$



91163501 丝锥
Tap, HA3.5




91320401 导针钻套
Drill Sleeve For Guide Pin, $\Phi 1.5\text{mm}$



91320501 钻套
Drill Sleeve, $\Phi 2.8\text{mm}$

Our company have final right to modify all the specification in this catalogue, If there is any inconformity/between contant, model, color, please refer to the object.20220304



* * 上肢锁定安装器械包 91329901

Upper Limb Locking Plate Instrument Set



91320601 钻套扳手
Drill Sleeve Wrench



91150402 限位器扳手
Limiter Wrench, SW2.5



91320801 加压孔钻套
Drill Sleeve For Compression Hole, Φ2.5mm



91320901 测深器
Depth Gauge




91133402 梅花起子
Torx Screwdriver, T3.4



91150901 持钉套
Screw Holding Sleeve



91133403 梅花起子(三角快装式)
Torx Screwdriver (Quick Coupling), T3.4



91321001 提拉钻
Pull Lifting Drill, Φ2.8mm



91321101 打滑取出器 (三角快装式)
Slippage Extractor



91150101 快装手柄 (三角快装式)
Quick Coupling Handle






91150201 转换接头 (三角快装式)
Exchanging Joint





* * 上肢锁定接骨板安装器械包 (PPSU) 91329902

Upper Limb Locking Plate Instrument Set



91121501 导针
Guide Pin, Φ1.5×150mm



91112504 钻头 (三角快装式)
Drill Bit(Quick Coupling), Φ2.5×125mm



91112802 钻头 (三角快装式)
Drill Bit(Quick Coupling), Φ2.8×150mm



91113006 钻头 (三角快装式)
Drill Bit(Quick Coupling), Φ3.0mm



91150813 限位器
Limiter, Φ2.8



91163502 锁定丝锥(三角快装式)
Tap(Quick Coupling), Φ3.5



91321201 断钉钻头
Tap



91320401 导针钻套
Drill Sleeve For Guide Pin, Φ1.5mm

91163503 螺钉丝锥 (三角快装式)
Nail break extractor, HA3.5



91164003 丝锥 (三角快装式)
Tap, HB4.0



91730409 断钉取出器 (三角快装式)
Drill Sleeve For Guide Pin, Φ4.5mm



91321301 埋头钻 (三角快装式)
Counter drill, (Quick Coupling)



91320501 钻套
Drill Sleeve, Φ2.8



91150402 限位器扳手
Limiter Wrench, SW2.5



91133403 梅花起子 (三角快装式)
Torx Screwdriver (Quick Coupling), T.3.4



* * 上肢锁定接骨板安装器械包 (PPSU) 91329902

Upper Limb Locking Plate Instrument Set



91133407 梅花起子
Torx Screwdriver ,T3.4



91150901 持钉套
Limiter Wrench



91320901 测深器
Depth Gauge



91321101 打滑取出器 (三角快装式)
Slippage Extractor



91321001 提拉钻
Pull Lifting Drill, $\Phi 2.8$



91150201 转换接头 (三角快装式)
Exchanging Joint



91150101 快装手柄 (三角快装式)
Quick Coupling Handle



91151102 扭力手柄 (三角快装式)
The torsion lever, 1.5N.M



91320801 导向器
Guider, $\Phi 3.5$



91220301 DCP 导向器
DCP Guider, $\Phi 2.5$



91150301 螺钉夹持钳
Screw clamp forceps



91221001 折弯器
Screw clamp forceps



91221702 尖头复位钳
Tip reduction forceps



91221802 弯型复位钳
Bent reduction forceps





91221901 中心化持骨钳
Centralizing bone holding forceps



* 上肢锁定安装器械包 91329903

Upper Limb Locking Plate Instrument Set



91121501 导针
Guide Pin, $\Phi 1.5 \times 150\text{mm}$



91112502 钻头
Drill Bit, $\Phi 2.5 \times 125\text{mm}$



91112801 钻头
Drill Bit, $\Phi 2.8 \times 150\text{mm}$



91150813 限位器
Limiter, $\Phi 2.8$



91112802 钻头 (三角快装式)
Drill Bit (Quick Coupling), $\Phi 2.8 \times 150\text{mm}$



91163502 丝锥 (三角快装式)
Tap (Quick Coupling), $\Phi 3.5\text{mm}$



91163501 丝锥
Tap, HA3.5



91320401 导针钻套
Drill Sleeve For Guide Pin, $\Phi 1.5\text{mm}$



91320501 钻套
Drill Sleeve, $\Phi 2.8\text{mm}$

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object.20220629



* 上肢锁定安装器械包 91329903

Upper Limb Locking Plate Instrument Set



91320601 钻套扳手
Drill Sleeve Wrench



91150402 限位器扳手
Limiter Wrench,SW2.5



91320801 加压孔钻套
Drill Sleeve For Compression Hole,Φ2.5mm



91320901 测深器
Depth Gauge



91133402 梅花起子
Torx Screwdriver,T3.4



91150901 持钉套
Screw Holding Sleeve



91133403 梅花起子 (三角快装式)
Torx Screwdriver (Quick Coupling),T3.4



91321001 提拉钻
Pull Lifting Drill,Φ2.8mm



91321101 打滑取出器 (三角快装式)
Slippage Extractor



91150101 快装手柄 (三角快装式)
Quick Coupling Handle




91150201 转换接头 (三角快装式)
Exchanging Joint




* * 下肢锁定安装器械包 91339901


Lower Limb Locking Plate Instrument Set



91105001 空心钻
hollow drill, $\Phi 5.0 \times 250\text{mm} / \Phi 3.0$



91122002 导针
Guide Pin, $\Phi 2.0 \times 250\text{mm}$



91123002 螺纹导针
Threaded Guide Pin, $\Phi 3.0 \times 300\text{mm}$



91164501 丝锥
Tap, HA4.5



91166502 丝锥 (三角快装式)
Tap (Quick Coupling), $\Phi 6.5\text{mm}$



91165001 丝锥 (三角快装式)
Tap (Quick Coupling), $\Phi 5.0\text{mm}$



91330401 导针钻套
Drill Sleeve For Guide Pin, $\Phi 2.0\text{mm}$



91330501 钻套
Drill Sleeve, $\Phi 4.3\text{mm}$



91330601 导针钻套
Drill Sleeve For Guide Pin, $\Phi 3.0\text{mm}$

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object. 20220629



* * 下肢锁定安装器械包 91339901

Lower Limb Locking Plate Instrument Set



91114302 钻头 (三角快装式)
Drill Bit (Quick Coupling), $\Phi 4.3 \times 280\text{mm}$



91113202 钻头
Drill Bit, $\Phi 3.2 \times 150\text{mm}$



91330701 钻套
Drill Sleeve, $\Phi 5.0\text{mm}$



91114303 钻头
Drill Bit, $\Phi 4.3 \times 280\text{mm}$



91115001 钻头
Drill Bit, $\Phi 5.0 \times 280\text{mm}$



91320601 钻套扳手
Drill Sleeve Wrench, $\Phi 4.5\text{mm}$



91150403 限位器扳手
Limiter Wrench, SW3.0



91134303 梅花起子(三角快装式)
Torx Screwdriver(Quick Coupling), T4.3



91134302 梅花起子
Torx Screwdriver, T4.3



91135302 梅花起子(三角快装式)
Torx Screwdriver(Quick Coupling), T5.3



91150902 持钉套
Screw Holding Sleeve



91145301 空心梅花起子
Hollow Torx Screwdriver, T5.3




91331001 测深器
Depth Gauge



91331101 加压孔钻套
Drill Sleeve For Compression Hole, $\Phi 3.2\text{mm}$



91331201 打滑取出器 (三角快装式)
Slippage Extractor



91331301 提拉钻
Pull Lifting Drill, $\Phi 4.0\text{mm}$



91150102 快装手柄 (三角快装式)
Quick Coupling Handle



91150201 转换接头 (三角快装式)
Exchanging Joint



91150814 限位器
Limiter, $\Phi 4.3\text{mm}$





91150812 限位器
Limiter, $\Phi 5.0\text{mm}$




下肢锁定接骨板安装器械包 (PPSU) 91339902

Lower Limb Locking Plate Instrument Set



91122002 导针
Guide Pin, $\Phi 2.0 \times 250\text{mm}$



91123002 螺纹导针
Threaded Guide Pin, $\Phi 3.0 \times 300\text{mm}$



91115001 钻头
Drill Bit, $5.0 \times 280\text{mm}$



91150812 限位器
Limiter, $\Phi 5.0\text{mm}/\text{SW}3.0$



91114302 钻头(三角快装式)
Drill Bit(Quick Coupling), $\Phi 4.3 \times 280\text{mm}$



91150814 限位器
Limiter, $\Phi 4.3\text{mm}/\text{SW}3.0$



91113206 钻头(三角快装式)
Drill Bit, $\Phi 3.2 \times 150\text{mm}$



91114503 钻头(三角快装式)
hollow drill(Quick Coupling), $\Phi 4.5\text{mm}$



91331401 埋头钻(三角快装式)
Counter drill



下肢锁定接骨板安装器械包 (PPSU) 91339902

Lower Limb Locking Plate Instrument Set



91105002 空心钻
hollow drill, $\Phi 5.0 \times 280\text{mm}$



91331501 合金钻头
Alloy bit, $\Phi 6.0\text{mm}$



91164503 丝锥 (三角快装式)
Tap, HA4.5



91166504 丝锥 (三角快装式)
Tap, HB6.5



91165001 丝锥 (三角快装式)
tap, $\Phi 5.0\text{mm}$



91330401 导针钻套
Drill Sleeve For Guide Pin, $\Phi 2.0\text{mm}$



91330501 钻套
Drill Sleeve, $\Phi 4.3\text{mm}$



91330601 导针钻套
Drill Sleeve For Guide Pin, $\Phi 3.0\text{mm}$



91330701 钻套
Drill Sleeve, $\Phi 5.0\text{mm}$



91150403 限位器扳手
Limiter Wrench, SW3.0



91134308 梅花起子
Torx Screwdriver, T4.3



91150902 持钉套
Limiter Wrench



91145302 空心梅花起子
Torx Screwdriver, T5.3



91135305 梅花起子
Torx Screwdriver, T5.3



91134303 梅花起子 (三角快装式)
Torx Screwdriver (Quick Coupling), T4.3



下肢锁定接骨板安装器械包 (PPSU) 91339902

Lower Limb Locking Plate Instrument Set



91135302 梅花起子 (三角快装式)
Torx Screwdriver(Quick Coupling),T5.3



91331001 测深器
Depth Gauge



91331201 打滑取出器 (三角快装式)
Slippage Extractor



91331301 提拉钻
Pull Lifting Drill,Φ4.0mm



91150201 转换接头 (三角快装式)
Exchanging Joint



91150102 快装手柄 (三角快装式)
Quick Coupling Handle



91151101 扭力手柄 (三角快装式)
The torsion lever,4.0N.M



91730408 断钉取出器 (三角快装式)
Nail break extractor, Φ5.5mm



91331101 导向器
Guider, Φ4.5mm



91230101 DCP 导向器
DCP Guider,Φ3.2mm



91150301 螺钉夹持钳
Screw clamp forceps



91230901 钢板折弯器
Bender(Left)



91231601 尖头复位钳
Tip reduction forceps




91231701 弯型复位钳
Bent reduction forceps



91231801 中心化持骨钳
Centralizing bone holding forceps

普通钢板螺钉系列





* * 上肢骨折安装器械包 91229901

Upper Limb Fracture Instrument Set



91112502 钻头
Drill Bit, $\Phi 2.5 \times 125\text{mm}$



91113001 钻头
Drill Bit, $\Phi 3.0 \times 125\text{mm}$



91113502 钻头
Drill Bit, $\Phi 3.5 \times 125\text{mm}$



91220101 导向器
Guider, HA3.5/ $\Phi 2.5\text{mm}$



91220201 导向器
Guider, HB4.0/ $\Phi 3.0\text{mm}$



91220301 双头导钻
Double-End Guide Drill, $\Phi 2.5\text{mm}$



91163501 丝锥
Tap, HA3.5



91164001 丝锥
Tap, HB4.0



91220601 测深器
Depth Gauge





* * 上肢骨折安装器械包 91229901

Upper Limb Fracture Instrument Set



91220701 模板(直型)
Template (Straight)



91220801 模板(重建)
Template (Reconstruction)



91221001 钢板折弯器
Plate Bender



91133402 梅花起子
Torx Screwdriver, T3.4



91150901 持钉套
Screw Holding Sleeve



91150302 螺钉夹持钳
Screw Holding Clamp



91221201 埋头钻
Countersink Drill



91221301 断钉取出器
Break-Off Screw Extractor



91221401 打滑取出器
Slippage Extractor



91221501 骨膜剥离器(平)
Periosteal Detacher (Flat), $\Phi 12\text{mm}$



91221502 骨膜剥离器(圆)
Periosteal Detacher (Round), $\Phi 6.5\text{mm}$



91221601 骨撬
Bone Lever



91221701 带尖复位钳
Reduction Plier With Points



91221801 弯型复位钳
Curved Reduction Plier



91221901 自动中心化持骨钳
Self Centering Bone Holding Plier



* * 下肢骨折安装器械包 91239901

Lower Limb Fracture Instrument Set



91113202 钻头
Drill Bit, Φ3.2×150mm



91114502 钻头
Drill Bit, Φ4.5×150mm



91230101 双头导钻
Double-End Guide Drill, Φ3.2mm



91230201 双头导向器
Double-End Guider, HA4.5/Φ3.2mm



91230301 双头导向器
Double-End Guider, HB6.5/Φ4.5mm



91230701 测深器
Depth Gauge



91231901 导向器
Guider, Φ3.2mm



91231902 导向器
Guider, Φ4.5mm




91230401 C型导向器
C Type Guider




*** * 下肢骨折安装器械包 91239901**


Lower Limb Fracture Instrument Set



91164501 丝锥
Tap, HA4.5



91166501 丝锥
Tap, HB6.5



91230801 模板(直型)
Template (Straight)



91134302 梅花起子
Torx Screwdriver, T4.3




91150902 持钉套
Screw Holding Sleeve



91231101 埋头钻
Countersink Drill



91231201 断钉取出器
Break-Off Screw Extractor



91231301 打滑取出器
Slippage Extractor



91231401 骨膜剥离器(平)
Periosteal Detacher (Flat), $\Phi 16\text{mm}$



91231402 骨膜剥离器(圆)
Periosteal Detacher (Round), $\Phi 8\text{mm}$



91231701 弯型复位钳
Curved Reduction Plier



91150302 螺钉夹持钳
Screw Holding Clamp



91231601 带尖复位钳
Reduction Plier With Points



91231501 骨撬
Bone Lever



91231801 自动中心化持骨钳
Self Centering Bone Holding Plier




91230901 钢板折弯器
Plate Bender



* * 四肢骨折器械包 (重建) 91259901

Limbs Fracture Instrument Set (Reconstruction)



91250101 顶棒
Bucking Bar, $\Phi 1.5 \times 8\text{mm}$



91250201 侧弯长复位钳
Lateral Curved Long Reduction Clamp



91250301 侧弯短复位钳
Lateral Curved Short Reduction Clamp




91250401 直复位钳
Straight Reduction Clamp



91250501 小持骨钳
Bone Holding Plier (Small)



91250601 大持骨钳
Bone Holding Plier (Large)



91250701 高低双爪复位钳
Asymmetric Reduction Clamp With Two Pointed Ball Tips



91250801 三爪复位钳
Reduction Clamp With Three Pointed Ball Tips



91250901 双爪复位钳
Reduction Clamp With Two Pointed Ball Tips

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object. 20220629



* * 四肢骨折器械包 (重建) 91259901

Limbs Fracture Instrument Set (Reconstruction)



91251001 活动钳头(直)
Adjustable Clamp Head(Flat)



91251101 活动钳头(弧形)
Adjustable Clamp Head(Curved)



91251201 髋臼复位钳
Acetabular Reduction Clamp



91251301 弯板钳
Plate Bending Clamp



91251401 拉钩
Pulling Hook



91251501 T型拉钩
T Pulling Hook



91251601 牵开器
Retractor



91251701 软钻
Soft Drill Bit,Φ2.5



91251801 软钻导向器
Guider For Soft Drill Bit,Φ2.5



91251901 万向梅花起子
Universal Torx Screwdriver,T3.4



91220801 模板(直型)
Template (Straight)



91252101 模板(弧形)
Template(Curved)



91112503 钻头
Drill Bit,Φ2.5×150mm



91113002 钻头
Drill Bit,Φ3.0×150mm



91220101 双头导向器
Double-End Guider,Φ2.5/HA3.5



* * 四肢骨折器械包 (重建) 91259901


Limbs Fracture Instrument Set (Reconstruction)




91220201 双头导向器
Double-End Guider, $\Phi 3.0$ /HB4.0



91220601 测深器
Depth Gauge



91163501 丝锥
Tap, HA3.5



91164001 丝锥
Tap, HB4.0



91221001 钢板折弯器
Plate Bender



91252201 螺钉夹持钳
Screw Holding Forceps



91133402 梅花起子
Torx screwdriver ,T3.4






91150901 持钉套
Screw Holding Sleeve





* * 鹅头钉安装器械包 91269901

DHS Fracture Plate Instrument Set



91260201 DHS组合铰刀 (三角快装式)
DHS Combination Reamer



91260301 DHS角度导向器 (三角快装式)
DHS Angled Guider, 130°



91260401 导针复位器
Reduction Device For Guide Pin, $\Phi 2.5$



91260601 丝锥扶持器
Support For Tap



91150509 导针测深器
Depth Gauge For Guide Pin



91122505 螺纹导针
Threaded Guide Pin, $\Phi 2.5 \times 250\text{mm}$



91122504 导针
Guide Pin, $\Phi 2.5 \times 250\text{mm}$



91261201 打入器
Impactor



91260501 拉力螺钉丝锥
Tap For Lag Screw, $\Phi 12.5\text{mm}$



91260901 拉力螺钉扳手
Wrench For Lag Screw, SW7.0



91261001 导杆
Guide Rod, M5



91260602 上钉扶持器
Support For Screw Holding



* * 鹅头钉安装器械包 91269901

DHS Fracture Plate Instrument Set



91150102 T型快装手柄 (三角快装式)
T Quick Coupling Handle



91113202 钻头
Drill Bit, Φ3.2×150mm



91164501 丝锥
Tap, HA4.5



91114502 钻头
Drill Bit, Φ4.5×150mm



91166501 丝锥
Tap, HB6.5



91230201 双头导向器
Double-End Guider, HA4.5/Φ3.2



91230301 双头导向器
Double-End Guider, HB6.5/Φ4.5



91230101 双头导钻
Double-End Guide Drill, Φ3.2mm



91230701 测深器
Depth Gauge



91134302 梅花起子
Torx Screwsriver, T4.3




91150902 持钉套
Screw Holding Sleeve



91331201 打滑取出器 (三角快装式)
Sliping Extractor, SW3.5



91261301 加压器
Compressor



91150706 清洁导针
Cleaning Guide Pin, Φ2.5



91260701 导套
Guide Sleeve, SW7.0




91260801 导杆
Guide Rod, M5/2.5



* 断钉取出器安装器械包 91739901

Break-Off Screw Extractor Instrument Set



91730101 断钉取出钳
Break-Off Screw Removal Clamp



91730201 打滑取出器
Slippage Extractor,SW2.5×0.5mm



91730202 打滑取出器
Slippage Extractor,SW2.5×1mm



91730203 打滑取出器
Slippage Extractor,SW3.5×0.7mm



91730204 打滑取出器
Slippage Extractor,SW4.0×1mm



91730301 环钻
Trephine, Φ3.5



91730302 环钻
Trephine, Φ4.5



91730303 环钻
Trephine, Φ6.0



91730304 环钻
Trephine, Φ6.5



* 断钉取出器安装器械包 91739901

Break-Off Screw Extractor Instrument Set



91131501 六角起子
Hexagonal Screwdriver,SW1.5



91132002 六角起子
Hexagonal Screwdriver,SW2.0



91132502 六角起子
Hexagonal Screwdriver,SW2.5



91133003 六角起子
Hexagonal Screwdriver,SW3.0



91133502 六角起子
Hexagonal Screwdriver,SW3.5



91134002 六角起子
Hexagonal Screwdriver,SW4.0



91134501 六角起子
Hexagonal Screwdriver,SW4.5



91135001 六角起子
Hexagonal Screwdriver,SW5.0



91135501 六角起子
Hexagonal Screwdriver,SW5.5



91132307 梅花起子
Torx Screwdriver,T2.3



91132701 梅花起子
Torx Screwdriver,T2.7



91132903 梅花起子
Torx Screwdriver,T2.9



91133201 梅花起子
Torx Screwdriver,T3.2



91133404 梅花起子
Torx Screwdriver,T3.4



91133801 梅花起子
Torx Screwdriver,T3.8



* 断钉取出器安装器械包 91739901

Break-Off Screw Extractor Instrument Set



91134305 梅花起子
Torx Screwdriver,T4.3



91135303 梅花起子
Torx Screwdriver,T5.3



91130601 四方起子
Square Screwdriver,□1.1



91130602 四方起子
Square Screwdriver,□1.2



91130603 四方起子
Square Screwdriver,□1.5



91730401 断钉取出器
Break-Off Screw Extractor,Φ2.7



91730402 断钉取出器
Break-Off Screw Extractor,Φ3.0



91730403 断钉取出器
Break-Off Screw Extractor,Φ4.0



91730404 断钉取出器
Break-Off Screw Extractor,Φ5.0



91730405 断钉取出器
Break-Off Screw Extractor,Φ6.5



91730406 断钉取出器
Break-Off Screw Extractor,Φ8.0



91730407 断钉取出器
Break-Off Screw Extractor,Φ10



91730501 一字起子(小号)
Straight Screwdriver (Small)



91730502 一字起子(大号)
Straight Screwdriver (Large)



91130204 十字起子
Crosshead Screwdriver,Φ2.3



* 断钉取出器安装器械包 91739901

Break-Off Screw Extractor Instrument Set



91730601 尖锥
Awl, 165mm



91150104 快换手柄(T型)
Quick Coupling Handle (T Type)



91150112 快换手柄(直型)
Quick Coupling Handle (Straight Type)



91730701 骨凿
Osteotome, 8mm



91730801 通用断钉取出器
Universal Break-Off Screw Extractor, $\Phi 3.5$



91730802 通用断钉取出器
Universal Break-Off Screw Extractor, $\Phi 4.5$



91730901 尖钩
Sharp Hook, $\Phi 2.2$



91731001 合金钻头
Alloy Drill Bit, $\Phi 4$



91731002 合金钻头
Alloy Drill Bit, $\Phi 6$



91731003 合金钻头
Alloy Drill Bit, $\Phi 8$



91731101 固定钉上钉器
Screw Screwdriver, $\Phi 4.0\text{mm}$



91731102 万向钉上钉器
Universal Screw Screwdriver, T3.4



91731103 万向钉上钉器
Universal Screw Screwdriver, SW3.0





91731104 万向钉上钉器
Universal Screw Screwdriver, SW3.5



*** * 断钉取出器安装器械包 91739902**

Break-Off Screw Extractor Instrument Set



91730101 断钉取出钳
Break-Off Screw Removal Clamp



91730201 打滑取出器
Slippage Extractor, SW2.5



91730202 打滑取出器
Slippage Extractor, SW2.5



91730203 打滑取出器
Slippage Extractor, SW3.5



91730204 打滑取出器
Slippage Extractor, SW4.0



91730301 环钻
Trepine, $\Phi 3.5$




91730303 环钻
Trepine, $\Phi 6.0$



91730304 环钻
Trepine, $\Phi 6.5$

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object.20220629



* * 断钉取出器安装器械包 91739902

Break-Off Screw Extractor Instrument Set



91131501 六角起子
Hexagonal Screwdriver,SW1.5



91132002 六角起子
Hexagonal Screwdriver,SW2.0



91132502 六角起子
Hexagonal Screwdriver,SW2.5



91133003 六角起子
Hexagonal Screwdriver,SW3.0



91133502 六角起子
Hexagonal Screwdriver,SW3.5



91134002 六角起子
Hexagonal Screwdriver,SW4.0



91134501 六角起子
Hexagonal Screwdriver,SW4.5



91135001 六角起子
Hexagonal Screwdriver,SW5.0



91135501 六角起子
Hexagonal Screwdriver,SW5.5



91132307 梅花起子
Torx Screwdriver,T2.3



91132701 梅花起子
Torx Screwdriver,T2.7



91132903 梅花起子
Torx Screwdriver,T2.9



91133201 梅花起子
Torx Screwdriver,T3.2



91133404 梅花起子
Torx Screwdriver,T3.4



91133801 梅花起子
Torx Screwdriver,T3.8



* * 断钉取出器安装器械包 91739902

Break-Off Screw Extractor Instrument Set



91134305 梅花起子
Torx Screwdriver,T4.3



91135303 梅花起子
Torx Screwdriver,T5.3



91130601 四方起子
Square Screwdriver,□1.1



91130602 四方起子
Square Screwdriver,□1.2



91130603 四方起子
Square Screwdriver,□1.5



91730401 断钉取出器
Break-Off Screw Extractor,Φ2.7



91730402 断钉取出器
Break-Off Screw Extractor,Φ3.0



91730403 断钉取出器
Break-Off Screw Extractor,Φ4.0



91730404 断钉取出器
Break-Off Screw Extractor,Φ5.0



91730405 断钉取出器
Break-Off Screw Extractor,Φ6.5



91730406 断钉取出器
Break-Off Screw Extractor,Φ8.0



91730407 断钉取出器
Break-Off Screw Extractor,Φ10



91730501 一字起子(小号)
Straight Screwdriver (Small)



91730502 一字起子(大号)
Straight Screwdriver (Large)



91130204 十字起子
Crosshead Screwdriver,Φ2.3



* * 断钉取出器安装器械包 91739902

Break-Off Screw Extractor Instrument Set



91730601 尖锥
Awl, 165mm



91150104 快换手柄(T型)
Quick Coupling Handle (T Type)



91150112 快换手柄(直型)
Quick Coupling Handle (Straight Type)



91730701 骨凿
Osteotome, 8mm



91730801 通用断钉取出器
Universal Break-Off Screw Extractor, $\Phi 3.5$



91730802 通用断钉取出器
Universal Break-Off Screw Extractor, $\Phi 4.5$



91730901 尖钩
Sharp Hook, $\Phi 2.2$



91731001 合金钻头
Alloy Drill Bit, $\Phi 4$



91731002 合金钻头
Alloy Drill Bit, $\Phi 6$



91731003 合金钻头
Alloy Drill Bit, $\Phi 8$



91731101 固定钉上钉器
Screw Screwdriver, $\Phi 4.0\text{mm}$



91731102 万向钉上钉器
Universal Screw Screwdriver, T3.4




91731103 万向钉上钉器
Universal Screw Screwdriver, SW3.0



91731104 万向钉上钉器
Universal Screw Screwdriver, SW3.5




空心钉、骨针、钛缆





* 钛缆固定系统安装器械包 91729901

Titanium Cable Fixation System Instrument Set



91720101 直引针
Straight Guide Pin



91720201 弯引针-1
Curve Guide Pin



91720202 弯引针-2
Curve Guide Pin



91720301 前压缆锁
Cable Compression Lock(Front), $\Phi 1.3/\Phi 1.8$



91720401 后压缆锁
Cable Compression Lock(Back), $\Phi 1.3/\Phi 1.8$



91720501 紧缆钳
Cable Tighten Plier



91720601 剪缆钳
Cable Cut Plier



91720701 骨针剪
Skeletal Pins Scissor



91720801 锁紧扣压紧钳
Compression Plier For Fasten Buckle

* 钛缆固定系统安装器械包 91729901

Titanium Cable Fixation System Instrument Set



91720901 骨针折弯器
Skeletal Pins Bender






91721001 引缆钩
Cable Guiding Hook





* * 钛缆固定系统安装器械包 91729902

Titanium Cable Fixation System Instrument Set



91720101 直引针
Straight Guide Pin



91720201 弯引针-1
Curve Guide Pin



91720202 弯引针-2
Curve Guide Pin



91720301 前压缆锁
Cable Compression Lock(Front), $\Phi 1.3/\Phi 1.8$



91720401 后压缆锁
Cable Compression Lock(Back), $\Phi 1.3/\Phi 1.8$



91720501 紧缆钳
Cable Tighten Plier



91720602 剪缆钳
Cable Cut Plier



91720701 骨针剪
Skeletal Pins Scissor



91720801 锁紧扣压紧钳
Compression Plier For Fasten Buckle



* * 钛缆固定系统安装器械包 91729902

Titanium Cable Fixation System Instrument Set



91720901 骨针折弯器
Skeletal Pins Bender



91721002 引缆钩
Cable Guiding Hook



91720302 前压缆锁(选配)
Cable Compression Lock(Front), Φ 1.3




91720303 前压缆锁 (选配)
Cable Compression Lock(Front), Φ 1.8



91720402 后压缆锁 (选配)
Cable Compression Lock(Back), Φ 1.3




91720403 后压缆锁 (选配)
Cable Compression Lock(Back), Φ 1.8



* * $\Phi 1.7$ 空心钉安装器械包 91459901


$\Phi 1.7$ Cannulated Screw Instrument Set




91101401 空心钻
Hollow Drill Bit, $\Phi 1.4 \times 70$ mm




91120601 导针
Guide Pin, $\Phi 0.6 \times 70$ mm




91133006 六角起子
Hexagonal Screwdriver, SW2.0




91143005 空心六角起子
Hollow Hexagonal Screwdriver, SW2.0



91150301 螺钉夹持钳
Screw Holding Forceps



91150401 限位器扳手
Limiter Wrench, SW2.0



91150511 导针测深器
Depth Gauge For Guide Pin



91150710 清洁导针
Cleaning Guide Pin, $\Phi 0.6 \times 70$ mm



91150823 钻头限位器
Limiter For Drill Bit, $\Phi 1.4$ mm




* * $\Phi 1.7$ 空心钉安装器械包 91459901


$\Phi 1.7$ Cannulated Screw Instrument Set



91450101 导向器
Guider




91150110 快换手柄
Two Holes Guider



91150204 快换接头
Opening Wrench



91450201 空心丝锥
Hollow Tap





91450301 空心埋头钻
Hollow Countersink Drill






* * $\Phi 7.3$ 空心钉安装器械包 91419901

$\Phi 7.3$ Cannulated Screw Instrument Set



91122504 导针
Guide Pin, $\Phi 2.5 \times 250\text{mm}$



91122505 螺纹导针
Threaded Guide Pin, $\Phi 2.5 \times 250\text{mm}$



91105002 空心钻
Hollow Drill Bit, $\Phi 5.0 \times 250\text{mm}$



91150812 钻头限位器
Limiter For Drill Bit, $\Phi 5.0\text{mm}$



91150506 导针测深器
Depth Gauge For Guide Pin, $\Phi 7.3\text{mm}$



91410201 导针钻套
Drill Sleeve For Guide Pin, $\Phi 2.5\text{mm}$



91410301 扩张器
Expander



91410401 钻套
Drill Sleeve, $\Phi 5.0\text{mm}$



91410501 扶持器
Support



* * $\Phi 7.3$ 空心钉安装器械包 91419901

$\Phi 7.3$ Cannulated Screw Instrument Set



91410601 多孔导向器
Multi Holes Guider



91410701 两孔导向器
Two Holes Guider



91410801 开口扳手
Opening Wrench, SW10



91150403 限位器扳手
Limiter Wrench, SW3.0



91410901 空心丝锥
Hollow Tap, $\Phi 7.3$ mm



91144001 空心六角起子
Hollow Hexagonal Screwdriver, SW4.0



91154001 持钉套
Screw Holding Sleeve



91134001 六角起子
Hexagonal Screwdriver, SW4.0




91411001 空心套筒
Hollow Sleeve, SW6.0



91411101 空心埋头钻
Hollow Countersink Drill



91411201 打滑取出器
Slippage Extractor, SW4.0




91150706 清洁导针
Cleaning Guide Pin, $\Phi 2.5 \times 250$ mm



*** * Φ4.5 空心钉安装器械包 91429901**

Φ4.5 Cannulated Screw Instrument Set



91150402 限位器扳手
Limiter Wrench, SW2.5



91103201 空心钻
Hollow Drill Bit, Φ3.2×150mm



91420301 钻套
Drill Sleeve, Φ3.2mm



91150802 限位器
Limiter, Φ3.2mm



91420401 空心丝锥
Hollow Tap, Φ4.5mm



91420201 导针钻套
Drill Sleeve For Guide Pin, Φ1.5mm



91121501 导针
Guide Pin, Φ1.5×150mm



91150907 持钉套
Screw Holding Sleeve



91143002 空心六角起子
Hollow Hexagonal Screwdriver, SW3.0



91121504 螺纹导针
Threaded Guide Pin, Φ1.5×150mm



91150302 螺钉夹持钳
Screw Holding Forceps



91133002 六角起子
Hexagonal Screwdriver, SW3.0



91420501 空心埋头钻
Hollow Countersink Drill



91420601 打滑取出器
Slippage Extractor, SW3.0



91150502 导针测深器
Depth Gauge For Guide Pin




91150705 清洁导针
Cleaning Guide Pin, Φ1.5mm



* $\Phi 4.5$ 空心钉安装器械包 PPSU 91429902

$\Phi 4.5$ Headless Compression Cannulated Screw Instrument Set



91420701 螺钉盒
Tray



91121501 导针
Guide Pin, $\Phi 1.5 \times 150$ mm



91121504 螺纹导针
Threaded Guide Pin, $\Phi 1.5 \times 150$ mm



91420401 空心丝锥
Hollow Tap, $\Phi 4.5$ mm



91420301 钻套
Drill Sleeve, $\Phi 3.2$ mm



91420201 导针钻套
Drill Sleeve For Guide Pin, $\Phi 1.5$ mm



91150402 限位器扳手
Limiter Wrench, SW2.5



91103201 空心钻头
Hollow Drill Bit, $\Phi 3.2 \times 150$ mm



91143004 空心六角起子
Hollow Hexagonal Screwdriver, SW3.0



91150802 限位器
Limiter, $\Phi 3.2$ mm/SW2.5



91420501 空心埋头钻
Hollow Countersink Drill



91133005 六角起子
Hexagonal Screwdriver, SW3.0



91420601 打滑取出器
Slippage Extractor



91150907 持钉套
Screw Holding Sleeve



91150502 导针测深器
Depth Gauge For Guide Pin



91150302 螺钉夹持钳
Screw Holding Forceps




91150705 清洁导针
Cleaning Guide Pin, $\Phi 1.5$ mm



* * $\Phi 4.0$ 空心钉安装器械包 91439901

$\Phi 4.0$ Cannulated Screw Instrument Set



91121201 导针
Guide Pin, $\Phi 1.2 \times 150\text{mm}$

91103001 空心钻
Hollow Drill Bit, $\Phi 3.0 \times 150\text{mm}$

91430301 钻套
Drill Sleeve, $\Phi 3.0\text{mm}$

91121202 螺纹导针
Threaded Guide Pin, $\Phi 1.2 \times 150\text{mm}$

91430401 空心丝锥
Hollow Tap, $\Phi 4.0\text{mm}$

91430201 导针钻套
Drill Sleeve For Guide Pin, $\Phi 1.2\text{mm}$

91150402 限位器扳手
Limiter Wrench, SW2.5

91150906 持钉套
Screw Holding Sleeve

91142501 空心六角起子
Hollow Hexagonal Screwdriver, SW2.5

91150811 钻头限位器
Limiter For Drill Bit, $\Phi 3.0\text{mm}$

91150302 螺钉夹持镊
Screw Holding Forceps

91132501 SW2.5六角起子
Hexagonal Screwdriver, SW2.5

91440501 空心埋头钻
Hollow Countersink Drill

91440601 打滑取出器
Slippage Extractor, SW2.5


91150502 导针测深器
Depth Gauge For Guide Pin

91150704 清洁导针
Cleaning Guide Pin, $\Phi 1.2\text{mm}$



* $\Phi 4.0$ 空心钉安装器械包 PPSU 91439902

$\Phi 4.0$ Headless Compression Cannulated Screw Instrument Set



91430501 螺钉盒
Tray



91150402 限位器扳手
Limiter Wrench,SW2.5



91150811 钻头限位器
Limiter For Drill Bit, $\Phi 3.0$ mm



91121201 导针
Guide Pin, $\Phi 1.2 \times 150$ mm

91121202 螺纹导针
Threaded Guide Pin, $\Phi 1.2 \times 150$ mm



91430301 钻套
Drill Sleeve, $\Phi 3.0$ mm



91430201 导针钻套
Drill Sleeve For Guide Pin, $\Phi 1.2$ mm



91440501 空心埋头钻
Hollow Countersink Drill



91103001 空心钻头
Hollow Drill Bit, $\Phi 3.0 \times 150$ mm



91430401 空心丝锥
Hollow Tap, $\Phi 4.0$ mm



91142502 空心六角起子
Hollow Hexagonal Screwdriver,SW2.5



91132503 六角起子
Hexagonal Screwdriver,SW2.5



91440601 打滑取出器
Slippage Extractor



91150906 持钉套
Screw Holding Sleeve



91150302 螺钉夹持钳
Screw Holding Forceps



91150502 导针测深器
Depth Gauge For Guide Pin




91150704 清洁导针
Cleaning Guide Pin, $\Phi 1.2$ mm



* * $\Phi 3.5$ 空心钉安装器械包 91449901

$\Phi 3.5$ Cannulated Screw Instrument Set



91440401 空心丝锥
Hollow Tap, $\Phi 3.5\text{mm}$



91121003 导针
Guide Pin, $\Phi 1.0 \times 150\text{mm}$



91150402 限位器扳手
Limiter Wrench, SW2.5



91150810 钻头限位器
Limiter For Drill Bit, $\Phi 2.5\text{mm}$



91440301 钻套
Drill Sleeve, $\Phi 2.5\text{mm}$



91440201 $\Phi 1.0$ 导针钻套
Drill Sleeve For Guide Pin, $\Phi 1.0\text{mm}$



91142501 空心六角起子
Hollow Hexagonal Screwdriver, SW2.5



91132501 SW2.5六角起子
Hexagonal Screwdriver, SW2.5



91440501 空心埋头钻
Hollow Countersink Drill



91440601 打滑取出器
Slippage Extractor, SW2.5



91150502 导针测深器
Depth Gauge For Guide Pin



91150703 清洁导针
Cleaning Guide Pin, $\Phi 1.0\text{mm}$



91150906 持钉套
Screw Holding Sleeve




91150302 螺钉夹持钳
Screw Holding Forceps



* $\Phi 3.5$ 空心钉安装器械包 PPSU 91449902

$\Phi 3.5$ Headless Compression Cannulated Screw Instrument Set



91440701 螺钉盒
Tray

91121003 导针
Guide Pin, $\Phi 1.0 \times 150\text{mm}$



91440501 空心埋头钻
Hollow Countersink Drill



91150810 钻头限位器
Limiter For Drill Bit, $\Phi 2.5\text{mm}$



91150402 限位器扳手
Limiter Wrench, SW2.5

91102501 空心钻头
Hollow Drill Bit, $\Phi 2.5 \times 150\text{mm}$



91440301 钻套
Depth Gauge For Guide Pin, $\Phi 2.5\text{mm}$



91142502 空心六角起子
Hollow Hexagonal Screwdriver, SW2.5



91440401 空心丝锥
Hollow Tap, $\Phi 3.5\text{mm}$



91440201 导针钻套
Threaded Guide Pin, $\Phi 2.5\text{mm}$



91132503 六角起子
Hexagonal Screwdriver, SW2.5



91150502 导针测深器
Double-End Guider



91150703 清洁导针
Cleaning Guide Pin, $\Phi 1.0\text{mm}$



91440601 打滑取出器
Slippage Extractor



91150906 持钉套
Screw Holding Sleeve




91150302 螺钉夹持钳
Screw Holding Forceps




* * $\Phi 2.5$ 双头加压空心钉安装器械包 91519901

$\Phi 2.5$ Double-End Compression Cannulated Screw Instrument Set



91510601 托盘
Tray



91121001 导针
Guide Pin, $\Phi 1.0 \times 120\text{mm}$



91150401 限位器扳手
Limiter Wrench, SW2.0



91142301 空心梅花起子(D型快装)
Hollow Torx Screwdriver, $\Phi 2.3$



91121002 螺纹导针
Threaded Guide Pin, $\Phi 1.0 \times 120\text{mm}$



91150807 限位器
Limiter, $\Phi 1.9\text{mm}$



91132303 梅花起子(D型快装)
Torx Screwdriver (Quick Coupling), $\Phi 2.3$



91510201 空心台阶钻 (D型快装)
Hollow Step Drill, $\Phi 1.9\text{mm}$



91510101 双头导向器
Double-End Guider, $\Phi 1.0/\Phi 1.9\text{mm}$



91150110 快换手柄(D型快装)
Quick Coupling Handle



91510401 空心埋头钻 (D型快装)
Hollow Countersink Drill



91150601 加压手柄
Compression Handle



91150501 导针测深器
Depth Gauge for Guide Pin



91150301 螺钉夹持钳
Screw Holding Forceps



91510501 加压套筒
Compression Sleeve



91150701 清洁导针
Cleaning Guide Pin, $\Phi 1.0\text{mm}$



91150204 空心快换接头 (D型快装)
Hollow Quick Coupling Joint



*** * Φ3.0 双头加压空心钉安装器械包 91529901**

Φ3.0 Double-End Compression Cannulated Screw Instrument Set



91520601 托盘
Tray



91121001 导针
Guide Pin,Φ1.0×120mm



91150401 限位器扳手
Limiter Wrench,SW2.0



91142001 空心六角起子 (D型快装)
Hollow Hexagonal Screwdriver,SW2.0

91121002 螺纹导针
Threaded Guide Pin,Φ1.0×120mm



91150808 限位器
Limiter,Φ2.4mm



91132006 六角起子 (D型快装)
Hexagonal Screwdriver (Quickn Coupling),SW2.0



91520201 空心台阶钻 (D型快装)
Hollow Step Drill,Φ2.4×120mm



91520101 双头导向器
Double-End Guider,Φ1.0/Φ2.4



91150110 快换手柄 (D型快装)
Quick Coupling Handle



91520401 空心埋头钻 (D型快装)
Hollow Countersink Drill,Φ3.0mm



91150601 加压手柄
Compression Handle



91150501 导针测深器
Depth Gauge for Guide Pin



91150301 螺钉夹持钳
Screw Holding Forceps



91520501 加压套筒
Compression Sleeve,Φ3.0mm



91150701 清洁导针
Cleaning Guide Pin,Φ1.0mm



91150204 空心快换接头 (D型快装)
Hollow Quick Coupling Joint



* * $\Phi 4.0$ 双头加压空心钉安装器械包 91539901

$\Phi 4.0$ Double-End Compression Cannulated Screw Instrument Set



91530601 托盘
Tray



91121502 导针
Guide Pin, $\Phi 1.5 \times 120\text{mm}$



91150401 限位器扳手
Limiter Wrench, SW2.0



91143001 空心六角起子 (D型快装)
Hollow Hexagonal Screwdriver, SW3.0



91121503 螺纹导针
Threaded Guide Pin, $\Phi 1.5 \times 120\text{mm}$



91150809 限位器
Limiter, $\Phi 3.1\text{mm}$



91133001 六角起子 (D型快装)
Hexagonal Screwdriver, SW3.0



91530201 空心钻头 (D型快装)
Hollow Drill Bit, $\Phi 3.1 \times 120\text{mm}$



91530101 双头导向器
Double-End Guide, $\Phi 3.1/\Phi 1.5$



91150110 快换手柄 (D型快装)
Quick Coupling Handle



91530401 空心埋头钻 (D型快装)
Hollow Countersink Drill



91150601 加压手柄
Compression Handle



91150501 导针测深器
Depth Gauge For Guide Pin



91150301 螺钉夹持钳
Screw Holding Forceps



91530501 加压套筒
Compression Sleeve, $\Phi 4.0\text{mm}$



91150702 清洁导针
Cleaning Guide Pin, $\Phi 1.5\text{mm}$



91150204 空心快换接头 (D型快装)
Hollow Quick Coupling Joint




* $\Phi 4.0$ 双头加压空心钉安装器械包 91539902

$\Phi 4.0$ Double-End Compression Cannulated Screw Instrument Set



91530601 托盘
Tray



91121502 导针
Guide Pin, $\Phi 1.5 \times 120\text{mm}$



91150401 限位器扳手
Limiter Wrench, SW2.0



91143003 空心六角起子 (三角快装)
Hollow Hexagonal Screwdriver, SW3.0



91121503 螺纹导针
Threaded Guide Pin, $\Phi 1.5 \times 120\text{mm}$



91150809 限位器
Limiter, $\Phi 3.1\text{mm}$



91133004 六角起子 (三角快装)
Hexagonal Screwdriver, SW3.0



91530202 空心钻头 (三角快装)
Hollow Drill Bit, $\Phi 3.1 \times 120\text{mm}$



91530101 双头导向器
Double-End Guide, $\Phi 3.1 / \Phi 1.5$



91150113 快换手柄 (三角快装)
Quick Coupling Handle



91530402 空心埋头钻 (三角快装)
Hollow Countersink Drill



91150602 加压手柄
Compression Handle



91150501 导针测深器
Depth Gauge For Guide Pin



91150301 螺钉夹持钳
Screw Holding Forceps



91530501 加压套筒
Compression Sleeve, $\Phi 4.0\text{mm}$



91150702 清洁导针
Cleaning Guide Pin, $\Phi 1.5 \times 120\text{mm}$




91150205 空心快换接头 (三角快装)
Hollow Quick Coupling Joint



* $\Phi 4.0$ 双头加压空心钉安装器械包 PPSU 91539903

$\Phi 4.0$ Double-End Compression Cannulated Screw Instrument Set



91531001 螺钉盒
Tray



91121501 导针
Guide Pin, $\Phi 1.5 \times 150\text{mm}$



91150811 限位器
Limiter, $\Phi 3.0\text{mm}$



91150402 L 型扳手
Limiter Wrench, SW2.5

91430301 钻套
Drill Sleeve, $\Phi 3.0\text{mm}$



91530701 导针钻套
Drill Sleeve For Guide Pin



91103002 空心钻头
Hollow Drill Bit, $\Phi 3.0 \times 150\text{mm} / \Phi 1.5\text{mm}$



91530901 空心丝锥
Hollow Tap, $\Phi 4.0\text{mm} / \Phi 1.5\text{mm} / \text{P1.5}$



91143004 空心六角起子
Hollow Hexagonal Screwdriver, SW3.0/ $\Phi 1.5\text{mm}$



91530403 沉孔空心钻
Hollow Countersink Drill



91133005 六角起子
Hexagonal Screwdriver, SW3.0



91420601 打滑取出器
Slippage Extractor



91150502 导针测深器
Depth Gauge For Guide Pin



91150912 持钉套
Screw Holding Sleeve



91150705 清洁导针
Cleaning Guide Pin, $\Phi 1.5\text{mm}$




91150302 螺钉夹持钳
Screw Holding Forceps



* $\Phi 4.0$ 双头加压空心钉安装器械包 91539904

$\Phi 4.0$ Double-End Compression Cannulated Screw Instrument Set



91121501 导针
Guide Pin, $\Phi 1.5 \times 150\text{mm}$



91150811 限位器
Limiter, $\Phi 3.0\text{mm}$



91150402 L 型扳手
Limiter Wrench, SW2.5



91430301 钻套
Drill Sleeve, $\Phi 3.0\text{mm}$



91530701 导针钻套
Drill Sleeve For Guide Pin



91103002 空心钻头
Hollow Drill Bit, $\Phi 3.0 \times 150\text{mm} / \Phi 1.5\text{mm}$



91530403 沉孔空心钻
Hollow Countersink Drill



91530901 空心丝锥
Hollow Tap, $\Phi 4.0\text{mm} / \Phi 1.5\text{mm} / P1.5$



91143004 空心六角起子
Hollow Hexagonal Screwdriver, SW3.0/ $\Phi 1.5\text{mm}$



91133005 六角起子
Hexagonal Screwdriver, SW3.0



91420601 打滑取出器
Slippage Extractor



91150502 导针测深器
Depth Gauge For Guide Pin



91150705 清洁导针
Cleaning Guide Pin, $\Phi 1.5\text{mm}$



91150912 持钉套
Screw Holding Sleeve




91150302 螺钉夹持钳
Screw Holding Forceps



* 2.5 无头加压空心钉安装器械包 91619901

Φ2.5 Headless Compression Cannulated Screw Instrument Set




91120803 导针
Guide Pin,Φ0.8×120mm



91140101 空心三角起子
Hollow Triangular Screwdriver,△2.0



91610501 导针钻套
Drill Sleeve For Guide Pin,Φ0.8mm



91120804 螺纹导针
Threaded Guide Pin,Φ0.8×120mm



91130103 三角起子
Triangular Screwdriver,△2.0



91610601 钻套
Drill Sleeve,Φ2.0mm



91610701 导向器
Guider



91610801 空心钻
Hollow Drill Bit,Φ1.7/Φ2.0×120mm



91150402 限位器扳手
Limiter Wrench,SW2.5



91610901 限位器
Limiter,Φ2.0mm



91611001 沉孔空心钻
Hollow Countersink Drill



91150501 导针测深器
Depth Gauge For Guide Pin



91150708 清洁导针
Cleaning Guide Pin,Φ0.8mm



91150301 螺钉夹持镊子
Screw Holding Forceps



* * $\Phi 2.5$ 无头加压空心钉安装器械包 91619902

$\Phi 2.5$ Headless Compression Cannulated Screw Instrument Set



91611101 托盘
Tray



91120801 导针
Guide Pin, $\Phi 0.8 \times 100\text{mm}$



91140102 空心三角起子(D型快装式)
Hollow Triangular Screwdriver(Quick Coupling), $\Delta 2.0$



91610201 空心钻(D型快装式)
Hollow Drill Bit (Quick Coupling), $\Phi 1.7\text{mm}$



91120802 螺纹导针
Threaded Guide Pin, $\Phi 0.8 \times 100\text{mm}$



91130102 三角起子(D型快装式)
Triangular Screwdriver(Quick Coupling), $\Delta 2.0$



91610401 沉孔空心钻(D型快装式)
Hollow Countersink Drill (Quick Coupling)



91610101 双头导向器
Double-End Guider, $\Phi 0.8 / \Phi 1.7\text{mm}$



91150110 快换手柄 (D型快装)
Quick Coupling Handle



91150507 导针测深器
Depth Gauge For Guide Pin



91150707 导清洁导针
Cleaning Guide Pin, $\Phi 0.8\text{mm}$



91150301 螺钉夹持镊子
Screw Holding Forceps



91150804 限位器
Limiter, $\Phi 1.7\text{mm}$



91150401 限位器扳手
Limiter Wrench, SW2.0



91150204 空心快换接头 (D型快装)
Hollow Quick Coupling Joint



* 3.5 无头加压空心钉安装器械包 91629901

Φ3.5 Headless Compression Cannulated Screw Instrument Set



91620601 钻套
Drill Sleeve,Φ2.0mm



91150402 限位器扳手
Limiter Wrench,SW2.5



91620701 导向器
Guider



91620501 导针钻套
Drill Sleeve For Guide Pin,Φ1.0mm



91150805 限位器
Limiter,Φ2.1mm



91150501 导针测深器
Depth Gauge For Guide Pin



91150701 清洁导针
Cleaning Guide Pin,Φ1.0mm



91150301 螺钉夹持镊子
Screw Holding Forceps

91621001 沉孔空心钻
Hollow Countersink Drill



* * $\Phi 3.5$ 无头加压空心钉安装器械包 91629902

$\Phi 3.5$ Headless Compression Cannulated Screw Instrument Set




91621101 托盘
Tray



91620101 双头导向器
Double-End Guider, $\Phi 1.0/\Phi 2.1$ mm



91121001 导针
Guide Pin, $\Phi 1.0 \times 120$ mm



91121002 螺纹导针
Threaded Guide Pin, $\Phi 1.0 \times 120$ mm



91620201 空心钻(D型快装式)
Hollow Drill Bit (Quick Coupling), $\Phi 2.1$ mm



91620401 沉孔空心钻(D型快装式)
Hollow Countersink Drill (Quick Coupling)



91150805 限位器
Limiter, $\Phi 2.1$ mm



91150401 限位器扳手
Limiter Wrench,SW2.0



91142301 空心梅花起子(D型快装式)
Hollow Torx Screwdriver(Quick Coupling),T2.3



91132303 梅花起子(D型快装式)
Torx Screwdriver(Quick Coupling),T2.3



91150110 快换手柄 (D型快装)
Quick Coupling Handle



91150501 导针测深器
Depth Gauge For Guide Pin



91150301 螺钉夹持镊子
Screw Holding Forceps



91150204 空心快换接头 (D型快装)
Hollow Quick Coupling Joint



91150701 清洁导针
Cleaning Guide Pin, $\Phi 1.0$ mm



* 4.5 无头加压空心钉安装器械包 91639901

Φ4.5 Headless Compression Cannulated Screw Instrument Set



91121001 导针
Guide Pin,Φ1.0×120mm



91630501 导针钻套
Drill Sleeve For Guide Pin,Φ1.0mm



91630701 导向器
Guider



91121002 螺纹导针
Threaded Guide Pin,Φ1.0×120mm



91630601 钻套
Drill Sleeve,Φ2.5mm



91630801 空心钻
Hollow Drill Bit,Φ2.5×120mm



91631001 沉孔空心钻
Hollow Countersink Drill



91150402 限位器扳手
Limiter Wrench,SW2.5



91150806 限位器
Limiter,Φ2.4mm



91142902 空心梅花起子
Hollow Plum Screwdriver,T2.9



91150501 导针测深器
Depth Gauge For Guide Pin



91132902 梅花起子
Plum Screwdriver,T2.9



91150701 清洁导针
Cleaning Guide Pin,Φ1.0mm

91150301 螺钉夹持镊子
Screw Holding Forceps



* * $\Phi 4.5$ 无头加压空心钉安装器械包 91639902

$\Phi 4.5$ Headless Compression Cannulated Screw Instrument Set



91631101 托盘
Tray



91630101 双头导向器
Double-End Guider, $\Phi 1.0/\Phi 2.7$ mm



91121001 导针
Guide Pin, $\Phi 1.0 \times 120$ mm

91121002 螺纹导针
Threaded Guide Pin, $\Phi 1.0 \times 120$ mm



91630201 空心钻(D型快装式)
Hollow Drill Bit (Quick Coupling), $\Phi 2.7 \times 120$ mm



91630401 沉孔空心钻(D型快装式)
Hollow Countersink Drill (Quick Coupling)



91150806 限位器
Limiter, $\Phi 2.7$ mm



91150401 限位器扳手
Limiter Wrench,SW2.0



91142901 空心梅花起子(D型快装式)
Hollow Torx Screwdriver(Quick Coupling),T2.9mm



91132901 梅花起子(D型快装式)
Torx Screwdriver(Quick Coupling),T2.9



91150110 快换手柄 (D型快装)
Quick Coupling Handle



91150501 导针测深器
Depth Gauge For Guide Pin, $\Phi 1.0$



91150701 清洁导针
Cleaning Guide Pin, $\Phi 1.0$ mm



91150301 螺钉夹持镊子
Screw Holding Forceps



91150204 空心快换接头 (D型快装)
Hollow Quick Coupling Joint



* $\Phi 4.5$ 无头加压空心钉安装器械包 91639903

$\Phi 4.5$ Headless Compression Cannulated Screw Instrument Set



91631101 托盘
Tray



91630101 双头导向器
Double-End Guider, $\Phi 1.0/\Phi 2.7$ mm



91630202 空心钻(三角快装式)
Hollow Drill Bit (Quick Coupling), $\Phi 2.7 \times 120$ mm



91150806 限位器
Limiter, $\Phi 2.7$ mm



91150401 限位器扳手
Limiter Wrench, SW2.0



91630402 沉孔空心钻(三角快装式)
Hollow Countersink Drill (Quick Coupling)



91121001 导针
Guide Pin, $\Phi 1.0 \times 120$ mm



91121002 螺纹导针
Threaded Guide Pin, $\Phi 1.0 \times 120$ mm



91142903 空心梅花起子(三角快装式)
Hollow Torx Screwdriver (Quick Coupling), T2.9mm



91132906 梅花起子(三角快装式)
Torx Screwdriver (Quick Coupling), T2.9



91150113 快换手柄 (三角快装)
Quick Coupling Handle



91150501 导针测深器
Depth Gauge For Guide Pin, $\Phi 1.0$



91150701 清洁导针
Cleaning Guide Pin, $\Phi 1.0 \times 120$ mm



91150301 螺钉夹持镊子
Screw Holding Forceps



* * ϕ 6.5无头加压空心钉安装器械包 91649901

ϕ 6.5 Headless Compression Cannulated Screw Instrument Set



91640501 托盘
Instrument box PPSU



91122507 导针
Guide Pin, ϕ 2.5 x 180



91122508 螺纹导针
Threaded Guide Pin, ϕ 2.5 x 180



91150102 空心T型快换手柄(三角快装)
Trial Template



91150202 空心快换接头(三角快装)
Hollow Quick Coupling Joint



91150125 空心快装手柄(三角快装)
Trial Template



91150301 螺钉夹持钳
Screw Holding Forceps



91640101 双头导向器
Double-End Guider, ϕ 4.5



91150403 限位器扳手
Limiter Wrench, SW3.0



91150819 限位器
Limiter, SW3.0



91150709 清洁导针
Clean the needle, ϕ 2.5 x 200



91134003 六角起子(三角快装)
Trial Template, SW4.0



91640301 空心埋头钻(三角快装)
Hollow Countersink Drill



91150510 导针测深器
Depth Gauge for Guide Pin, ϕ 2.5 x 180



91144002 空心六角起子(三角快装)
Trial Template, ϕ 2.5/SW4.0



91650401 空心钻头
Drill Bit, ϕ 4.5/ ϕ 2.5



* * $\Phi 2.0 / \Phi 2.5$ 自攻自钻螺钉安装器械包 91249901

$\Phi 2.0 / \Phi 2.5$ Seft-Tapping & Drilling Screw Instrument Set



91111201 骨钻
Bone Drill Bit, $\Phi 1.2 \times 50\text{mm}$



91240101 埋头钻(D型快装)
Countersink Drill



91240201 一字起子 (D型快装)
Straight Screwdriver, $\Phi 2.0\text{mm}$



91240301 一字起子 (D型快装)
Straight Screwdriver, $\Phi 2.5\text{mm}$



91240401 三角套筒 (D型快装)
Triangular Sleeve, $\Phi 2.0\text{mm}$



91240501 三角套筒 (D型快装)
Triangular Sleeve, $\Phi 2.5\text{mm}$




91150110 快装手柄 (D型快装)
Quick Coupling Handle



91150303 螺钉夹持钳
Screw Holding Forceps




髓内钉系列安装器械包





* * QML 多维股骨倒打髓内钉安装器械包 92279901

QML Multi-Dimensional Femoral Reverse Intramedullary Nail Instrument Set



92270101 显影量尺
Developing Measure Gauge



91123204 导针
Guide Pin, $\Phi 3.2 \times 280\text{mm}$



92100501 转换接头
Exchanging Joint



92100401 快换手柄
Quick Coupling Handle



92270201 开髓钻
Tray, $\Phi 13.5\text{mm}$




92100709 中空开孔器
Hollow Awl



92100303 复位棒 (导针引入器)
Reduction Rod (Guide Pin Guider)



92100101 导针把持器
Guide Pin Holder



92100201 球头导针
Ball Head Guide Pin, $\Phi 2.5 \times 800\text{mm}$

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between contant, model, color, please refer to the object.20220629



*** * QML 多维股骨倒打髓内钉安装器械包 92279901**

QML Multi-Dimensional Femoral Reverse Intramedullary Nail Instrument Set



91109801 软扩刀头
Soft Expansion Tool Bit, Φ8.5mm



91109803 软扩刀头
Soft Expansion Tool Bit, Φ9.5mm



91109805 软扩刀头
Soft Expansion Tool Bit, Φ10.5mm



91109802 软扩刀头
Soft Expansion Tool Bit, Φ9mm



91109804 软扩刀头
Soft Expansion Tool Bit, Φ10mm



91109806 软扩刀头
Soft Expansion Tool Bit, Φ11mm



91109807 软扩刀头
Soft Expansion Tool Bit, Φ11.5mm



91109808 软扩刀头
Soft Expansion Tool Bit, Φ12mm



92100901 软扩刀头盒
Soft Expansion Tool Box



92100801 软扩刀杆
Soft Expansion Tool Bar



92270301 手柄
Handle



92270401 连接螺杆
Connecting Rod



92270501 连接螺杆扳手
Connecting Rod Wrench, SW6.5



92270601 主钉打拔连接杆
Insert-Extractor Connecting Rod For Main Nail



92101405 开口扳手
Opening Wrench, SW11




92270701 主钉打拔杆
T Universal Wrench For Connecting Rod



92100601 滑锤
Slide Hammer



92270801 踝部端导杆
Distal Guide rod



* * QML 多维股骨倒打髓内钉安装器械包 92279901

QML Multi-Dimensional Femoral Reverse Intramedullary Nail Instrument Set



92270901 髌部端连接导杆锁轮
Locking Nut For Connecting Handle



92271001 短型钉粗隆端导杆
Proximal Guide Rod For Short Screw



92271101 粗隆端导杆连接锁轮
Locking Nut Connecting Guide Rod



92271201 定位架
Proximal Located Bracket



92271301 前后位定位块
Anteroposterior Position Located Block



92271401 前后定位块锁轮
Locking Nut For Anteroposterior Position Located Block



92271501 定位杆套管
Located Rod Sleeve



92271601 定位杆套管针
Located Rod Trocar, M8



92271701 定位杆钻套
Drill Sleeve For Located Rod, Φ5.2



92271801 定位杆骨钻
Bone Drill For Located Rod, Φ5.2mm



92271901 定位杆平头钻
Flat Drill Bit For Located Rod



92272001 定位杆
Located Rod




92272101 定位卡块
Located Rod Plug



92272201 锁钉套管
Locking Screw Sleeve



92272301 粗隆端锁钉钻套
Proximal Locking Sleeve, Φ4.0



* * QML 多维股骨倒打髓内钉安装器械包 92279901

QML Multi-Dimensional Femoral Reverse Intramedullary Nail Instrument Set



92272401 粗隆端锁钉套管针
Proximal Locking Trocar ,Φ4.0mm



91124001 导针
Proximal Guide Pin,Φ4.0×300mm



91114001 骨钻
Proximal Bone Drill,Φ4.0×300mm



91150801 限位器
Proximal Bone Drill Limit,,Φ4.0/SW3.5



92272501 粗隆端临时定位杆
Proximal Temporary Located Rod



92272601 锁钉测深器
Depth Gauge For Locking Screw



92272701 T型锁钉扳手
T Locking Screw Wrench,,SW4.0



92272801 锁钉扳手
Locking Screw Wrench,SW4.0



92272901 锁钉扳手杆
Locking Screw Wrench Rod,SW3.5



92273001 髁部端锁钉钻套
Distal Locking Drill Sleeve ,Φ3.2mm



92273101 髁部端锁钉套管针
Distal Locking Screw Trocar , Φ3.2



91113205 骨钻
Distal Locking Screw Trocar ,Φ3.2×280mm



91150821 限位器
Ball Head Guide Pin,Φ3.2/SW3.5



92273201 尾帽六角持杆
Guide Pin Holder , SW4.0






92273401 长型钉粗隆端导杆
Guide



* * KTO - PFNA安装器械包加强型 短款 92219906

KTO-PFNA Gamma Intramedullary Nail Instrument Set



92100708 中空开孔器
Hollow Awl




92220601 近端空心钻
Proximal Hollow Drill Bit, $\Phi 16.5\text{mm}$



92233501 T型锁钉扳手
Quick Coupling Handle



92100501 转换接头
Exchanging Joint



92100201 球头导针
Ball Head Guide Pin, $\Phi 2.5 \times 800\text{mm}$



91109823 软扩刀头
Soft Expansion Tool Bit, $\Phi 8.5\text{mm}$



91109824 软扩刀头
Soft Expansion Tool Bit, $\Phi 9.0\text{mm}$



91109826 软扩刀头
Soft Expansion Tool Bit, $\Phi 10\text{mm}$



91109828 软扩刀头
Soft Expansion Tool Bit, $\Phi 11\text{mm}$



91109825 软扩刀头
Soft Expansion Tool Bit, $\Phi 9.5\text{mm}$



91109827 软扩刀头
Soft Expansion Tool Bit, $\Phi 10.5\text{mm}$



91109829 软扩刀头
Soft Expansion Tool Bit, $\Phi 11.5\text{mm}$

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object. 20220629



* * KTO - PFNA安装器械包加强型 短款 92219906

KTO-PFNA Gamma Intramedullary Nail Instrument Set



91109830 软扩刀头
Soft Expansion Tool Bit, Φ12mm



91109831 软扩刀头
Soft Expansion Tool Bit, Φ12.5mm



91109832 软扩刀头
Soft Expansion Tool Bit, Φ13mm



92100904 软扩刀盒
Soft Expansion Tool Box



92210701 手把
Handle



92210801 连接螺栓
Connecting Bolt, SW7.0



92231001 连接螺栓扳手
Connecting Bolt Wrench, SW7.0



92101401 组合扳手
Combination Wrench



92211201 近端导杆
Proximal Guide Rod

曾用编码 92210901




92211501 支持加压螺母
Supporting Compression Nut



92211601 钻套
Drill Sleeve, Φ3.2mm



92211701 套管针
Trocar Φ3.2mm



91123201 螺纹导针
Threaded Guide Pin, Φ3.2×400mm



92230101 导针测深器
Depth Gauge Of Guide Pin



92222301 近端皮质骨空心钻
Hollow Drill Bit For Proximal Cortical Bone



* * KTO - PFNA安装器械包加强型 短款 92219906

KTO-PFNA Gamma Intramedullary Nail Instrument Set



92222401 近端限位空心钻
Proximal Limit Hollow Drill Bit, Φ10.5mm



92222501 限位器
Limiter



92212101 刀片钉打入锁定起子
Locking Screwdriver For Inserting Blade Nail



92212201 刀片钉退出连接杆
Connecting Rod For Extracting Blade Screw



92212302 远端导杆
Distal Guide Rod



92100601 滑锤
Slide Hammer



92211401 刀片钉套管
Blade Nail Sleeve



92223301 远端锁钉套管
Distal Locking Screw Sleeve



92223401 钻套
Drill Sleeve, Φ4.0mm



92223501 套管针
Trocar, Φ4.0mm



91114001 骨钻
Bone Drill, Φ4.0×300mm



91150801 限位器
Limiter, Φ4.0/SW3.5



91124001 导针
Guide Pin, Φ4.0×300mm



92233301 测深器
Depth Gauge



92233501 远端锁钉扳手
Depth Gauge, SW3.5




* * KTO - PFNA安装器械包加强型 短款 92219906

KTO-PFNA Gamma Intramedullary Nail Instrument Set



92233801 T型连接螺杆万向扳手
T Universal Wrench For Connecting Rod,SW7.0



92213101 尾帽把持扳手
End Cap Holding Wrench,SW7.0



92213201 刀片钉取出器
Blade Nail Extractor



92235401 清洁导针
Cleaning Guide Pin ,3.0mm



92100301 复位棒(导针引入器)
Reduction Rod(Guide Pin Guider)



92231201 导杆锁轮
Locking Nut For Guide Rod



92224301 打入取出连接杆
Inserter For Main Nail



92100101 导针把持器
Guide Pin Holder




92100401 快装手柄
Quick Coupling Handle



92210801 连接螺栓
Connecting Bolt, SW7.0



92224401 主钉取出杆
Extracting Rod For Main Nail,M12



92213101 尾帽把持扳手
End Cap Holding Wrench,SW7.0



92220401 保护套筒
Protective Sleeve, Φ16.5mm



92221201 连接螺杆
Connecting Rod,SW7.0






92100807 股骨软杆
Soft Expansion Tool Bar



* * KTO - PFNA安装器械包加强型 长款 92219905

KTO-PFNA Gamma Intramedullary Nail Instrument Set



92100708 中空开孔器
Hollow Awl



92101001 保护套筒
Protective Sleeve




92220601 近端空心钻
Proximal Hollow Drill Bit, $\Phi 16.5\text{mm}$



92233501 T型锁钉扳手
Quick Coupling Handle



92100501 转换接头
Exchanging Joint



92100201 球头导针
Ball Head Guide Pin, $\Phi 2.5 \times 800\text{mm}$

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object. 20220629



* * KTO - PFNA安装器械包加强型 长款 92219905

KTO-PFNA Gamma Intramedullary Nail Instrument Set



91109823 软扩刀头
Soft Expansion Tool Bit, Φ8.5mm



91109824 软扩刀头
Soft Expansion Tool Bit, Φ9.0mm



91109825 软扩刀头
Soft Expansion Tool Bit, Φ9.5mm



91109826 软扩刀头
Soft Expansion Tool Bit, Φ10mm



91109827 软扩刀头
Soft Expansion Tool Bit, Φ10.5mm



91109828 软扩刀头
Soft Expansion Tool Bit, Φ11mm



91109829 软扩刀头
Soft Expansion Tool Bit, Φ11.5mm



91109830 软扩刀头
Soft Expansion Tool Bit, Φ12mm



91109831 软扩刀头
Soft Expansion Tool Bit, Φ12.5mm



91109832 软扩刀头
Soft Expansion Tool Bit, Φ13mm



92100904 软扩刀盒
Soft Expansion Tool Box



92210701 手把
Handle



92210801 连接螺栓
Connecting Bolt, SW7.0



92231001 连接螺栓扳手
Connecting Bolt Wrench, SW7.0



92101401 组合扳手
Combination Wrench

曾用编码 92210901



* * KTO - PFNA安装器械包加强型 长款 92219905

KTO-PFNA Gamma Intramedullary Nail Instrument Set



92211201 近端导杆
Proximal Guide Rod




92211501 支持加压螺母
Supporting Compression Nut



92211601 钻套
Drill Sleeve, $\Phi 3.2\text{mm}$



92211701 套管针
Trocar $\Phi 3.2\text{mm}$



91123201 螺纹导针
Threaded Guide Pin, $\Phi 3.2 \times 400\text{mm}$



92230101 导针测深器
Depth Gauge Of Guide Pin



92222301 近端皮质骨空心钻
Hollow Drill Bit For Proximal Cortical Bone



92222401 近端限位空心钻
Proximal Limit Hollow Drill Bit, $\Phi 10.5\text{mm}$



92222501 限位器
Limiter



92212101 刀片钉打入锁定起子
Locking Screwdriver For Inserting Blade Nail



92212201 刀片钉退出连接杆
Connecting Rod For Extracting Blade Screw



92234401 定位杆套管针
Trocar For Located Rod



92100601 滑锤
Slide Hammer



92211401 刀片钉套管
Blade Nail Sleeve



92223301 远端锁钉套管
Distal Locking Screw Sleeve



* * KTO - PFNA安装器械包加强型 长款 92219905

KTO-PFNA Gamma Intramedullary Nail Instrument Set



92223401 钻套
Drill Sleeve ,Φ4.0mm



92223501 套管针
Trocars,Φ4.0mm



91114001 骨钻
Bone Drill,Φ4.0×300mm



91150801 限位器
Limiter,Φ4.0/SW3.5



91124001 导针
Guide Pin,Φ4.0×300mm



92233301 测深器
Depth Gauge



92101201 套筒手柄
Sleeve Handle




92233801 T型连接螺杆万向扳手
T Universal Wrench For Connecting Rod,SW7.0



92213101 尾帽把持扳手
End Cap Holding Wrench,SW7.0



92213201 刀片钉取出器
Blade Nail Extractor



92235401 清洁导针
Cleaning Guide Pin ,3.0mm



92224401 主钉取出杆
Extracting Rod For Main Nail,M12



92100301 复位棒(导针引入器)
Reduction Rod(Guide Pin Guider)



92218301 长钉远端连接导杆
Distal Connecting Guide Rod For Long Screw



92218401 长钉远端导杆
Distal Guide Rod For Long Screw



* * KTO - PFNA安装器械包加强型 长款 92219905

KTO-PFNA Gamma Intramedullary Nail Instrument Set



92218501 长钉远端定位架 (左)
Distal Located Bracket For Long Screw



92218601 长钉远端定位架 (右)
Distal Located Bracket For Long Screw



92218701 定位杆
Located Rod



92218801 定位卡块
Located Piu



92231201 导杆锁轮
Locking Nut For Guide Rod



92234801 远端连接锁轮
Locking Nut For Connecting Guide Rod



91115201 定位杆骨钻
Bone Drill For Located Rod, $\Phi 5.2 \times 300\text{mm}$



92251301 定位杆平头钻
Flat Drill Bit For Located Rod, $\Phi 5.2\text{mm}$



92250501 定位杆钻套
Drill Sleeve For Located Rod, $\Phi 5.2\text{mm}$



92250601 定位杆套管
Located Rod Sleeve



92224301 打入取出连接杆
Inserter For Main Nail



92100101 导针把持器
Guide Pin Holder




92100401 快装手柄
Quick Coupling Handle



92100806 软杆刀杆
Soft Expansion Tool Bar






92101101 蜂窝式导针套筒
Honeycomb Guide Pin Sleeve



* * KTO-PFNA伽玛带锁髓内钉安装器械包 92219901

KTO-PFNA Gamma Intramedullary Nail Instrument Set



92100708 中空开孔器
Hollow Awl



92220401 保护套筒
Protective Sleeve, $\Phi 16.5\text{mm}$



92220501 导针套管
Guide Pin Sleeve, $\Phi 3.2\text{mm}$




92220601 近端空心钻
Proximal Hollow Drill Bit, $\Phi 16.5\text{mm}$



92100401 快装手柄
Quick Coupling Handle



92100501 转换接头
Exchanging Joint



92100201 球头导针
Ball Head Guide Pin, $\Phi 2.5 \times 800\text{mm}$




92210201 软扩
Soft Expansion, $\Phi 9.0\text{mm}$




92210301 软扩
Soft Expansion, $\Phi 10\text{mm}$

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object.20220629



* * KTO-PFNA伽玛带锁髓内钉安装器械包 92219901

KTO-PFNA Gamma Intramedullary Nail Instrument Set



92210401 软扩
Soft Expansion, $\Phi 11\text{mm}$



92210501 软扩
Soft Expansion, $\Phi 12\text{mm}$



92210601 软扩
Soft Expansion, $\Phi 13\text{mm}$



92210701 手把
Handle



92210801 连接螺栓
Connecting Bolt, SW7.0



92231001 连接螺栓扳手
Connecting Bolt Wrench, SW7.0
曾用编码 92210901



92224301 打入取出连接杆
Connecting Rod For Inserting And Extracting, M12
曾用编码 92211001



92101401 组合扳手
Combination Wrench, SW11



92100601 滑锤
Slide Hammer



92211201 近端导杆
Proximal Guide Rod



92231201 导杆锁轮
Locking Nut For Guide Rod



92211401 刀片钉套管
Blade Nail Sleeve



92211501 支持加压螺母
Supporting Compression Nut



92211601 钻套
Drill Sleeve, $\Phi 3.2\text{mm}$




92211701 套管针
Trocar



* * KTO-PFNA伽玛带锁髓内钉安装器械包 92219901

KTO-PFNA Gamma Intramedullary Nail Instrument Set



91123201 螺纹导针
Threaded Guide Pin,Φ3.2×400mm



92230101 导针测深器
Depth Gauge Of Guide Pin
曾用编码 92220101



92222301 近端皮质骨空心钻
Hollow Drill Bit For Proximal Cortical Bone
曾用编码 92211801



92222401 近端限位空心钻
Proximal Limit Hollow Drill Bit,Φ10.5mm
曾用编码 92211901



92222501 限位器
Limiter
曾用编码 92212001



92212101 刀片钉打入锁定起子
Locking Screwdriver For Inserting Blade Nail



92212201 刀片钉退出连接杆
Connecting Rod For Extracting Blade Screw



92212301 远端导杆
Distal Guide Rod



92223301 远端锁钉套管
Distal Locking Screw Sleeve
曾用编码 92212401



92223401 钻套
Drill Sleeve ,Φ4.0mm
曾用编码 92212501



92223501 套管针
Trocac,Φ4.0mm
曾用编码 92212601



91114001 骨钻
Bone Drill,Φ4.0×300mm



91124001 导针
Guide Pin,Φ4.0×300mm



91150801 限位器
Limiter,Φ4.0/SW3.5



92233301 测深器
Depth Gauge
曾用编码 92212801



* * KTO-PFNA伽玛带锁髓内钉安装器械包 92219901


KTO-PFNA Gamma Intramedullary Nail Instrument Set



92233501 T型锁钉扳手
T Locking Wrench, SW3.5
曾用编码 92212901




92233801 T型连接螺杆万向扳手
T Universal Wrench For Connecting Rod,SW7.0
曾用编码 92213001



92213101 尾帽把持扳手
End Cap Holding Wrench,SW7.0



92213201 刀片钉取出器
Blade Nail Extractor



92224401 主钉取出杆
Extracting Rod For Main Nail,M12
曾用编码 92213301



92235401 清洁导针
Cleaning Guide Pin ,Φ3.0mm
曾用编码 92213401



92100101 导针把持器
Guide Pin Holder






92100301 复位棒(导针引入器)
Reduction Rod(Guide Pin Guider)





* * KTO - PFNA伽玛带锁髓内钉安装器械包(加长版) 92219902

KTO-PFNA Gamma Intramedullary Nail Instrument Set



92100708 中空开孔器
Hollow Awl



92220401 保护套筒
Protective Sleeve, $\Phi 16.5\text{mm}$



92220501 导针套管
Guide Pin Sleeve, $\Phi 3.2\text{mm}$



92220601 近端空心钻
Proximal Hollow Drill Bit, $\Phi 16.5\text{mm}$



92213501 通用T型扳手
Quick Coupling Handle



92213502 通用小扳手
Quick Coupling Handle

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object. 20220629




* * KTO - PFNA伽玛带锁髓内钉安装器械包(加长版) 92219902

KTO-PFNA Gamma Intramedullary Nail Instrument Set



92213601 转换接头
Exchanging Joint



92100201 球头导针
Ball Head Guide Pin, $\Phi 2.5 \times 800\text{mm}$



92213701 软扩
Soft Expansion, $\Phi 9.0\text{mm}$



92213702 软扩
Soft Expansion, $\Phi 10\text{mm}$



92213703 软扩
Soft Expansion, $\Phi 11\text{mm}$



92213704 软扩
Soft Expansion, $\Phi 12\text{mm}$



92213705 软扩
Soft Expansion, $\Phi 13\text{mm}$



92213706 软扩
Soft Expansion, $\Phi 14\text{mm}$



92213707 软扩
Soft Expansion, $\Phi 15\text{mm}$



92213708 软扩
Soft Expansion, $\Phi 16\text{mm}$



92213801 手把
Handle



92210801 连接螺栓
Connecting Bolt, SW7.0



92231001 连接螺栓扳手
Connecting Bolt Wrench, SW7.0



92101401 组合扳手
Combination Wrench



92214001 滑锤
Slide Hammer

曾用编码 92210901



* * KTO - PFNA伽玛带锁髓内钉安装器械包(加长版) 92219902

KTO-PFNA Gamma Intramedullary Nail Instrument Set



92214101 打入锤
Insert Hammer



92214201 近端导杆
Proximal Guide Rod



92211401 刀片钉套管
Blade Nail Sleeve




92211501 支持加压螺母
Supporting Compression Nut



92211601 钻套
Drill Sleeve, $\Phi 3.2\text{mm}$



92211701 套管针
Trocars



91123201 螺纹导针
Threaded Guide Pin, $\Phi 3.2 \times 400\text{mm}$



92214301 导针测深器
Depth Gauge Of Guide Pin



92222301 近端皮质骨空心钻
Hollow Drill Bit For Proximal Cortical Bone



92222401 近端限位空心钻
Proximal Limit Hollow Drill Bit, $\Phi 10.5\text{mm}$



92222501 限位器
Limiter



92212101 刀片钉打入锁定起子
Locking Screwdriver For Inserting Blade Nail



92212201 刀片钉退出连接杆
Connecting Rod For Extracting Blade Screw



92214401 远端导杆
Distal Guide Rod, 180~240



92214501 远端导杆
Distal Guide Rod, 180、220



* * KTO - PFNA伽玛带锁髓内钉安装器械包(加长版) 92219902

KTO-PFNA Gamma Intramedullary Nail Instrument Set



92214601 远端导杆
Distal Guide Rod, 200、240



92223301 远端锁钉套管
Distal Locking Screw Sleeve



92223401 钻套
Drill Sleeve , Φ 4.0mm



92223501 套管针
Trocar, Φ 4.0mm



91114001 骨钻
Bone Drill, Φ 4.0 $\times$ 300mm




91150801 限位器
Limiter, Φ 4.0/SW3.5



91124001 导针
Guide Pin, Φ 4.0 $\times$ 300mm



92233301 测深器
Depth Gauge




92233501 远端锁钉扳手
Depth Gauge,SW3.5



92215101 远端锁钉扳手
Locking Screw Wrench,SW4.0




92233801 T型连接螺杆万向扳手
T Universal Wrench For Connecting Rod,SW7.0




92213101 尾帽把持扳手
End Cap Holding Wrench,SW7.0



92213201 刀片钉取出器
Blade Nail Extractor



92215201 清洁导针
Cleaning Guide Pin , Φ 2.0mm



92215301 清洁导针
Cleaning Guide Pin , Φ 2.8mm



* * KTO - PFNA伽玛带锁髓内钉安装器械包(加长版) 92219902

KTO-PFNA Gamma Intramedullary Nail Instrument Set



93334401 主钉取出杆
Extraction Rod For Main Nail, M12



92215401 显影模板
Developing Measure Gauge



92100301 复位棒(导针引入器)
Reduction Rod (Guide Pin Guider)



92215501 长钉远端连接导杆
Distal Connecting Guide Rod For Long Screw



92215601 长钉远端导杆
Distal Guide Rod For Long Screw



92215701 长钉远端定位架
Distal Located Bracket For Long Screw



92215801 长钉远端连接导杆锁轮
Distal Connecting Guide Rod Locking Nut For Long Screw



92215901 定位杆
Located Rod



92216001 定位卡块
Located Pin



92216101 长钉远端导杆锁轮
Distal Guide Rod Locking Nut For Long Screw



92216201 长钉远端定位架锁轮
Distal Located Bracket Locking Nut For Long Screw



91115201 定位杆骨钻
Bone Drill For Located Rod, Φ5.2×300mm



92216301 定位杆平头钻
Flat Drill Bit For Located Rod, Φ5.2mm



92216401 定位杆钻套
Drill Sleeve For Located Rod, Φ5.2mm



92216501 定位杆套管
Located Rod Sleeve



* * KTO - PFNA伽玛带锁髓内钉安装器械包(加长版) 92219902

KTO-PFNA Gamma Intramedullary Nail Instrument Set



92213901 打入器连接杆
Inserter For Main Nail





92213902 打入器连接块
Inserter For Main Nail





* KTO - PFNA伽玛带锁髓内钉安装器械包（碳纤维） 92219904

KTO-PFNA Gamma Intramedullary Nail Instrument Set



92216601 中空开孔器
Hollow Awl



92220401 保护套筒
Protective Sleeve, $\Phi 16.5\text{mm}$



92220501 导针套管
Guide Pin Sleeve, $\Phi 3.2\text{mm}$




92220601 近端空心钻
Proximal Hollow Drill Bit, $\Phi 16.5\text{mm}$



92100401 快装手柄
Quick Coupling Handle



92100501 转换接头
Exchanging Joint



92100201 球头导针
Ball Head Guide Pin, $\Phi 2.5 \times 800\text{mm}$



92100801 软扩刀杆
Soft Expansion Tool Bar



92100901 软扩刀盒
Soft Expansion Tool Box

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object. 20220629



* KTO - PFNA伽玛带锁髓内钉安装器械包（碳纤维） 92219904

KTO-PFNA Gamma Intramedullary Nail Instrument Set



91109801 软扩刀头
Soft Expansion Tool Bit,Φ8.5mm



91109802 软扩刀头
Soft Expansion Tool Bit,Φ9.0mm



91109803 软扩刀头
Soft Expansion Tool Bit,Φ9.5mm



91109804 软扩刀头
Soft Expansion Tool Bit,Φ10mm



91109805 软扩刀头
Soft Expansion Tool Bit,Φ10.5mm



91109806 软扩刀头
Soft Expansion Tool Bit,Φ11mm



91109807 软扩刀头
Soft Expansion Tool Bit,Φ11.5mm



91109808 软扩刀头
Soft Expansion Tool Bit,Φ12mm



91109809 软扩刀头
Soft Expansion Tool Bit,Φ12.5mm



91109810 软扩刀头
Soft Expansion Tool Bit,Φ13mm



92216701 手把
Handle



92210801 连接螺栓
Connecting Bolt, SW7.0



92231001 连接螺栓扳手
Connecting Bolt Wrench,SW7.0
曾用编码 92210901



92211001 打入取出连接杆
Connecting Rod For Inserting And Extracting,M12



92101401 组合扳手
Combination Wrench



* KTO - PFNA伽玛带锁髓内钉安装器械包（碳纤维） 92219904

KTO-PFNA Gamma Intramedullary Nail Instrument Set



92100601 滑锤
Slide Hammer



92216801 近端导杆
Proximal Guide Rod



92216901 导杆锁轮
Locking Nut For Guide Rod



92217001 刀片钉套管
Blade Nail Sleeve




92211501 支持加压螺母
Supporting Compression Nut



92217101 钻套
Drill Sleeve, Φ3.2mm



92217201 套管针
Trocár, Φ3.2mm



91123201 螺纹导针
Threaded Guide Pin, Φ3.2×400mm



92217301 导针测深器
Depth Gauge Of Guide Pin



92222301 近端皮质骨空心钻
Hollow Drill Bit For Proximal Cortical Bone



92222401 近端限位空心钻
Proximal Limit Hollow Drill Bit, Φ10.5mm



92222501 限位器
Limiter



92212101 刀片钉打入锁定起子
Locking Screwdriver For Inserting Blade Nail



92212201 刀片钉退出连接杆
Connecting Rod For Extracting Blade Screw



92217601 远端导杆
Distal Guide Rod



* KTO - PFNA伽玛带锁髓内钉安装器械包（碳纤维） 92219904

KTO-PFNA Gamma Intramedullary Nail Instrument Set



92217701 远端锁钉套管
Distal Locking Screw Sleeve



92217801 钻套
Drill Sleeve , $\Phi 4.0\text{mm}$



92217901 套管针
Trocar, $\Phi 4.0\text{mm}$



91114001 骨钻
Bone Drill, $\Phi 4.0 \times 300\text{mm}$



91124001 导针
Guide Pin, $\Phi 4.0 \times 300\text{mm}$



91150801 限位器
Limiter, $\Phi 4.0/\text{SW}3.5$



92218001 测深器
Depth Gauge



92218101 T型锁钉扳手
T Locking Wrench, SW3.5



92233801 T型连接螺杆万向扳手
T Universal Wrench For Connecting Rod,SW7.0



92218201 尾帽把持扳手
End Cap Holding Wrench,SW7.0



92213201 刀片钉取出器
Blade Nail Extractor



92213301 主钉取出杆
Extracting Rod For Main Nail,M12



92235401 清洁导针
Cleaning Guide Pin , $\Phi 3.0\text{mm}$



92100101 导针把持器
Guide Pin Holder







92100301 复位棒(导针引入器)
Reduction Rod(Guide Pin Guider)



* * KTO 滑动加压伽玛钉安装器械包 92229901

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set



91123201 螺纹导针
Threaded Guide Pin, $\Phi 3.2 \times 400\text{mm}$



91124001 导针
Guide Pin, $\Phi 4.0 \times 300\text{mm}$



92230101 导针测深器
Depth Gauge Of Guide Pin



92220201 显影量尺
Developing Measure Gauge



92100701 开孔器
Awl



92220401 保护套筒
Protective Sleeve, $\Phi 16.5\text{mm}$



92220501 导针套管
Guide Pin Sleeve, $\Phi 3.2\text{mm}$



92220601 近端空心钻
Proximal Hollow Drill Bit, $\Phi 16.5\text{mm}$




92100101 导针把持器
Guide Pin Holder

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object. 20220620



*** * KTO 滑动加压伽玛钉安装器械包 92229901**

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set



92100201 球头导针
Ball Head Guide Pin, $\Phi 2.5 \times 800\text{mm}$



92100301 复位棒(导针引入器)
Reduction Rod (Guide Pin Guider)



92221001 手把
Handle



92221101 近端导杆
Proximal Guide Rod



92221201 连接螺杆
Connecting Rod, SW7.0



92231001 连接螺杆扳手
Connecting Rod Wrench, SW7.0
曾用编码 92210901



92221401 远端导杆
Distal Guide Rod



92231201 导杆锁轮
Locking Nut For Guide Rod



92221601 近端套管
Proximal Sleeve



92221701 近端皮质套管
Proximal Cortical Sleeve



92231501 近端套管针
Proximal Trocar



92221901 支持加压螺母
Supporting Compression Nut



92222001 子钉开口钻
Opening Drill Bit For Sub Screw, $\Phi 7.3\text{mm}$



92222101 子钉钻头
Drill Bit For Sub Screw, $\Phi 6.5\text{mm}$



92222201 子钉孔定位器
Locator For Sub Screw Hole, 65mm



*** * KTO 滑动加压伽玛钉安装器械包 92229901**

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set



92222301 母钉皮质骨空心钻
Cortical Bone Hollow Drill Bit For Parent Screw, Φ10.5mm



92222401 母钉限位空心钻
Limit Hollow Drill Bit For Parent Screw, Φ10.3mm



92222501 限位器
Limiter



92222601 母钉空心丝锥
Hollow Tap For Parent Screw



92232401 母钉持杆
Holding Rod For Parent Screw



92232501 母钉中空持杆
Hollow Holding Rod For Parent Screw



92222901 母钉上钉器
Parent Screw Screwdriver



92232701 母钉加压螺母
Compression Nut For Parent Screw



92233101 母钉导针套管
Guiding Pin Sleeve For Parent Screw, Φ3.2mm



92232901 子钉上钉器
Sub Screw Screwdriver, SW4.0/M2.5



91114001 骨钻
Bone Drill Bit, Φ4.0×300mm



91150801 限位器
Limiter, Φ4.0/SW3.5



92233001 远端锁钉套管
Distal Locking Screw Sleeve



92233101 钻套
Drill Bit Sleeve, Φ4.0mm



92233201 套管针
Trocator, Φ4.0mm



*** * KTO 滑动加压伽玛钉安装器械包 92229901**

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set



92233301 测深器
Depth Gauge



92233401 锁钉自持杆
Self Holding Rod For Locking Screw,SW3.5



92233501 T型六角扳手
T Hexagonal Wrench , SW3.5



92233601 六角起子(快装式)
Hexagonal Screwdriver(Quick Coupling),SW4.0



92233701 锁钉扳手
Locking Screw Wrench,SW4.0/M2.5



92233801 连接螺栓/尾帽扳手
Wrench For Connecting Bolt/End Cap,SW7.0



92233901 软杆尾帽万向扳手
Universal Wrench For Soft Rod End Cap,SW7.0



92100401 快装手柄
Quick Coupling Handle



92100501 转换接头
Exchanging Joint



92224301 打入取出连接杆
Connecting Rod For Inserting And Extracting, M12



92224401 主钉取出杆
Extracting Rod For Main Nail,M12



92100601 滑锤
Inserter-Extractor



92101401 组合扳手
Combination Wrench



92235401 清洁导针
Cleaning Guide Pin,Φ3.0mm
曾用编码 92213401



92100801 软扩刀杆
Soft Expansion Tool Bar



* * KTO 滑动加压伽玛钉安装器械包 92229901

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set



92100901 软扩刀盒
Soft Expansion Tool Box



91109801 软扩刀头
Soft Expansion Tool Bit,Φ8.5mm



91109802 软扩刀头
Soft Expansion Tool Bit,Φ9.0mm



91109803 软扩刀头
Soft Expansion Tool Bit,Φ9.5mm



91109804 软扩刀头
Soft Expansion Tool Bit,Φ10mm



91109805 软扩刀头
Soft Expansion Tool Bit,Φ10.5mm



91109806 软扩刀头
Soft Expansion Tool Bit,Φ11mm



91109807 软扩刀头
Soft Expansion Tool Bit,Φ11.5mm



91109808 软扩刀头
Soft Expansion Tool Bit,Φ12mm



91109809 软扩刀头
Soft Expansion Tool Bit,Φ12.5mm



91109810 软扩刀头
Soft Expansion Tool Bit,Φ13mm








92222202 子钉孔定位器
Locator For Sub Screw Hole, 95mm



*** * KTO 滑动加压伽玛钉安装器械包(加长版) 92229902**

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)



91123201 螺纹导针
Threaded Guide Pin, $\Phi 3.2 \times 400\text{mm}$



91124001 导针
Guide Pin, $\Phi 4.0 \times 300\text{mm}$



92230101 导针测深器
Depth Gauge Of Guide Pin



92220201 显影量尺
Developing Measure Gauge




92100701 开孔器
Awl



92220401 保护套筒
Protective Sleeve, $\Phi 16.5\text{mm}$

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object.20220629



*** * KTO 滑动加压伽玛钉安装器械包(加长版) 92229902**

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)




92220501 导针套管
Guide Pin Sleeve, $\Phi 3.2\text{mm}$



92220601 近端空心钻
Proximal Hollow Drill Bit, $\Phi 16.5\text{mm}$



92100101 导针把持器
Guide Pin Holder



92100201 球头导针
Ball Head Guide Pin, $\Phi 2.5 \times 800\text{mm}$



92100301 复位棒(导针引入器)
Reduction Rod (Guide Pin Guider)



92221001 手把
Handle



92221101 近端导杆
Proximal Guide Rod



92221201 连接螺杆
Connecting Rod, SW7.0



92231001 连接螺杆扳手
Connecting Rod Wrench, SW7.0
曾用编码 92210901



92221401 远端导杆
Distal Guide Rod



92231201 导杆锁轮
Locking Nut For Guide Rod



92221601 近端套管
Proximal Sleeve



92221701 近端皮质套管
Proximal Cortical Sleeve



92231501 近端套管针
Proximal Trocar



92221901 支持加压螺母
Supporting Compression Nut



* * KTO 滑动加压伽玛钉安装器械包(加长版) 92229902

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)



92222001 子钉开口钻
Opening Drill Bit For Sub Screw, Φ7.3mm



92222101 子钉钻头
Drill Bit For Sub Screw, Φ6.5mm



92222201 子钉孔定位器
Locator For Sub Screw Hole, 65mm



92222301 母钉皮质骨空心钻
Cortical Bone Hollow Drill Bit For Parent Screw, Φ10.5mm



92222401 母钉限位空心钻
Cortical Bone Hollow Drill Bit For Parent Screw, Φ10.3mm



92222501 限位器
Limiter



92222601 母钉空心丝锥
Hollow Tap For Parent Screw



92222701 母钉持杆
Holding Rod For Parent Screw



92222801 母钉中空持杆
Hollow Holding Rod For Parent Screw



92222901 母钉上钉器
Parent Screw Screwdriver



92232701 母钉加压螺母
Compression Nut For Parent Screw



92223101 母钉导针套管
Guiding Pin Sleeve For Parent Screw, Φ3.2mm




92232901 子钉上钉器
Sub Screw Screwdriver, SW4.0/M2.5



91114001 骨钻
Bone Drill Bit, Φ4.0×300mm



91150801 限位器
Limiter, Φ4.0/SW3.5



* * KTO 滑动加压伽玛钉安装器械包(加长版) 92229902

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)



92233001 远端锁钉套管
Distal Locking Screw Sleeve



92233101 钻套
Drill Sleeve, $\Phi 4.0\text{mm}$



92233201 套管针
Trocar, $\Phi 4.0\text{mm}$



92233301 测深器
Depth Gauge



92233401 锁钉自持杆
Self Holding Rod For Locking Screw, SW3.5



92233501 T型六角扳手3.5
T Hexagonal Wrench 3.5



92233601 六角起子(快装式)
Hexagonal Screwdriver(Quick Coupling), SW4.0



92233701 锁钉扳手
Locking Screw Wrench, SW4.0/M2.5



92233801 连接螺栓/尾帽扳手
Wrench For Connecting Bolt/End Cap, SW7.0



92233901 软杆尾帽万向扳手
Universal Wrench For Soft Rod End Cap, SW7.0



92100401 快装手柄
Quick Coupling Handle



92100501 转换接头
Exchanging Joint




92224301 打入取出连接杆
Connecting Rod For Inserting And Extracting, M12



92224401 主钉取出杆
Extracting Rod For Main Nail, M12



92100601 滑锤
Slider-Extractor



*** * KTO 滑动加压伽玛钉安装器械包(加长版) 92229902**

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)



92101401 组合扳手
Combination Wrench, SW11



92235401 清洁导针
Cleaning Guide Pin, Φ3.0mm
曾用编码 92213401



92100801 软扩刀杆
Soft Expansion Tool Bar



92100901 软扩刀盒
Soft Expansion Tool Box



91109801 软扩刀头
Soft Expansion Tool Bit, Φ8.5mm



91109802 软扩刀头
Soft Expansion Tool Bit, Φ9.0mm



91109803 软扩刀头
Soft Expansion Tool Bit, Φ9.5mm



91109804 软扩刀头
Soft Expansion Tool Bit, Φ10mm



91109805 软扩刀头
Soft Expansion Tool Bit, Φ10.5mm



91109806 软扩刀头
Soft Expansion Tool Bit, Φ11mm



91109807 软扩刀头
Soft Expansion Tool Bit, Φ11.5mm



91109808 软扩刀头
Soft Expansion Tool Bit, Φ12mm




91109809 软扩刀头
Soft Expansion Tool Bit, Φ12.5mm



91109810 软扩刀头
Soft Expansion Tool Bit, Φ13mm



91155201 定位杆钻头
Located Rod Drill Bit, Φ5.2×300mm



*** * KTO 滑动加压伽玛钉安装器械包(加长版) 92229902**

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)



92234201 定位杆平头钻
Flat Drill Bit For Located Rod ,Φ5.2mm



92224801 定位杆
Located Rod



92234401 定位杆套管针
Located Rod Trocar



92234501 定位杆套管
Located Rod Sleeve



92234601 定位杆钻套
Drill Sleeve For Located Rod,Φ5.2mm



92225201 定位卡块
Located Plug



92225301 远端连接锁轮
Distal Connecting Locking Nut



92225401 股骨近端瞄准架
Proximal Femoral Diopter



92225501 长钉远端连接导杆
Distal Connecting Guide Rod For Long Screw



92225601 长钉远端导杆
Distal Guide Rod For Long Screw



92225701 长钉远端定位架左
Distal Located Bracket For Long Screw(Left)




92225801 长钉远端定位架右
Distal Located Bracket For Long Screw(Right)







92222202 子钉孔定位器
Locator For Sub Screw Hole, 95mm





* * KTO 滑动加压伽玛钉安装器械包 加强型(小号) 92239903

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set(Small)



91123201 螺纹导针
Threaded Guide Pin,Φ3.2×400mm



91124001 导针
Guide Pin,Φ4.0×300mm



92230101 导针测深器
Depth Gauge For Guide Pin



92230201 显影量尺
Developing Measure Gauge



92100702 开孔器
Awl



92230401 保护套筒
Protective Sleeve,Φ15mm



92230501 导针套管
Guide Pin Sleeve,Φ15mm



92230601 近端空心钻
Proximal Hollow Drill Bit,Φ15mm




92100101 导针把持器
Guide Pin Holder

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object.20220629



*** * KTO 滑动加压伽玛钉安装器械包 加强型(小号) 92239903**

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set(Small)



92100201 球头导针
Ball Head Guide Pin, $\Phi 2.5 \times 800\text{mm}$



92100301 复位棒(导针引入器)
Reduction Rod (Guide Pin Guider)



92230701 手把
Handle



92230804 近端导杆
Proximal Guide Rod



92230901 连接螺杆
Connecting Rod, SW7.0



92231001 连接螺杆扳手
Connecting Rod Wrench, SW7.0
曾用编码 92210901



92231104 远端导杆
Distal Guide Rod



92231201 导杆锁轮
Locking Nut For Guide Rod



92231301 近端套管
Proximal Sleeve



92231401 近端皮质套管
Proximal Cortical Sleeve



92231501 近端套管针
Proximal Trocar



92231601 支持加压螺母
Supporting Compression Nut




92231701 子钉开口钻
Opening Drill Bit For Sub Screw, $\Phi 7.1\text{mm}$



92231801 子钉钻头
Drill Bit For Sub Screw, $\Phi 6.5\text{mm}$



92231901 子钉孔定位器
Locator For Sub Screw Hole, 65mm



* * KTO 滑动加压伽玛钉安装器械包 加强型(小号) 92239903

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set(Small)



92232001 母钉皮质骨空心钻
Cortical Bone Hollow Drill Bit For Parent Screw, Φ9.6mm



92232101 母钉限位空心钻
Limit Hollow Drill Bit For Parent Screw, Φ9.0mm



92232201 限位器
Limiter



92232301 母钉空心丝锥
Hollow Tap For Parent Screw



92232401 母钉持杆
Holding Rod For Parent Screw



92232501 母钉中空持杆
Hollow Holding Rod For Parent Screw



92232601 母钉上钉器
Parent Screw Screwdriver



92232701 母钉加压螺母
Compression Nut For Parent Screw



92232801 母钉导针套管
Guide Pin Sleeve For Parent Screw, Φ3.2mm



92232901 子钉上钉器
Sub Screw Screwdriver, SW4.0/M2.5



91114001 骨钻
Bone Drill Bit, Φ4.0×300mm



91150801 限位器
Limiter, Φ4.0/SW3.5



92233001 远端锁钉套管
Distal Locking Screw Sleeve, Φ8.2mm



92233101 钻套
Drill Sleeve, Φ4.2mm



92233201 套管针
Trocator, Φ4.0mm



*** * KTO 滑动加压伽玛钉安装器械包 加强型(小号) 92239903**

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set(Small)



92233301 测深器
Depth Gauge



92233401 锁钉自持杆
Self Holding Rod For Locking Screw,SW3.5



92233501 T型六角扳手
T Hexagonal Wrench ,SW3.5



92233601 六角起子(快装式)
Hexagonal Screwdriver(Quick Coupling),SW4.0



92233701 锁钉扳手
Locking Screw Wrench,SW4.0/M2.5



92233801 连接螺栓/尾帽扳手
Wrench For Connecting Bolt/End Cap,SW7.0



92233901 软杆尾帽万向扳手
Universal Wrench For Soft Rod End Cap,SW7.0



92100401 快装手柄
Quick Coupling Handle



92100501 转换接头
Exchanging Joint



92234001 打入取出连接杆
Connecting Rod For Inserting And Extracting



92234101 主钉取出杆
Extracting Rod For Main Nail,M10



92100601 滑锤
Inserter-Extractor



92101401 组合扳手
Combination Wrench, SW11



92235401 清洁导针
Cleaning Guide Pin,Φ3.0mm
曾用编码 92213401



92100801 软扩刀杆
Soft Expansion Tool Bar



* * KTO 滑动加压伽玛钉安装器械包 加强型(小号) 92239903

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set(Small)



92100902 软扩刀盒
Soft Expansion Tool Box



91109801 软扩刀头
Soft Expansion Tool Bit,Φ8.5mm



91109802 软扩刀头
Soft Expansion Tool Bit,Φ9.0mm



91109803 软扩刀头
Soft Expansion Tool Bit,Φ9.5mm



91109804 软扩刀头
Soft Expansion Tool Bit,Φ10mm



91109805 软扩刀头
Soft Expansion Tool Bit,Φ10.5mm



91109806 软扩刀头
Soft Expansion Tool Bit,Φ11mm



91109807 软扩刀头
Soft Expansion Tool Bit,Φ11.5mm



91109808 软扩刀头
Soft Expansion Tool Bit,Φ12mm








92231902 子钉孔定位器
Locator For Sub Screw Hole,95mm





* * KTO 滑动加压伽玛钉安装器械包小号加强型(加长版) 92239904

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)



91123201 螺纹导针
Threaded Guide Pin, $\Phi 3.2 \times 400\text{mm}$



91124001 导针
Guide Pin, $\Phi 4.0 \times 300\text{mm}$



92230101 导针测深器
Depth Gauge For Guide Pin



92230201 显影量尺
Developing Measure Gauge



92100702 开孔器
Awl



92230401 保护套筒
Protective Sleeve, $\Phi 15\text{mm}$



*** * KTO 滑动加压伽玛钉安装器械包小号加强型(加长版) 92239904**

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)




92230501 导针套管
Guide Pin Sleeve ,Φ15mm



92230601 近端空心钻
Proximal Hollow Drill Bit,Φ15mm



92100101 导针把持器
Guide Pin Holder



92100201 球头导针
Ball Head Guide Pin,Φ2.5×800



92100301 复位棒(导针引入器)
Reduction Rod(Guide Pin Guider)



92230701 手把
Handle



92230804 近端导杆
Proximal Guide Rod



92230901 连接螺杆
Connecting Rod,SW7.0



92231001 连接螺杆扳手
Connecting Rod Wrench,SW7.0



92231104 远端导杆
Distal Guide Rod



92231201 导杆锁轮
Locking Nut For Guide Rod



92231301 近端套管
Proximal Sleeve



92231401 近端皮质套管
Proximal Cortical Sleeve



92231501 近端套管针
Proximal Trocar



92231601 支持加压螺母
Supporting Compression Nut



*** * KTO 滑动加压伽玛钉安装器械包小号加强型(加长版) 92239904**

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)



92231701 子钉开口钻
Opening Drill Bit For Sub Screw, $\Phi 7.1\text{mm}$



92231801 子钉钻头
Drill Bit For Sub Screw, $\Phi 6.5\text{mm}$



92231901 子钉孔定位器
Locator For Sub Screw Hole, 65mm



92232001 母钉皮质骨空心钻
Cortical Bone Hollow Drill Bit For Parent Screw, $\Phi 9.6\text{mm}$



92232101 母钉限位空心钻
Limit Hollow Drill Bit For Parent Screw, $\Phi 9.0\text{mm}$



92232201 限位器
Limiter



92232301 母钉空心丝锥
Hollow Tap For Parent Screw, $\Phi 9\text{mm}$



92232401 母钉持杆
Holding Rod For Parent Screw, M5



92232501 母钉中空持杆
Hollow Holding Rod For Parent Screw, M5



92232601 母钉上钉器
Parent Screw Screwdriver



92232701 母钉加压螺母
Compression Nut For Parent Screw



92232801 母钉导针套管
Guide Pin Sleeve For Parent Screw, $\Phi 3.2\text{mm}$



92232901 子钉上钉器
Sub Screw Screwdriver, SW4.0/M2.5



91114001 骨钻
Bone Drill Bit, $\Phi 4.0 \times 300\text{mm}$



91150801 限位器
Limiter, $\Phi 4.0/\text{SW}3.5$



* * KTO 滑动加压伽玛钉安装器械包小号加强型(加长版) 92239904

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)



92233001 远端锁钉套管
Distal Locking Screw Sleeve, $\Phi 8.2\text{mm}$



92233101 钻套
Drill Sleeve, $\Phi 4.2\text{mm}$



92233201 套管针
Trocar, $\Phi 4\text{mm}$



92233301 测深器
Depth Gauge



92233401 锁钉自持杆
Self Holding Rod For Locking Screw, SW3.5



92233501 T型六角扳手
T Hexagonal Wrench, SW3.5



92233601 六角起子(快装式)
Hexagonal Screwdriver(Quick Coupling), SW4.0



92233701 锁钉扳手
Locking Screw Wrench, SW4.0/M2.5



92233801 连接螺栓/尾帽扳手
Wrench For Connecting Bolt/End Cap, SW7.0



92233901 软杆尾帽万向扳手
Universal Wrench For Soft Rod End Cap, SW7.0



92100401 快装手柄
Quick Coupling Handle



92100501 转换接头
Exchanging Joint



92234001 打入取出连接杆
Connecting Rod For Inserting And Extracting, M10



92234101 主钉取出杆
Extracting Rod For Main Nail, M10



92100601 滑锤
Inserter-Extractor



*** * KTO 滑动加压伽玛钉安装器械包小号加强型(加长版) 92239904**

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)



92101401 组合扳手
Combination Wrench, SW11



92235401 清洁导针
Cleaning Guide Pin, $\Phi 3.0\text{mm}$
曾用编码92213401



92100801 软扩刀杆
Soft Expansion Tool Bar



92100902 软扩刀盒
Soft Expansion Tool Box



91109801 软扩刀头
Soft Expansion Tool Bit, $\Phi 8.5\text{mm}$



91109802 软扩刀头
Soft Expansion Tool Bit, $\Phi 9.0\text{mm}$



91109803 软扩刀头
Soft Expansion Tool Bit, $\Phi 9.5\text{mm}$



91109804 软扩刀头
Soft Expansion Tool Bit, $\Phi 10\text{mm}$



91109805 软扩刀头
Soft Expansion Tool Bit, $\Phi 10.5\text{mm}$



91109806 软扩刀头
Soft Expansion Tool Bit, $\Phi 11\text{mm}$



91109807 软扩刀头
Soft Expansion Tool Bit, $\Phi 11.5\text{mm}$



91109808 软扩刀头
Soft Expansion Tool Bit, $\Phi 12\text{mm}$



91115201 定位杆钻头
Located Rod Drill Bit, $\Phi 5.2 \times 300\text{mm}$



92234201 定位杆平头钻
Flat Drill Bit For Located Rod, $\Phi 5.2\text{mm}$



92234301 定位杆
Located Rod



* * KTO 滑动加压伽玛钉安装器械包小号加强型(加长版) 92239904

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)



92234401 定位杆套管针
Trocar For Located Rod



92234501 定位杆套管
Located Rod Sleeve



92234601 定位杆钻套
Drill Sleeve For Located Rod, $\phi 5.2\text{mm}$



92234701 定位卡块
Located Plug



92234801 远端连接锁轮
Distal Connecting Nut



92234901 股骨近端瞄准架
Proximal Femoral Diopter



92235001 长钉远端连接导杆
Distal Connecting Guide Rod For Long Screw



92235101 长钉远端导杆
Distal Guide Rod For Long Screw



92235201 长钉远端定位架左
Distal Located Bracket For Long Screw(Left)



92235301 长钉远端定位架右
Distal Located Bracket For Long Screw(Right)



92231902 子钉孔定位器
Locator For Sub Screw Hole, 95mm




92235601 近端瞄准架
Locator For Sub Screw Hole



*** * KTO滑动加压伽玛钉安装器械包小号加强型(加长)(加强型)(碳纤维款)92239905**

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)



91123201 螺纹导针
Threaded Guide Pin, Φ3.2×400mm



91124001 导针
Guide Pin, Φ4.0×300mm



92230101 导针测深器
Depth Gauge For Guide Pin



92230201 显影量尺
Developing Measure Gauge



92100702 开孔器
Awl



92230401 保护套筒
Protective Sleeve, Φ15mm

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between contant, model, color, please refer to the object.20220629



*** * KTO滑动加压伽玛钉安装器械包小号加强型(加长)(加强型)(碳纤维款)92239905**

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)




92230501 导针套管
Guide Pin Sleeve ,Φ15mm



92230601 近端空心钻
Proximal Hollow Drill Bit,Φ15mm



92100101 导针把持器
Guide Pin Holder



92100201 球头导针
Ball Head Guide Pin,Φ2.5×800



92100301 复位棒(导针引入器)
Reduction Rod(Guide Pin Guider)



92230701 手把
Handle



92230801 近端导杆
Proximal Guide Rod



92230901 连接螺杆
Connecting Rod,SW7.0



92231001 连接螺杆扳手
Connecting Rod Wrench,SW7.0



92231101 远端导杆
Distal Guide Rod



92231201 导杆锁轮
Locking Nut For Guide Rod



92231301 近端套管
Proximal Sleeve



92231401 近端皮质套管
Proximal Cortical Sleeve



92231501 近端套管针
Proximal Trocar



92231601 支持加压螺母
Supporting Compression Nut



* * KTO滑动加压伽玛钉安装器械包小号加强型(加长)(加强型)(碳纤维款)92239905

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)



92231701 子钉开口钻
Opening Drill Bit For Sub Screw, Φ7.1mm



92231801 子钉钻头
Drill Bit For Sub Screw, Φ6.5mm



92231901 子钉孔定位器
Locator For Sub Screw Hole, 65mm



92232001 母钉皮质骨空心钻
Cortical Bone Hollow Drill Bit For Parent Screw, Φ9.6mm



92232101 母钉限位空心钻
Limit Hollow Drill Bit For Parent Screw, Φ9.0mm



92232201 限位器
Limiter



92232301 母钉空心丝锥
Hollow Tap For Parent Screw, Φ9mm



92232401 母钉持杆
Holding Rod For Parent Screw, M5



92232501 母钉中空持杆
Hollow Holding Rod For Parent Screw, M5



92232601 母钉上钉器
Parent Screw Screwdriver



92232701 母钉加压螺母
Compression Nut For Parent Screw



92232801 母钉导针套管
Guide Pin Sleeve For Parent Screw, Φ3.2mm



92232901 子钉上钉器
Sub Screw Screwdriver, SW4.0/M2.5



91114001 骨钻
Bone Drill Bit, Φ4.0×300mm

91150822 限位器
Limiter, Φ4.0



* * KTO滑动加压伽玛钉安装器械包小号加强型(加长)(加强型)(碳纤维款)92239905

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)



92233001 远端锁钉套管
Distal Locking Screw Sleeve, $\Phi 8.2\text{mm}$



92233101 钻套
Drill Sleeve, $\Phi 4.2\text{mm}$



92233201 套管针
Trocar, $\Phi 4\text{mm}$



92233301 测深器
Depth Gauge



92233401 锁钉自持杆
Self Holding Rod For Locking Screw, SW3.5



92233501 T型六角扳手
T Hexagonal Wrench, SW3.5



92233601 六角起子(快装式)
Hexagonal Screwdriver (Quick Coupling), SW4.0



92233701 锁钉扳手
Locking Screw Wrench, SW4.0/M2.5



92233801 连接螺栓/尾帽扳手
Wrench For Connecting Bolt/End Cap, SW7.0



92233901 软杆尾帽万向扳手
Universal Wrench For Soft Rod End Cap, SW7.0



92100401 快装手柄
Quick Coupling Handle



92100501 转换接头
Exchanging Joint



92234001 打入取出连接杆
Connecting Rod For Inserting And Extracting, M10



92234101 主钉取出杆
Extracting Rod For Main Nail, M10



92100601 滑锤
Slider-Extractor



* * KTO滑动加压伽玛钉安装器械包小号加强型(加长)(加强型)(碳纤维款)92239905

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)



92101401 组合扳手
Combination Wrench, SW11



92238401 清洁导针
Cleaning Guide Pin, Φ3.0mm



92100801 软扩刀杆
Soft Expansion Tool Bar



92100902 软扩刀盒
Soft Expansion Tool Box



91109801 软扩刀头
Soft Expansion Tool Bit, Φ8.5mm



91109802 软扩刀头
Soft Expansion Tool Bit, Φ9.0mm



91109803 软扩刀头
Soft Expansion Tool Bit, Φ9.5mm



91109804 软扩刀头
Soft Expansion Tool Bit, Φ10mm



91109805 软扩刀头
Soft Expansion Tool Bit, Φ10.5mm



91109806 软扩刀头
Soft Expansion Tool Bit, Φ11mm



91109807 软扩刀头
Soft Expansion Tool Bit, Φ11.5mm



91109808 软扩刀头
Soft Expansion Tool Bit, Φ12mm



92231902 子钉孔定位器
Locator For Sub Screw Hole, 95mm

92235601 近端瞄准架
Flat

股骨近端瞄准架锁轮 M8×1
Located



* * KTO滑动加压伽玛钉安装器械包小号加强型(加长)(加强型)(碳纤维款)92239905

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)

92235701 近端瞄准架-左
Trocar For Located Rod

92235801 近端瞄准架-右
Located Rod Sleeve

92235901 远端瞄准架
Drill Sleeve For Located Rod, $\Phi 5.2\text{mm}$




92236001 显影套针 $\Phi 11$
Located Plug

92234801 远端连接锁轮
Distal Connecting Nut




92234901 股骨近端瞄准架
Proximal Femoral Diopter





*** * KTO 滑动加压伽玛钉安装器械包(小号) 92239901**

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set(Small)



91123201 螺纹导针
Threaded Guide Pin, $\Phi 3.2 \times 400\text{mm}$



91124001 导针
Guide Pin, $\Phi 4.0 \times 300\text{mm}$



92230101 导针测深器
Depth Gauge For Guide Pin



92230201 显影量尺
Developing Measure Gauge



92100702 开孔器
Awl



92230401 保护套筒
Protective Sleeve, $\Phi 15\text{mm}$



92230501 导针套管
Guide Pin Sleeve, $\Phi 15\text{mm}$



92230601 近端空心钻
Proximal Hollow Drill Bit, $\Phi 15\text{mm}$




92100101 导针把持器
Guide Pin Holder

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object. 20220629



* * KTO 滑动加压伽玛钉安装器械包(小号) 92239901

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set(Small)



92100201 球头导针
Ball Head Guide Pin, $\Phi 2.5 \times 800\text{mm}$



92100301 复位棒(导针引入器)
Reduction Rod(Guide Pin Guider)



92230701 手把
Handle



92230801 近端导杆
Proximal Guide Rod



92230901 连接螺杆
Connecting Rod, SW7.0



92231001 连接螺杆扳手
Connecting Rod Wrench, SW7.0
曾用编码 92210901



92231101 远端导杆
Distal Guide Rod



92231201 导杆锁轮
Locking Nut For Guide Rod



92231301 近端套管
Proximal Sleeve



92231401 近端皮质套管
Proximal Cortical Sleeve



92231501 近端套管针
Proximal Trocar



92231601 支持加压螺母
Supporting Compression Nut



92231701 子钉开口钻
Opening Drill Bit For Sub Screw, $\Phi 7.1\text{mm}$



92231801 子钉钻头
Drill Bit For Sub Screw, $\Phi 6.5\text{mm}$



92231901 子钉孔定位器
Locator For Sub Screw Hole, 65mm



* * KTO 滑动加压伽玛钉安装器械包(小号) 92239901

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set(Small)



92232001 母钉皮质骨空心钻
Cortical Bone Hollow Drill Bit For Parent Screw, Φ9.6mm



92232101 母钉限位空心钻
Limit Hollow Drill Bit For Parent Screw, Φ9.0mm



92232201 限位器
Limiter



92232301 母钉空心丝锥
Hollow Tap For Parent Screw



92232401 母钉持杆
Holding Rod For Parent Screw



92232501 母钉中空持杆
Hollow Holding Rod For Parent Screw



92232601 母钉上钉器
Parent Screw Screwdriver



92232701 母钉加压螺母
Compression Nut For Parent Screw



92232801 母钉导针套管
Guide Pin Sleeve For Parent Screw, Φ3.2mm



92232901 子钉上钉器
Sub Screw Screwdriver, SW4.0/M2.5



91114001 骨钻
Bone Drill Bit, Φ4.0×300mm



91150801 限位器
Limiter, Φ4.0/SW3.5



92233001 远端锁钉套管
Distal Locking Screw Sleeve, Φ8.2mm



92233101 钻套
Drill Sleeve, Φ4.2mm



92233201 套管针
Trocar, Φ4.0mm



* * KTO 滑动加压伽玛钉安装器械包(小号) 92239901

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set(Small)



92233301 测深器
Depth Gauge



92233401 锁钉自持杆
Self Holding Rod For Locking Screw,SW3.5



92233501 T型六角扳手
T Hexagonal Wrench ,SW3.5



92233601 六角起子(快装式)
Hexagonal Screwdriver(Quick Coupling),SW4.0



92233701 锁钉扳手
Locking Screw Wrench,SW4.0/M2.5



92233801 连接螺栓/尾帽扳手
Wrench For Connecting Bolt/End Cap,SW7.0



92233901 软杆尾帽万向扳手
Universal Wrench For Soft Rod End Cap,SW7.0



92100401 快装手柄
Quick Coupling Handle



92100501 转换接头
Exchanging Joint



92234001 打入取出连接杆
Connecting Rod For Inserting And Extracting



92234101 主钉取出杆
Extracting Rod For Main Nail,M10



92100601 滑锤
Inserter-Extractor



92101401 组合扳手
Combination Wrench, SW11



92235401 清洁导针
Cleaning Guide Pin,Φ3.0mm
曾用编码 92213401



92100801 软扩刀杆
Soft Expansion Tool Bar



* * KTO 滑动加压伽玛钉安装器械包(小号) 92239901

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set(Small)



92100902 软扩刀盒
Soft Expansion Tool Box



91109801 软扩刀头
Soft Expansion Tool Bit,Φ8.5mm



91109802 软扩刀头
Soft Expansion Tool Bit,Φ9.0mm



91109803 软扩刀头
Soft Expansion Tool Bit,Φ9.5mm



91109804 软扩刀头
Soft Expansion Tool Bit,Φ10mm



91109805 软扩刀头
Soft Expansion Tool Bit,Φ10.5mm



91109806 软扩刀头
Soft Expansion Tool Bit,Φ11mm



91109807 软扩刀头
Soft Expansion Tool Bit,Φ11.5mm



91109808 软扩刀头
Soft Expansion Tool Bit,Φ12mm





92231902 子钉孔定位器
Locator For Sub Screw Hole,95mm





*** * KTO 滑动加压伽玛钉安装器械包小号(加长版) 92239902**

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)



91123201 螺纹导针
Threaded Guide Pin, $\Phi 3.2 \times 400\text{mm}$



91124001 导针
Guide Pin, $\Phi 4.0 \times 300\text{mm}$



92230101 导针测深器
Depth Gauge For Guide Pin



92230201 显影量尺
Developing Measure Gauge



92100702 开孔器
Awl



92230401 保护套筒
Protective Sleeve, $\Phi 15\text{mm}$

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between contant, model, color, please refer to the object.20220629



*** * KTO 滑动加压伽玛钉安装器械包小号(加长版) 92239902**

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)




92230501 导针套管
Guide Pin Sleeve ,Φ15mm



92230601 近端空心钻
Proximal Hollow Drill Bit,Φ15mm



92100101 导针把持器
Guide Pin Holder



92100201 球头导针
Ball Head Guide Pin,Φ2.5×800



92100301 复位棒(导针引入器)
Reduction Rod(Guide Pin Guider)



92230701 手把
Handle



92230801 近端导杆
Proximal Guide Rod



92230901 连接螺杆
Connecting Rod,SW7.0



92231001 连接螺杆扳手
Connecting Rod Wrench,SW7.0



92231101 远端导杆
Distal Guide Rod



92231201 导杆锁轮
Locking Nut For Guide Rod



92231301 近端套管
Proximal Sleeve



92231401 近端皮质套管
Proximal Cortical Sleeve



92231501 近端套管针
Proximal Trocar



92231601 支持加压螺母
Supporting Compression Nut



* * KTO 滑动加压伽玛钉安装器械包小号(加长版) 92239902

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)



92231701 子钉开口钻
Opening Drill Bit For Sub Screw, Φ7.1mm



92231801 子钉钻头
Drill Bit For Sub Screw, Φ6.0mm



92231901 子钉孔定位器
Locator For Sub Screw Hole, 65mm



92232001 母钉皮质骨空心钻
Cortical Bone Hollow Drill Bit For Parent Screw, Φ9.6mm



92232101 母钉限位空心钻
Limit Hollow Drill Bit For Parent Screw, Φ9.0mm



92232201 限位器
Limiter



92232301 母钉空心丝锥
Hollow Tap For Parent Screw, Φ9mm



92232401 母钉持杆
Holding Rod For Parent Screw, M5



92232501 母钉中空持杆
Hollow Holding Rod For Parent Screw, M5



92232601 母钉上钉器
Parent Screw Screwdriver



92232701 母钉加压螺母
Compression Nut For Parent Screw



92232801 母钉导针套管
Guide Pin Sleeve For Parent Screw, Φ3.2mm



92232901 子钉上钉器
Sub Screw Screwdriver, SW4.0/M2.5



91114001 骨钻
Bone Drill Bit, Φ4.0×300mm



91150801 限位器
Limiter, Φ4.0/SW3.5



*** * KTO 滑动加压伽玛钉安装器械包小号(加长版) 92239902**

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)



92233001 远端锁钉套管
Distal Locking Screw Sleeve, $\Phi 8.2\text{mm}$



92233101 钻套
Drill Sleeve, $\Phi 4.2\text{mm}$



92233201 套管针
Trocar, $\Phi 4\text{mm}$



92233301 测深器
Depth Gauge



92233401 锁钉自持杆
Self Holding Rod For Locking Screw, SW3.5



92233501 T型六角扳手
T Hexagonal Wrench, SW3.5



92233601 六角起子(快装式)
Hexagonal Screwdriver(Quick Coupling), SW4.0



92233701 锁钉扳手
Locking Screw Wrench, SW4.0/M2.5



92233801 连接螺栓/尾帽扳手
Wrench For Connecting Bolt/End Cap, SW7.0



92233901 软杆尾帽万向扳手
Universal Wrench For Soft Rod End Cap, SW7.0



92100401 快装手柄
Quick Coupling Handle



92100501 转换接头
Exchanging Joint



92234001 打入取出连接杆
Connecting Rod For Inserting And Extracting, M10



92234101 主钉取出杆
Extracting Rod For Main Nail, M10



92100601 滑锤
Inserter-Extractor



*** * KTO 滑动加压伽玛钉安装器械包小号(加长版) 92239902**

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)



92101401 组合扳手
Combination Wrench, SW11



92235401 清洁导针
Cleaning Guide Pin, $\Phi 3.0\text{mm}$
曾用编码 92213401



92100801 软扩刀杆
Soft Expansion Tool Bar



92100902 软扩刀盒
Soft Expansion Tool Box



91109801 软扩刀头
Soft Expansion Tool Bit, $\Phi 8.5\text{mm}$



91109802 软扩刀头
Soft Expansion Tool Bit, $\Phi 9.0\text{mm}$



91109803 软扩刀头
Soft Expansion Tool Bit, $\Phi 9.5\text{mm}$



91109804 软扩刀头
Soft Expansion Tool Bit, $\Phi 10\text{mm}$



91109805 软扩刀头
Soft Expansion Tool Bit, $\Phi 10.5\text{mm}$



91109806 软扩刀头
Soft Expansion Tool Bit, $\Phi 11\text{mm}$



91109807 软扩刀头
Soft Expansion Tool Bit, $\Phi 11.5\text{mm}$



91109808 软扩刀头
Soft Expansion Tool Bit, $\Phi 12\text{mm}$



91115201 定位杆钻头
Located Rod Drill Bit, $\Phi 5.2 \times 300\text{mm}$



92234201 定位杆平头钻
Flat Drill Bit For Located Rod, $\Phi 5.2\text{mm}$



92234301 定位杆
Located Rod



* * KTO 滑动加压伽玛钉安装器械包小号(加长版) 92239902

KTO Sliding Compression Gamma Intramedullary Nail Instrument Set (Extended Edition)



92234401 定位杆套管针
Trocar For Located Rod



92234501 定位杆套管
Located Rod Sleeve



92234601 定位杆钻套
Drill Sleeve For Located Rod,Φ5.2mm



92234701 定位卡块
Located Plug



92234801 远端连接锁轮
Distal Connecting Nut



92234901 股骨近端瞄准架
Proximal Femoral Diopter



92235001 长钉远端连接导杆
Distal Connecting Guide Rod For Long Screw



92235101 长钉远端导杆
Distal Guide Rod For Long Screw



92235201 长钉远端定位架左
Distal Located Bracket For Long Screw(Left)



92235301 长钉远端定位架右
Distal Located Bracket For Long Screw(Right)





92231902 子钉孔定位器
Locator For Sub Screw Hole,95mm



*** * 股骨倒打带锁髓内钉安装器械包 92249901**

Femoral Intramedullary Nail Instrument Set



92240101 显影量尺
Developing Measure Gauge



91123203 螺纹导针
Threaded Guide Pin, $\Phi 3.2 \times 300\text{mm}$



92240201 保护套筒
Protective Sleeve



92240301 近端空心钻
Proximal Hollow Drill Bit, $\Phi 14\text{mm}$



92100401 快装手柄
Quick Coupling Handle



92100501 转换接头
Exchanging Joint



92100710 中空开孔器
Hollow Awl



92100101 导针把持器
Guide Pin Holder



92100304 复位棒 (导针引入器)
Reduction Rod (Guide Pin Guider), $\Phi 2.5 \times 320\text{mm}$

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object. 20220629



*** * 股骨倒打带锁髓内钉安装器械包 92249901**

Femoral Intramedullary Nail Instrument Set

92100201 球头导针 Ball Head Guide Pin,Φ2.5×800	91109801 软扩刀头 Soft Expansion Tool Bit,Φ8.5mm	91109802 软扩刀头 Soft Expansion Tool Bit,Φ9.0mm
91109803 软扩刀头 Soft Expansion Tool Bit,Φ9.5mm	91109805 软扩刀头 Soft Expansion Tool Bit,Φ10.5mm	91109807 软扩刀头 Soft Expansion Tool Bit,Φ11.5mm
91109804 软扩刀头 Soft Expansion Tool Bit,Φ10mm	91109806 软扩刀头 Soft Expansion Tool Bit,Φ11mm	91109808 软扩刀头 Soft Expansion Tool Bit,Φ12mm
91109809 软扩刀头 Soft Expansion Tool Bit,Φ12.5mm	91109810 软扩刀头 Soft Expansion Tool Bit,Φ13mm	92100803 软扩刀杆 Soft Expansion Tool Bar
92240401 手把 Handle	92240501 连接螺杆 Connecting Rod,SW6.5	92240601 连接螺杆扳手 Connecting Rod Wrench,SW6.5
92240701 延长连接螺栓 Extended Connecting Bolt	92101405 开口扳手 Opening Wrench,SW11	92240801 主钉打拔杆 Insert-Extractor Rod For Main Nail , M10



* * 股骨倒打带锁髓内钉安装器械包 92249901

Femoral Intramedullary Nail Instrument Set



92100601 滑锤
Inserter Extraactor



92240901 近端导杆
Proximal Guide Rod



92241001 手把连接锁轮
Connecting Nut For Handle



91150404 锁轮扳手
Locking Nut Wrench, SW5.0



92241101 远端导杆
Distal Guide Rod



92241201 远端连接锁轮
Distal Locking Nut, M8



92241301 定位架
Located Bracket



92241401 外套管/定位杆套管
Outer Sleeve



92241501 软组织分离器
Soft Tissue Separator



92241601 远端钻套/定位杆钻套
Distal Drill Sleeve/Drill Sleeve For Located Rod, Φ5.2mm



92241701 T型平头钻
T Flat Drill Bit, Φ5.2mm



91115202 定位杆骨钻
Bone Drill For Located Rod, Φ5.2×250mm



92241801 定位杆
Located Rod



92241901 定位卡块
Located Rod Plug, Φ9mm, Φ10mm



92242001 定位卡块
Located Rod Plug, Φ11mm, Φ12mm



* * 股骨倒打带锁髓内钉安装器械包 92249901

Femoral Intramedullary Nail Instrument Set



92242101 近端钻套
Proximal Drill Sleeve, Φ4.0mm



91114001 近端骨钻
Proximal Bone Drill Bit, Φ4.0×300mm



91150801 限位器
Limiter, Φ4.0/SW3.5



92242201 直角定位杆
Right-Angled Located Rod, Φ4.0mm



92242301 测深器
Depth Gauge



92242401 T型锁钉扳手
T Locking Screw Wrench, SW3.5



92242501 锁钉扳手
Locking Screw Wrench, SW3.5



91150803 限位器
Depth Gauge, Φ5.2/SW3.5



92242601 远端锁定丝锥
Distal Locking Tap



92242701 尾帽球头扳手
Ball Head Wrench For End Cap, SW3.5



92242801 尾帽持杆
Holding Rod For End Cap, SW3.5



92242901 主钉拔出连接杆
Main Nail Extractor Connecting Rod For Main Nail



91124001 导针
Guide Pin, Φ4.0×300mm



92100901 软扩刀盒
Soft Expansion Tool Box



91115201 远端骨钻
Bone Drill For Located Rod, Φ5.2×300mm



* * KTO 股骨带锁髓内钉安装器械包 92259901

KTO Femoral Intramedullary Nail Instrument Set



92250101 锁钉钻套
Drill Sleeve For Locking Screw, Φ4.0mm



92250201 导针钻套
Drill Sleeve For Guide Pin, Φ2.5mm



92250301 锁钉套管
Locking Screw Sleeve



92250401 软组织分离器(长)
Soft Tissue Separator (Long)



92250501 定位杆钻套
Drill Sleeve For Located Rod, Φ5.2mm




92250601 定位杆套管
Located Rod Sleeve



92250701 软组织分离器(短)
Soft Tissue Separator, (Short)



91122502 螺纹导针
Threaded Guide Pin, Φ2.5×320mm



91122503 导针
Guide Pin, Φ2.5×320mm

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between contant, model, color, please refer to the object.20220629



* * KTO 股骨带锁髓内钉安装器械包 92259901

KTO Femoral Intramedullary Nail Instrument Set



91123201 螺纹导针
Threaded Guide Pin, Φ3.2×400mm



92250801 导针测深器
Depth Gauge Of Guide Pin, Φ2.5×320mm



92250901 测深器
Depth Gauge



91114001 远端骨钻
Distal Bone Drill, Φ4.0×300mm



91150801 限位器
Limiter, Φ4.0/SW3.5



91124001 导针
Guide Pin, Φ4.0×300mm



92251101 空心台阶钻
Hollow Step Drill Bit, Φ6.4mm



92251201 限位器
Limiter



91115201 定位杆骨钻
Bone Drill For Located Rod, Φ5.2×300mm



92251301 定位杆平头钻
Flat Drill Bit For Located Rod, Φ5.2mm



92223801 T型六角起子
T Hexagonal Screwdriver, SW3.5



92251501 T型股骨颈钉/尾帽扳手
T Femoral Neck Screw / End Cap Wrench, SW4.0/M2.5



92251601 连接螺杆扳手
Connecting Rod Wrench, SW6.5



92251701 六角起子(快装式)
Hexagonal Screwdriver (Quick Coupling), SW3.5



92251801 股骨颈钉/尾帽扳手(快装式)
Femoris Neck Screw / End Cap Wrench, SW4.0



* * KTO 股骨带锁髓内钉安装器械包 92259901

KTO Femoral Intramedullary Nail Instrument Set



92251901 T型连接螺杆万向扳手

T Universal Wrench For Connecting Rod,SW6.5



92252001 软杆尾帽万向扳手

Universal Wrench For Soft Rod End Cap,SW4.0



92252101 扳手

Wrench,Φ4.0mm



92252201 近端导杆

Proximal Guide Rod



92252301 远端导杆

Distal Guide Rod



92252401 股骨近端钉远端导杆

Distal Guide Rod For Proximal Femoral Screw



92252501 远端定位架

Distal Located Bracket



92252601 导杆锁轮

Locking Nut For Guide Rod



92252602 导杆连接锁轮

Locking Nut For Connecting Guide Rod



92252701 定位卡块

Located Plug



92252801 手把

Handle



92252901 主钉拔出连接杆

Connecting Rod For Extracting Nail,M8



92253001 连接螺杆

Connecting Rod



92253101 T型定位杆

T Located Rod



92253201 显影量尺

Developing Measure Gauge



* * KTO 股骨带锁髓内钉安装器械包 92259901

KTO Femoral Intramedullary Nail Instrument Set



92253301 直角定位杆
Right-Angled Located Rod



92253401 近端空心钻
Proximal Hollow Drill, Φ13.8mm



92100801 软扩刀杆
Soft Expansion Tool Bar



92100901 软扩刀盒
Soft Expansion Tool Box



91109801 软扩刀头
Soft Expansion Tool Bit, Φ8.5mm



91109802 软扩刀头
Soft Expansion Tool Bit, Φ9.0mm



91109803 软扩刀头
Soft Expansion Tool Bit, Φ9.5mm



91109804 软扩刀头
Soft Expansion Tool Bit, Φ10mm



91109805 软扩刀头
Soft Expansion Tool Bit, Φ10.5mm



91109806 软扩刀头
Soft Expansion Tool Bit, Φ11mm



91109807 软扩刀头
Soft Expansion Tool Bit, Φ11.5mm



91109808 软扩刀头
Soft Expansion Tool Bit, Φ12mm



91109809 软扩刀头
Soft Expansion Tool Bit, Φ12.5mm



91109810 软扩刀头
Soft Expansion Tool Bit, Φ13mm



92100702 开孔器
Awl




* * KTO 股骨带锁髓内钉安装器械包 92259901

KTO Femoral Intramedullary Nail Instrument Set



92100301 复位棒(导针引入器)
Reduction Rod(Guide Pin Guider)



92100201 球头导针
Ball Head Guide Pin, $\Phi 2.5 \times 800$



92100101 导针把持器
Guide Pin Holder



92253501 打拔器螺帽
Insert-Extractor Nut



92253601 打拔锤
Insert-Extractor Hammer



92253701 打拔器杆
Insert-Extractor Rod



92101001 保护套筒
Protective Sleeve



92101101 蜂窝式导针套筒
Honeycomb Guide Pin Sleeve, $\Phi 3.2\text{mm}$



92101201 套筒手柄
Sleeve Handle



92100401 快装手柄
Quick Coupling Handle






92100501 转换接头
Exchanging Joint





*** * 股骨带锁髓内钉安装器械包 92269901**

Femoral Intramedullary Nail Instrument Set



92260101 髓腔铰刀
Medullary Canal Reamer, Φ9mm



92260201 髓腔铰刀
Medullary Canal Reamer, Φ10.4mm



92320401 髓腔铰刀
Medullary Canal Reamer, Φ11mm



92260401 髓腔铰刀
Medullary Canal Reamer, Φ12mm



92260501 髓腔铰刀
Medullary Canal Reamer, Φ12.5mm



92320501 钻套
Drill Sleeve, Φ4.0mm



92320601 钻套
Drill Sleeve, Φ5.2mm



92320701 外套管
Outer Sleeve



92320801 软组织分离器
Soft Tissue Separator

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object.20220629



* * 股骨带锁髓内钉安装器械包 92269901

Femoral Intramedullary Nail Instrument Set



91114001 钻头
Drill Bit, $\Phi 4.0 \times 300\text{mm}$



91150818 限位器
Limiter, $\Phi 4.0/\text{SW}3.0$



91115201 钻头
Drill Bit, $\Phi 5.2 \times 300\text{mm}$



92321001 T型平头钻
T Flat Drill Bit, $\Phi 5.2 \times 220\text{mm}$



92322401 T型六角起子
T Hexagonal Wrench, SW3.5



91133504 封帽起子
Sealing Cap Screwdriver, SW3.5



91150908 持钉套
Screw Holding Sleeve



92261501 测深器
Depth Gauge



92261601 定位架锁轮
Locking Nut For Located Bracket



92261701 导杆锁轮
Locking Nut For Guide Rod



92261801 U型卡块
U Plug



91150404 锁轮扳手
Locking Nut Wrench, SW5.0



91150403 限位器扳手
Limiter Wrench, SW3.0



92261901 瞄准器手把
Diopter Handle



92262001 连接螺杆
Connecting Rod



* * 股骨带锁髓内钉安装器械包 92269901

Femoral Intramedullary Nail Instrument Set



92262101 远端定位架
Distal Located Bracket



92262201 导杆
Guide Rod



92262301 定位杆
Located Rod



92262401 打入器手柄
Impactor Handle



92262501 打拔器
Insert-Extractor



92262601 万向接头
Universal Joint, M8




92101402 开口扳手
Opening Wrench, SW12



92100706 开孔器
Awl



92322201 直角定位杆
Right-Angled Located Rod



92100202 导针
Guide Pin, Φ2.5×660mm




92100402 快装手柄
Quick Coupling Handle




91124001 导针
Guide Pin, Φ4.0×300mm





* * QML 多维胫骨髌上入路髓内钉安装器械包 92339901

QML Multi-Dimensional Tibial Suprapatellar approach Intramedullary Nail Instrument Set




92330101 显影量尺
Developing Measure Gauge

91123203 螺纹导针
Threaded Guide Pin, $\Phi 3.2 \times 300\text{mm}$

91123202 螺纹导针
Threaded Guide Pin, $\Phi 3.2 \times 400\text{mm}$



92101201 套筒手柄
Sleeve Handle



92101002 保护套筒
Protective Sleeve, $\Phi 12.8 \times 125\text{mm}$



92101102 蜂窝式导针套筒
Honeycomb Guide Pin Sleeve, $\Phi 3.2 \times 150\text{mm}$



92100501 转换接头
Exchanging Joint




92100401 快换手柄
Quick Coupling Handle



92230201 近端扩孔钻
Proximal Expanding Drill, $\Phi 12.8\text{mm}$

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object. 20220629




*** * QML 多维胫骨髌上入路髓内钉安装器械包 92339901**

QML Multi-Dimensional Tibial Suprapatellar approach Intramedullary Nail Instrument Set




92230301 软管近端空心钻
Proximal Hollow Drill Bit for soft rod



92330401 开孔器内芯
Inner core for awl, $\Phi 3.2 \times 410\text{mm}$



92330501 中空开孔器
Hollow awl, $\Phi 9.5\text{mm}$



92100903 软扩刀盒
Soft Expansion Tool Box



91109811 软扩刀头
Soft Expansion Tool Bit, $\Phi 7.5\text{mm}$



91109812 软扩刀头
Soft Expansion Tool Bit, $\Phi 8.0\text{mm}$



91109801 软扩刀头
Soft Expansion Tool Bit, $\Phi 8.5\text{mm}$



91109802 软扩刀头
Soft Expansion Tool Bit, $\Phi 9.0\text{mm}$



91109803 软扩刀头
Soft Expansion Tool Bit, $\Phi 9.5\text{mm}$



91109804 软扩刀头
Soft Expansion Tool Bit, $\Phi 10\text{mm}$



91109805 软扩刀头
Soft Expansion Tool Bit, $\Phi 10.5\text{mm}$



91109806 软扩刀头
Soft Expansion Tool Bit, $\Phi 11\text{mm}$



91109807 软扩刀头
Soft Expansion Tool Bit, $\Phi 11.5\text{mm}$



91109808 软扩刀头
Soft Expansion Tool Bit, $\Phi 12\text{mm}$



92100802 软扩刀杆
Soft Expansion Tool Bar



* * QML 多维胫骨髌上入路髓内钉安装器械包 92339901


QML Multi-Dimensional Tibial Suprapatellar approach Intramedullary Nail Instrument Set



92334101 复位棒
Reduction Rod



92100101 导针把持器
Guide Pin Holder



92100201 球头导针
Ball Head Guide Pin, $\Phi 2.5 \times 800\text{mm}$



92330601 球头导针测量器
Ball head guide pin gauge



92330701 手把
Handle



92330801 连接螺栓 (长)
Connecting bolt (Long)



92330901 手把 (短)
Handle (Short)



92331001 连接螺栓 (短)
Connecting bolt (Short)



92331101 连接螺杆扳手
Connecting Rod Wrench, SW6.5



92331201 主钉打拔连接杆
Inserter-Extractor connector for main nail, M8mm



92101401 开口扳手
Opening Wrench, SW11



92331301 主钉取出杆
Inserter-Extractor connector for main nail, M8mm




92100601 滑锤
Slide Hammer



92331401 导杆
Guide rod



92331501 导杆锁轮
Locking Nut For Guide Rod



* * QML 多维胫骨髌上入路髓内钉安装器械包 92339901

QML Multi-Dimensional Tibial Suprapatellar approach Intramedullary Nail Instrument Set



92331601 远端瞄准架
Distal diopther



92331701 远端瞄准架锁轮
Locking Nut For Distal diopther



92331801 超远端瞄准架
Super Distal diopther



92331901 超远端瞄准架锁轮
Locking Nut For Super Distal diopther



92332001 锁钉套管
Locking Screw Sleeve



92332101 软组织分离器
Soft Tissue Separator



92332201 钻套
Drill Sleeve, Φ5.2mm



91115202 定位杆骨钻
Bone Drill For Located Rod, Φ5.2×250mm



92332301 定位杆平头钻
Flat Drill For Located Rod, Φ5.2mm



92332401 定位杆
Located Rod



92332501 前后位定位杆
Anteroposterior position located rod



92332601 定位卡块
Located Plug




92332701 钻套
Drill Sleeve, Φ4.0mm



91124001 导针
Guide Pin, Φ4.0×300mm



91114001 远端骨钻
Distal Bone Drill, Φ4.0×300mm



* * QML 多维胫骨髌上入路髓内钉安装器械包 92339901

QML Multi-Dimensional Tibial Suprapatellar approach Intramedullary Nail Instrument Set



91150818 限位器
Limiter, $\Phi 4.0\text{mm}$ /SW3.0



91150403 L型六角起子
L hexagon Screwdriver, SW3.0mm



92232801 临时定位杆
Temporary Located Rod, $\Phi 4.0\text{mm}$



92332901 测深器
Depth gauge



92333001 锁钉扳手
Tr, SW4.0mm



92333101 快装六角起子
Hexagon screwdriver, SW4.0mm



92333201 近端瞄准架 (铝)
Proximal diopter (aluminium)



92333301 近端导杆锁轮
Locking Nut For Proximal guide rod



92333401 锁钉套管
Locking Screw Sleeve



92333501 导针钻套
Trocar, $\Phi 3.2\text{mm}$



92333601 套管针
Trocar, $\Phi 3.2\text{mm}$



91113204 骨钻
Bone Drill, $\Phi 3.2 \times 250\text{mm}$



91150820 限位器
Limiter, $\Phi 3.2\text{mm}$ /SW3.0



92333701 加压螺栓 (长)
Compression bolt (Long)



92333801 加压螺栓 (短)
Compression bolt (Short)

* * QML 多维胫骨髌上入路髓内钉安装器械包 92339901

QML Multi-Dimensional Tibial Suprapatellar approach Intramedullary Nail Instrument Set



92333901 连接螺杆万向扳手

Universal Wrench For Connecting Rod, SW6.5




92334001 尾帽六角持杆

Hexagonal Holding Rod For End Cap, SW3.5




92100704 中空开孔器

Hollow Awl



* * KTO 胫骨带锁髓内钉安装器械包 92319901

KTO Tibial Intramedullary Nail Instrument Set



92310101 显影量尺
Developing Measure Gauge

91123203 螺纹导针
Threaded Guide Pin, $\Phi 3.2 \times 300\text{mm}$

92101001 保护套筒
Protective Sleeve, $\Phi 17\text{mm}$



92101101 蜂窝式导针套筒
Honeycomb Guide Pin Sleeve, $\Phi 3.2\text{mm}$



92101201 套筒手柄
Sleeve Handle



92310201 近端扩孔钻
Proximal Expanding Drill, $\Phi 12\text{mm}$



92100401 快换手柄
Quick Coupling Handle



92100501 转换接头
Exchanging Joint



92100704 中空开孔器
Hollow Awl

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object. 20220629



*** * KTO 胫骨带锁髓内钉安装器械包 92319901**


KTO Tibial Intramedullary Nail Instrument Set



92100301 复位棒(导针引入器)
Reduction Rod(Guide Pin Guider)



92100101 导针把持器
Guide Pin Holder



92100201 球头导针
Ball Head Guide Pin,Φ2.5×800mm



92100903 软扩刀盒
Soft Expansion Tool Box



91109811 软扩刀头
Soft Expansion Tool Bit,Φ7.5mm



91109812 软扩刀头
Soft Expansion Tool Bit,Φ8.0mm



91109801 软扩刀头
Soft Expansion Tool Bit,Φ8.5mm



91109802 软扩刀头
Soft Expansion Tool Bit,Φ9.0mm



91109803 软扩刀头
Soft Expansion Tool Bit,Φ9.5mm



91109804 软扩刀头
Soft Expansion Tool Bit,Φ10mm



91109805 软扩刀头
Soft Expansion Tool Bit,Φ10.5mm



91109806 软扩刀头
Soft Expansion Tool Bit,Φ11mm



91109807 软扩刀头
Soft Expansion Tool Bit,Φ11.5mm



91109808 软扩刀头
Soft Expansion Tool Bit,Φ12mm



92100802 软扩刀杆
Soft Expansion Tool Bar



* * KTO 胫骨带锁髓内钉安装器械包 92319901

KTO Tibial Intramedullary Nail Instrument Set



92310301 手把
Handle



92310401 连接螺栓
Connecting Bolt, SW6.5



92310501 连接螺杆扳手
Connecting Rod Wrench, SW6.5



92310601 打入器连接块
Impactor Connecting Block



92310701 主钉打拔连接器
Connector For Inserting And Extrating Nail



92101405 开口扳手
Opening Wrench, SW11



92310901 扳手
Wrench, Φ4.0mm



92311001 滑锤导向杆
Slide Hammer Guide Rod, M8



92311101 滑锤导向杆螺帽
Slide Hammer Guide Rod Nut, M10



92100601 滑锤
Slide Hammer



92311201 导杆
Guide Rod



92311301 连接导杆锁轮
Locking Nut For Connecting Guide Rod, M8



92311401 定位杆钻套
Drill Sleeve For Located Rod, Φ5.2mm



91115201 定位杆骨钻
Bone Drill For Located Rod, Φ5.2×300mm



92311501 定位杆平头钻
Flat Drill For Located Rod, Φ5.2mm



* * KTO 胫骨带锁髓内钉安装器械包 92319901

KTO Tibial Intramedullary Nail Instrument Set



92311601 定位杆
Located Rod



92311701 定位杆卡块
Located Rod Plug



92311801 远端瞄准架
Distal Diopter



92311901 远端瞄准架锁轮
Locking Nut For Distal Diopter



92312001 外套管
Outer Sleeve



92312101 软组织分离器
Soft Tissue Separator



92312201 钻套
Drill Sleeve, Φ4.0mm



91124001 导针
Guide Pin, Φ4.0×300mm



91114001 远端骨钻
Distal Bone Drill, Φ4.0×300mm



91150801 限位器
Limiter, Φ4.0/SW3.5



92312401 临时定位杆
Temporary Located Rod, Φ4.0mm



92312501 测深器
Depth Gauge




92312601 T型锁钉扳手
T Locking Screw Wrench, SW3.5



92312701 锁钉扳手
Locking Screw Wrench, SW3.5



92312801 近端瞄准架
Proximal Diopter



* * KTO 胫骨带锁髓内钉安装器械包 92319901

KTO Tibial Intramedullary Nail Instrument Set



92312901 近端松质骨瞄准架
Proximal Cancellous Bone Diopter



92313001 近端瞄准架锁轮
Locking Nut For Proximal Diopter, M6



92313101 锁钉套管
Locking Screw Sleeve, Φ11



92313201 钻套
Drill Sleeve, Φ3.2mm



92313301 套管针
Trocar, Φ3.2mm



91113204 骨钻
Bone Drill, Φ3.2×250mm



91150817 限位器
Limiter, Φ3.2mm/SW3.5



92313501 加压螺栓
Compression Bolt



92313601 连接螺杆万向扳手
Universal Wrench For Connecting Rod, SW6.5



92313701 尾帽球头扳手
Ball Head Wrench For End Cap, SW3.5



92313801 尾帽六角持杆
Hexagonal Holding Rod For End Cap, SW3.5





92314001 近端空心钻保护套筒
Protective Sleeve For Proximal Hollow Drill Bit



* * 胫骨带锁髓内钉安装器械包 92329901

Tibial Intramedullary Nail Instrument Set



92320101 硬扩
Hard Expansion, $\Phi 8\text{mm}$



92320201 硬扩
Hard Expansion, $\Phi 9.4\text{mm}$



92320301 硬扩
Hard Expansion, $\Phi 10\text{mm}$



92320401 硬扩
Hard Expansion, $\Phi 11\text{mm}$



92320501 钻套
Drill Sleeve, $\Phi 4.0\text{mm}$



92320601 钻套
Drill Sleeve, $\Phi 5.2\text{mm}$



92320701 套管
Sleeve



92320801 软组织分离器
Soft Tissue Separator



91114001 钻头
Drill Bit, $\Phi 4.0 \times 300\text{mm}$

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object.20220629



* * 胫骨带锁髓内钉安装器械包 92329901

Tibial Intramedullary Nail Instrument Set



91150818 限位器
Limiter For Drill Bit, $\Phi 4.0$ /SW3.0



91115201 钻头
Drill Bit, $\Phi 5.2 \times 300$ mm



92321001 T型平头钻
T Flat Drill Bit, $\Phi 5.2$ mm



92321101 T型定位杆
T Located Rod



92321201 定位架锁轮
Locking Nut For Located Bracket



92321301 导杆锁轮
Locking Nut For Guide Rod



92321401 U型卡块
U Plug



92321501 连接螺杆
Connecting Rod



92321601 瞄准器手把
Diopter Handle



92321701 近端瞄准架
Proximal Diopter



92321801 远端瞄准架
Distal Diopter



92321901 导杆
Guide Rod



91150403 限位器扳手
Limiter Wrench, SW3.0



91150404 锁轮扳手
Locking Nut Wrench, SW5.0



92100707 开孔器
Awl



*** * 胫骨带锁髓内钉安装器械包 92329901**

Tibial Intramedullary Nail Instrument Set



92100402 快装手柄
Quick Coupling Handle



92322201 直角定位杆
Right-Angled Located Rod



92322301 测深器
Depth Gauge



92322401 T型六角起子
T Hexagonal Screwdriver,SW3.5



91133504 六角起子
Hexagonal Screwdriver,SW3.5



91150908 持钉套
Screw Holding Sleeve



92322701 打拔器
Inserter-Extractor



92322801 万向接头
Universal Joint,M6



92101402 开口扳手
Opening Wrench,SW12







91124001 导针
Guide Pin,Φ4.0×300mm




* * QML 多维肱骨带锁髓内钉安装器械包 92529901

QML Multi-Dimensional Humeral Intramedullary Nail Instrument Set



92100705 中空开孔器
Hollow Awl



92100207 球头导针
Ball Head Guide Pin, $\Phi 2.5 \times 680\text{mm}$



92100501 转换接头
Exchanging Joint



92100401 T型快装手柄
T Quick Coupling Handle



92520101 手把
Handle



92520201 连接螺栓
Connecting Bolt

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object.20220629



* * QML 多维肱骨带锁髓内钉安装器械包 92529901

QML Multi-Dmensional Humeral Intramedullary Nail Instrument Set



92520301 加压螺杆
Compression Bolt



92100601 滑锤
Slide Hammer



92520401 软扩
Soft Expansion, $\Phi 7\text{mm}$



92520501 软扩
Soft Expansion, $\Phi 8\text{mm}$



92520601 软扩
Soft Expansion, $\Phi 9\text{mm}$



92520701 软扩
Soft Expansion, $\Phi 10\text{mm}$



92101404 开口扳手
Opening Wrench, SW11



91121502 导针
Guide Pin, $\Phi 1.5 \times 120\text{mm}$



91122504 导针
Guide Pin, $\Phi 2.5 \times 250\text{mm}$



92520801 止深导针
H, $\Phi 2.5 \times 200\text{mm}$



92520901 测深器
Depth Gauge



92521001 测深器
Depth Gauge



91113801 骨钻
Bone Drill, $\Phi 3.8 \times 270\text{mm}$



92521101 指示块
Indication Block, $\Phi 3.8\text{mm}$



91112804 骨钻
Bone Drill, $\Phi 2.8 \times 200\text{mm}$



* * QML 多维肱骨带锁髓内钉安装器械包 92529901

QML Multi-Dimensional Humeral Intramedullary Nail Instrument Set



92521201 指示块
Indication Block, $\Phi 2.8\text{mm}$



91113004 骨钻
Bone Drill, $\Phi 3.0 \times 300\text{mm}$



91150815 限位器
Limiter, $\Phi 3.0/\text{SW}3.0$



92100302 复位棒
Reduction Rod



92522701 导针把持器
Guide Pin Holder



92522801 近端导杆
Proximal Screw Guide Rod



92522901 近端导杆锁轮
Locking Nut For Proximal guide rod



92523001 多维锁钉套管
Multi-Dimensional Locking Sleeve



92523101 多维锁钉钻套
Multi-Dimensional Locking Drill Sleeve, $\Phi 3.8\text{mm}$



92523201 套管针
Distal Located Bracket, $\Phi 3.8\text{mm}$



92523301 锁钉套管
Locking Sleeve



92523401 锁钉钻套
Locking Drill Sleeve, $\Phi 3.0\text{mm}$



92523501 套管针
Distal Located Bracket, $\Phi 3.0\text{mm}$



92523601 环钻
Trephine, $\Phi 10\text{mm}$



92523602 环钻
Trephine, $\Phi 11.5\text{mm}$



* * QML 多维肱骨带锁髓内钉安装器械包 92529901

QML Multi-Dmensional Humeral Intramedullary Nail Instrument Set



92523701 球头导针测量器
Ball head guide pin gauge



92523801 保护套筒
Protective Sleeve, $\Phi 11.5 \times 70\text{mm}$



92523901 梅花起子杆
Trox Screwdriver, T3.4mm



91151104 1.5N.m扭力起子
Torque Screwdriver, T3.4mm



92521401 梅花起子 (短)
Torx Screwdriver (Short), T3.4mm



92521501 梅花起子 (长)
Torx Screwdriver (Long), T3.4mm



92521601 近端前位导杆
Proximal Anteposition Guide Rod



92521701 远端导杆
Distal Guide Rod



92521801 远端导杆锁轮
Distal Guide Rod Locking Screw



92521901 定位杆平头钻
Flat Drill Bit For Located Rod, $\Phi 3.8 \times 220\text{mm}$



92522001 定位卡快



92522101 定位杆
Located Rod



92522201 远端定位架 (右)
Distal Located Bracket (Right)



92522202 远端定位架 (左)
Distal Located Bracket (Left)



92522301 连接锁轮
Connecting Locking Nut



* * QML 多维肱骨带锁髓内钉安装器械包 92529901

QML Multi-Dimensional Humeral Intramedullary Nail Instrument Set



92522401 临时定位杆
Temporay Located Rod



92522501 尾帽测量器
End Cap Measurer



92522601 梅花起子持杆
Torx Screwdriver Holding Rod, T3.4mm



92524001 取出杆
Extraction Rod



92524101 埋头钻
Countersink Drill Bit



92524201 显影量尺
Deveiping Measure Gauge



92524301 多维锁钉扳手
Multi-Dimensional Locking Screw Wrench



92524401 锁定钉钻套
Drill Sleeve For Locking Screw, $\Phi 2.8\text{mm}$



92524501 锁定钉定位套
Located Sleeve For Lccking Screw



91150403 限位器扳手
Limiter Wrench, SW3.0



91150404 锁轮扳手
Locking Nut Wrench, SW5.0






92524601 滑锤导向杆
Guide Rod For Slide Hammer, M6





KTO 肱骨带锁髓内钉安装器械包 92519901

KTO Humeral Intramedullary Nail Instrument Set



92510201 导针把持器
Guide Pin Holder



92510301 显影量尺
Developing Measure Gauge



92100705 中空开孔器
Hollow Awl



92510401 保护套筒
Protective Sleeve, $\Phi 9.5$



92510501 近端空心钻
Proximal Hollow Drill, $\Phi 9.5$



92100403 T型快装手柄
T Quick Coupling Handle



92100502 转换接头
Exchanging Joint



91122506 螺纹导针
Threaded Guide Pin, $\Phi 2.5 \times 200$



92100206 球头导针
Ball Head Guide Pin, $\Phi 2.0 \times 660\text{mm}$

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object. 20220629



KTO 肱骨带锁髓内钉安装器械包 92519901

KTO Humeral Intramedullary Nail Instrument Set



92510801 埋头锥
Countersink Tap



92100302 复位棒
Reduction Rod



92510101 钉盒
Screw Box



91114002 骨钻
Bone Drill Bit, Φ4.0×120mm



92511001 软扩
Soft Expansion, Φ7.0mm



92511201 软扩
Soft Expansion, Φ8.5mm



92510701 埋头钻
Countersink Drill Bit



92511101 软扩
Soft Expansion, Φ8.0mm



92511301 软扩
Soft Expansion, Φ9.0mm



92511401 手把
Handle



92511501 导杆
Guide Rod



92511601 近端钉导杆
Proximal Screw Guide Rod



92511701 远端定位架A
Distal Located Bracket



92511801 远端定位架B
Distal Located Bracket



92511901 连接锁轮
Connecting Locking Nut, M6/SW5



92512001 定位卡快
Located Plug



92512101 连接螺栓
Connecting Bolt



92512201 加压螺栓
Compression Bolt



KTO 肱骨带锁髓内钉安装器械包 92519901

KTO Humeral Intramedullary Nail Instrument Set



92512301 辅助定位杆
Auxiliary Located Rod



92101403 组合扳手
Combination Wrench, SW10



92512501 临时定位杆
Temporary Located Rod, $\Phi 3.0\text{mm}$



91150403 限位器扳手
Limiter Wrench, SW3.0



91150404 锁轮扳手
Locking Nut Wrench, SW5.0



92512601 定位杆
Located Rod, $\Phi 3.6/\Phi 7.1\text{mm}$



92512701 外套管
Outer Sleeve, $\Phi 7.1/\Phi 9\text{mm}$



92512801 套管针
Trocar, $\Phi 7.1\text{mm}$



92512901 定位杆钻套
Drill Sleeve For Located Rod, $\Phi 3.6\text{mm}$



91113601 定位杆骨钻
Bone Drill For Located Rod, $\Phi 3.6 \times 180\text{mm}$



92513101 平头钻
Flat Drill Bit, $\Phi 3.6\text{mm}$



92513201 梅花起子持杆
Torx Screwdriver Holding Rod, T3.4



92513301 梅花起子
Torx Screwdriver, T3.4



92513401 梅花起子 (快装式)
Torx Screwdriver (Quick Coupling), T3.4



91122003 导针
Guide Pin, $\Phi 2.0 \times 200\text{mm}$



KTO 肱骨带锁髓内钉安装器械包 92519901

KTO Humeral Intramedullary Nail Instrument Set



91122004 螺纹导针
Threaded Guide Pin, Φ2.0×200mm



92513501 导针测深器
Depth Gauge For Guide Pin, Φ2.0mm



92513601 空心钻
Hollow Drill Bit



92513701 叶片插入器
Blade Insert



92513801 叶片连接杆
Blade Connecting Rod, Φ2.0/M5



92513901 锁钉套管
Locking Sleeve, Φ7.1/Φ13



92514001 导针钻套
Drill Sleeve For Guide Pin, Φ2.0/Φ7.1



92514101 套管针
Trocars, Φ2.0mm



92100602 滑锤
Slide Hammer



91123003 导针
Guide Pin, Φ3.0×220



91113003 骨钻
Bone Drill, Φ3.0×220mm



91150815 限位器
Limiter, Φ3.0/SW3.0



92514401 钻套
Drill Sleeve, Φ3.0/Φ7.1mm



92514501 测深器
Depth Gauge




92514601 滑锤导向杆
Slide Hammer For Guide Rod



* * 金属弹性髓内针安装器械包 92619901

Metal Elastic Intramedullary Pins Instrument Set



91112702 钻头
Drill Bit, $\Phi 2.7 \times 115\text{mm}$



91113202 钻头
Drill Bit, $\Phi 3.2 \times 150\text{mm}$
曾用编码 91113203



91114502 钻头
Drill Bit, $\Phi 4.5 \times 150\text{mm}$



92610101 螺纹杆
Threaded Rod



92610201 三角扳手(小)
Triangle Wrench(Small)



92610202 三角扳手(大)
Triangle Wrench(Large)



92610301 直型打击器
Straight Impactor



92610302 斜型打击器
Oblique Impactor



92610401 取出钳
Removal Clamp

Our company have final right to modify all the specification in this catalogue, If there is any inconformity between content, model, color, please refer to the object. 20220629



*** * 金属弹性髓内针安装器械包 92619901**

Metal Elastic Intramedullary Pins Instrument Set



92610501 断钉器
Screw Break Off Device



92100601 滑锤
Slide Hammer



92610701 开孔锥
Opening Awl



92610801 插入器
Insert Device



92610901 双头导钻
Double-End Guide Drill, $\Phi 4.5/\Phi 3.2$




92611001 快装手柄
Quick Coupling Handle





*** * 尺桡骨髓内钉安装器械包系统 92629901**

Instrument kit for installing radioulnar internal bone marrow nail Instrument Set



91121601 导针 (镍钛合金)
Guide Pin, $\Phi 1.6 \times 180\text{mm}$



91121602 导针(不锈钢)
Guide Pin, $\Phi 1.6 \times 180\text{mm}$



91121603 螺纹导针
Threaded Guide Pin, $\Phi 1.6 \times 180\text{mm}$



92620101 中空开孔器
Hollow Awl



92620201 骨科复位钳
Orthopedic reduction forceps



92620301 保护套筒
Protective sleeve, $\Phi 6.5 \times 45\text{mm}$



92620401 导针套筒
Needle sleeve, $\Phi 6.5/1.6 \times 60\text{mm}$



92620501 近端空心钻
Proximal hollow drill, $\Phi 6.3/1.6\text{mm}$



91150101 T型快装手柄
Ttype quick mount handle



*** * 尺桡骨髓内钉安装器械包系统 92629901**

Instrument kit for installing radioulnar internal bone marrow nail Instrument Set



92620601 导针把持器
Guide pin holder, $\Phi 1.5/2.5\text{mm}$



91121604 导针(镍钛合金)
Guide pin, $\Phi 1.6 \times 350\text{mm}$



91121605 导针(不锈钢)
Guide pin, $\Phi 1.6 \times 350\text{mm}$



92620701 导针测深器
Guide needle sounder, $\Phi 1.6 \times 350\text{mm}$



92620801 扩髓钻
Expanding medullary, $\Phi 3.2\text{mm}$



92620901 扩髓钻
Expanding medullary, $\Phi 3.8\text{mm}$



92621001 复位杆
Release link



92621101 复位套筒
Reset the sleeve, $\Phi 1.6\text{mm}$



92621201 手把
Handlebar



92621301 连接螺栓
Connecting bolt, M4



92621401 小扳手
The little wrench, $\Phi 4.0\text{mm}$



92621501 尺骨钉导杆
Ulna nail guide rod



92621601 桡骨钉导杆
Radial screw guide rod



92621701 导杆锁轮
Guide wheel lock, M5



92621801 显影杆连接块
Developing rod connection block



* * 尺桡骨髓内钉安装器械包系统 92629901

Instrument kit for installing radioulnar internal bone marrow nail Instrument Set



92621901 显影杆
Developing pole



92622001 锁钉套管
Locking sleeve, $\Phi 8.0/5.2 \times 75\text{mm}$



92622101 锁钉钻套
Locking drill sets, $\Phi 5.2/2.5 \times 88\text{mm}$



92622201 套管针
Trocars, $\Phi 2.5\text{mm}$



92622301 骨钻
Bone drill, $\Phi 2.5\text{mm}$



91150810 限位器
Limiter, $\Phi 2.5/\text{SW}2.5$



91150402 限位器扳手
Limiter wrench, SW2.5



92622401 测深器 (直型)
Depth finder (Straight)



92622501 测深器 (钩型)
Depth finder (Hook type)



92622601 锁钉丝锥
Locking screw tap, $\Phi 3.0\text{mm}$



92622701 斜孔钻套
Incliner drill sets, $\Phi 2.5\text{mm}$



92623101 梅花起子
The plum blossom screwdriver, T3.4



92623201 梅花起子杆
Club driver, T3.4




92622801 主钉取出杆
Main nail extraction lever, M4/ $\Phi 6\text{mm}$



92100603 滑锤
Slide Hammer

外固定支架系列





* * 外固定支架器械包（上肢） 91759901

External Fixation Instrument Set (Upper Limb)



91150402 L型扳手
L Wrench, SW2.5




91750101 T型扳手
T Wrench



91750201 定位器
Locator, φ5x4x3



91750301 弓形钻架
Arch Drill Frame, φ4



91750401 拇指轮
Speedy Screwdriver



91750501 上钉器
Screwdriver, φ4



91113002 钻头 (D型)
L, φ3.0x150




91750601 稳定复位扳手
Reduction Wrench For Stabilization



* * 外固定支架器械包（下肢） 91769901

External Fixation Instrument Set (Lower Limb)



91760101 扭力扳手
Torque Wrench




91760201 T型扳手
T Wrench



91760301 定位器
Locator, $\phi 6 \times 5 \times 4$



91760401 弓形钻架
Arch Drill Frame, $\phi 5$



91760501 拇指轮
Speedy Screwdriver



91760601 上钉器
Screwdriver, $\phi 5$




91114001 钻头（D型）
L, $\phi 4.0 \times 300$



91760701 稳定复位扳手
Reduction Wrench For Stabilization





肋骨钢板系列





肋骨安装器械包 91719901


Rib Plate Instrument Set




91710101 单边施夹钳
Single Edge Holding Clamp



91710201 弯形施夹钳
Curved Holding Clamp



91710301 双边施夹钳
Double Edges Holding Clamp



91710401 枪式施夹钳
Gun Type Holding Clamp



91710501 骨膜剥离子 (小号)
Periosteal Detacher(Small)



91710601 骨膜剥离子
Periosteal Detacher



91710701 剪板钳
Plate Cutting Clamp






91710801 测量器 (选配)
Meter



* * 肋骨安装器械包 91719902

Rib Plate Instrument Set




91710101 单边施夹钳
Single Edge Holding Clamp



91710201 弯形施夹钳
Curved Holding Clamp



91710301 双边施夹钳
Double Edges Holding Clamp



91710401 枪式施夹钳
Gun Type Holding Clamp



91710502 骨膜剥离子 (小号)
Periosteal Detacher(Small)



91710601 骨膜剥离子
Periosteal Detacher




91710701 剪板钳
Plate Cutting Clamp



91710801 测量器 (选配)
Plate Cutting Clamp

足踝安装器械包






* * 足踝安装器械包

Veal bone fusion device installation instrument Set



***** 撑开钳（窄）
Trial Template,φ



***** 撑开钳（宽）
Trial Template,φ




***** 导针撑开钳
Trial Template,φ




***** 普通复位钳
Trial Template,φ




***** 带尖复位钳
Trial Template,φ



***** 骨刀
Trial Template,φ10mm




***** 骨刀
Trial Template,φ15mm



***** 骨刀
Trial Template,φ20mm



***** 骨撬
Trial Template,φ6mm



* * 足踝安装器械包

Veal bone fusion device installation instrument Set



***** 深部拉钩
Trial Template, $\phi 6 \times 20\text{mm}$



***** 深部拉钩
Trial Template, $\phi 9 \times 22\text{mm}$



***** 骨膜剥离器
Trial Template, $\phi 19\text{mm}$



***** 骨膜剥离器
Trial Template, $\phi 14\text{mm}$



***** 骨膜剥离器
Trial Template, $\phi 7\text{mm}$




***** 测深器
Trial Template



***** 钢板折弯器
Trial Template



***** KFM 加压导向器
Trial Template



***** 转换接头
Trial Template



***** KFM 梅花起子杆
Trial Template



***** KFM 2.7 两孔导向器
Trial Template,



***** KFM 3.5 两孔导向器
Trial Template,



***** 2.7 锁定导向器
Trial Template,



***** 3.5 锁定导向器
Trial Template,



***** KFM 带球克氏针 (螺纹)
Trial Template,



* * 足踝安装器械包

Veal bone fusion device installation instrument Set



***** KFM 带球克氏针（光杆）
Trial Template,



***** D字快装钻头（2.7螺钉用）
Trial Template, $\phi 2.8\text{mm}$



***** D字快装钻头
Trial Template, $\phi 2.1\text{mm}$



***** D字快装钻头
Trial Template, $\phi 3.5\text{mm}$



***** D字快装钻头（3.5螺钉用）
Trial Template, $\phi 2.8\text{mm}$




***** 2.7 万向导向器
Trial Template



***** 3.5万向导向器
Trial Template


其他系列安装器械包





*** * 骨卡环骑缝钉安装器械包 97119902**

Sawtooth-arm block hoop straddle nail Instrument Set



91351401 转换接头
Exchange joint



97110101 滑锤
Slide hammer



97110601 钻头
Drill bit, ϕ 2.8mm



97110201 导向器
Guider, 10mm



97110202 导向器
Guider, 12mm



97110203 导向器
Guider, 14mm



97110204 导向器
Guider, 16mm



97110205 导向器
Guider, 19mm



97110206 导向器
Guider, 21mm



*** * 骨卡环骑缝钉安装器械包 97119902**

Sawtooth-arm block hoop straddle nail Instrument Set



97110207 导向器
Guider, 26mm



97110301 撑开器
Distraction device, 10mm



97110302 撑开器
Distraction device, 12mm



97110303 撑开器
Distraction device, 14mm



97110304 撑开器
Distraction device, 16mm



97110305 撑开器
Distraction device, 19mm



97110306 撑开器
Distraction device, 21mm



97110307 撑开器
Distraction device, 26mm



97110401 打入器
Impactor



97110703 尺寸轮
Size wheel (Small)



97110901 限位套筒
Limit the sleeve, ϕ 2.1mm



97110902 限位套筒
Limit the sleeve, ϕ 2.8mm



97111301 钉盒
Nail box



97110501 钻头
Drill bit, ϕ 2.1mm



97110801 定位棒
Located rod, ϕ 2.0mm




97110802 定位棒
Located rod, ϕ 2.6mm



* * 恒天腰椎椎间融合器安装器械包 95219901

Veal bone fusion device installation instrument Set



95210101 试模
Trial Template, ϕ 7mm



95210103 试模
Trial Template, ϕ 9mm



95210105 试模
Trial Template, ϕ 11mm



95210102 试模
Trial Template, ϕ 8mm



95210104 试模
Trial Template, ϕ 10mm



95210106 试模
Trial Template, ϕ 12mm



95210202 侧方刮匙
Lateral curet, ϕ 8mm



95210203 侧方刮匙
Lateral curet, ϕ 9mm



95210204 侧方刮匙
Lateral curet, ϕ 10mm



95210201 侧方刮匙
Lateral curet, ϕ 7mm



95210205 侧方刮匙
Lateral curet, ϕ 11mm



95210206 侧方刮匙
Lateral curet, ϕ 12mm



95210301 带齿刮匙
Toothed curet, ϕ 6mm



* * 恒天腰椎椎间融合器安装器械包 95219901

Veal bone fusion device installation instrument Set



95210501 水滴形刮匙 (左)
Water drop curette (Left)



95210502 水滴形刮匙 (右)
Water drop curette (Right)



95210601 植骨推送器
Bone graft push device



95210701 植骨漏斗
Bone graft funnel



95210801 神经拉钩 (大)
Nerve retractor (Large)



95210802 神经拉钩 (小)
Nerve retractor (Small)



95210302 带齿刮匙
Toothed curette, $\phi 8\text{mm}$



95210401 方头拉钩
Coach retractor, $\phi 6\text{mm}$



95211101 持取器
A take device



95210901 剥离器 (大)
Detacher (Large)



95211001 剥离器 (小)
Detacher (Small)



95211201 六角扳手
Hexagon wrench



95211501 植入物推送器
Implant pusher



95211301 植入物推送器
Implant pusher



95211401 T型手柄
T handle





* * BES 2.5螺钉安装器械包 95139901

BES 2.5 screw instrument set



95130106 BES 阶梯钻头
Trial Template, $\phi 2.5 \times 10\text{mm}$



95130105 BES 阶梯钻头
Trial Template, $\phi 2.5 \times 12\text{mm}$



95130104 BES 阶梯钻头
Trial Template, $\phi 2.5 \times 14\text{mm}$



95130103 BES 阶梯钻头
Trial Template, $\phi 2.5 \times 16\text{mm}$



95130102 BES 阶梯钻头
Trial Template, $\phi 2.5 \times 18\text{mm}$



95130201 上钉器
Trial Template



95130401 丝锥
Lateral curet



91150124 快装手柄
Lateral curet





* * BES 3.5螺钉安装器械包 95109901

BES 3.5 screw instrument set



95100105 BES 阶梯钻头
Trial Template, $\phi 3.5 \times 20\text{mm}$



95100104 BES 阶梯钻头
Trial Template, $\phi 3.5 \times 25\text{mm}$



95100103 BES 阶梯钻头
Trial Template, $\phi 3.5 \times 30\text{mm}$



95100102 BES 阶梯钻头
Trial Template, $\phi 3.5 \times 35\text{mm}$



95100101 BES 阶梯钻头
Trial Template, $\phi 3.5 \times 40\text{mm}$



95100201 上钉器
Trial Template



95100301 丝锥
Lateral curet



* * BES 3.5空心钉螺钉安装器械包 95119901

BES 3.5 screw instrument set



95100505 BES 阶梯钻头
Trial Template, $\phi 3.5 \times 20\text{mm}$



95100504 BES 阶梯钻头
Trial Template, $\phi 3.5 \times 25\text{mm}$



95100503 BES 阶梯钻头
Trial Template, $\phi 3.5 \times 30\text{mm}$



95100502 BES 阶梯钻头
Trial Template, $\phi 3.5 \times 35\text{mm}$



95100501 BES 阶梯钻头
Trial Template, $\phi 3.5 \times 40\text{mm}$




95100202 上钉器
Trial Template



95100302 丝锥
Lateral curet





康力·器械产品目录

Kangli · Instrument Product Catalogue



20250530版