

Locking Connecting Plate Series


CE 0197 ISO 13485 FDA



KangLi

KANGLI ORTHOPAEDICS
INSTRUMENT

Locking connecting plate series



Character Of The System

- Suitable for comminuted fracture or fracture around the prosthesis.
- 2.7、 3.5、 5.0 Three series are optional for use, flexible use, wide indications.
- Locking connecting hole is fixed with main plate through locking connecting screw, and the side holes can be reinforced with locking screw.

Product Specification

JY 3.5 Locking Connecting Plate (Anodized) 23016



Product Name	Material	Specification	Product No.
JY 3.5 Locking Connecting Plate (Anodized)	Pure Titanium	4 x 28	23016428
		4 x 31	23016431
		4 x 34	23016434
Note: Use ϕ 3.5 Locking Connecting Screw (Anodized), ϕ 3.5 Locking Screw (Anodized).			

Locking connecting plate series

**JY 5.0 Locking Connecting Plate (Anodized)
23017**



Product Name	Material	Specification	Product No.
JY 5.0 Locking Connecting Plate (Anodized)	Pure Titanium	4 x 35	23017435
		4 x 38	23017438
		4 x 41	23017441

Note: Use ϕ 5.0 Locking Connecting Screw (Anodized), ϕ 5.0 Locking Screw (Anodized).

**Locking Connecting Screw (Anodized)
33004**



Product Name	Material	Specification	Product No.
3.5 Locking Connecting Screw (Anodized)	Titanium Alloy	ϕ 3.5 x 10	33004035
5.0 Locking Connecting Screw (Anodized)		ϕ 5.0 x 14	33004050

**2.7 Locking Connecting Plate
23017**



Product Name	Material	Specification	Product No.
2.7 Locking Connecting Plate	Pure Titanium	4 x 22	23017422
		4 x 25	23017425
		4 x 28	23017428

Note: Use ϕ 2.7 Locking Connecting Screw, ϕ 2.7 Locking Screw.

**3.5 Locking Connecting Plate
23018**



Product Name	材质 Material	规格 Specification	产品编号 Product No.
3.5 Locking Connecting Plate	Pure Titanium	4 x 28	23018428
		4 x 31	23018431
		4 x 34	23018434

Note: Use ϕ 3.5 Locking Connecting Screw, ϕ 3.5 Locking Screw .

**5.0 Locking Connecting Plate
23019**



Product Name	Material	Specification	Product No.
5.0 Locking Connecting Plate	Pure Titanium	4 x 35	23019435
		4 x 38	23019438
		4 x 41	23019441

Note: Use ϕ 5.0 Locking Connecting Screw , ϕ 5.0 Locking Screw .

**Locking Connecting Screw
33070**



Product Name	Material	Specification	Product No.
2.7 Locking Connecting Screw	Titanium Alloy	ϕ 2.7 x 8	33072027
3.5 Locking Connecting Screw		ϕ 3.5 x 10	33072035
5.0 Locking Connecting Screw		ϕ 5.0 x 14	33072050