

# MCF Multifunction craniomaxillofacial fixation system


CE 0197 ISO 13485 FDA




KangLi

KANGLI ORTHOPAEDICS  
INSTRUMENT

# MCF Multifunction craniomaxillofacial fixation system



**MCF straight Plate**  
T705A、T707A



Product Name	Material	Thickness(mm)	Bar Length(mm)	Product No.
MCF straight Plate	Pure Titanium	0.6	2 × 15	T705A0200
			4 × 20	T705A0400
		0.8	2 × 16	T707A0200
			4 × 21	T707A0400


**MCF straight Plate**  
T705A、T707A



Product Name	Material	Thickness(mm)	Bar Length(mm)	Product No.
MCF straight Plate	Pure Titanium	0.6	12 × 53	T705A1200
		0.8	12 × 54	T707A1200


## MCF Multifunction craniomaxillofacial fixation system

**MCF Titanium Oblique L Plate (left,right)**  
T705B、T707B




Product Name	Material	Thickness(mm)	Bar Length(mm)	Product No.
MCF Titanium Oblique L Plate (left,right)	Pure Titanium	0.6	4 × 17L	T705B0401
			4 × 17R	T705B0402
		0.8	4 × 17L	T707B0401
			4 × 17R	T707B0402

**MCF Titanium Oblique L Plate (left,right)**  
T705B、T707B




Product Name	Material	Thickness(mm)	Bar Length(mm)	Product No.
MCF Titanium Oblique L Plate (left,right)	Pure Titanium	0.6	5 × 21L	T705B0501
			5 × 21R	T705B0502
		0.8	5 × 22L	T707B0501
			5 × 22R	T707B0502

**MCF Titanium Oblique L Plate (left,right)**  
T705B、T707B



Product Name	Material	Thickness(mm)	Bar Length(mm)	Product No.
MCF Titanium Oblique L Plate (left,right)	Pure Titanium	0.6	7 × 26L	T705B0701
			7 × 26R	T705B0702
		0.8	7 × 27L	T707B0701
			7 × 27R	T707B0702


**MCF T Plate**  
T705C、T707C



Product Name	Material	Thickness(mm)	Bar Length(mm)	Product No.
MCF T Plate	Pure Titanium	0.6	7 × 25	T705C0700
		0.8	7 × 26	T707C0700


# MCF Multifunction craniomaxillofacial fixation system

**MCF Y Plate**  
T705D、T707D




Product Name	Material	Thickness(mm)	Bar Length(mm)	Product No.
MCF Y Plate	Pure Titanium	0.6	5 × 20	T705D0500
		0.8	5 × 21	T707D0500

**MCF Arc-Shaped Plate**  
T705E、T707E




Product Name	Material	Thickness(mm)	Bar Length(mm)	Product No.
MCF Arc-Shaped Plate	Pure Titanium	0.6	6 × 30	T705E0600
			8 × 39	T705E0800
		0.8	6 × 30	T707E0600
			8 × 39	T707E0800

**MCF Double Y Plate**  
T705F、T707F



Product Name	Material	Thickness(mm)	Bar Length(mm)	Product No.
MCF Double Y Plate	Pure Titanium	0.6	6 × 20	T705F0600
		0.8	6 × 21	T707F0600


**MCF X plate**  
T705G、T707G




Product Name	Material	Thickness(mm)	Bar Length(mm)	Product No.
MCF X plate	Pure Titanium	0.6	4 × 10	T705G0400
		0.8	4 × 11	T707G0400

# MCF Multifunction craniomaxillofacial fixation system


## MCF anatomic L Plate orthognathic T707H




## MCF oblique L plate orthognathic T707J



## MCF double curved iaw plate orthognathic T710C



## MCF straight sagittal split plate orthognathic T710B



## MCF curved sagittal split plate orthognathic T710A



## MCF Orbital floor anatomical reconstruction plate T705H




Product Name	Material	Thickness(mm)	Bar Length(mm)	Long bridge	Product No.
MCF anatomic L Plate orthognathic	Pure Titanium	0.8	6L	4	T707H6041
			6R	4	T707H6042
			6L	8	T707H6081
			6R	8	T707H6082
			6L	12	T707H6121
			6R	12	T707H6122
MCF oblique L plate orthognathic	Pure Titanium	0.8	4L	2	T707J4021
			4R	2	T707J4022
			4L	6	T707J4061
			4R	6	T707J4062
			4L	11	T707J4111
			4R	11	T707J4112
			6L	4	T707J6041
			6R	4	T707J6042
			6L	8	T707J6081
			6R	8	T707J6082
			6L	13	T707J6131
			6R	13	T707J6132
MCF double curved iaw plate orthognathic	Pure Titanium	0.8	5	4	T710C0504
				6	T710C0506
				8	T710C0508
				10	T710C0510

Product Name	Material	Thickness(mm)	Bar Length(mm)	Long bridge	Product No.
MCF straight sagittal split plate orthognathic	Pure Titanium	1.0	4	6	T710B0406
				8	T710B0408
				10	T710B0410
				12	T710B0412
MCF curved sagittal split plate orthognathic	Pure Titanium	1.0	6	6	T710A0660
				8	T710A0860
				10	T710A1060
				12	T710A1260

Product Name	Material	Thickness(mm)	Bar Length(mm)	Material	Product No.
MCF Orbital floor anatomical reconstruction plate	Pure Titanium	0.6		Big Left	42 x 37 T705H001L
				Big Right	42 x 37 T705H001R
				Small Left	42 x 36 T705H002L
				Small Right	42 x 36 T705H002R

# MCF Multifunction craniomaxillofacial fixation system

## MCF 0.6 Thin Maxillofacial Plate Series



### Matching screws

MCF  $\Phi 1.5$  Self drilling and tapping screw




MCF  $\Phi 1.5$  Seft tapping Screw



MCF  $\Phi 1.7$  Correction Screw



## MCF 0.8 Thin Maxillofacial Plate Series



### Matching screws

MCF  $\Phi 2.0$  Self drilling and tapping screw



MCF  $\Phi 2.0$  Seft tapping Screw



## MCF 0.8 Standard Orthognathic Maxillofacial Plate Series



## MCF 1.0 Standard Orthognathic Maxillofacial Plate Series



## MCF Orbital floor anatomical reconstruction plate



### Matching screws

MCF  $\Phi 2.0$  Self drilling and tapping screw



MCF  $\Phi 2.0$  Seft tapping Screw



## MCF Multifunction craniomaxillofacial fixation system

### MCF $\Phi$ 1.5 Self drilling and tapping screw 32325



Product Name	Material	Bar Length(mm)	Thickness(mm)	Product No.
MCF $\Phi$ 1.5 Self drilling and tapping screw	Titanium Alloy	$\Phi$ 1.5	4	32325004
			5	32325005
			6	32325006
			7	32325007

### MCF $\Phi$ 1.5 Seft tapping Screw 32321



Product Name	Material	Bar Length(mm)	Thickness(mm)	Product No.
MCF $\Phi$ 1,5 Seft tapping Screw	Titanium Alloy	$\Phi$ 1.5	4	32321004
			5	32321005
			6	32321006
			7	32321007

### MCF $\Phi$ 1.7 Correction Screw 32322



Product Name	Material	Bar Length(mm)	Thickness(mm)	Product No.
MCF $\Phi$ 1.7 Correction Screw	Titanium Alloy	$\Phi$ 1.7	5	32322005
			6	32322006
			7	32322007
			8	32322008

### MCF $\Phi$ 2.0 Self drilling and tapping screw 32326



Product Name	Material	Bar Length(mm)	Thickness(mm)	Product No.
MCF $\Phi$ 2.0 Self drilling and tapping screw	Titanium Alloy	$\Phi$ 2.0	5	32326005
			6	32326006
			7	32326007
			8	32326008



### MCF $\Phi$ 2.0 Seft tapping Screw 32323



Product Name	Material	Bar Length(mm)	Thickness(mm)	Product No.
MCF $\Phi$ 2.0 Seft tapping Screw	Titanium Alloy	$\Phi$ 2.0	5	32323005
			6	32323006
			7	32323007
			8	32323008

# MCF Multifunction craniomaxillofacial fixation system

## Matched Instrument Set



91130206  
MCF Filler Cross Screwdriver Rod,40mm



91130207  
MCF Filler Cross Screwdriver Rod,65mm



91270101  
MCF Awl,  $\phi 1.4 \times 68\text{mm}$



91270102  
MCF Awl,  $\phi 1.4 \times 71\text{mm}$



91270103  
MCF Awl,  $\phi 1.4 \times 73\text{mm}$



91150124  
MCF Quick Coupling Handle



91270201  
MCF Plate Bending Plier(Left)



91270202  
MCF Plate Bending Plier(Right)



91270301  
MCF Screw Box,  $\phi 1.5$



91270302  
MCF Screw Box,  $\phi 2.0$



91270605  
MCF Screw Lenght Colum,4



91270601  
MCF Screw Lenght Colum,5



91270602  
MCF Screw Lenght Colum, 6



91270603  
MCF Screw Lenght Colum,7



91270604  
MCF Screw Lenght Colum,8



91270401  
MCF Steel Plate Box, 0.6



91270402  
MCF Steel Plate Box,0.8



91270501  
MCF Profiling Plate Box,0.6



91270502  
MCF Profiling Plate Box, 0.8