

Prue Titanium Rib Plate (Sterile Packaging)

CE 0197 ISO 13485 FDA




KangLi

KANGLI ORTHOPAEDICS
INSTRUMENT

Prue Titanium Rib Plate (Sterile Packaging)



About Patent



About Registration Certificate



医疗器械注册登记表	
注册号: 国食药监械(注)字2013第3461819号	
生产企业名称	苏州市康力骨科器械有限公司
企业注册地址	张家港市沙洲东路
生产地址	张家港市沙洲东路、张家港市锦丰镇建设路
产品名称	金属解剖型接骨板
型号、规格	见附表
产品标准	YZB/T 039-2013 安全金属解剖型接骨板
产品性能指标及结构	“Ⅱ类”产品符合GB4218.1-2008标准规定 GB18311-1963(不锈钢或符合GB7124)中规定142、143 纯钛材料制成,纯钛产品在原包装中,非灭菌包装。
产品适用范围	该产品用于骨科手术材料作固定于髌骨、锁骨、肋骨、骨盆骨折内固定用。
产品禁忌症	见说明书
备注	

Rib Plate (Sterile Packaging)



Prue Titanium Rib Plate (Sterile Packaging)

Character Of The System

- Minimally invasive implant, each part is fixed with clamp, easy operation
- The pure titanium has a light weight and good biocompatibility. It doesn't contain ferric ion and won't affect CT & MRI examination. Besides, this flexible immobilization doesn't oppress intercostals nerve and the implant needn't be removed out after the operation.


Indications

- It is used to the fixation for floating chest wall and obvious dislocation of single or multiple ribs fractures.
- It can be customized to be used for the reconstruction of ribs after the rib tumor resection.
- It can also be used to the bone-setting after the thoracotomy.

Matters Deeding Attention

- Common sterilization on non-sterile packaging rib plate and clip applicator prior to use.
- Don't need cut the rib periosteum or remove the intercostals nerve.
- Hemostasis and drainage should be done completely during operative process.
- Adopting the closed thoracic drainage.

Rib Plate T132



Product Name	Packaging From	Material	Width	Length	Product No.
Rib Plate	Clean Packaging	Pure Titanium	11	49	T1321149
			15	34	T1321534
			15	49	T1321545
			15	59	T1321559
			19	34	T1321934
			19	49	T1321949
	Sterile Packaging		19	59	T1321959
			11	49	T1321149-S
			15	34	T1321534-S
			15	49	T1321545-S
			15	59	T1321559-S
			19	34	T1321934-S
			19	49	T1321949-S
			19	59	T1321959-S

Supporting Instruments

